FOCUS *CR,*
FOCUS *CR*e*,*
FOCUS *CR TTL*

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 **Ki** Mobility

I. INTRODUCTION

Thank you for purchasing a Focus CR wheelchair!

Please do not use this wheelchair without first reading this entire manual. BEFORE riding, you should be trained in the safe use of this chair by an Assistive Technology Practitioner (ATP) or clinical professional.

If you have any questions or concerns about any aspect of this wheelchair, this manual, or the service provided by us or your retail supplier, please do not hesitate to contact us by telephone at:

715-254-0991

In writing at:

Ki Mobility
5201 Woodward Drive
Stevens Point, WI 54481
U.S.A

Via email at:

sales@kimobility.com

Or via our Authorized EU Representative:

James Leckey Design
19C Ballinderry Road
Lisburn
BT28 2SA
Phone: 0800 318265 (UK) or 1800 626020 (ROI)
www.leckey.com

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I.	INTRODUCTION	
II.	TABLE OF CONTENTS	
III.	NOTICE - READ BEFORE USE	
	A. Your Safety and Stability	4
IV.	WARNINGS	
	A. Signal Words.....	4
	B. General Warnings.....	5
	C. Positioning Belts	6
	D. Riding Your Wheelchair	7
	E. Power Drives.....	8
	F. Ascending Stairs	8
	G. Descending Stairs	8
	H. Transfers.....	9
	I. Your Wheelchair and the Environment.....	9
	J. Modifying your Wheelchair	10
	K. Wheelchair Stability	10
V.	SET UP & USE OF YOUR WHEELCHAIR	
	A. Your Focus CR It's Parts.....	13
	B. Transit Use.....	14
	C. Height Adjustable T-Arms	17
	D. Dual Post Height Adjustable Armrest.....	19
	E. Angle Adjustable Locking Flip Up Extendable Armrest.....	20
	F. Armrest Warnings	20
	G. Swing Away Hangers	21
	H. Swing Away Hangers with 4-Way Latch.....	22
	I. Extension Tubes	23
	J. Elevated Leg Rest	24
	K. Pro Elevated Leg Rest.....	25
	L. Pediatric Pro Elevating Leg Rest.....	26
	M. Contracture Footrest.....	27
	N. Multi-Angle Footrest	28
	O. Backrest	29
	P. Reclining Backrest	31
	Q. Tilt Mechanism	34
	R. Seat Frame Set Up.....	36
	S. Casters	37
	T. Axle Plate.....	38
	U. Rear Wheels.....	40
	V. Wheel Locks	41

II. TABLE OF CONTENTS

W.	Drum Brake.....	43
X.	Frame Width	44
Y.	Anti-Tips	45
Z.	Seat Pan	46
AA.	Bearings.....	47
BB.	Vent and Battery Tray.....	48
CC.	Oxygen Tank Holder	48
DD.	Power Tilt Interface Control Box	49
EE.	Canopy	55
VI.	MAINTENANCE	
A.	Inspecting Your Wheelchair	56
B.	Routine Maintenance	57
C.	Cleaning.....	58
D.	Storage.....	58
VII.	TROUBLESHOOTING.....	59
VIII.	WARRANTY	61

III. NOTICE - READ BEFORE USE

A. Your Safety and Stability

Ki Mobility manufactures many different wheelchairs that might meet your needs. You should consult an Assistive Technology Professional when selecting which model would best meet your particular requirements and how the wheelchair should be set up and adjusted. Final selection of the type of wheelchair, options and adjustments rests solely with you and your medical professional. The options you choose and the set-up and adjustment of the wheelchair have a direct impact on its stability. Factors to consider that affect your safety and stability are:

- a. Your personal abilities and capabilities including strength, balance and coordination.
- b. The types of hazards and obstacles you might encounter during your day.
- c. The specific dimensions, options and set up. In particular, the seat height, seat depth, seat angle, back angle, size and position of the rear wheels and size and position of the front casters. Any change to any of these items will change the stability of your wheelchair. You should only make changes after consulting with a qualified professional.

IV. WARNINGS

A. Signal Words

Within this manual you will find what are referred to as “Signal” words. These words are used to identify and convey the severity of varying hazards. Before using this chair you, and each person who may assist you, should read this entire manual. Please note the Signal word and consider any warnings, cautions or dangers. Make sure to follow all instructions and use your chair safely. The Signal word refers to a hazard or unsafe practice that may cause severe injury or death to you or to other persons. The “Warnings” are in three main categories, as follows:

DANGER – Danger indicates an imminently hazardous situation which, if not avoided, will result in serious injury or death.

WARNING – Warning indicates a potentially hazardous situation which, if not avoided, could result in serious injury or death.

CAUTION – Caution indicates a potentially hazardous situation which, if not avoided, could result in injury or damage to your wheelchair.

These signal words will be placed throughout the manual, where appropriate to highlight the hazardous situation. Refer to the following list for hazardous situations that will apply to the general use of this wheelchair.

IV. WARNINGS

B. General Warnings

WARNING

WARNING: Do not exceed the max weight capacity of 300 lbs (136 kg) on the standard frame or 400 lbs (181 kg) on a heavy duty frame. Reference serial number label for chair max capacity. This is the combined weight of user and all items carried. Exceeding the weight limit can cause damage to your chair or increase the likelihood of a fall or tip back resulting in severe injury or death to the user or others.

DANGER

DANGER: Do not use this chair for weight training. The movement of the additional weight will alter the center of gravity of the wheelchair increasing the likelihood of a tip-over which can cause damage to your chair or cause severe injury or death to the user or others.

WARNING

WARNING: If your wheelchair is equipped with inflatable tires, make sure the tires have been inflated to the correct tire pressure as indicated on the side wall of the tire. Your wheelchair provider can determine if you have inflatable tires. In many cases the wheel locks will not function properly with under-inflated tires. Using your wheelchair without properly inflated tires can have an affect on the stability of the wheelchair causing it to tip over resulting in death or injury to the user.

DANGER

DANGER: Do not attempt to push your wheelchair up or down ramps or traverse across a slope of greater than 9 degrees. This is dangerous and increases the likelihood of a fall or tip back resulting in severe injury or death to the user or others.

WARNING

WARNING: Do not attempt to push your wheelchair up an incline that is slick or coated with ice, oil or water. This can cause an unstable situation resulting in death or injury to the user or others.

WARNING

WARNING: Do not lean over the side or back of the wheelchair to extend your reach. This may cause you to fall out of the wheelchair or the wheelchair to tip over resulting in injury or death.

DANGER

DANGER: Do not attempt to lift the wheelchair by holding on to removable parts such as the arms or footrests. Only lift the wheelchair by holding on to the frame. This may cause a fall or loss of control and result in serious injury or death.

CAUTION

CAUTION: Do not overtighten the bolts and hardware that attaches various components together on the frame. This could cause serious damage and affect the safety and durability of the wheelchair.

IV. WARNINGS

C. Positioning Belts

Positioning belts are designed to assist with proper positioning within the wheelchair. They are not designed as seat belts. Use positioning belts ONLY to help support the user's posture. Misuse of positioning belts may cause severe injury to or death of the user.

- Ensure the user does not slide underneath the positioning belt in the wheelchair seat. If this occurs, the user's breathing may be hampered causing death or serious injury.
- The positioning belt should have a snug fit; tight enough to hold their position, but not so tight as to restrict breathing. You should be able to slide your hand between the positioning belt and the user.
- **NEVER** Use Positioning Belts:
 - a. As a restraint. A restraint requires a doctor's order.
 - b. On a user who is unconscious.
 - c. As an occupant restraint in a motor vehicle. A positioning belt is not designed to replace a seat belt that is attached to the frame of a vehicle, which would be required of an effective seat belt. During a sudden stop, with the force of the stop, the user would be thrown forward. Wheelchair seat belts will not prevent this, and further injury may result from the belts or straps. See Transit Use (V. Set Up & Use of Your Wheelchair - Section B).

DANGER: Failure to comply with the instruction above could result in serious injury or death.

IV. WARNINGS

D. Riding Your Wheelchair

Your chair is designed for use on solid, flat surfaces such as concrete, asphalt and flooring. Use caution if you push your wheelchair on a wet or slick surface.

WARNING

WARNING: Do not push your chair in sand, loose soil or over rough terrain. This may cause a loss of stability and result in a fall or loss of control and cause serious injury or death.

DANGER

DANGER: In most states, wheelchairs are not legal for use on public roads. If you find you must push on a public road, be alert to the danger of motor vehicles. Use of a wheelchair on a public road can cause serious injury or death.

WARNING

WARNING: Obstacles and road hazards (such as potholes and broken pavement) can damage your chair and may cause a fall, tip-over or loss of control. Failure to comply with this instruction could result in serious injury or death.

DANGER

DANGER: Do not ride your wheelchair on an escalator. Use of a wheelchair on an escalator can cause serious injury or death.

To minimize these risks:

- 1) Keep a lookout for danger-scan the area well ahead of your chair as you ride.
- 2) Make sure the floor areas where you live and work are level and free of obstacles.
- 3) Remove or cover threshold strips between rooms.
- 4) Install a ramp at entry or exit doors. Make sure there is not a drop off at the bottom of the ramp.
- 5) To Help Correct Your Center of Balance:
 - a. Lean your upper body **FORWARD** slightly as you go **UP** over an obstacle.
 - b. Press your upper body **BACKWARD** as you go **DOWN** from a higher to a lower level.
- 6) If your chair has anti-tip tubes, lock them in place before you go **UP** over an obstacle.
- 7) Keep both of your hands on the handrims as you go over an obstacle.
- 8) Never push or pull on an object (such as furniture or a doorjamb) to propel your chair.
- 9) Do not operate your wheelchair on roads, streets or highways.
- 10) Do not attempt to push over obstacles without assistance.

IV. WARNINGS

E. Power Drives

Ki Mobility does not recommend the installation of power drive systems on any Focus CR wheelchair.

Focus CR wheelchairs have not been designed or tested as power wheelchairs. If you add a power drive system to a Focus CR wheelchair, be sure the manufacturer of the power drive system has validated and approved the combination of the power drive system and Focus CR wheelchair as safe and effective.

WARNING

WARNING: Use of a power drive system that has not been properly validated could result in serious injury or death.

F. Ascending Stairs

- Have at least two people, who have sufficient strength and skill to handle the weight of the user and wheelchair, assist when trying to go up a set of stairs in this wheelchair.
- Move the wheelchair and user backwards up the stairs.
- Position one person behind the user, one person in front. The person in front must hold onto a non-removable part of the wheelchair.
- The rear attendant tilts the chair back and they both lift together. Take one step at a time.
- This may require the anti-tips be flipped up or removed. Make sure the anti-tips are reattached or flipped back down before using the wheelchair.

DANGER

DANGER: Failure to comply with the instructions above could result in serious injury or death.

G. Descending Stairs

- When descending a set of stairs the user should be facing forward.
- A person behind the user, who has sufficient strength and skill to handle the weight of the user and the wheelchair, should tilt the chair backward and let the chair down the stairs one step at a time on the rear wheels.
- This may require the anti-tips be flipped up or removed. Make sure the anti-tips are reattached or flipped back down before using the wheelchair.

DANGER

DANGER: Failure to comply with the instructions above could result in serious injury or death.

IV. WARNINGS

H. Transfers

A transfer requires good balance and stability. You should receive training from your therapist before attempting to do a transfer on your own.

- Before transferring out of your wheelchair every caution should be taken to reduce the gap between the two surfaces.
- Engage the wheel locks to lock the rear wheels.
- Rotate the casters forward to increase the wheelbase of the wheelchair.
- Remove or swing away the footrests.
- Have someone assist you unless you are well experienced in transfers.

It is dangerous to transfer on your own. It requires good balance and agility. Be aware there is a point during every transfer when the wheelchair seat is not below you.

WARNING

WARNING: Failure to perform a transfer properly can result in a fall and can cause severe injury or death.

I. Your Wheelchair and the Environment

- Your wheelchair is made of many different materials including metal and fabric. Exposure to water or excessive moisture may cause the metal in the wheelchair to rust or corrode and the fabric to tear. Dry your chair as soon as possible if exposed to water.
- **DO NOT USE YOUR WHEELCHAIR IN A SHOWER, POOL OR BODY OF WATER.** This will cause your wheelchair to rust or corrode and eventually fail.
- Do not operate your wheelchair in sand. Sand can get into the wheel bearings and moving parts. This will cause damage and eventually will cause the wheelchair to fail.
- Make sure any ramp, slope or curb cut you may attempt to ride on is compliant with ADA guidelines. Riding across, up or down any slope that is too great may cause a loss of stability.

ADA Guidelines and more information about accessible design are available at:
<http://www.ada.gov/>

WARNING

WARNING: Failure to comply with the instructions above may cause a fall or loss of control, which may result in serious injury or death.

IV. WARNINGS

J. Modifying your Wheelchair

Your wheelchair was engineered and manufactured under strict design controls. An integral part of this process is ensuring the various components work together correctly; they have been tested to various standards to ensure quality and are approved to work together.

NO ONE SHOULD MODIFY THIS WHEELCHAIR EXCEPT BY ADJUSTING IT ACCORDING TO THIS MANUAL OR BY ADDING KI MOBILITY APPROVED OPTIONS. THERE ARE NO APPROVED OPTIONS THAT INVOLVE DRILLING OR CUTTING THE FRAME BY ANYONE OTHER THAN A TRAINED KI MOBILITY ASSOCIATE. Contact Ki Mobility or an authorized Ki Mobility supplier before adding any accessories or components not provided by Ki Mobility.

DANGER: Failure to comply to these instructions may cause the wheelchair to fail and result in serious injury or death.

K. Wheelchair Stability

To ensure proper stability of your wheelchair, you must make sure the center of gravity and the wheelchairs base of support is correct for your balance and abilities. Many factors can affect these two elements:

- Seat height
- Seat depth
- Back angle
- Seat angle
- Size and position of rear wheels
- Size and position of front casters
- Any seating system components
- Tilt position

Generally, the most important factor is the position of the rear wheels for rearward stability. There are other actions than can have an adverse effect on your stability. You should consult with your wheelchair provider and clinicians familiar with your needs and capabilities in determining how this affects your use.

WARNING: This chair is equipped with a system to vary the tilt angle of the seat frame. Make sure the chair is stable throughout the range of tilt angle change. Failure to comply to the instruction above could result in serious injury or death.

WARNING: Moving the rear wheels forward increases the likelihood of the wheelchair tipping backwards. Make small adjustments and proceed slowly until you learn the new balance point of your wheelchair. Failure to comply with the instruction above could result in serious injury or death.

IV. WARNINGS

K. Wheelchair Stability

WARNING

WARNING: The farther rearward you place the front casters the greater the likelihood of the wheelchair tipping forwards. If possible, have your casters mounted forward and whenever doing a static activity which involves shifting your weight, rotate the casters forward to increase your wheel base. Failure to comply with this instruction above could result in serious injury or death.

WARNING

WARNING: Always have a qualified technician set up your wheelchair with the accessories you plan to use daily.

WARNING

WARNING: Changes to how you sit or changes in your weight require your chair to be readjusted by a qualified technician. Always use anti-tips while you acclimate to any changes in your chair set up. Failure to comply with the instruction above could result in serious injury or death.

WARNING

WARNING: Changes to your Center of Gravity during your daily activities may occur many times a day, changing and affecting the stability of your wheelchair. You should be aware of these activities and take precautions to minimize the risk of a fall. Failure to comply with the instruction above could result in serious injury or death.

WARNING

WARNING: Dressing in your wheelchair produces movements and momentary positions that can reduce stability. Ensure that your anti-tips are in place and rotate your casters forward. Failure to comply with the instruction above could result in serious injury or death.

WARNING

WARNING: Be very careful when reaching for objects if this movement requires you to shift in your seat. This changes your center of gravity. Ensure that your anti-tips are in place. Failure to comply with the instruction above could result in serious injury or death.

WARNING

WARNING: Pushing up an incline shifts your center of gravity rearward and can reduce stability. Ensure your anti-tips are in place. Failure to comply with the instruction above could result in serious injury or death.

IV. WARNINGS

K. Wheelchair Stability

WARNING

WARNING: If attempting a wheelie to get over a curb or obstacle, ensure your anti-tips are in place and lean forward. Do not attempt a wheelie unless you have been trained and always have an attendant behind you to provide assistance if needed. Failure to comply with the instruction above could result in serious injury or death.

WARNING

WARNING: Placing items on the back or front of your wheelchair, such as a backpack or briefcase, alters the balance and center of gravity of the wheelchair. Since the weight of these items can vary greatly at each use do not assume you are accustomed to the balance point. Failure to comply with the instruction above could result in serious injury or death.

BE AWARE THAT CARRYING HEAVY OBJECTS ON YOUR WHEELCHAIR CAN HAVE AN ADVERSE EFFECT ON THE BALANCE WHICH MAY CAUSE A TIP-OVER WHICH MAY RESULT IN SERIOUS INJURY OR DEATH TO THE USER.

WARNING

WARNING: Ensure your anti-tips are in place. You should discuss how you plan to use your wheelchair or any changes you are planning with your clinician. Failure to comply with this instruction may create a potential hazardous situation which, if not avoided, could result in serious injury or death.


V. SET UP & USE OF YOUR WHEELCHAIR

A. Your Focus CR It's Parts

1. Inspect and maintain this chair. See *MAINTENANCE* on page 56.
2. If you detect a problem, make sure to service or repair the chair before use.
3. Have a complete inspection, safety check and service of your chair performed by an authorized supplier annually.

⚠ WARNING ⚠

WARNING: Failure to read or comply with these instructions may result in damage to your wheelchair, a fall, or loss of control causing severe injury to the user or others.

- | | |
|------------------------|-----------------------|
| 1. Base Tube | 12. Anti-Tip |
| 2. Rotary Frame | 13. Wheel Lock |
| 3. Seat Tube | 14. Caster Fork |
| 4. Caster Arm | 15. Caster Wheel |
| 5. Back Cane | 16. Swing Away Hanger |
| 6. Push Handle | 17. Footrest |
| 7. Tilt Lock Trigger | 18. Swing Away Lever |
| 8. Tilt Cable | 19. Axle Plate |
| 9. Trigger Lock Pin | 20. Seat Pan |
| 10. Rear Wheel | 21. Dual Post Armrest |
| 11. Quick Release Axle | |

V. SET UP & USE OF YOUR WHEELCHAIR

B. Transit Use

It is always safest to transfer out of your wheelchair onto a seat in a motor vehicle with appropriate seat and shoulder belts. **Never use this wheelchair as a seat in a motor vehicle unless it has been equipped with the Transit Option.**

The Focus CR wheelchair equipped with the Transit Option has been tested to and passed the RESNA WC-4:2012, Section 19: Wheelchairs used as seats in motor vehicles and ISO 7176-19:2008 Wheelchairs -- Part 19: Wheeled mobility devices for use as seats in motor vehicles. RESNA and ISO standards are designed to test the structural integrity of the wheelchair as a seat for use in a motor vehicle. These standards are also designed to create compatibility with Wheelchair Tie-down and Occupant Restraint Systems (WTORS).

Not all configurations of the Focus CR wheelchairs are compatible with the Transit Option. Ki Mobility manages the configuration and does not offer the Focus CR wheelchair except in compatible configurations. If you make changes to your Focus CR wheelchair after you receive it, you should contact your wheelchair provider or Ki Mobility to make sure it is appropriate to continue to use your wheelchair as a seat in a motor vehicle.

Aftermarket seating may have replaced the original equipment seat and back support designed and tested as part of the Transit Option. Your wheelchair provider should tell you if the seating they provided is original equipment or replacement aftermarket seating. A complete system of wheelchair frame, seating, Wheelchair Tie-down and Occupant Restraint Systems and a properly equipped motor vehicle, that have all complied with the standards mentioned in this section, should be in place before using a Focus CR wheelchair equipped with the Transit Option as a seat in a motor vehicle.


When using your wheelchair as a seat in a motor vehicle you should always observe the following instructions:

- The rider must be in a forward-facing position.
- The rider and all items carried must not weigh more than 300 lbs (136 kg) (400 lbs (181 kg) for Heavy Duty option).
- Backpacks and pouches should be removed and secured separately in the motor vehicle. In the event of an accident these items can become dangerous projectiles, which may injure or kill you or other occupants of the motor vehicle.

V. SET UP & USE OF YOUR WHEELCHAIR

B. Transit Use

- The rider must use a Wheelchair Tie-down and Occupant Restraint System that complies with RESNA WC-4:2012, Section 18: Wheelchair tie-down and occupant restraint systems for use in motor vehicles or ISO 10542-1:2012 Technical systems and aids for disabled or handicapped persons -- Wheelchair tie-down and occupant-restraint systems -- Part 1: Requirements and test methods for all systems.

- Attach wheelchair tie-downs to the four securement points (two front, two rear) on the Focus CR wheelchair with the Transit Option (Fig. 1) in accordance with the wheelchair tie-down manufacturer's instructions and RESNA WC-4:2012, Section 18 or ISO 10542-1:2012 - Part 1.
- Attach occupant restraints in accordance with the occupant restraint manufacturer's instructions and RESNA WC-4:2012, Section 18 or ISO 10542-1:2012, Part 1.
- Use of lap belts, chest straps, shoulder harnesses, any other positioning strap system or positioning accessory should not be used, or relied on as an occupant restraint, unless it is marked as such by the manufacturer in accordance with RESNA WC-4:2012, Section 18 or ISO 10542-1:2012, Part 1.
- Use of headrests, lateral supports or other positioning accessories should not be used, or relied on as an occupant restraint, unless it is marked as such by the manufacturer in accordance with RESNA WC-4:2012, Section 18 or ISO 10542-1:2012, Part 1 or RESNA WC-4:2012, Section 20: Wheelchair seating systems for use in motor vehicles or ISO 16840-4:2009 Wheelchair seating -- Part 4: Seating systems for use in motor vehicles.
- Any aftermarket seating should be tested to comply with RESNA WC-4:2012, Section 20 or ISO 16840-4:2009 - Part 4.
- Attach the seating to the wheelchair frame in accordance with the seating manufacturer's instructions and RESNA WC-4:2012, Section 20 or ISO 16840-4:2009 - Part 4.
- Aftermarket accessories such as trays, oxygen tank holders, oxygen tanks, IV poles, back packs, pouches and other items not manufactured by Ki Mobility should be removed and secured separately in the motor vehicle. In the event of an accident, these items can become dangerous projectiles which may injure or kill you or other occupants of the motor vehicle.
- If the wheelchair has been involved in an accident, you should not continue to use it, as it may have suffered fatigue that may not be visible.

DANGER

DANGER: Failure to comply with any of these instructions could result in severe injury or death!

V. SET UP & USE OF YOUR WHEELCHAIR

B. Transit Use

- When using this wheelchair as a seat in a motor vehicle, you must remove any items attached to the vent tray, battery tray or oxygen tank holder and properly secure them separately.

DANGER: Failure to remove items attached to the vent tray, battery tray or oxygen tank holder and properly secure them separately in the motor vehicle can result in these items becoming dangerous projectiles in the event of an accident. Oxygen tanks contain a highly pressurized gas that vigorously accelerates combustion. Failure to heed this warning can result in severe injury or death.

- If the chair is equipped with an oxygen tank holder never use the wheelchair as a seat in a motor vehicle.

DANGER: Using a Focus wheelchair equipped with an oxygen tank holder as a seat in a motor vehicle can result in severe injury or death.

NOTE: To obtain copies of RESNA or ISO standards please contact the standards organizations below:

RESNA

1700 N Moore St Ste 1540
Arlington, VA 22209
Phone: 703-524-6686
Fax: 703-524-6630
Email: technicalstandards@resna.org

ANSI/RESNA Standards:

RESNA WC-4:2012, Section 18:
Wheelchair tie-down and occupant restraint systems for use in motor vehicles.

RESNA WC-4:2012, Section 19:
Wheelchairs used as seats in motor vehicles.

RESNA WC-4:2012, Section 20:
Wheelchair seating systems for use in motor vehicles.

International Organization for Standardization (ISO)

BIBC II
Chemin de Blandonnet 8
CP 401
1214 Vernier, Geneva
Switzerland
Phone: +41 22 749 01 11
Fax: +41 22 733 34 30
Email: central@iso.org

ISO Standards:

ISO 10542-1:2012 Technical systems and aids for disabled or handicapped persons -- Wheelchair tie-down and occupant-restraint systems -- Part 1:
Requirements and test methods for all systems.

ISO 16840-4:2009 Wheelchair seating – Part 4:
Seating systems for use in motor vehicles.

ISO 7176-19:2008 Wheelchairs – Part 19:
Wheeled mobility devices for use as seats in motor vehicles.

V. SET UP & USE OF YOUR WHEELCHAIR

C. Height Adjustable T-Arms

1. Installation


- a. Slide the outer armpost into the receiver mounted to the wheelchair frame.
- b. The armrest will automatically lock into place. Check to make sure the locking lever is as shown (Fig. 2:B).

2. Height Adjustment

- a. Rotate release lever (Fig. 2:A).
- b. Slide armrest pad up or down to desired height.
- c. Return lever to locked position against arm post.
- d. Push arm pad until upper arm locks firmly into place. Check to make sure the locking lever is as shown (Fig. 2:A).

3. Removing Armrest

- a. Squeeze release lever (Fig. 2:B) and remove the armrest.


V. SET UP & USE OF YOUR WHEELCHAIR

C. Height Adjustable T-Arms

4. Adjusting Position

- a. Remove both screws (Fig. 3:A) from each armrest receiver (Fig. 3:B) retaining the spacers (Fig. 3:C) if installed. If not installed at factory, they are included in a separate bag in the packaging.
- b. Move receiver to desired location on seat frame and reinstall screws with washers and nuts and tighten.
- c. If screws for new armrest receiver location coincide with seat frame screws (Fig. 3:D and E), remove seat frame screws and replace with armrest receiver screws. Retain (do not discard) these screws for potential future use, or optionally, they may be retained by reinstalling them in an alternate location on the seat frame.
- d. If one of the armrest receiver screws lies on the center growth section of the seat frame, be sure to use the spacers (as described in step 1) around growth tube in location of mounting hole.

⚠ DANGER ⚠

DANGER: Failure to comply with the instructions above may result in the armrest accidentally disconnecting from the wheelchair and result in a fall or loss of control and may cause serious injury or death.

⚠ DANGER ⚠

DANGER: Never attempt to lift the chair by the armrests; they may break or disconnect resulting in a fall or loss of control and may cause serious injury or death.

V. SET UP & USE OF YOUR WHEELCHAIR

D. Dual Post Height Adjustable Armrest

1. Installation

- a. Ensure both levers (Fig. 4:A and B) are flipped upward.
- b. Set both posts of armrest into receivers until they are seated against the bottoms of the receivers.
- c. Flip levers back downward. Ensure locking pins are seated by lifting on armrest.

2. Height Adjustment

- a. Flip lever (Fig. 4:C) forward at the top of the front post of the lower frame of the armrest.
- b. Adjust height of armrest by pushing down or pulling upwards on arm pad near the center between the two posts of the armrest.
- c. When proper height is achieved, flip lever rearward again and move armrest up or down to allow pin to seat in nearest hole.


3. Flipping Back Armrest

- a. Flip lever (Fig. 4:A) on front armrest receiver.
- b. Lift upwards on arm pad and cross tube, allowing the armrest to flip back.

NOTE: Do not push outward or backward on armrest when it is flipped back or damage to rear receiver may occur.

4. Removal

- a. Flip the levers (Fig. 4:A and B) on both receivers upward.
- b. Lift armrest out of receivers by pulling up on arm pad and cross tube.

V. SET UP & USE OF YOUR WHEELCHAIR


E. Angle Adjustable Locking Flip Up Extendable Armrest

Adjustments

1. Set the angle of the armrest based on your preferences. There are five holes that can be used to set the angle (Fig. 5:A). Tighten the bolt once angle is set.
2. Set the length of the armrest based on your preferences. To adjust the length, remove the bolts and spacer (Fig. 5:B) on the tube and the screw (Fig. 5:C) closest to the back of the chair. Slide the armrest to desired length available by the predrilled holes and reinstall the screw and bolts.
3. Set the height of the armrest on the chair based on your preferences. There are four holes on the armrest that allow for two different height settings for each set of holes on the back tube. Use the holes that provide the correct height setting for the user. The two bolts (Fig. 5:D) pass through the spacer (Fig. 5:E), sleeves (Fig. 5:F), back posts and into the armrest.

Use

1. Press the release lever (Fig. 5:G) up to release the armrest and swing it upwards. To return the armrest to the operating position, push the armrest down until the lever clicks and locks.

F. Armrest Warnings

- All Ki Mobility armrests are designed to detach from the chair and will not bear even the weight of this chair.
- NEVER lift this chair by its armrests. The armrests will release and the user may fall.
- Lift this chair only by non-detachable parts of the main frame.

⚠ DANGER ⚠

DANGER: Always make sure the arms are locked in place before using them for repositioning yourself, failure to do so may result in a fall or loss of control and may cause serious injury or death.

⚠ WARNING ⚠

WARNING: Failure to heed these instructions may result in a fall, tip-over or loss of control causing severe injury to the user or others.

V. SET UP & USE OF YOUR WHEELCHAIR

G. Swing Away Hangers

1. Installation



- a. Place swing away pivot saddle into the receiver on front frame tube with the footrest facing either inward or outward from the frame (Fig. 6).
- b. Rotate the footrest so that it aligns with the frame until it locks into place in the latch block (Fig. 7).

2. Swinging the Footrest Away

- a. Push release latch toward the frame.
- b. Rotate footrest outward or inward as desired.

3. Removal

- a. To remove footrest, push release latch toward the frame.
- b. Lift the footrest straight upward to remove. You may also swing the footrest inward or outward before lifting it off.

⚠ WARNING ⚠

WARNING: Always make sure the hangers are locked in place before using them or riding the wheelchair. Failure to do so may result in a fall or loss of control and may cause serious injury or death.

⚠ DANGER ⚠

DANGER: Never attempt to lift the chair by the hangers; they may break or disconnect resulting in a fall or loss of control and may cause serious injury or death.

V. SET UP & USE OF YOUR WHEELCHAIR

H. Swing Away Hangers with 4-Way Latch

1. Installation

- Place swing away pivot saddle into the receiver on front frame tube (Fig. 8:A).
- Rotate the footrest so that it aligns with the frame and locks into place in the latch block (Fig. 8:B).

2. Swinging the Footrest Away

- Push or pull on release latch.
- Rotate footrest outward or inward as desired.

3. Removal

- To remove footrest, push or pull release latch.
- Lift the footrest straight upward to remove. You may also swing the footrest inward or outward before lifting it off.

4. Rotating (Fig. 9)




The 4-way latch has eight possible configurations, four with the curve of the lever facing outward and four with the curve of the lever facing inward. See Fig. 10 for possible configurations.

- To rotate the 4-way latch, remove the screw using a 3mm Allen wrench while the hanger is still on the chair (spring must be engaged to remove and reinsert screw and keeping the hanger on the chair keeps the spring engaged). Ensure nut does not fall out.
- Rotate the 4-way latch to desired orientation and reinsert screw with 3mm Allen wrench. Ensure that nut stays in position while tightening the screw. Do not over-tighten screw or mechanism will bind.

NOTE: To reverse the 4-way latch, the same screw is removed, but the hanger has to be removed from the latch block. Once removed,

slide the latch off, flip over and reinstall. Ensure spring is engaged, by pushing and holding the latch button in, and nut stays in position while reinstalling the screw.

NOTE: In-line position is not achievable with the Pro ELR Footrest option.

V. SET UP & USE OF YOUR WHEELCHAIR

I. Extension Tubes



1. Adjustment

- a. Remove mounting fasteners from each side of the hanger tube.
- b. Slide footrest extension tube to the desired height.
- c. Line up holes and reassemble fasteners in desired hole through hanger and extension tube.
- d. Follow same procedure on opposite side (Fig. 11).

At the lowest point, footrests should be **AT LEAST 2 ½ INCHES (6.35 cm)** off the ground. If set too low, they may "catch" on obstacles you would expect to find in normal use. This may cause the chair to stop suddenly and tip forward (Fig. 12).

To Avoid a Trip or Fall When You Transfer:

- Make sure your feet do not get caught in the space between the footrests.
- Avoid putting weight on the footrests as the chair may tip forward.

V. SET UP & USE OF YOUR WHEELCHAIR

J. Elevated Leg Rest

1. Installation

- Place swing away pivot saddle into the receiver on front frame tube with the footrest facing either inward or outward from the frame. Similar to swing away footrest pictured in Fig. 6 and 7.
- Rotate the footrest so it aligns with the frame until it locks into place in the latch block.

2. Removal


- To remove footrest, push release latch toward the frame.
- Rotate footrest outward and lift.

3. Extension Tube Adjustment

- Remove mounting bolt that is threaded into the leg rest tube in the center of the pad bracket (Fig. 13:C). In some circumstances this bolt may be inserted in the tube below the pivot bracket.
- Slide footrest extension tube to the desired height.
- Line up holes and replace the bolt through the leg rest tube and extension tube.
- Follow same procedure on opposite side.

4. Angle Adjustment

- To raise, lift up on the rear of the extension tube (Fig. 13:A). The rod will slide through the ratchet in this direction. Stop at desired position.
- To lower, firmly hold the leg from behind the extension tube (Fig. 13:A). Pull forward on the lever (Fig. 13:B) and while holding the lever, lift the leg rest. Releasing the lever will cause the leg rest to lock into position.

To Avoid a Trip or Fall When You Transfer:

- Make sure your feet do not get caught in the space between the footrests.
- Avoid putting weight on the footrests as the chair may tip forward.

⚠ WARNING ⚠

WARNING: At the lowest point, footrests should be AT LEAST 2 ½ INCHES (6.35 cm) off the ground. If set too low, they may "catch" on obstacles you would expect to find in normal use. This may cause the chair to stop suddenly and tip forward and could result in a fall or loss of control and may cause serious injury or death.

⚠ DANGER ⚠

DANGER: Never attempt to lift the chair by the hangers; they may break or disconnect resulting in a fall or loss of control and may cause serious injury or death.

V. SET UP & USE OF YOUR WHEELCHAIR

K. Pro Elevated Leg Rest

Installation

NOTE: The Pro ELR is mounted onto the chair the same way as a swing away hanger.


1. Place swing away pivot saddle into the receiver on front frame tube with the leg rest facing either inward or outward from the frame.
2. Rotate the leg rest so that it aligns with the frame and locks into place in the latch block.

Adjusting Height of Calf Pad

1. Loosen nut (Fig. 14:A) using a 10mm wrench.
2. Slide calf pad arm up or down to desired location. Retighten nut.


Adjusting Depth of Calf Pad

1. Remove screw (Fig. 14:B) and nut (Fig. 14:C) on calf pad arm using a 5mm Allen wrench and a 13mm wrench.
2. Pick the desired location based on the four predrilled holes and reinstall screw and nut.


Adjusting Length of Footrest

1. Loosen set screw (Fig. 15:A) with a 4mm Allen wrench.
2. Slide extension tube in or out to get to desired length and secure by tightening set screw.

Adjusting Knee Height

1. Using two 10mm socket wrenches, loosen the two nuts (Fig. 16:A) on the cover.
2. Adjust knee height to desired setting.
3. Retighten the two nuts (Fig. 16:A) after desired height is attained.

Use


1. To raise the Pro ELR, lift the leg rest tube (Fig. 17:A) to desired angle of elevation.
2. To lower the Pro ELR, press and hold the lever lock (Fig. 17:B) while pushing the leg rest tube (Fig. 17:A) down.

NOTE: Remove the Pro ELR from the chair or remove weight from the Pro ELR while lowering to avoid a sudden drop when the lever lock is pushed.

NOTE: The calf pad can swing outward to clear the front of the chair for transfers.

Removal

1. To remove leg rest, push or pull release latch.
2. Lift the leg rest straight upward to remove. You may also swing the leg rest inward or outward before lifting it off.

V. SET UP & USE OF YOUR WHEELCHAIR

L. Pediatric Pro Elevating Leg Rest

Installation

NOTE: The Pediatric Pro ELR is mounted onto the chair the same way as a swing away hanger.

1. Place swing away pivot saddle into the receiver on front frame tube with the leg rest facing either inward or outward from the frame.
2. Rotate the leg rest so that it aligns with the frame and locks into place in the latch block.


Adjusting Height of Calf Pad

1. Loosen nut (Fig. 18:D) using a 10mm wrench.
2. Slide calf pad arm up or down to desired location. Retighten nut.

Adjusting Depth of Calf Pad

1. Remove screw (Fig. 18:C), washer (Fig. 18:B) and nut (Fig. 18:A) on calf pad arm using a 5mm Allen wrench and a 13mm wrench.
2. Pick the desired location based on the four predrilled holes and reinstall screw and nut.

Fig. 18


Adjusting Knee Height

1. Using two 10mm socket wrenches, loosen the two nuts (Fig. 19:A) on the cover.
2. Adjust knee height to desired setting.
3. Retighten the two nuts (Fig. 19:A) after desired height is attained.

Use

1. To raise the Pediatric Pro ELR, lift the leg rest tube (Fig. 20:A) to desired angle of elevation.
2. To lower the Pediatric Pro ELR, press and hold the lever lock (Fig. 20:B) while pushing the leg rest tube (Fig. 20:A) down.

Fig. 19


NOTE: Remove the Pediatric Pro ELR from the chair or remove weight from the Pediatric Pro ELR while lowering to avoid a sudden drop when the lever lock is pushed.

NOTE: The calf pad can swing outward to clear the front of the chair for transfers.

Removal

1. To remove leg rest, push or pull release latch.
2. Lift the leg rest straight upward to remove. You may also swing the leg rest inward or outward before lifting it off.

Fig. 20


V. SET UP & USE OF YOUR WHEELCHAIR

M. Contracture Footrest

1. Changing Footrest and Footplate Angle


NOTE: Contracture Footrest is shown in the outside mount configuration. A center mount can be achieved by moving the leg rests to the opposite side.

- The footrests on a contracture footrest can be angled toward the chair and away from the chair by simply adjusting to upper pivot point (Fig. 21:A).
- The footplates on a contracture footrest can also be moved front to rear (Fig. 21:B), are angle adjustable (Fig. 21:C) and can be flipped up for transfer (Fig. 21:D).

2. Changing Footrest Height

- Remove screw (Fig. 22:B) and nut (Fig. 22:C) from the long pivot (Fig. 22:A) and slide extension tube to desired height. See Chart in Fig. 22 for footrest heights that correspond with each hole.


V. SET UP & USE OF YOUR WHEELCHAIR

N. Multi-Angle Footrest


1. Changing Footrest Angle (Fig. 23)

The Multi-Angle Footrest utilizes a ball joint which allows the footrest to be set in a variety of angles just by applying pressure to different areas on the top of the footplate.

2. Changing Footrest Height

- Remove footplate assembly from hanger by removing screw (Fig. 24:G), washer (Fig. 24:F), washer (Fig. 24:E), saddle spacer (Fig. 24:D), washer (Fig. 24:C), washer (Fig. 24:B) and nut (Fig. 24:A) using a 5mm Allen wrench and a 10mm wrench.
- Slide the footplate assembly to desired height and reinstall hardware in nearest hole.

V. SET UP & USE OF YOUR WHEELCHAIR




O. Backrest

1. Folding Backrest Down

- If chair is equipped with a fold down backrest, push levers (Fig. 25:A) toward the front of the chair to disengage the latches. If backrest seating is attached, both latches must be disengaged simultaneously in order to allow the backrest to fold.
- Push forward on the backrest to rotate it downward. Once it begins to rotate forward, the latch levers may be released.

2. Angle Adjustment

- Using a 4mm Allen wrench and a 10mm wrench, remove button head screws (Fig. 26:A) with their washers and nuts.
- Rotate backrest to desired angle. Each hole represents 5° of rotation. The upper row of holes correspond to -5°, 5°, 15° and 25° and the lower row of holes correspond to 0° (vertical), 10° and 20° (Fig. 27).
- Replace button head screw and tighten with washer and nut.
- Repeat for opposite side of chair.

V. SET UP & USE OF YOUR WHEELCHAIR

O. Backrest

3. Adjusting Height



- If chair is equipped with a height adjustable backrest, remove both screws (Fig. 28:A) of the telescoping backrest tubes using a 4mm Allen wrench and a 10mm wrench.
- Move upper tubes up or down to set height as desired.
- Replace screws with washers and nuts to set height.

4. Push Handle Angle Adjustment

- If chair is equipped with a height adjustable backrest, push buttons of locking hinges (Fig. 28:B). Both buttons must be pressed simultaneously to adjust angle.
- While holding buttons down, rotate push handle to desired position.
- Release buttons when desired angle is achieved.
- Verify that hinge is locked by pulling upwards on push handles.

5. Removable Stroller Handle

- Depress buttons (Fig. 29:A) on both sides of chair toward lower end of stroller tubes to disengage locking tabs.
- Pull handle out of receivers to remove or push stroller handle in to receivers to install. For best results, insert both sides at the same time.
- Verify snap buttons are engaged in receivers.

⚠ DANGER ⚠

DANGER: Do not attempt to lift or tilt the wheelchair by holding the center strut that joins the Adjustment Height Back push handles (Fig. 28:C) or the center strut that joins Removable Stroller Handle push handles (Fig. 29:B). These struts may break resulting in a fall or loss of control and may cause serious injury or death.

6. Dynamic Back

- See instructions for this back provided separately with chair.

V. SET UP & USE OF YOUR WHEELCHAIR

P. Reclining Backrest

WARNING

WARNING: Back height, back angle, seat depth, tilt angle, recline angle, seat height, size/position of rear wheels, size/position of front casters and any added weight to the back of the chair such as accessories, backpacks and oxygen tanks all relate to and affect the stability of the wheelchair. Any adjustments or change to one or more of the items listed above can decrease the stability of the wheelchair. Take caution when making changes to the above items. Consult a qualified technician when adjusting back height, back angle, seat depth, seat height, size/position of rear wheels and size/position of front casters. Inappropriate adjustments can result in serious injury.

Reclining Back Notes for Operation

1. After making any adjustments to the system make sure all attaching hardware is securely tightened.
2. Always make sure the wheelchair is stable in the fully reclined position. Depending on the angle of the backrest, the full recline position can be 85°, 75° or 65° and is relative to the full upright position. Full upright position is measured as 0°.
3. In order to ensure maximum safety and security for the patient in the chair, it is important to make sure the patient is properly positioned before reclining or inclining (sitting up) the backrest.
4. Always engage the wheel locks while reclining or inclining to ensure maximum safety of the patient.
5. Before using the reclining backrest, make sure that the anti-tippers are adjusted properly to ensure the wheelchair will not become unstable when in use. See anti-tip instructions on page 45 of this owner's manual.

DANGER

DANGER: Tilting or reclining the chair beyond level greatly increases the risk of falling out of the chair, which can result in serious injury or death.



V. SET UP & USE OF YOUR WHEELCHAIR

P. Reclining Backrest

Instructions for Use:

Recline / Incline:

1. To unlock the system prior to recline or incline adjustment by removing the locking pin from locking hole on each recline trigger (Fig. 30).
2. Make sure wheel locks are engaged and wheelchair is on a level surface.
3. Grip the handles securely and slowly depress both reclining triggers (Fig. 31) at the same time to release the gas cylinders. This will allow you to adjust the backrest position as follows:
 - a. **Reclining** - to open the angle of the backrest relative to the seat frame, or recline, by applying downward force on gas springs.
 - b. **Inclining** - return to upright position by applying force upward and pushing the backrest towards the front of the wheelchair.
4. Once desired backrest recline angle has been achieved, slowly release both reclining triggers.
5. To lock the system, insert the locking pin into the locking pin hole (Fig. 30). It is recommended that you always lock both reclining backrest triggers when you are done positioning the patient.

V. SET UP & USE OF YOUR WHEELCHAIR


P. Reclining Backrest

Angle Adjustment of Recline Backrest Plates:

⚠ WARNING ⚠

WARNING: Remove patient from wheelchair prior to making this adjustment. Failure to do so could result in death or serious injury.

1. The reclining option comes with three different backrest angle adjustment options. See Fig. 32 for range information prior to making the adjustment.
2. Using a 4mm Allen wrench and a 10mm wrench, remove button head screws (Fig. 32) with their washers and nuts.
3. Rotate backrest to desired angle. Each hole represents 65° of rotation (Fig. 32).
4. Replace button head screw and tighten with washer and nut until securely fastened.
5. Repeat for opposite side of chair.

V. SET UP & USE OF YOUR WHEELCHAIR

Q. Tilt Mechanism


1. Hand Operated (Trigger) Tilt

a. Tilting Seat Frame

- i. Squeeze trigger (Fig. 33:A).
- ii. While keeping triggers depressed, rotate seat frame to desired position.
- iii. Release trigger to lock into place.

b. Cable Adjustment


- i. Unlock jam nut (Fig. 33:B) from cable adjuster body (Fig. 33:C) using a 10mm and 3/8" open ended wrench.
- ii. Thread cable adjuster (Fig. 33:D) in or out until slack in cable is eliminated but also such that the cogs are still fully seated in the teeth of the rotary frames.
- iii. Verify function of the cable. Tilt mechanism should lock at all locations and cogs should not drag along the teeth at any location. If the cogs drag or if there is excessive play in the triggers, the male adjuster must be unthreaded more. If the cog is disengaged at all from the teeth, the adjuster must be threaded into the body further.
- iv. Once cable is adjusted and function is verified, re-lock the jam nut against the cable adjuster body by tightening the nut only, not the male adjuster or the adjuster body.

c. Locking Triggers


i. To secure the triggers in the locked position:

- a. Determine which style back canes and trigger were provided with the chair by checking against Fig. 34 and 35.
- b. Make sure the chair is locked by pushing on the back handles.
- c. Insert the pin attached to the trigger by a lanyard into hole number 1.
- d. Make sure the trigger no longer releases the tilt locks by gently squeezing the trigger and pushing on the back canes to tilt the chair.

ii. To secure the triggers in the unlocked position:

- a. Determine which style back canes and trigger were provided with the chair by checking against Fig. 34 and 35.
- b. Squeeze the triggers to unlock the tilt mechanism while holding the back canes and controlling the tilt position.
- c. Allow the seat frame to tilt until it comes to rest in a stable position.
- d. Insert the pin attached to the trigger by a lanyard into hole number 2 (Fig. 34:A and B).

NOTE: This feature is intended to be used solely for the purposes of adjusting CG which is covered in Section V: R: CG (Center of Gravity) Optimization.


WARNING: Never leave someone seated in the chair unattended with the triggers secured in the unlocked position. This may lead to serious injury to the person seated in the wheelchair.

V. SET UP & USE OF YOUR WHEELCHAIR

Q. Tilt Mechanism


2. Foot Operated Tilt

- a. Tilting Seat Frame
 - i. Depress foot tilt pedal (Fig. 36:A).
 - ii. Rotate seat frame to desired position while keeping pedal depressed.
 - iii. Release pedal to lock into place.
- b. Linkage Adjustment
 - i. Using a 3mm Allen wrench, loosen locking collar (Fig. 36:B) by loosening screw in collar.
 - ii. With cog (Fig. 36:C) fully engaged in rotary frame teeth and with tilt pedal released, reposition locking collar against end of tube (Fig. 36:D) and tighten clamp screw.

3. Tilt Stops


- a. Place tilt stop pads (Fig. 37:A) on both sides of each rotary frame (Fig. 37:B) at desired location. Pads should nest together.
- b. Place thrust plates (Fig. 37:C) into pockets on pads.
- c. Insert and tighten screw (Fig. 37:D) with nut and washer to lock into place.
- d. To change, remove tilt stops and reassemble according to the instructions above at the new desired location.

NOTE: Cable actuated mechanisms require periodic maintenance (See Fig. 38). Verify that cables are:

- Adjusted properly.
- Routed correctly.
- Bends in the cable are minimized to reduce friction.
- Replace the cables if they have any kinks or cannot be adjusted.

Fig. 38

V. SET UP & USE OF YOUR WHEELCHAIR


R. Seat Frame Set Up

1. Setting Seat Depth

- Using a 5mm Allen wrench and a 13mm wrench, remove depth adjustment screws from seat pan, or remove seat seat pan (See Section Z).
- Using a 4mm Allen wrench and a 10mm wrench, remove button head screw (Fig. 39:A) with nut and washer from front seat frame (Fig. 39:B) on each side of chair.
- Slide rear seat frame in or out as desired to set the seat depth.
- Replace button head screws in new location, placing screws in rearmost hole in front seat frame.

2. CG (Center of Gravity) Optimization

- Place occupant in chair, with all equipment that attaches to upper frame, such as foot boxes, backpacks, etc.
- Lock wheel locks (if equipped).
- Squeeze triggers (if hand operated tilt) or depress foot pedal (if foot operated tilt). Lock triggers in this position using locking pins (See Section Q. Step 1 -C.), or if foot operated tilt, use a zip tie or other means of holding foot pedal down.
- Allow rotary frame to find its resting place. A slight rocking motion applied to the push handles can be used to determine where the center of the resting place lies.
- Observe the tilt angle by viewing the angular scale on the side of the rotary frame. The angle is indicated by the rear edge of the bearing plate (Fig. 39:C). Resting angle for optimum performance is 17-23°. If frame rests outside of this range, remove the seat frame position locking screws (Fig. 39:D) using a 5mm Allen wrench.
- If resting angle is larger than 23°, pull seat frame toward the front of the chair until approximate 20° resting angle is achieved. If frame rests at an angle smaller than 17°, push seat frame toward the rear until approximately 20° is achieved.
- Different occupants will respond differently, but a good starting point is to move the seat frame ½" for every 5° that the angle is off from 20°.
- Once optimum position is achieved, reinstall and tighten positioning screws in holes closest to optimum position. Screws may be placed in front or rear hole of seat frame saddle (Fig. 39:E).



V. SET UP & USE OF YOUR WHEELCHAIR

S. Casters

1. Adjusting Caster Arm Position

- Using 13mm wrench, remove 8mm hex bolt (Fig. 40:A).
- Slide caster arm (Fig. 40:B) to desired location, aligning holes in caster arm with hole in base frame.
- Replace hex bolt thru aligned holes. Nut for bolt is captured in pocket inside of base frame.
- Ensure both sides are positioned identically.
- The back of hanger and the receiver tube (Fig. 40:C and 41:C) should not be positioned beyond the front of the caster stem bearing tube (Fig. 40:D and 41:D) as shown in Fig. 41. If it does, follow steps A. through D. until it appears as shown in Fig. 40.

NOTE: The hanger receiver may be positioned behind the caster stem bearing tube.


⚠ DANGER ⚠

DANGER: Failure to comply with the instructions above may cause the wheelchair to tip forward during use and result in a fall or loss of control and may cause serious injury or death.

2. Adjusting Caster Height

- Remove hex bolt (Fig. 42:A) and nut using 13mm wrenches.
- Remove caster wheel (Fig. 42:B) and spacers (Fig. 42:C).
- Feed hex bolt thru desired hole in one side of the caster fork.
- Slide one spacer over hex bolt.
- Slide caster wheel over hex bolt and push hex bolt thru caster wheel, flush with the opposite side of the caster wheel bearing.
- Slide second spacer between caster wheel and caster fork, aligning spacer with bearing hole.
- Push hex bolt thru spacer and hole in opposite side of caster fork.
- Retighten hex bolt with nut.
- Repeat for opposite side of chair. Ensure both casters are set up identically.


V. SET UP & USE OF YOUR WHEELCHAIR

T. Axle Plate

1. Adjusting Axle Plate Position

- a. Using two 13mm wrenches, remove axle plate screws (Fig. 43:A) from both sides of chair, noting orientation of anti-tip receiver (Fig. 43:B).
- b. Slide axle plates (Fig. 43:C) simultaneously forward or rearward to desired location.
- c. Replace and tighten axle plate screws.
- d. If equipped with foot operated tilt, the tilt lock linkage will need to be readjusted (See Section Q, Step 2:B).
- e. Ensure both axle plates are in the same locations on the base frames.

⚠ DANGER ⚠

DANGER: Adjusting the axle plate too far forward can result in a chair that may tip over backward resulting in a fall or loss of control that could cause severe injury or death.

⚠ WARNING ⚠

WARNING: Increasing the seat to back angle will change the center of gravity of the wheelchair. Make sure your axle position is adjusted properly in order to maintain constant stability of the wheelchair. Reference instruction for axle plate position at top of page and Fig. 43.

V. SET UP & USE OF YOUR WHEELCHAIR

T. Axle Plate

2. Adjusting Axle Sleeve Height

- a. Remove rear wheels.
- b. Using two 30mm wrenches, remove axle sleeve nut (Fig. 44:A) and washer (Fig. 44:B) from outside of axle plate.
- c. Remove axle sleeve (Fig. 44:C) from axle plate and reposition in desired height location.
- d. Reinstall and tighten axle sleeve nut.
- e. Additionally, the axle plate screws may be removed (See step 1) and axle plate may be moved up or down by 1" (2.54cm) by using alternate set of holes in axle plate. Using the lower set of holes (Fig. 44:D) yields a seat height 1" (2.54 cm) lower than using the upper set of holes (Fig. 44:E) for a given axle sleeve location. If moving axle plate up or down, ensure that the anti-tip receiver (Fig. 44:F) moves up or down (respectively) with the axle plate. All four holes of the anti-tip receiver should line up with all four holes in the axle plate. See Section V for further instructions on the anti-tip receiver orientation.
- f. Ensure that both axle sleeves and axle plates are set up identically.

3. Adjusting Axle Spacing

- a. Loosen axle sleeve nut (Fig. 44:A) using two 30mm wrenches.
- b. Thread axle sleeve (Fig. 44:C) in and out of nuts to achieve desired wheel spacing.
- c. Retighten axle sleeve nuts.
- d. Ensure both sides of the chair have the same axle spacing.
- e. Changing the axle spacing may require spacing the wheel locks in or out (See Section V).




Fig. 44

V. SET UP & USE OF YOUR WHEELCHAIR

U. Rear Wheels

1. Installing Wheels

- Push in the quick release button on the axle to allow the locking balls to retract. Make note of the difference between the extended and depressed position of the axle release button and its effect on the locking balls on the other end of the axle (Fig. 45).
- Insert the axle into the bearing housing on the wheel if it's separate.
- Push on quick release button again and slide axle into axle sleeve (Fig. 45).
- Release the button to lock axle in sleeve. If release button does not fully extend and the locking balls do not move into the locked position after releasing the button, the axle length needs to be adjusted (Fig. 46).

2. Adjusting Axles

- To adjust the axle you will need a 19mm wrench to turn the outside axle nuts.
- You will also need an 11mm wrench to hold the ball bearings on the opposite end of the axle and prevent the axle from turning.
- When the wheel is installed, if the axle does not lock, loosen the axle nut until it locks and axle can move in and out slightly when locked. If the axle locks, but moves in and out slightly, proceed to next step.
- Turn the axle nut counterclockwise until it is tight, then turn the axle nut clockwise 1/3 turn at a time. Check to see if it is properly adjusted after each 1/3 turn. When properly adjusted, there should be no perceptible movement of the axle in and out in the sleeve, but the axle release button should be easy to push in.

3. Removing Wheels

- Hold the wheel close to the hub and push in the button on the outside end of the axle.
- While still holding the button, pull the wheel and axle out of the axle sleeve.

DANGER: Make sure the push button is completely extended and the locking balls on the inside of the chair are fully engaged before operating the wheelchair. Failure to do so may result in the wheel falling off and may cause serious injury or death.

Review and understand Section U. Rear Wheel Installation and Removal before attempting an axle adjustment.

V. SET UP & USE OF YOUR WHEELCHAIR


V. Wheel Locks

1. Adjusting Wheel Locks

- a. Use a 10mm wrench or socket to loosen (do not remove) hex bolt (Fig. 47:A) beneath mounting block of base frame.
- b. Slide mount plate (Fig. 47:B) forward or rearward to achieve proper locking.
- c. Retighten hex bolt.

2. Reversing Wheel Lock Mount

- a. Using a 10mm wrench, remove the nuts (Fig. 47:C) from the wheel lock studs and remove wheel lock mechanism (Fig. 47:D).
- b. Remove hex bolt (Fig. 47:E) and nut (Fig. 47:F) using 10mm wrenches.
- c. Flip wheel lock mount (Fig. 47:G) forward or reverse.
- d. Reinstall and tighten hex bolt (Fig. 47:E) and nut (Fig. 47:F).
- e. Reinstall wheel lock mechanism (Fig. 47:D) and fasten with nuts (Fig. 47:C). Note that large diameter washer (Fig. 47:H) goes beneath the nut that covers the scalloped hole (Fig. 47:I). Always install the wheel lock using the scalloped hole that allows the front stud of the wheel lock mechanism to lie below the rear stud, i.e., if the single hole is toward the front of the chair, install wheel lock to that the rearmost stud of the wheel lock mechanism is in the upper hole of the scalloped hole. If the single hole is toward the rear of the chair, place the front stud of the mechanism in the lower hole of the scalloped hole.
- f. Note that if the chair is equipped with attendant foot lock and you are changing the wheel size, you will need to obtain a new linkage rod. Please contact Ki Mobility Customer Service to order.


V. SET UP & USE OF YOUR WHEELCHAIR

V. Wheel Locks

3. Reversing Wheel Lock Mounting Plate


- Remove wheel lock mounts according to steps 2:A and 2:B on previous page.
- Using a 10mm wrench, loosen hex bolt (Fig. 48:A) and slide wheel lock mounting plate (Fig. 48:B) off of mounting block on base frame. Repeat for opposite side of chair.
- Note lateral offset of wheel lock mounting plates. Choose mounting configuration of wheel lock mounting plates that best suits the axle spacing of the chair. Standard or wide mounting (Fig. 49) may be achieved by switching the wheel lock mounting plates to the opposite side of the chair. For foot lock and/or dual post armrests, wide configuration should always be used.

- Once width configuration is chosen, slide wheel lock mounting plates back onto mount blocks, indexing T-nuts back into slots of mounting plates. Do not tighten hex bolt yet.
- Reattach wheel lock mounts in proper configuration according to steps 2:C and 2:E on previous page.
- Adjust wheel lock mounting plate position and tighten hex nut when proper mount plate location is achieved.

4. Adjusting Foot Pedal Angle (Attendant Foot Lock Only)

- Using 10mm wrenches, remove hex bolt (Fig. 50:A), nut (Fig. 50:B), plastic washer (Fig. 50:C) and bushing (Fig. 50:D) from wheel lock index plate (Fig. 50:E) and linkage rod (Fig. 50:F) on both sides of chair.
- Choose new hole according to preference to mount hardware.
- Retighten hardware, ensuring plastic washer (Fig. 50:C) is installed between index plate (Fig. 50:E) and linkage rod (Fig. 50:F). Note that the head of the hex bolt is to be installed toward the outside of the chair.

NOTE: Always loosen and tighten wheel hardware by alternating between the two bolts while loosening/tightening a little at a time. This prevents overclamping on one set of hardware which can lead to binding of the fasteners and increased difficulty in removal.

V. SET UP & USE OF YOUR WHEELCHAIR

W. Drum Brake

1. Brake Operation

- a. Squeeze brake levers (Fig. 51:A) on both sides of chair to activate brakes. Each side is operated independently and may be squeezed with more or less force to get more or less braking action.
- b. Let go of brake levers to release brake.
- c. For parking brake function, squeeze each brake lever and, with trigger squeezed, use index finger to flip locking lever (Fig. 51:B) down into teeth on lever mount. While holding locking lever down, release brake handles and parking brake will remain engaged.
- d. To release parking brake, squeeze brake levers until spring loaded locking levers snap back. Release brake levers to deactivate the brake.




2. Brake Adjustment

- a. Using a 10mm wrench, unlock jam nut (Fig. 52:A) on adjuster (Fig. 52:B). Unthread adjuster to get tighter braking action or thread adjuster in to reduce braking action.
- b. It is recommended to adjust the brake until there is a slight bit of drag on the wheel when spun. Then thread the adjuster inward about $\frac{1}{2}$ turn, or whatever is necessary to remove all drag from the spinning wheel.
- c. Retighten jam nut to lock the adjustment in place. Repeat for opposite side of chair.

3. Axle Sleeve Position Adjustment

- a. Remove both wheels. Using a 24mm wrench, loosen axle sleeve nut (Fig. 53:A) and washer. Repeat for opposite side of chair.
- b. Move both brake assemblies to desired position.
- c. Replace axle sleeve washer and nut and tighten nut.
- d. Drum brake cables may need to be adjusted after repositioning axle sleeve.
- e. Caster wheel position will also need to be adjusted the same amount as the axle sleeve. Ensure that, when seat frame is at 0° , the seat height, as measured at the front of the chair, is within $\frac{1}{4}$ " of the seat height as measured at the rear of the chair. Adjust caster or axle sleeve as necessary.

NOTE: Axle sleeve spacing may not be adjusted with drum brake option.


V. SET UP & USE OF YOUR WHEELCHAIR

X. Frame Width

1. Frame Width Adjustment


- a. Remove backrest seating.
- b. If equipped with a seat pan, remove seat pan and adjust to new width (See Section Z).
- c. Remove screws in growth tubes of all cross struts - Base Frame (Fig. 54:A), Axle Plate (Fig. 54:B), Seat Frame (Fig. 54:C), Rotary Frame (Fig. 54:D), Tilt Lock Mechanism (Fig. 54:E) and Angle Adjustable Push Handle / Removable Stroller Handle (Fig. 54:F) (if equipped). Two 10mm wrenches are needed for the screws on items 1, 2 and 3, two 8mm wrenches for the screws on items 4 and 5, and a 4mm Allen wrench for screws in item 6.
- d. Set chair to desired width, aligning holes in growth tubes/struts.
- e. Reinstall but **DO NOT TIGHTEN** all screws/nuts/washers into growth tube holes.
- f. Once all screws are reinstalled, place chair on a flat surface with wheels installed. Ensure all growth tubes are centered on the chair (push handle/removable stroller handle will be offset by ½" on odd chair widths).
- g. Retighten all strut fasteners. Ensure chair tilts easily.
- h. If binding occurs when tilting, loosen axle plate/base frame and rotary frame, tilt rotary frame 2 - 4 times to allow chair to align itself and retighten fasteners, ending with rotary frame struts.

V. SET UP & USE OF YOUR WHEELCHAIR

Y. Anti-Tips

1. Flipping Anti-Tips Up
 - a. Pull anti-tip (Fig. 55:A) rearward to disengage rotary lock studs (Fig. 55:B).
 - b. Rotate anti-tip upwards.
 - c. Reset anti-tip back onto rotary lock studs, ensuring studs are nested into alignment slots.
2. Removing Anti-Tips
 - a. Depress buttons (Fig. 55:C) on anti-tips.
 - b. Pull anti-tip straight back to remove anti-tip from receiver.
 - c. To reinstall, press buttons in and push anti-tips back into receivers, aligning holes in receiver with front detent button.
3. Setting for Seat Height
 - a. Using a 4mm Allen wrench and a 10mm wrench, remove cross bolt (Fig. 55:D) from height adjuster (Fig. 55:E).
 - b. Slide height adjuster up or down as necessary, aligning holes in adjuster with hole in anti-tip. The distance between the bottom of the anti-tip wheel and the floor should be between 1.5" and 2" (3.81 and 5.08 cm).
 - c. Reinstall cross bolt through height adjuster and tighten.
 - d. Repeat for opposite side of chair. Ensure anti-tips on both sides of chair are set up identically.

V. SET UP & USE OF YOUR WHEELCHAIR

Z. Seat Pan

1. Removal

- Using a 5mm Allen wrench, remove button head screws (Fig. 56:A) from corners of seat pan.
- Remove seat pan assembly (Fig. 56:B), seat pan brackets (Fig. 56:C), seat pan clamps (Fig. 56:D) and nuts (Fig. 56:E) from seat tubes.

2. Installation

- Set seat pan onto seat frame.
- Slide one seat pan bracket (Fig. 56:C) onto seat tube and link seat pan clamp (Fig. 56:D) onto bracket, wrapping around seat tube. Hold nut (Fig. 56:E) in slot of clamp as shown and thread button head screw (Fig. 56:A) through seat pan into clamp. Do not tighten. Repeat for all four corners of seat pan.


Fig. 56

- Adjust seat pan in so it is square to the chair and in the desired position, with the front edge of the seat pan approximately 3/8" behind vertical footrest hanger tube.
- Tighten all four button head screws using a 5mm Allen wrench when proper position is obtained.

3. Length Adjustment

- Using a 5mm Allen wrench and a 13mm wrench, remove all four button head screws (Fig. 56:F) from center area of seat pan.
- With 5mm Allen wrench, loosen all four button head screws (Fig. 56:A) from seat pan clamps (Fig. 56:C).
- Slide seat pans forward or rearward, aligning proper holes to achieve desired seat pan length.
- Reinstall and tighten button head screws (Fig. 56:F) with nuts.
- Reposition seat pan on seat frame as necessary and tighten button head screws (Fig. 56:A) in clamps.

4. Width Adjustment

NOTE: Seat pan width may be adjusted up to 2" (5.08 cm) wider for even widths or 1" (2.54 cm) wider for odd widths. For greater adjustment than this, new seat pans must be obtained. Contact Ki Mobility Customer Service.


- Remove seat pan (See step 1 above) and adjust width of chair to desired width (See Section X).
- Align seat pan clamp holes with desired set of holes in seat pan to achieve larger frame width.
- Reinstall and tighten button head screws and nuts through seat pan clamps (See step 2 above).

V. SET UP & USE OF YOUR WHEELCHAIR

AA. Bearings

1. Bearing Replacement (it is recommended this procedure be performed by a qualified technician)

- a. Using a 17mm wrench and a 6mm Allen wrench, loosen and remove all four bearing locknuts (Fig. 57:A).
- b. Remove flat head screws (Fig. 57:B and C) from bearing subassemblies.
- c. Remove all bearing assembly components from chair.
- d. Using new bearing/liners (Fig. 57:D), slide bearing shaft (Fig. 57:E) through one bearing/liner subassembly as shown.
- e. Place shaft with bearing/liner through slot in rotary frame from the outside and slide spacer (Fig. 57:F) and second bearing/liner over shaft. "Noses" of bearing liners should face each other.
- f. Align bearing shaft hole with hole in front bearing plate of base frame (Fig. 57:G) and slide 10mm x 50mm (Fig. 57:B) flat head screw through shaft.
- g. Install washers (Fig. 57:H) and locknut (Fig. 57:A) back onto flat head screw, finger tight.
- h. Repeat steps d through g for opposite side of chair.
- i. Install bearing shaft through bores of two more bearings with spacer (Fig. 57:F) in between, as in step e above.
- j. Slide onto 10mm x 60mm (Fig. 57:C) flat head screw through rear bearing plate holes (Fig. 57:I) as shown. Repeat on opposite side of chair.
- k. Set rotary frame on both sets of bearings, placing ridge in bottom of rotary frame into groove between bearing liners. Install swivel bushing (Fig. 57:J) with flange toward bearing and slide swivel plate (Fig. 57:K) over swivel bushing with retainer bearing nested in inside groove of rotary frame. Secure with washer (Fig. 57:E) and lock nut (Fig. 57:A). Repeat for opposite side of chair.

V. SET UP & USE OF YOUR WHEELCHAIR

BB. Vent and Battery Tray

1. Instructions for Use:

Consult your authorized Ki Mobility supplier to ensure the vent and battery tray are appropriate for the devices being secured to them.

Adjustments to the vent and battery tray should only be done by an authorized Ki Mobility supplier.

2. Instructions for Transit Use:

When using this wheelchair as a seat in a motor vehicle you must remove any items attached to the vent and battery tray and properly secure them separately. Only authorized dealers should install, remove or adjust the vent and battery tray components.

DANGER: Failure to remove items attached to the vent and battery tray and properly secure them separately in the motor vehicle can result in these items becoming dangerous projectiles in the event of an accident. This can cause severe injury or death.

CC. Oxygen Tank Holder

1. Instructions for Use:

This oxygen tank holder has been designed to secure D cylinders to your wheelchair. Make sure that the bottle is well seated in the bottom of the holder and that the clamp is tightened to the point that the cylinder cannot be pulled out.

Oxygen tanks can be dangerous, make sure you follow all instructions for use from the tank's provider.


DANGER: Failure to follow these instructions and the instructions for use of the cylinder from the cylinder's provider can result in severe injury or death.

V. SET UP & USE OF YOUR WHEELCHAIR

DD. Power Tilt Interface Control Box

The Focus Power Tilt uses a Reac RCB10 battery and interface control box when operating with adaptive switches. The Reac Control Box - RCB10 is a compact and flexible 4 channel control box with a 2 Ah battery pack and an external charger. The RCB10 is designed to comply with standards for powered wheelchairs (EN 12184:2014) and assistive products (EN 12182:2012). See Fig. 58.

Interface Control Box Mounting

Mount the interface control box and cable in a location that the cable will not be pinched, cut, trapped by moving components or obstruct tilt movement. The interface control box must be close enough to the battery so that the cable is not stretched or pulled free. Ensure cable is not exposed to being caught or snagged by any other objects.

Connecting Adaptive Switches to Interface Control Box (See image above)

For a single switch, connect the single switch interface cable to port three on the interface control box. For a dual switch, connect the two dual switch interface cables to port one and two of the interface box. Route the adaptive switch cable(s) so it will not be pinched, cut, trapped by moving components or obstruct tilt movement. Ensure there is enough slack in the cable(s) so they are not stretched or pulled free from the interface control box or mounting location.

NOTE: The interface control box cable accepts a 3.5mm / 1/8" mono phono jack.

WARNING


WARNING: Ensure hand control box and cables are secure and out of the way of any moving parts to prevent possible damage and/or a possible fall which can result in injury or death.

V. SET UP & USE OF YOUR WHEELCHAIR

DD. Power Tilt Interface Control Box

Initial Charge

1. Identify the box with the charger.
2. Connect the charger cable to the extension cable (Fig. 59:A).
3. Plug the other end of the charger into a wall outlet for an initial charge.
4. Fully charge battery before first use.

Charging


For more information on charging and the battery charger, refer to the battery charger owner's manual that is included in the battery charger box.

- The battery is charged through the external charger.
- Charging is only allowed with charger provided.
- Charging time for a completely discharged system is approximately 6 hours.
- Charging must take place in a room with good air ventilation and a temperature of 0-40°C (32-104°F).
- Operation is prevented by the system during charging.
- If stored, the battery must be charged before use.

Center of Gravity Adjustment

NOTE: A Center of Gravity adjustment is required for initial use and after any adjustment is made to the chair.


1. With the occupant in the chair, ensure two people are present to assist in the Center of Gravity adjustment.
2. One person must support the back of the chair to prevent a sudden drop while the other person assisting pulls the motor disconnect pin (Fig. 60:A) on the power tilt assembly.

V. SET UP & USE OF YOUR WHEELCHAIR

DD. Power Tilt Interface Control Box

6. Allow the chair to settle into position and read the angle measurement on the rotary frame. See Fig. 61.


7. The chair should settle around 20°, if not, the seat frame set up needs adjustment. See the Seat Frame Set Up section in your owner's manual.
8. If the chair does settle around 20°, reinsert the motor disconnect pin.

⚠ WARNING ⚠

WARNING: Ensure the motor disconnect pin is fully engaged before attempting a transfer. Failure to comply with the instruction above could result in serious injury or death.

Toggle Switch Hand Control Use

1. Ensure hand control toggle switch and the motor are connected to the battery before use. See Fig. 62.


V. SET UP & USE OF YOUR WHEELCHAIR

DD. Power Tilt Interface Control Box


2. Pulling the hand control toggle switch up will tilt the chair (Fig. 63).
3. Pushing the hand control toggle switch down will undo the tilt (Fig. 63).

NOTE: The power tilt motor has limits set to protect the system if too much energy is required to change the tilt angle. If a limit is triggered, the motor will reset to an active state within a few seconds. These same limits are responsible for the motor discontinuing when the rotary frame reaches the maximum tilt range or a tilt stop. Power tilt chairs with recline are more susceptible to a motor reaching a limit and discontinuing when the client is reclined because the center of gravity has changed. It may be necessary to reduce the recline angle of the client for the motor to operate normally.

Rotating Toggle Switch Hand Control (Fig. 64)

1. Loosen bolt holding clamp and hardware to the toggle switch until the toggle control can be detached using a 10mm wrench.
2. Rotate the toggle switch 90° so toggle switch is perpendicular with the armrest.
3. Reinstall the toggle switch onto the clamp hardware assembly and secure by tightening bolt with a 10mm wrench.
4. Now, pulling the hand control toggle switch back will tilt the chair and pushing the switch forward will undo the tilt.


V. SET UP & USE OF YOUR WHEELCHAIR

DD. Power Tilt Interface Control Box

Pendant Hand Control Use


1. Ensure hand control toggle switch and the motor are connected to the battery before use. See Fig. 65.

2. Push the down arrow button to tilt the chair (See Fig. 66).
3. Push the up arrow button to undo the tilt (See Fig. 66).
4. Push the power button to turn on or turn off the battery level lights (See Fig. 66).
5. Push and hold the up and down arrow buttons for three seconds to lock the key functions. Repeat to unlock key functions (See Fig. 66).

NOTE: The pendant hand control has a hook on the end of the control. This hook can be used to mount, or hang, the control on the wheelchair. Wherever the hand control is mounted, ensure cable is not in a position where it will be pinched, cut, trapped by moving components or obstruct tilt movement.

NOTE: The power tilt motor has limits set to protect the system if too much energy is required to change the tilt angle. If a limit is triggered, the motor will reset to an active state within a few seconds. These same limits are responsible for the motor discontinuing when the rotary frame reaches the maximum tilt range or a tilt stop. Power tilt chairs with recline are more susceptible to a motor reaching a limit and discontinuing when the client is reclined because the center of gravity has changed. It may be necessary to reduce the recline angle of the client for the motor to operate normally.

V. SET UP & USE OF YOUR WHEELCHAIR

DD. Power Tilt Interface Control Box

WARNING

WARNING: In the unlikely event of battery acid leakage and you get acid in your eyes or skin - rinse thoroughly with water and contact your local healthcare professional for further consultation.

WARNING

WARNING: CHOKING HAZARD - The system is equipped with a charger connector plug to seal the connector hole, be aware that this plug can be removed and swallowed by a small child.

Transfers

WARNING

WARNING: Avoid the switch system and other adaptive switches when transferring to prevent unintentional chair tilting. Failure to comply with the instruction above could result in serious injury or death.

Cleaning and Disinfection

The RCS system is designed to fulfill IPx4 (or IPx6 as an option), which allows cleaning with water. For disinfection, it is recommended to use soap or equal and clean with a paper towel/rag.

Standard Features
Controls: 1-4 channel individual actuators
Current Limit Control: Default 8 A
System Power: 24VDC
Batteries: <ul style="list-style-type: none">•2 Ah (sealed lead accumulators)•Ventilated housing
Charger: Power supply 120 VAC - 230 VAC
Protection Class: IPX4
Color: Black and orange
Weight: 1.9 kg (4.19 lbs) (including batteries)
EMC Conformity: <ul style="list-style-type: none">•EN 61000-6-2•EN 61000-6-3•EN 60601-1-2•ISO 7176-21
Operating Temperature: -15°C to +50°C (5°F to 122°F)
Storage Temperature: -40°C to +60°C (-40°F to 140°F)
Soft Start / Stop
Flammability Rating: UL94V-0
Single Fault Protection

V. SET UP & USE OF YOUR WHEELCHAIR

EE. Canopy

1. Folding (Fig. 67)


- Close - Lift the front end of the canopy and take the front end all the to the back end to fold the canopy out of the way.
- Open - Take the front end and pull out and down until the canopy is fully stretched out.

2. Quick Removal and Quick Install (Fig. 67)

- Push in the latches on both sides of the canopy and lift up and out of receivers.
- To reinstall canopy, install both ends into the receivers until latches catch the bottoms of the tube mount clamps and "lock" into place.

3. Rain Catch (Fig. 67)

- The edge of the canopy will fold up to catch and direct rain away from the chair user.

VI. MAINTENANCE

A. Inspecting Your Wheelchair

Regular and routine maintenance will extend the life of your wheelchair while improving its performance. Wheelchair repairs and the replacement of parts should be done by a qualified technician of an authorized Ki Mobility supplier.

1. General Inspections:

- a. Clean your chair at least once per month. You may need to clean your chair more frequently if you operate it in dirty environments, such as a worksite.
- b. Check to be sure that all fasteners are tight. Unless otherwise noted, fasteners should be tightened to 40 in./lbs (4.52 N*m).
- c. Check tires and casters:
 - Check the tire for tread wear. Replace the tires if they have flat spots, visible cracks or if the tread is worn off.
 - If you have inflatable tires with a valve stem, check the pressure and set to the pressure listed on the tire sidewall.

DANGER

DANGER: Replace worn tires. The wheel locks will not grip properly if you fail to maintain the air pressure shown on tire sidewall. This could result in a fall or loss of control and cause severe injury or death.

- d. Check spoke wheels for loose spokes.
- e. Check your wheel locks. As tires wear the wheel locks should be adjusted. See Section V. Wheel Locks.

VI. MAINTENANCE

B. Routine Maintenance

1. Weekly:

- a. Check wheel locks to be sure they are adjusted correctly.
- b. Check axle sleeves to ensure the axle sleeve nuts are tight.
- c. Check for broken, bent or loose spokes.
- d. Check that casters spin freely.
- e. Inspect tires and casters for wear spots.
- f. Check pneumatic tires for proper inflation.
- g. Wipe rotary frame glide surfaces with a damp cloth.
- h. Adjust tilt cable to ensure proper tilt lock function (See Section V., Part Q Step 1:B or 2:B).
- i. Inspect roller bearings to ensure proper function and wear.
- j. Ensure hand grips do not rotate or pull off.

2. Monthly:

- a. Inspect rear wheel axles and tighten if necessary.
- b. Inspect caster housing bearings for hair build up and remove if necessary.
- c. Inspect wheel locks to be sure assembly is tight. Make sure wheel locks properly engage the tires.
- d. Check that all fasteners are tight and secure.
- e. Inspect frame for any deformities, defects, cracks, dimples or bends. These could be signs of fatigue in the frame which could result in a failure of the chair. Discontinue use of the wheelchair immediately and contact your authorized Ki Mobility dealer.

3. Annually:

- a. Have wheelchair checked and adjusted by a qualified technician.

WARNING

WARNING: After adjustments and before using this wheelchair, make sure that all of the fasteners are tight and secure or injury or damage may occur.

CAUTION

CAUTION: Do not overtighten fasteners as this could damage the frame tubing.

DANGER

DANGER: Failure to read and comply with these instructions may result in a fall or loss of control causing severe injury or death to the user or others.

VI. MAINTENANCE

C. Cleaning

1. Axles and Wheels:

- a. Clean around the axles and wheels WEEKLY with a damp rag.
- b. Hair and lint will lodge in the caster housing. Disassemble the caster housing every six months to remove entangled hair.

NOTE: Do not use WD-40 or any other penetrating oil on this wheelchair. This will destroy the sealed bearings.

NOTE: Do not use any chemical cleaning agents on casters or tires.

WARNING

WARNING: Failure to read and comply with these instructions may result in damage to your wheelchair, a fall or loss of control causing severe injury to the user or others.

D. Storage

1. When not in use, keep your chair in a clean, dry area. Failure to do so may result in your chair rusting and/or corroding.
2. If your chair has been in storage for more than a few weeks you should make sure it is working properly. You should inspect and service, if necessary, all items in Section A. Inspecting Your Wheelchair.
3. If your chair has been in storage for more than two months, it should be serviced and inspected by your authorized supplier before you use it.

WARNING

WARNING: Failure to read and comply with these instructions may result in damage to your wheelchair, a fall, or loss of control causing severe injury to the user or others.

VII. TROUBLESHOOTING

A. Symptoms and Solutions

For optimum performance, your Focus CR wheelchair will need to be adjusted occasionally. On the next page is a list of potential symptoms that may occur when adjusting your chair, as well as one or more suggested solutions for each symptom. For best results, when troubleshooting, try one solution at a time before proceeding to a second solution. If symptoms are not relieved, it is recommended to take the chair to an authorized service center for repair.

Please see chart on next page.

Symptoms										Solutions
Chair does not track straight	Looseness in chair	Sluggishness in turning	One caster does not contact the floor	Caster flutters	Squeaking or rattling	Will not tilt/ Rattles when tilting	Difficulty in tilting	Chair tilts too quickly and is difficult to recover	Chair is difficult to tilt and recovers too quickly	
✓				✓						Add spacer washer between caster stem and bottom caster housing bearing on left or right side of chair.
✓		✓	✓							Ensure all tire pressures are correct and equal on both sides of chair.
✓	✓			✓			✓			Ensure all fasteners are tightened.
					✓					Use Teflon® based lubricant (such as Tri-Flow®) between modular frame components.
						✓				Ensure cable/foot tilt mechanism is adjusted and cogs disengage teeth (see section V. N. 1. b or V. N. 2. b).
						✓	✓			Inspect roller bearings for wear. Replace if necessary (see section V. X).
						✓	✓			Inspect rotary frame bearing surfaces for damage. Clean as necessary.
								✓	✓	Optimize CG setting (Section V. O. 2).

VIII. WARRANTY

A. Warranty Information

Ki Mobility warrants the frame, hangers, tilt assembly (including bearings) and quick-release axles of this wheelchair against defects in materials and workmanship for the life of the original purchaser. Power tilt components are warranted for two years except the batteries. All other Ki Mobility-made parts and components of this wheelchair are warranted against defects in materials and workmanship for one year from the date of first consumer purchase.

The expected life of the frame is five years.

Limitations to the Warranty


1. We do not warrant:
 - a. Wear items: Upholstery, tires, armrest pads, tubes, armrests and push-handle grips.
 - b. Damage resulting from neglect, misuse or from improper installation or repair.
 - c. Damage from exceeding weight limit.
2. This warranty is VOID if the original chair serial number tag is removed or altered.
3. This warranty is VOID if the original chair has been modified from its original condition and it is determined the modification resulted in failure.
4. This warranty applies in the USA only. Check with your supplier to find out if international warranties apply.

Ki Mobility's Responsibility

Ki Mobility's only liability is to replace or repair, at our discretion, the covered parts. There are no other remedies, expressed or implied.

Your Responsibility

- a. Notify Ki Mobility, via an authorized supplier, prior to the end of the warranty period and get a return authorization (RA) for the return or repair of the covered parts.
- b. Have the supplier send the authorized return, freight pre-paid, to:
Ki Mobility
5201 Woodward Drive
Stevens Point, WI 54481
- c. Pay any charges for labor to repair or install parts.

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