



Kimba neo

EN Instructions for Use (Qualified Personnel) 3

Table of contents

EN

| | | |
|----------|--|----------|
| 1 | Foreword | 6 |
| 2 | Application | 6 |
| 2.1 | Intended Use | 6 |
| 2.2 | Indications | 6 |
| 2.3 | Contraindications | 6 |
| 2.4 | Necessary qualifications | 6 |
| 2.5 | Service life | 6 |
| 3 | Safety | 7 |
| 3.1 | Explanation of Warning Symbols | 7 |
| 3.2 | General safety instructions | 7 |
| 3.3 | Safety instructions for assembly | 7 |
| 3.4 | Name plate and warning labels | 8 |
| 3.4.1 | Name plate | 8 |
| 3.4.2 | Warning Labels | 8 |
| 4 | Delivery | 9 |
| 4.1 | Scope of Delivery | 9 |
| 4.2 | Options | 9 |
| 4.3 | Storage | 9 |
| 5 | Preparation for Use | 9 |
| 5.1 | Assembly | 9 |
| 5.2 | Installing options/accessories | 10 |
| 5.2.1 | Headrest supports (multifunctional seating unit) | 10 |
| 5.2.1.1 | Headrest supports | 10 |
| 5.2.1.2 | Integrated headrests | 11 |
| 5.2.2 | Safety belts | 11 |
| 5.2.2.1 | Strap guide | 11 |
| 5.2.2.2 | Lap belt | 12 |
| 5.2.2.3 | Four-point belt | 14 |
| 5.2.2.4 | Five-point belt | 15 |
| 5.2.3 | Positioning aids | 17 |
| 5.2.3.1 | Thoracic supports | 17 |
| 5.2.3.2 | Abductor | 17 |
| 5.2.3.3 | Pelvic band | 17 |
| 5.2.3.4 | Chest and shoulder support | 18 |
| 5.2.3.5 | Fixation vest | 19 |
| 5.2.3.6 | Groin strap | 20 |
| 5.2.3.7 | Neck Support | 21 |
| 5.2.3.8 | Lumbar Padding | 21 |
| 5.2.3.9 | Hook & loop foot straps | 22 |
| 5.2.4 | Armrests | 22 |
| 5.2.5 | Tray | 22 |
| 5.2.6 | Canopy adapter | 22 |
| 5.2.7 | Rain cover | 23 |
| 5.2.8 | Summer/winter slip sack | 23 |
| 5.2.9 | Frame padding | 23 |
| 5.2.10 | Storage basket | 23 |
| 5.2.11 | Platform for Respirator | 24 |
| 5.2.12 | Buggy Board | 24 |
| 5.2.13 | Cup Holder | 24 |
| 5.2.14 | Footrest Padding | 24 |
| 5.2.15 | Anti-Slip Upholstery | 25 |
| 5.2.16 | Foot block for shortening the lower leg length | 25 |

| | | |
|-----------|--|-----------|
| 5.2.17 | Incontinence Cover | 25 |
| 6 | Adjustments..... | 26 |
| 6.1 | Prerequisites | 26 |
| 6.2 | Multifunctional seating unit (standard seat)..... | 26 |
| 6.2.1 | Adjusting the Seat Depth | 26 |
| 6.2.1.1 | Removing the seat depth reducer..... | 27 |
| 6.2.2 | Adjusting the seat width | 27 |
| 6.2.3 | Adjusting the backrest height | 28 |
| 6.2.4 | Further settings | 28 |
| 6.3 | Headrest supports (multifunctional seating unit)..... | 28 |
| 6.3.1 | Adjusting the headrest support | 28 |
| 6.3.2 | Fitting the integrated headrests..... | 28 |
| 6.4 | Positioning aids | 29 |
| 6.4.1 | Adjusting the thoracic supports..... | 29 |
| 6.4.2 | Adjusting the abductor..... | 29 |
| 6.4.3 | Adjusting the chest and shoulder support | 29 |
| 6.4.4 | Adjusting the fixation vest | 30 |
| 6.4.5 | Adjusting the groin strap | 30 |
| 6.4.6 | Adjusting the pelvic band | 30 |
| 6.4.7 | Adjusting the Neck Support..... | 30 |
| 6.4.8 | Adjusting the Lumbar Padding..... | 30 |
| 6.4.9 | Adjusting the hook & loop foot straps..... | 30 |
| 6.5 | Footplate | 30 |
| 6.5.1 | Adjusting the footplate height (lower leg length)..... | 30 |
| 6.5.2 | Adjusting the footplate angle (option)..... | 31 |
| 6.5.3 | Further settings | 31 |
| 6.6 | Safety belts | 31 |
| 6.7 | Spring tension | 31 |
| 6.8 | Pushbar | 31 |
| 7 | Use in a Wheelchair Accessible Vehicle..... | 31 |
| 7.1 | Instruction of Users | 31 |
| 7.2 | General Safety Instructions | 32 |
| 7.3 | Permitted Use..... | 32 |
| 7.4 | Required Options..... | 32 |
| 7.4.1 | Mygo Seating System - Specifics | 33 |
| 7.5 | Adjusting the Seating Unit..... | 33 |
| 7.6 | Adjusting and Securing the Seat Angle..... | 33 |
| 7.6.1 | Multifunctional Seating Unit..... | 33 |
| 7.6.2 | Mygo Seating System | 34 |
| 7.7 | Adjusting and securing the back angle | 34 |
| 7.7.1 | Multifunctional Seating Unit..... | 34 |
| 7.7.2 | Mygo Seating System | 35 |
| 7.8 | Securing the Product in the Vehicle..... | 35 |
| 7.9 | Securing the User in the Product | 36 |
| 8 | Handover..... | 37 |
| 8.1 | Final inspection | 37 |
| 8.2 | Transport to the customer | 37 |
| 8.3 | Delivery | 37 |
| 9 | Maintenance/Repair | 38 |
| 10 | Disposal | 38 |
| 10.1 | Disposal information..... | 38 |
| 10.2 | Information on re-use..... | 38 |
| 11 | Legal Information..... | 38 |
| 11.1 | Liability | 38 |

11.2 CE Conformity.....38
11.3 Warranty Terms and Conditions.....39
12 Appendices.....39
12.1 Technical Data.....39

1 Foreword

INFORMATION

Date of the last update: 2014-01-28

- ▶ Please read this document carefully.
- ▶ Follow the safety instructions.

INFORMATION

Missing instructions for use or service manuals can be downloaded from the homepage at www.ottobock.com. Missing instructions for use can also be ordered directly from the manufacturer's service department (see back cover or inside back cover for addresses).

The product was delivered as specified on the order form. Settings must be adjusted and the product must be fine-tuned according to the medical indications or the needs of the user; this must be done in the presence of the user. These instructions for use provide the information necessary for adjusting the settings. Adjustments should be made in close cooperation with the user or attendant.

Please note the following:

- The Instructions for Use (Qualified Personnel) are intended only for qualified personnel and remain with them.
- The operation of the product is described in the Instructions for Use (User). These instructions must be given to the user or attendant.
- We recommend checking the product settings regularly in order to assure optimal fit over the long term. A fitting is required every six months for children and youth in particular.
- The installation of all optional components is generally described in the service manual.
- The manufacturer assumes no liability for combinations with seating systems that have not been approved by the manufacturer (see Page 38).

2 Application

2.1 Intended Use

The rehab buggy is intended solely for transporting children with walking disabilities or who are unable to walk.

The rehab buggy has been approved for transporting only one child aged 1 to 10 years (depending on the size and weight of the child). It is operated by an attendant (conveyance by another person).

The rehab buggy can be used indoors and outdoors. It is fully adjustable to provide an optimised, individual fitting.

The product may be used only with the options listed on the product order form.

The Otto Bock warranty applies only when the product is used according to the specified conditions and for the intended purposes.

2.2 Indications

The broad range of settings allow it to be used for a disability or walking impairment, for example due to:

- Paralysis (paraplegia/tetraplegia)
- Loss of limbs (leg amputation)
- Limb defect/deformity
- Joint contractures/defects
- Neurological and muscular diseases

2.3 Contraindications

None

2.4 Necessary qualifications

The installation and adjustments described below may only be carried out by qualified personnel.

2.5 Service life

Based on market observations and the current state of technology, the manufacturer has calculated that the product can be used for a period of **4 years**, provided that it is used properly and that the service and maintenance instructions are observed. Storage times at the dealer or with paying parties are not included in this period.

However, we would like to emphasise that the product is reliable far beyond this defined period of time, provided that it is cared for and maintained properly.

If the service life is reached, the user or a responsible attendant should contact the qualified personnel who fitted the product or the manufacturer's service (see inside or outside of rear cover for address).

3 Safety

3.1 Explanation of Warning Symbols

| | |
|------------------|---|
| ⚠ WARNING | Warnings regarding possible risks of severe accident or injury. |
| ⚠ CAUTION | Warnings regarding possible risks of accident or injury. |
| NOTICE | Warnings regarding possible technical damage. |

3.2 General safety instructions

| |
|---|
| ⚠ WARNING |
| Lack of instruction |
| Risk of falling, tipping over due to lack of knowledge |
| ▶ Instruct the user or the attendant in the proper use of the product when handing it over. |

| |
|---|
| ⚠ CAUTION |
| Use of unsuitable tools |
| Pinching, crushing or damaging the product |
| ▶ When working, use only tools that are suitable for the conditions at the workplace and whose proper use ensures safety and health protection. |

| |
|--|
| NOTICE |
| Tipping or falling of the product |
| Damage to the product |
| ▶ When you work on the product, secure it so that it cannot tip over or fall over. |
| ▶ Use a clamping device to secure the product for all work at a workbench. |

| |
|--|
| NOTICE |
| Use of unsuitable packaging |
| Damage caused by transport in the wrong packaging |
| ▶ Use only the original packaging for delivery of the product. |

3.3 Safety instructions for assembly

| |
|--|
| ⚠ CAUTION |
| Unsecured screw connections |
| Risk of pinching, crushing, tipping over, falling of user |
| ▶ After all adjusting/readjusting work, retighten the mounting screws/nuts firmly. Observe any torque settings which may be specified. |
| ▶ Any time you loosen a screw connection with thread lock, replace it with a new screw connection or secure the old screw connection with medium strength thread locker (e.g. Loctite® 241). |

| |
|--|
| ⚠ CAUTION |
| Improper use of self-locking nuts |
| Risk of the user falling or tipping over |
| ▶ Always replace self-locking nuts with new self-locking nuts after disassembly. |

3.4 Name plate and warning labels

3.4.1 Name plate

| Label | Meaning |
|---|---|
| <p>The name plate "Outdoor mobility base" is located on the back connecting tube of the base frame.</p> <p>The name plate "Seating unit" is located on the bottom back plate.</p> | A Type designation |
| | B Manufacturer article number |
| | C Maximum load capacity (see the chapter "Technical Data") |
| | D Manufacturer information/address/country of origin |
| | E Serial number |
| | F European Article Number/Global Trade Item Number |
| | G Read the instructions for use before using the product. |
| | H CE marking – product safety according to EU guidelines |

3.4.2 Warning Labels

| Label | Meaning |
|---|--|
| <p>Achtung: Bitte beachten Sie die maximale Zuladung des Untergestells!</p> <p>Attention: Please note the permissible vehicle payload!</p> | <p>Attention: Please observe the maximum load capacity of the mobility base! (See the section "Technical Data")</p> |
| | <p>Fixation point for fastening the product in wheelchair accessible vehicles</p> |
| <p>WARNING</p> <p>ACHTUNG: Lassen sie Ihr Kind nicht unbeaufsichtigt!</p> <p>ATTENTION: Never leave your child unattended!</p> <p>ACHTUNG: Dieser Sitz ist für Kinder unter 6 Monaten nicht geeignet!</p> <p>ATTENTION: This seating unit is unsuitable for children aged below 6 month!</p> <p>ottobock.</p> | <p>WARNUNG/WARNING</p> <p>ACHTUNG! Lassen Sie Ihr Kind nicht unbeaufsichtigt! ATTENTION! Never leave your child unattended!</p> <p>ACHTUNG! Dieser Sitz ist für Kinder unter 6 Monaten nicht geeignet! ATTENTION! This seating unit is unsuitable for children aged below 6 months!</p> |

4 Delivery

4.1 Scope of Delivery

The product is delivered pre-assembled, with the rear wheels and seating unit removed, and packed in a crate.

The scope of delivery includes:

- Outdoor mobility base, folded
- Multifunctional seating unit with padding (standard option) or an alternative seating system
- Further options, depending on the order
- Tools required
- Instructions for Use (Qualified Personnel), Instructions for Use (User)

INFORMATION

Only provide the Instructions for Use (User) to the user.

4.2 Options

The functionality and operation of the options are described in more detail in the instructions for use (user).

All of the available options/accessories are listed on the order form.

4.3 Storage

The product must be stored in a dry place.

For transportation and storage, an ambient temperature between 10 °C and +40 °C must be maintained.

Note regarding the tyres

- Tyres contain chemical substances that can react with other chemical substances (such as cleaning agents, acids, etc.).
- Direct exposure to sunlight or UV radiation causes the tyres to age more quickly. As a result, the tread surface hardens and corner pieces break out of the tread.
- Avoid unnecessary parking outdoors. The tyres should be replaced every 2 years regardless of wear and tear.

5 Preparation for Use

5.1 Assembly

⚠ CAUTION

Exposed pinch points

Risk of crushing or pinching the fingers

- ▶ Do not reach into the folding mechanism when folding/unfolding.

⚠ CAUTION

Improper assembly

Risk of the user falling or tipping over

- ▶ After each assembly, verify the proper fit of the removable wheels. The quick-release axles must be firmly locked in the wheel attachment.
- ▶ Before using the product, the folding mechanism must be firmly locked in place.
- ▶ When attaching the seating unit, make sure it clicks audibly onto the mobility base.
- ▶ Secure the seating unit using the seat lock.

The following steps are needed to make the rehab buggy ready for use:

Attaching the wheels (if wheels are detached)

- 1) Press the lock on the rear wheels and on the fixed front wheels.
- 2) Push the wheels onto the axle until the lock audibly clicks into place.

Unfolding the rehab buggy (when the seating unit is detached)



- 1) Unfold the pushbar (see Fig. 1, item 1).
 - 2) Unfold the rehab buggy by pulling the pushbar (see Fig. 1, item 2).
- The release handles of the folding mechanism click into place automatically after unfolding (see Fig. 1, item 3).

Attaching the seat (multifunctional seating unit)

INFORMATION

The product can also be equipped with alternative seating systems. In this case, please observe the information in the separate instructions for use for the seating system.

- 1) Set the seat onto the rear tube of the seat adapter at an angle of approximately 45°.
- 2) Press the front edge of the seat onto the mobility base until the clutch audibly snaps into the front tube of the seat adapter.
- 3) Slide the seat lock to the left to engage it (see Fig. 2).



Attaching the padding (multifunctional seating unit)

INFORMATION

- ▶ Mount the pads before installing the standard padding.
- ▶ Prior to the assembly of the padding, adjust the seating unit (see chapter "Adjustments").

- 1) Open the zip on the backrest padding.
- 2) Pull the back padding with its cap over the back insert and close the zip (see Fig. 3).
- 3) Pull the excess back padding down between the lower edge of the back and the seat and attach it to the twist locks on the outside of the back.
- 4) Place the seat cushion on the seat in such a way that the corresponding velcro/fleece surfaces can be connected (see Fig. 4).




5.2 Installing options/accessories

5.2.1 Headrest supports (multifunctional seating unit)

5.2.1.1 Headrest supports

- 1) Pull the proper sides of the headrest supports through the padding (see Fig. 5).
- 2) Open the zipper of the padding cap (see Fig. 3) and push the padding behind the seat to the side.

- 3) Put the mounting screws through the slotted hole in the headrest support and through the hole in the back plate (see Fig. 6).
INFORMATION: If the headrest supports are to be mounted to the top holes (= mounting points of the back plate), the longer mounting screws, which come with the device, will have to be used instead of the regular mounting screws.
- 4) Fasten the headrest supports on the back of the back plate using the knurled nuts.
- 5) If necessary, adjust the width. Hand-tighten the knurled nuts.
- 6) Close the padding again.



5.2.1.2 Integrated headrests

Assembly is done at the factory following receipt of an order.

To change the distance of the headrests: see Page 28.



5.2.2 Safety belts

5.2.2.1 Strap guide


For improved belt guidance in the pelvic area, strap guides must be installed in the following cases:

- Lap belt; lap-belt part of the four/five-point belt: installation position see Fig. 7, item 1
- Bottom mounting strap of the chest and shoulder support: installation position see Fig. 7, item 2

- 1) Put on the strap guide from the inside, connecting the flat side to the fastening point of the seating unit (see Fig. 8, item 1).
- 2) Put the mounting screw through the fastening point from the outside (see Fig. 8, item 2) and screw it into the strap guide.
- 3) Pre-mount the strap guide with the washer/cap nut.
- 4) Twist the strap guide in the direction of use and finish assembly (see Fig. 7).




5.2.2.2 Lap belt




| | | | |
|---|--------------------|---|----------------------------|
| 1 | Strap guide | 4 | Adjustment buckle |
| 2 | Attachment fitting | 5 | Belt pads |
| 3 | Belt restraint | 6 | Release button/belt buckle |

Mounting the lap belt

- 1) Mount the strap guide according to the chapter "Strap guide" (see Page 11).
 - 2) Pull the belt out of the belt restraint (3) (see Fig. 10, item 1).
 - 3) Press the adjustment buckle (4) together and pull out the belt (see Fig. 10, item 2).
INFORMATION: Make sure not to push the plastic part beyond the guidance in the metal part.
 - 4) Pull out a belt pad (5).
 - 5) Push the belt between the seat and the backrest through the half-moon-shaped connecting piece to the back and through the first strap guide (see Fig. 12, item 1).
 - 6) Put the belt behind the backrest and through the second strap guide (see Fig. 13, item 1).
INFORMATION: If necessary, you can run the belt inside the seat cushion (no illustration).
 - 7) Run the belt between the backrest and seat to the front (see Fig. 13, item 2).
 - 8) String up the belt pad (5).
 - 9) Run the end of the belt through the adjustment buckle (4) and, to complete, through the belt restraint (3) (see Fig. 10, item 1/2).
- If the belt buckle (6) can be closed while being in the centre of the seat, the lap belt has been installed properly.




5.2.2.3 Four-point belt



| | | | |
|---|-------------------------------------|---|---------------------------------|
| 1 | Clip with rivet | 5 | Lap-belt part |
| 2 | Shoulder belt, left/right | 6 | Belt restraint |
| 3 | Belt pads | 7 | Adjustment buckle |
| 4 | Strap guide with attachment fitting | 8 | Belt buckle with release button |

Installing the four-point belt

- 1) Mount the strap guide (4) according to the chapter "Strap guide" (see Page 11, see Fig. 15, item 1).
 - 2) Rivet 2 clips (1) to the holes of the bottom back plate (see Fig. 15, item 2).
 - 3) Run the shoulder belts (2) through the padding/back plate (see Fig. 16) and secure them in the clips (see Fig. 15, item 3).
 - 4) If necessary, cut off protruding belt ends of the shoulder belts and melt the free ends by applying heat (such as from a lighter).
 - 5) Mount the lap-belt part (5) according to the chapter "Lap belt" (see Page 12).
- The four-point belt is installed properly if the latches can be clipped into the buckle (8).



5.2.2.4 Five-point belt



| | | | |
|---|-------------------------------------|----|---|
| 1 | Clip with rivet | 6 | Strap guide, fixed |
| 2 | Shoulder belt, left/right | 7 | Belt restraint |
| 3 | Belt pads | 8 | Adjustment buckle |
| 4 | Strap guide with attachment fitting | 9 | Crotch belt with buckle and adjustment buckle |
| 5 | Lap-belt part, left/right | 10 | Stop plate |

Installation start

- 1) Mount the loose strap guide on the right (4) according to the chapter "Strap guide" (see Page 11).
- 2) Rivet 2 clips (1) to the holes of the bottom back plate (see Fig. 15, item 2).
- 3) Run the shoulder belts (2) through the padding/back plate (see Fig. 16) and secure them in the clips (see Fig. 15, item 3).


Mounting the lap-belt part (5)

- 1) Run the lap belt on the right (5) between the seat and the backrest through the half-moon-shaped connecting piece to the back, through the strap guide and let it hang (see Fig. 18).
 - 2) Mount the fixed strap guide (6) according to the chapter "Strap guide" (see Page 11).
 - 3) Run the lap belt on the left (5) between the seat and the backrest through the half-moon-shaped connecting piece to the back (see Fig. 12).
 - 4) String up the belt restraint (7).
 - 5) Run the lap belt on the right (5) behind the seat through the adjustment buckle (8) and, to complete, through the belt restraint (7) (see Fig. 19 - see Fig. 21).
- For pressing together/threading through the adjustment buckle: see Page 12.



Mounting the crotch belt (9)

- 1) Thread the crotch belt, with the proper side, through the front slot of the seat plate (see Fig. 22, Pos. 1).
- 2) Thread the crotch belt twice through the stop plate (10) and back up through the seat plate (see Fig. 22, Pos. 2).
- 3) Thread the crotch belt a second time below the crossbar of the adjustment buckle (see Fig. 23, Pos. 1).



Completing the installation


► If necessary, cut off protruding belt ends of the shoulder belts and melt the free ends by applying heat (such as from a lighter).

→ The five-point belt is installed properly if the latches can be clipped into the buckle (9).

5.2.3 Positioning aids

5.2.3.1 Thoracic supports

- 1) Run the proper sides of the thoracic supports through the padding.
- 2) Put the mounting screws through the slotted hole in the thoracic support and the back insert (see Fig. 24, item 1).
- 3) Fasten the thoracic supports on the back of the back insert using knurled nuts.
- 4) If necessary, adjust the width/height. Hand-tighten the knurled nuts.





5.2.3.2 Abductor

- 1) Place the abductor in the desired position on the seat plate (see Fig. 25, item 1) and fasten it with screws from below loosely (see Fig. 26, item 1).

INFORMATION: When using a slip sack, make a cut in it in the specified spot, insert the slip sack and fasten the abductor with screws loosely.


- 2) Adjust the depth of the abductor in the slotted hole if necessary.
- 3) Now fasten the abductor with screws tightly.



5.2.3.3 Pelvic band

For installation information, consult the instructions that come with the product.



5.2.3.4 Chest and shoulder support




| | | | |
|---|---------------------------------------|---|-------------------------------------|
| 1 | Pair of clips with attachment fitting | 4 | Belt end of bottom attachment strap |
| 2 | Shoulder belt | 5 | Sewn-on clip |
| 3 | Bottom attachment strap | 6 | Pad |

Mounting the chest and shoulder support

- 1) Mount the strap guide according to the chapter "Strap guide" (see Page 11).
- 2) Mount the pairs of clips (1) using mounting screws/knurled nuts in the slotted holes of the bottom back plate (see Fig. 28, item 1/see Fig. 29, item 1).
- 3) Place the chest and shoulder support on the seat bottom – the sewn-on clips (5) must point towards the seat bottom.
- 4) Run the shoulder belt (2) through the padding/back plate (see Fig. 16) and secure them in the top clips (see Fig. 29, item 2).
- 5) Run the bottom attachment straps (4) between the seat bottom and backrest to the back through the strap guide and secure them in the bottom clips (see Fig. 29, item 3).
INFORMATION: The bottom padded ends of the chest and shoulder support are to rest subsequently on the user's spinae.
- 6) If necessary, cut off protruding belt ends of the attachment straps and melt the free ends by applying heat (such as from a lighter).





5.2.3.5 Fixation vest




| | | | |
|---|-------------------------|---|-------------------------|
| 1 | Clip with rivet | 4 | Bottom attachment strap |
| 2 | Shoulder belt | 5 | Buckle |
| 3 | Middle attachment strap | 6 | Pad |

Mounting the fixation vest

- 1) Rivet 6 clips (1) to the holes of the bottom back plate (see Fig. 31, item 1/see Fig. 32, item 1).
- 2) Place the fixation vest on the seat bottom - the buckles (5) must point towards the seat bottom.
- 3) Run the shoulder belts (2) through the padding/back plate (see Fig. 16) and secure them in the clips (see Fig. 32, item 2).
- 4) Run the middle/bottom attachment straps (3/4) around the backrest tubes to the back and secure them in the clips (1) (see Fig. 32, item 3).
INFORMATION: For improved guidance of the middle/bottom attachment straps (3/4), you can additionally install strap guides.
- 5) If necessary, cut off protruding belt ends of the attachment straps and melt the free ends by applying heat (such as from a lighter).





5.2.3.6 Groin strap



| | | | |
|---|-------------------------------------|---|----------------------------------|
| 1 | Underside groin strap | 4 | Long attachment strap |
| 2 | Attachment strap with loop | 5 | Strap adjustment/plug connection |
| 3 | Free belt end for length adjustment | 6 | Leg strap |

Mounting the groin strap

- 1) Place the groin strap, underside (1) down, on the seat cushion; the attachment strap with the loop (2) must face backwards.
- 2) Run the attachment strap with the loop (2) between the back support and seat bottom and secure it under the seat bottom (see Fig. 35). Use the screws/knurled nuts to adjust the depth/width of the seat (see Fig. 35, item 1).
- 3) Rivet the clips to the holes of the bottom back plate (see Fig. 34, Pos. 1; see Fig. 35, item 2).
- 4) Run the belt ends of the long attachment strap (4) around the back support and secure them in the clips (see Fig. 35, item 3).
- 5) If necessary, cut off protruding belt ends of the long attachment strap (4) and melt the free ends by applying heat (such as from a lighter).



5.2.3.7 Neck Support

- 1) Determine the medically required position of the neck support.
- 2) Thread the neck support through the back padding and the back plate.
- 3) Adhere the two hook-and-loop straps for securing the neck support to the back plate. Do this flexibly, according to the installation situation.
- 4) Secure the neck support firmly.



5.2.3.8 Lumbar Padding

- 1) Open the zip on the backrest padding and fold the backrest pad forward.
- 2) Determine the medically required position of the lumbar padding.
- 3) Adhere the supplied hook-and-loop straps for securing the lumbar padding to the back plate (see Fig. 40) or use the existing hook-and-loop straps.
- 4) Secure the lumbar padding.
- 5) Pull the back padding with its cap over the back insert and close the zip.



5.2.3.9 Hook & loop foot straps

Attach the hook-and-loop foot straps to the heel part of the footplate with two rivets or screws positioned beside one another. The holes must be inserted in the footplate accordingly.

In addition, hook-and-loop foot straps can also be used for the forefoot area.



5.2.4 Armrests

For mounting information, see the Instructions for Use (User) - reference number 647G770.

5.2.5 Tray



For mounting information, see the Instructions for Use (User) - reference number 647G770.

5.2.6 Canopy adapter

INFORMATION

The canopy adapter is used for the subsequent installation of the canopy to a multifunctional seating unit. For information on this, see the Instructions for Use (User) - reference number 647G770.

- 1) Loosen the set screws at the back tube (see Fig. 42, item 1) and pull the back insert up.
- 2) Slide the canopy adapter with the proper side onto the backrest bar (see Fig. 43).
- 3) Pre-position the canopy adapter in a level manner (see Fig. 43).
- 4) Tighten the set screws of the canopy adapter again.
- 5) Put on the back insert and tighten again the set screws at the back tube.



5.2.7 Rain cover

For mounting information, see the Instructions for Use (User) - reference number 647G770.

5.2.8 Summer/winter slip sack

NOTICE

Incorrect attachment

Damage to the product

- ▶ During the attachment, make sure that the product cannot grind against the wheels.

- 1) When needed, make a cut in the slip sack for putting through the lap belt (see Fig. 44, item 1), the abductor/five-point belt (see Fig. 44, item 2) or the groin strap (see Fig. 44, item 3).
- 2) Fasten the attachment cords depending on the attachment situation:
 - Back: e.g., at the seat plate (see Fig. 45, item 1)
 - Front: e.g., at the footrest tube (see Fig. 45, item 2)



5.2.9 Frame padding

For mounting information, see the Instructions for Use (User) - reference number 647G770.

5.2.10 Storage basket



NOTICE

Incorrect attachment

Damage to the product

- ▶ During the attachment, make sure that the product cannot grind against the wheels.

- ▶ Mount the storage basket, with the proper side, to the frame using the snap fasteners and the hook-and-loop strap (see Fig. 46, item 1/see Fig. 47, item 1).



5.2.11 Platform for Respirator

⚠ WARNING

Use in a wheelchair accessible vehicle is forbidden

Severe injury in case of accidents

- Use of the product in wheelchair accessible vehicles in combination with the artificial respiration platform option is not permitted.

INFORMATION

The respiration platform cannot be used on a product with fixed front wheels.

The platform for a respirator, if ordered, is installed when delivered. The platform only needs to be moved if the seat is to be turned around. More information on the installation can be taken from the Kimba neo service manual – article number 647G796.

Mounting instructions

- The respiration equipment is secured with belts using the slots on the respiration platform.
- It is recommended to use an anti-slip underlay.



5.2.12 Buggy Board

The installation is usually carried out by the user's attendant.

For mounting information, see the Instructions for Use (User) - reference number 647G770.

5.2.13 Cup Holder

The installation is usually carried out by the user's attendant.

For mounting information, see the Instructions for Use (User) - reference number 647G770.

5.2.14 Footrest Padding

- 1) Cut the self-adhesive hook-and-loop straps to length and secure them to the outer sides of the footrest.
- 2) Secure the padding to the hook-and-loop straps on the footrest.



5.2.15 Anti-Slip Upholstery

- 1) Clean the supporting surface of the footrest.
- 2) Adhere the self-adhesive anti-slip upholstery to the footrest, covering the whole surface.



5.2.16 Foot block for shortening the lower leg length

- 1) Clean the supporting surface of the footrest.
- 2) Glue on the self-adhesive hook-and-loop strap at the correct distance (see the underside of the foot block for the distance).
- 3) Press the foot block on to the hook-and-loop strap.



5.2.17 Incontinence Cover

The installation is usually carried out by the user's attendant.

For mounting information, see the Instructions for Use (User) - reference number 647G770.

6 Adjustments

6.1 Prerequisites

⚠ WARNING

Incorrect settings

Risk of falling, tipping over or improper user posture due to incorrect settings

- ▶ Adjustments may be carried out only by qualified personnel.
- ▶ Only the settings described in these instructions for use may be carried out.
- ▶ Settings may only be changed within the allowable adjustment ranges; otherwise, the stability of the product may be impaired (see this chapter and the chapter "Appendixes"). If you have questions, contact the manufacturer's service (see back cover for addresses).
- ▶ Only conduct tests in the presence of an assistant.
- ▶ Unless expressly described, you may not change any settings with a person sitting in the product.
- ▶ Secure the user against falling out during all tests.
- ▶ Before testing setting changes with the user seated, firmly tighten all screw connections.
- ▶ Check for safe function before delivering the product.

⚠ CAUTION

Unsecured screw connections

Risk of pinching, crushing, tipping over, falling of user

- ▶ After all adjusting/readjusting work, retighten the mounting screws/nuts firmly. Observe any torque settings which may be specified.
- ▶ Any time you loosen a screw connection with thread lock, replace it with a new screw connection or secure the old screw connection with medium strength thread locker (e.g. Loctite® 241).

⚠ CAUTION

Improper use of self-locking nuts

Risk of the user falling or tipping over

- ▶ Always replace self-locking nuts with new self-locking nuts after disassembly.

Fine adjustments and setting changes to match the concrete physical and mental condition of the user should always be made in the user's presence.

All parts of the product should be cleaned thoroughly before adjustments are made.



6.2 Multifunctional seating unit (standard seat)

INFORMATION

The product can also be equipped with alternative seating systems. In this case, please observe the information in the separate instructions for use for the seating system.



6.2.1 Adjusting the Seat Depth

- 1) Loosen the 2 knurled nuts on the underside of the seat (see Fig. 55, item 1).
- 2) Adjust the seat depth by moving the footrest (see Fig. 56).
- 3) Hand-tighten the knurled nuts.





6.2.1.1 Removing the seat depth reducer

- 1) Remove the covers.
- 2) Remove the hip pads (see Fig. 57).
- 3) Loosen and remove the four mounting screws on each seat depth reducer (see Fig. 58, item 1/2). Remove the bearing bushings.
→ The back frame is now separated from the base frame.
- 4) Attach the back frame with the bearing bushings inserted (see Fig. 59, item 3) to the base frame (see Fig. 60, item 1) using two mounting screws (see Fig. 59, item 1/2).
- 5) Attach the hip pads and pull over the covers.



6.2.2 Adjusting the seat width



- 1) Loosen the 4 knurled nuts on the underside of the seat (see Fig. 61, item 1).
- 2) Adjust the seat width by moving the hip supports (see Fig. 62). This also allows for abducted settings.
- 3) Hand-tighten the knurled nuts.



6.2.3 Adjusting the backrest height

- 1) Loosen the 2 set screws at the backrest tubes (see Fig. 63, item 1).
- 2) Adjust the back insert to the desired height.
- 3) Tighten the set screws again.

INFORMATION: You can achieve a greater back height by additionally shifting the top back plate at the back insert. To do so, loosen/remove the knurled nuts, shift the back plate and hand-tighten the knurled nuts again.



6.2.4 Further settings

The following fitting tasks are described in more detail in the included Instructions for Use (User) – reference number 647G770:

- Adjusting the back angle
- Adjusting the seat angle (seat tilt)

6.3 Headrest supports (multifunctional seating unit)

6.3.1 Adjusting the headrest support

- 1) Open the zipper of the padding cap (see Fig. 3) and push the padding behind the seat to the side.
- 2) Loosen the knurled nuts on the back of the back plate.
- 3) Adjust the width by moving the supports in the slotted holes (see Fig. 64, item 1).
- 4) Adjust the height by shifting the supports in the holes of the back plate (see Fig. 64, item 2).

INFORMATION: A further height adjustment of the headrest supports is possible by adjusting the back height.

- 5) Hand-tighten the knurled nuts.
- 6) Close the padding again.

6.3.2 Fitting the integrated headrests

- 1) Open the zipper of the padding cap (see Fig. 3) and remove the padding from the integrated headrests.
- 2) To fit the free space in the head region, glue up to 3 foam pads, each 1 cm thick, on the headrests (see Fig. 65, item 1/see Fig. 66, item 1).

- 3) Put back/close the padding again.



6.4 Positioning aids

⚠ CAUTION


Improper adjustments

States of panic, falling of the user

- ▶ If the product is adjusted too tightly, it may lead to unnecessary pain or suffering of the user. Immediately modify the settings if you detect signs of discomfort or fear in the user when using the product.
- ▶ Adjusting the belt straps too loosely can cause them to slip upwards/downwards so that the user can slide into a position that is dangerous to him/her. In addition, the fastening snaps could open unintentionally if they slide against hard parts of clothing (e.g., buttons).

6.4.1 Adjusting the thoracic supports

- 1) Loosen the knurled nuts on the back of the back insert.
- 2) Adjust the width/height by moving the supports in the slotted holes (see Fig. 67).
INFORMATION: If necessary, you can adjust the height of the supports by shifting the mounting screws in the slotted hole of the back plate.
- 3) Hand-tighten the knurled nuts.



6.4.2 Adjusting the abductor

The abductor can be adjusted as follows:

- Shifting of the mounting screw to other holes/slotted holes of the seat plate (see Fig. 25).
- Moving the abductor in the slotted hole (see Fig. 26)

6.4.3 Adjusting the chest and shoulder support

INFORMATION

The application of the product is described in the Instructions for Use (User).

The chest and shoulder support can be adjusted as follows:

- Moving the belt in the 2 pairs of clips on the back plate (see Fig. 29, item 2/3)
- Moving the belt in the sewn-on clips of the chest and shoulder support (see Fig. 27, item 5).

6.4.4 Adjusting the fixation vest

INFORMATION

The application of the product is described in the Instructions for Use (User).

The fixation vest can be adjusted as follows:

- Moving the belt in the 6 clips on the back plate (see Fig. 32, item 1-3)

6.4.5 Adjusting the groin strap

INFORMATION

The application of the product is described in the Instructions for Use (User).

The groin strap can be adjusted as follows:

- Moving the belts in the 2 clips on the back plate (see Fig. 35, item 2/3)
- Moving the belts in the strap adjustment at the leg straps (see Fig. 33, item 5/6)

6.4.6 Adjusting the pelvic band

INFORMATION

The application of the product is described in the Instructions for Use (User).

The pelvic band can be adjusted as follows:

- Moving the supports on the crossbar (see included Instructions for Use 647G129)
- Sliding the band onto the plastic mounting rail (see included Instructions for Use 647G129)

6.4.7 Adjusting the Neck Support

The adjustments are made when the option is installed (see section "Preparation for Use" > "Installing options/accessories").

6.4.8 Adjusting the Lumbar Padding

The adjustments are made when the option is installed (see section "Preparation for Use" > "Installing options/accessories").

6.4.9 Adjusting the hook & loop foot straps

INFORMATION

The application of the product is described in the Instructions for Use (User).


The hook-and-loop foot strap can be adjusted as follows:

- Moving the hook-and-loop foot straps against each other (no illustration)

6.5 Footplate

6.5.1 Adjusting the footplate height (lower leg length)

- 1) Loosen the Allen head screws (see Fig. 68, item 1).
- 2) Adjust the footplate height (lower leg length).
- 3) Tighten the Allen head screws again.



6.5.2 Adjusting the footplate angle (option)

- 1) Loosen the fastening nuts of the locking screws in the slotted holes (see Fig. 69, item 1).
- 2) Move the locking screws. Pay attention to the setting, so that a parallel adjustment is assured.
- 3) Align the wide slots of the locking screws perpendicular to the slotted holes in the footrest, so as to ensure wear-free locking.
- 4) Tighten the fastening nuts of the locking screws again.

6.5.3 Further settings

The following fitting tasks are described in more detail in the included Instructions for Use (User) – reference number 647G770:

- Adjusting the knee angle

6.6 Safety belts

⚠ CAUTION

Improper adjustments

States of panic, falling of the user

- ▶ If the product is adjusted too tightly, it may lead to unnecessary pain or suffering of the user. Immediately modify the settings if you detect signs of discomfort or fear in the user when using the product.
- ▶ Adjusting the belt straps too loosely can cause them to slip upwards/downwards so that the user can slide into a position that is dangerous to him/her. In addition, the fastening snaps could open unintentionally if they slide against hard parts of clothing (e.g., buttons).

The following fitting tasks are described in more detail in the included Instructions for Use (User) – reference number 647G770:

- Adjusting the lap belt
- Adjusting the four-point belt
- Adjusting the five-point belt

6.7 Spring tension

The settings are described in more detail in the included Instructions for Use (User) - reference number 647G770.

6.8 Pushbar

The settings are described in more detail in the included Instructions for Use (User) - reference number 647G770.

7 Use in a Wheelchair Accessible Vehicle

7.1 Instruction of Users

Based on the information that follows, the users and attendants must be instructed in the safe use of the product in wheelchair accessible vehicles.

This information is also found in the Instructions for Use (User).

7.2 General Safety Instructions

WARNING

Use as a seat in a wheelchair accessible vehicle

Severe injury in case of accident

- ▶ Always use the seats and restraint systems in the wheelchair accessible vehicle first. This is the only way to provide optimal protection for the occupants in case of accident.
- ▶ The product was tested for use in a wheelchair accessible vehicle in combination with a multifunctional seating unit (standard seat) or with a Mygo seating system.
- ▶ Please note that the use of the Mygo seating system in a wheelchair accessible vehicle is permitted only if a second locking pin for securing the seat angle is fitted on the Kimba outdoor mobility base (see the section "Adjusting and Securing the Seat Angle").
- ▶ Use of the product in combination with the artificial respiration platform option is **not** permitted.
- ▶ The product may only be used as a seat in a wheelchair accessible vehicle with the securing elements offered by the manufacturer and with suitable restraint systems.
- ▶ Please note that you must always remove certain options before using the product in a wheelchair accessible vehicle (see the section "Permitted Use").
- ▶ In general, you should only use the product in wheelchair accessible vehicles with the seat facing forwards (see the section "Adjusting the Seating Unit").
- ▶ In general, you should only use the product in wheelchair accessible vehicles with the seating unit secured horizontally and the backrest locked in the upright position (see the section "Adjusting and securing the seat angle" and "Adjusting and securing the back angle").

WARNING

Incorrect transport weight when using in a wheelchair accessible vehicle

Severe injury in case of accident

- ▶ When using the rehab buggy as a seat in a wheelchair accessible vehicle, the product may only be loaded as follows:
 - Kimba neo, **size 1/size 2**: max. user weight **40 kg (88 lbs)**
- ▶ The storage basket must be emptied prior to use in a wheelchair accessible vehicle.

WARNING

Use of the belt system as a restraint system in a wheelchair accessible vehicle is forbidden

The user may fall out of the product or sustain severe injuries if there is an accident

- ▶ Under no circumstances may the safety belts and positioning aids offered with the product be used as part of a restraint system for transportation in a wheelchair accessible vehicle.
- ▶ Note that the safety belts and positioning aids offered with the product are only intended to help support the user sitting in the product.

7.3 Permitted Use

The product may be used for seating in a wheelchair accessible vehicle if a suitable restraint system is used. The safety belts and positioning aids offered by the manufacturer only serve to further stabilise the person using the product.


The following options must be detached from the rehab buggy prior to transport in the wheelchair accessible vehicle:

- Sun canopy
- Grab rail
- Tray
- Diaper bag
- Backpack
- Buggy board

7.4 Required Options

A belt set is required to secure the rehab buggy for use as a seat in a wheelchair accessible vehicle.

- The lashing points are integrated in the base frame:
 - Lashing points, front: see Fig. 70, item 1
 - Lashing points, rear: see Fig. 71, item 1
- The belt set can be obtained from a third-party manufacturer. For more information on suppliers, please refer to our brochure on using our products in wheelchair accessible vehicles, order number 646D158.



7.4.1 Mygo Seating System - Specifics

- Use of the Mygo seating system in a wheelchair accessible vehicle is permitted only if a second locking pin for securing the seat angle is fitted on the Kimba outdoor mobility base (see Page 33).
- Use of the Mygo seating system in a wheelchair accessible vehicle is permitted only when a chest and shoulder belt or a chest and shoulder support is used.

7.5 Adjusting the Seating Unit

Before using the product in a wheelchair accessible vehicle, the seating unit must be adjusted in the driving direction (so the user faces forward in the direction of travel).

To remove/rotate and mount the seating unit: see the Instructions for Use (User).



7.6 Adjusting and Securing the Seat Angle

7.6.1 Multifunctional Seating Unit

Activating the Locking Mechanism

Before the product can be used in a wheelchair accessible vehicle, the seat angle must be adjusted horizontally and secured.



- 1) Turn the top rotary knob to the right and adjust the seat bottom so it is horizontal.
- 2) Pull out the locking pin, rotate it to the left and reinsert it (see Fig. 72, item 1).
→ The locking pin is now engaged (see Fig. 73, item 1).
- 3) Tilt the seat adapter (see Fig. 73, item 2) against the locking pin.
→ The seat adapter must contact the locking pin without play.



Deactivating the Locking Mechanism

After use in a wheelchair accessible vehicle, the seat angle adjustment must be enabled again:

- 1) Pull out the locking pin and secure it in the disengaged position by turning it (see Fig. 74).
- 2) The locking pin is disengaged (see Fig. 75, item 1). The seat angle can now once again be adjusted by turning the top rotary knob.



7.6.2 Mygo Seating System

Activating the Locking Mechanism

Before the product can be used in a wheelchair accessible vehicle, the seat angle must be adjusted horizontally and secured.

- 1) Twist the knob on the upper front part of the frame and adjust the seat bottom so it is horizontal.
- 2) Pull out the first locking pin, rotate it to the left and reinsert it (see Fig. 72, item 1).
→ The locking pin is now engaged (see Fig. 73, item 1).
- 3) Tilt the seat adapter (see Fig. 73, item 2) against the locking pin.
→ The seat adapter must contact the locking pin without play.
- 4) For additional stability, pull out the second locking pin on the opposite side, rotate it to the left and reinsert it (no picture available).

Deactivating the Locking Mechanism

After use in a wheelchair accessible vehicle, the seat angle adjustment must be enabled again:

- 1) Pull out the locking pins on both sides and secure them in the disengaged position by turning them (see Fig. 74).
- 2) The locking pins are disengaged (see Fig. 75, item 1). The seat angle can now once again be adjusted by twisting the knob on the upper front part of the frame.

7.7 Adjusting and securing the back angle

7.7.1 Multifunctional Seating Unit

Activating the Locking Mechanism


Before using the product in a wheelchair accessible vehicle, the back angle must be locked in the vertical position.

- 1) Turn the bottom rotary knob to the left and adjust the back angle to vertical (see Fig. 76).
- 2) Turn the bottom rotary knob approx. **5 mm** further to the left in order to activate the back securing bolt.
- 3) Pull out the back securing bolt and turn it 90° (see Fig. 77). The protruding locking aids on the back securing bolt (see Fig. 77, item 1) have to engage in the horizontal slotted holes (see Fig. 78, item 1).
- 4) Turn the bottom rotary knob to the right to secure the back angle.

Deactivating the Locking Mechanism

After use in a wheelchair accessible vehicle, the back securing bolt must be deactivated:

- 1) Turn the bottom rotary knob to the left until the back securing bolt is released.
- 2) Pull out the back securing bolt and turn it 90° (see Fig. 77). The protruding locking aids on the back securing bolt must engage in the vertical slotted holes (see Fig. 78, item 2).
- 3) Turn the bottom rotary knob to the right to secure the back angle without the back securing bolt (see Fig. 76).



7.7.2 Mygo Seating System

Before using the product in a wheelchair accessible vehicle, the back angle must be tilted back gently to an angle of **approx. 95°**.

INFORMATION: Please note the information in the separate instructions for use provided for the seating system.



7.8 Securing the Product in the Vehicle

INFORMATION

Please refer to the additional information in the brochure on using our products in wheelchair accessible vehicles, order number 646D158.

For transportation in a wheelchair accessible vehicle, the rehab buggy must be secured with tie-down straps as follows:

- 1) In the direction of travel, slide the product approximately in the middle between or over the perforated rails anchored on the bottom of the vehicle and secure it with the wheel lock.
 - 2) Check that the belt system in the vehicle matches the lashing points installed on the product.
 - 3) Attach the tie-down straps mounted on the bottom of the vehicle to the lashing points on the product (see Fig. 70; front: item 1, rear: item 2):
 - **Front:** engage the hooks in the lashing points from the inside to the outside (see Fig. 79, item 1).
 - **Rear:** engage the hooks in the lashing points from the inside or the outside (see Fig. 80, item 1).
 - 4) Firmly tighten the tie-down straps in accordance with the manufacturer's specifications. The maximum angle of the front tie-down-straps is **45°** (see Fig. 81).
 - 5) Fold the push handle down to provide additional stability for the seating system (see Fig. 81, item 1).
 - 6) **For "swiveling" front wheel option only:** turn the caster wheels to face forwards for improved stability (see Fig. 81).
 - 7) Remove all prohibited components from the product (see Page 32).
 - 8) Adjust the suspension to "firm" (see the Instructions for Use (User)).
- The tie-down straps are correctly attached (see Fig. 81).



7.9 Securing the User in the Product

⚠ WARNING

Lack of additional belt system on the Mygo seating system

Severe injury in case of accident

- ▶ Please note that the use of the Mygo seating system in a wheelchair accessible vehicle is permitted only when a chest and shoulder belt or a chest and shoulder support is used.

⚠ CAUTION

Incorrect belt settings

User fear, strangulation due to incorrectly positioned safety belts

- ▶ The safety belts must be positioned close to the body. They must never be held away from the user's body by components of the product or other parts.
- ▶ Firmly tighten the safety belts after they are applied. However, the comfort of the user must be taken into consideration:
 - Excessively tight adjustment may lead to unnecessary pain or suffering for the user.
 - Adjusting the belt straps too loosely can cause them to shift up or down so that the user can slide into a position that is dangerous to him/her.
 - The belt straps must not be twisted but must lie flat on the user's body.
- ▶ The shoulder belt must run over the user's shoulder. If this is not possible, the user must not be transported in this position in a motor vehicle.



INFORMATION

Please refer to the additional information in the brochure on using our products in wheelchair accessible vehicles, order number 646D158.

For transportation in a wheelchair accessible vehicle, the user must be secured in the rehab buggy with safety belts as follows:

- 1) Place the user in a nearly upright seated position.

- 2) Attach the lap belt to the anchor points (pins) provided:
 - Pull one end of the lap belt from each side of the seat through to the bottom (see Fig. 82, see Fig. 83).
 - Engage the ends of the lap belt on the anchor points (pins) provided (see Fig. 84, item 1).
 - **INFORMATION: The lap belt must be worn at an angle of 45° to 75° to the horizontal. When this is not possible in exceptional cases, an angle between 30° and 45° to the horizontal can also be chosen** (see Fig. 85).
- 3) Engage the shoulder belt in the lap belt (see Fig. 86).



8 Handover

8.1 Final inspection

Carry out a final inspection before the product is delivered:

- Are all options installed according to the order form?
- Do the wheels turn easily?
- Has the seat been correctly fitted?
- Is the padding placed properly?
- Is the footrest adjusted properly?
- Can the quick-release axles be easily removed?
- Have all screws been properly tightened?

8.2 Transport to the customer

The rehab buggy should be transported to the user in disassembled state using the outer packaging.

8.3 Delivery

The following steps must be performed for the safe delivery of the product:

- Test out the seat with the user. Pay special attention to proper positioning according to medical considerations (seating position, footrest position).
- The attendants must be instructed in safe use of the product. In particular, the enclosed Instructions for Use (User) are to be used.
- The Instructions for Use (User) must be given to the attendant upon delivery of the product.

9 Maintenance/Repair

The manufacturer recommends regular maintenance of the product every **12 months**.

More information on cleaning/disinfection and maintenance/repair can be found in the instructions for use (user).

The service manual contains detailed information on repairs.

10 Disposal

10.1 Disposal information

In case of disposal, all components of the product must be disposed of in accordance with the respective national environmental laws.

10.2 Information on re-use

⚠ CAUTION

Use of used seat upholstery

Functional and hygiene risks for the skin

- ▶ Replace the seat upholstery if the wheelchair is to be re-used.

The product is suitable for re-use.

Similar to second-hand machines or cars, products that are being re-used are subject to increased strain. The characteristics and specifications must not change to an extent that would endanger the safety of patients and, if applicable, third parties during the lifespan of the product.

Based on market observations and the current state of technology, the manufacturer has calculated that the product can be used for a period of **4 years**, provided that it is used properly and that the service and maintenance instructions are observed. Storage times at the dealer or with paying parties are not included in this period. It should be emphasised, however, that the product is reliable far beyond this defined period of time if it is cared for and maintained properly.

The product must first be thoroughly cleaned and disinfected before it can be re-used. Then the product must be examined by an authorised specialist to check the condition and to look for wear and possible damages. Any worn and damaged components as well as components that do not fit or are not suitable for the new user must be replaced.

Detailed information on the replacement of parts and details of tools required and prescribed service intervals can be found in the service manual.

11 Legal Information

11.1 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purpose described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

The manufacturer is not responsible for damages caused by components and spare parts not approved by the manufacturer. Repairs must be carried out exclusively by authorised dealers or by the manufacturer.

11.2 CE Conformity

This device meets the requirements of the 93/42/EEC guidelines for medical devices. This product has been classified as a class I device according to the classification criteria outlined in appendix IX of the directive. The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

11.3 Warranty Terms and Conditions

Further information on the warranty terms and conditions is available from the manufacturer's service (see inside back cover for addresses).

12 Appendices

12.1 Technical Data

Dimensions (cm) and weights (kg)

| Outdoor mobility base | Size 1 | Size 2 ¹⁾ |
|--|------------------------|------------------------|
| Overall width | 60 | 70 |
| Front wheel diameter, "swiveling" front wheels | 17 | 17 |
| Front wheel diameter, "fixed" front wheels | 28 | --- |
| Rear wheel diameter | 28 | 28 |
| Push handle height (min./max.) | 71/116 | 71/116 |
| Seat tilt (seating unit angle (min./max.)) | -35°/+10° or -10°/+35° | -35°/+10° or -10°/+35° |
| Max. load capacity ²⁾ | 55 | 55 |
| Max. load storage basket ²⁾ | 7 | 7 |
| Weight of mobility base | 10.5 | 11.5 |
| Folded size without seat (LxWxH) | 79x60x47 | 79x67x47 |
| Min. folded size with seat facing forwards; seat tilted backwards 10°, back angle 100° (LxWxH) | 90x60x53 ³⁾ | 97x67x59 |

| Multifunctional seating unit | Size 1 | Size 2 |
|---|--------------------------------|------------|
| Seat width | 20 – 30 | 26 – 35 |
| Seat depth with seat depth reducer (option) | 19 – 31 possible from 14 cm | 24 – 40 |
| Backrest height | 41 – 61 | 56 – 70 |
| Backrest angle | 80° – 180° | 80° – 180° |
| Height of shoulder strap slots (min./max.) | 30/45 | 30/45 |
| Width of shoulder strap slots | 5 | 5 |
| Lower leg length | 19 – 31 | 20 – 37 |
| Max. load capacity ²⁾ | 40 | 40 |
| Max. load capacity for use in a wheelchair accessible vehicle ⁴⁾ | 40 | 40 |
| Weight | 7 | 8.5 |
| Folded size with lateral padding (LxWxH) | 58x37x42 | 59x43x58 |

¹⁾ The outdoor mobility base with fixed front wheels in size 2 is available only as a custom fabrication.

²⁾ The load capacity differs between the outdoor mobility base and the seating unit:

- Seating unit (max. 40 kg) + own weight of the seating unit, size 1 (7.5 kg) = 47.5 kg
- Outdoor mobility base (max. 55 kg) – fully loaded seating unit (47.5 kg) = 7.5 kg for other accessories + weight in storage basket

³⁾ Folded dimensions of the outdoor mobility base with fixed front wheels: 97x60x53 cm.

⁴⁾ Observe the information in the section "Use in a Wheelchair Accessible Vehicle".

Kundenservice/Customer Service

Europe

Otto Bock HealthCare Deutschland GmbH
Max-Näder-Str. 15 · 37115 Duderstadt · Germany
T +49 5527 848-3433 · F +49 5527 848-1460
healthcare@ottobock.de · www.ottobock.de

Otto Bock Healthcare Products GmbH
Kaiserstraße 39 · 1070 Wien · Austria
T +43 1 5269548 · F +43 1 5267985
vertrieb.austria@ottobock.com · www.ottobock.at

Otto Bock Adria Sarajevo D.O.O.
Omladinskih radnih brigada 5
71000 Sarajevo · Bosnia-Herzegovina
T +387 33 766200 · F +387 33 766201
obadria@bih.net.ba · www.ottobockadria.com.ba

Otto Bock Bulgaria Ltd.
41 Tzar Boris III Blvd. · 1612 Sofia · Bulgaria
T +359 2 80 57 980 · F +359 2 80 57 982
info@ottobock.bg · www.ottobock.bg

Otto Bock Suisse AG
Pilatusstrasse 2 · CH-6036 Dierikon
T +41 41 455 61 71 · F +41 41 455 61 70
suisse@ottobock.com · www.ottobock.ch

Otto Bock ČR s.r.o.
Protetická 460 · 33008 Zruč-Senec · Czech Republic
T +420 377825044 · F +420 377825036
email@ottobock.cz · www.ottobock.cz

Otto Bock Iberica S.A.
C/Majada, 1 · 28760 Tres Cantos (Madrid) · Spain
T +34 91 8063000 · F +34 91 8060415
info@ottobock.es · www.ottobock.es

Otto Bock France SNC
4 rue de la Réunion · CS 90011
91978 Courtaboeuf Cedex · France
T +33 1 69188830 · F +33 1 69071802
information@ottobock.fr · www.ottobock.fr

Otto Bock Healthcare plc
32, Parsonage Road · Englefield Green
Egham, Surrey TW20 0LD · United Kingdom
T +44 1784 744900 · F +44 1784 744901
bockuk@ottobock.com · www.ottobock.co.uk

Otto Bock Hungária Kft.
Tatai út 74. · 1135 Budapest · Hungary
T +36 1 4511020 · F +36 1 4511021
info@ottobock.hu · www.ottobock.hu

Otto Bock Adria d.o.o.
Dr. Franje Tuđmana 14 · 10431 Sveta Nedelja · Croatia
T +385 1 3361 544 · F +385 1 3365 986
ottobockadria@ottobock.hr · www.ottobock.hr

Otto Bock Italia Srl Us
Via Filippo Turati 5/7 · 40054 Budrio (BO) · Italy
T +39 051 692-4711 · F +39 051 692-4720
info.italia@ottobock.com · www.ottobock.it

Otto Bock Benelux B.V.
Ekkersrijt 1412 · 5692 AK
Son en Breugel · The Netherlands
T +31 499 474585 · F +31 499 476250
info.benelux@ottobock.com · www.ottobock.nl

Industria Ortopédica Otto Bock Unip. Lda.
Av. Miguel Bombarda, 21 - 2º Esq.
1050-161 Lisboa · Portugal
T +351 21 3535587 · F +351 21 3535590
ottobockportugal@mail.telepac.pt

Otto Bock Polska Sp. z o. o.
Ulica Korolowa 3 · 61-029 Poznań · Poland
T +48 61 6538250 · F +48 61 6538031
ottobock@ottobock.pl · www.ottobock.pl

Otto Bock Romania srl
Șos de Centura Chitila - Mogoșoia Nr. 3
077405 Chitila, Jud. Ilfov · Romania
T +40 21 4363110 · F +40 21 4363023
info@ottobock.ro · www.ottobock.ro

OOO Otto Bock Service
p/o Pultikovo, Business Park „Greenwood”,
Building 7, 69 km MKAD
143441 Moscow Region/Krasnogorskiy Rayon
Russian Federation
T +7 495 564 8360 · F +7 495 564 8363
info@ottobock.ru · www.ottobock.ru

Otto Bock Scandinavia AB
Koppargatan 3 · Box 623 · 60114 Norrköping · Sweden
T +46 11 280600 · F +46 11 312005
info@ottobock.se · www.ottobock.se

Otto Bock Slovakia s.r.o.
Röntgenova 26 · 851 01 Bratislava 5 · Slovak Republic
T +421 2 32 78 20 70 · F +421 2 32 78 20 89
info@ottobock.sk · www.ottobock.sk

Otto Bock Sava d.o.o.
Industrijska bb · 34000 Kragujevac · Republika Srbija
T +381 34 351 671 · F +381 34 351 671
info@ottobock.rs · www.ottobock.rs

Otto Bock Ortopedi ve
Rehabilitasyon Tekniği Ltd. Şti.
Ali Dursun Bey Caddesi · Lati Lokum Sokak
Meriç Sitesi B Blok No: 6/1
34387 Mecidiyeköy-Istanbul · Turkey
T +90 212 3565040 · F +90 212 3566688
info@ottobock.com.tr · www.ottobock.com.tr

Africa

Otto Bock Algérie E.U.R.L.
32, rue Ahcène Outaleb - Coopérative les Mimosas
Mackle-Ben Aknoun · Alger · DZ Algérie
T +213 21 913863 · F +213 21 913863
information@ottobock.fr · www.ottobock.fr

Otto Bock Egypt S.A.E.
28 Soliman Abaza St. Mohandessein - Giza · Egypt
T +202 330 24 390 · F +202 330 24 380
info@ottobock.com.eg · www.ottobock.com.eg

Otto Bock South Africa (Pty) Ltd
Building 3 Thornhill Office Park · 94 Bekker Road
Midrand · Johannesburg · South Africa
T +27 11 312 1255
info-southafrica@ottobock.co.za
www.ottobock.co.za

Americas

Otto Bock Argentina S.A.
Av. Belgrano 1477 · CP 1093
Ciudad Autónoma de Buenos Aires · Argentina
T +54 11 5032-8201 / 5032-8202
atencionclientes@ottobock.com.ar
www.ottobock.com.ar

Otto Bock do Brasil Tecnica Ortopédica Ltda
Alameda Maria Tereza 4036 - Bairro Dois Córregos
13278-181 - Valinhos-São Paulo · Brasil
T +55 19 3729 3500 · F +55 19 3269 6061
ottobock@ottobock.com.br · www.ottobock.com.br

Otto Bock HealthCare Canada
5470 Harvester Road
Burlington, Ontario, L7L 5N5, Canada
T +1 289 288-4848 · F +1 289 288-4837
infocanada@ottobock.com · www.ottobock.ca

Otto Bock HealthCare Andina Ltda.
Calle 138 No 53-38 · Bogotá · Colombia
T +57 1 8619988 · F +57 1 8619977
info@ottobock.com.co · www.ottobock.com.co

Otto Bock de Mexico S.A. de C.V.
Prolongación Calle 18 No. 178-A
Col. San Pedro de los Pinos
C.P. 01180 México, D.F. · Mexico
T +52 55 5575 0290 · F +52 55 5575 0234
info@ottobock.com.mx · www.ottobock.com.mx

Otto Bock HealthCare
Two Carlson Parkway North, Suite 100
Minneapolis, MN 55447 · USA
T +1 763 553 9464 · F +1 763 519 6153
usa.customerservice@ottobockus.com
www.ottobockus.com

Asia/Pacific

Otto Bock Australia Pty. Ltd.
Suite 1.01, Century Corporate Centre
62 Norwest Boulevard
Baulkham Hills NSW 2153 · Australia
T +61 2 8818 2800 · F +61 2 8814 4500
healthcare@ottobock.com.au · www.ottobock.com.au

Beijing Otto Bock Orthopaedic Industries Co., Ltd.
B12E, Universal Business Park
10 Jiuxianqiao Road, Chao Yang District
Beijing, 100015, P.R. China
T +8610 8598 6880 · F +8610 8598 0040
news-service@ottobock.com.cn
www.ottobock.com.cn

Otto Bock Asia Pacific Ltd.
Unit 1004, 10/F, Greenfield Tower, Concordia Plaza
1 Science Museum Road, Tsim Sha Tsui
Kowloon, Hong Kong · China
T +852 2598 9772 · F +852 2598 7886
info@ottobock.com.hk · www.ottobock.com

Otto Bock HealthCare India
Behind FairLawn Housing Society
St. Gregorios Lane, Sion Trombay Road
Chembur, Mumbai, 400071 · India
T +91 22 2520 1268 · F +91 22 2520 1267
information@indiaottobock.com · www.ottobock.in

Otto Bock Japan K. K.
Yokogawa Building 8F, 4-4-44 Shibaura
Minato-ku, Tokyo, 108-0023 · Japan
T +81 3 3798-2111 · F +81 3 3798-2112
ottobock@ottobock.co.jp · www.ottobock.co.jp

Otto Bock Korea HealthCare Inc.
4F Agaworld Building · 1357-74, Seocho-dong
Seocho-ku, 137-070 Seoul · Korea
T +82 2 577-3831 · F +82 2 577-3828
info@ottobockkorea.com · www.ottobockkorea.com

Otto Bock South East Asia Co., Ltd.
1741 Phaholyothin Road
Kwaeng Chatuchark · Khet Chatuchark
Bangkok 10900 · Thailand
T +66 2 930 3030 · F +66 2 930 3311
obsea@ottobock.co.th · www.ottobock.co.th

Other countries

Otto Bock HealthCare GmbH
Max-Näder-Straße 15 · 37115 Duderstadt · Germany
T +49 5527 848-1590 · F +49 5527 848-1676
reha-export@ottobock.de · www.ottobock.com

Ihr Fachhändler | Your specialist dealer



Versandanschrift für Rücksendungen/Adress for Returns:
Otto Bock Manufacturing Königsee GmbH
Lindenstraße 13 · 07426 Königsee-Rottenbach/Germany



Otto Bock Mobility Solutions GmbH
Lindenstraße 13 · 07426 Königsee-Rottenbach/Germany
www.ottobock.com

Ottobock has a certified Quality Management System in accordance with ISO 13485.