

CE



Illustrative photo



mouse

**Operation and maintenance manual**

Edition 8 date 2023-06-21

**Special stroller Mouse**

A modern stroller with excellent driving characteristics

**This is a medical device.  
Use it according to the instruction manual or label.**

ATTENTION! Users with visual, reading or cognitive impairment are requested to report to the Manufacturer.

## Purpose of the Manual

The operation and maintenance manual contains basic information necessary for the correct use of the Mouse stroller, including the principles of safe use, description of its construction, operation and adjustment, preparation of the stroller for operation, how to maintain and clean the stroller, and warranty conditions.

It is designed for people who directly care of children with disabilities (childcare assistants, parents, guardians), which will be transported in a stroller. It is also intended for doctors and physiotherapists operating the stroller.

ATTENTION! Before using the stroller, each User is obliged to read and act in accordance with its guidelines

This manual is an integral part of the product and must always be attached to the product wherever being sold. **Keep this manual for future reference!** Keep the manual in an easily accessible place. The manual in the electronic version is available on the website [www.mywam.eu](http://www.mywam.eu).

ATTENTION! In the event of any medical incidents related to the presence of the Mouse stroller, the Manufacturer must be immediately notified in writing by contacting it via registered letter or e-mail with acknowledgment of receipt.

Manufacturer's contact details:

- Postal address: MyWam sp.j., ul. Szczecińska 10, 41-516 Chorzów, Poland
- contact phone number: +48 32 733 11 31
- e-mail address: [mywam@mywam.pl](mailto:mywam@mywam.pl)



User and maintenance manual version: EN **Edition: 8 Dated 2023-06-21**

MyWam Kupiec, Bartold, Angres sp.j. reserves the right to make technical and commercial changes to the content of the Operation and Maintenance Manual without prior notice. Each change will be determined by the date the instructions were last updated.

**ATTENTION! The photos contained in the Instructions for Use are indicative and may differ from the purchased product.**

## Table of contents

<b>1.</b>	<b>PURPOSE AND INDICATIONS FOR THE USE OF THE MOUSE STROLLER</b>	<b>4</b>
<b>2.</b>	<b>BASIC INFORMATION ABOUT THE MOUSE STROLLER</b>	<b>4</b>
<b>3.</b>	<b>QUICK START</b>	<b>5</b>
<b>4.</b>	<b>IMPORTANT INSTRUCTIONS AND SAFETY RULES</b>	<b>6</b>
<b>5.</b>	<b>NAMEPLATE</b>	<b>7</b>
<b>6.</b>	<b>WARNING LABELS</b>	<b>7</b>
<b>7.</b>	<b>SYMBOLS DESIGNATION</b>	<b>7</b>
<b>8.</b>	<b>GENERAL CONSTRUCTION OF THE STROLLER</b>	<b>8</b>
<b>9.</b>	<b>CONFIGURATION AND USE – STROLLER FRAME</b>	<b>9</b>
9.1	<i>Frame folding mechanism</i>	9
9.2	<i>Unfolding the stroller frame. Transport lock.</i>	9
9.3	<i>Folding the frame</i>	9
9.4	<i>Swivel front wheels, solid. Locking the direction of travel of the front wheels</i>	10
9.5	<i>Rear wheels, solid. Installation and removal of rear wheels on the axle</i>	10
9.6	<i>Rear wheel amortization system</i>	11
9.7	<i>Ergonomic handle for pushing the stroller. Mechanism and angle adjustment of the handle</i>	11
9.8	<i>Central brake. Applying and deactivating the brake</i>	11
9.9	<i>Tilting foot</i>	12
9.10	<i>The transport handle</i>	12
9.11	<i>Footrest. Adjustment and use.</i>	12
<b>10.</b>	<b>ASSEMBLY AND DISASSEMBLY OF THE SEAT TO THE FRAME</b>	<b>12</b>
10.1.	<i>Assembly and disassembly of the seat to the frame</i>	12
10.2.	<i>Adjustable seat, „reclining back”</i>	13
10.3.	<i>Adjustable backrest</i>	13
<b>11.</b>	<b>BASIC EQUIPMENT</b>	<b>13</b>
11.1.	<i>Five-point belts</i>	13
11.2.	<i>Safety barrier</i>	14
11.3.	<i>Adjustable canopy with window</i>	14
11.4.	<i>Shopping basket</i>	14
11.5.	<i>Abduction and stabilizing belts</i>	15
11.6.	<i>Wedge</i>	15
11.7.	<i>Lininig pads</i>	15
11.8.	<i>Upholstery</i>	15
<b>12.</b>	<b>ACCESSORIES</b>	<b>15</b>
12.1.	<i>Adjustable head supports - ↔ and ↕</i>	15
12.2.	<i>Adjustable chest and loin pads - ↔ and ↕</i>	15
12.3.	<i>Therapeutic table</i>	16
12.4.	<i>Seat minimizing kit</i>	16
12.5.	<i>Travel bag</i>	17
12.6.	<i>Sun umbrella</i>	17
12.7.	<i>Winter gloves</i>	17
12.8.	<i>Winter sleeping bag</i>	17
12.9.	<i>Mosquito net</i>	17
12.10.	<i>Raincover</i>	17
<b>13.</b>	<b>TECHNICAL SPECIFICATION OF THE MOUSE STROLLER</b>	<b>18</b>
13.1.	<i>Material composition of the stroller</i>	18

---

<b>14.</b>	<b>DELIVERY AND TOOLS</b>	<b>18</b>
<b>15.</b>	<b>ASSEMBLY AND DISASSEMBLY OF THE STROLLER</b>	<b>19</b>
<b>16.</b>	<b>PREPARATION FOR USE</b>	<b>19</b>
<b>17.</b>	<b>USING THE STROLLER</b>	<b>20</b>
17.1.	<i>Carrying the stroller</i>	20
17.2.	<i>Getting on and off the stroller</i>	20
17.3.	<i>Correct position in a stroller</i>	20
17.4.	<i>Riding in a stroller</i>	21
17.5.	<i>Moving on slopes</i>	21
17.6.	<i>Overcoming thresholds or curbs</i>	21
17.7.	<i>Transport in a car and other vehicles</i>	21
17.8.	<i>Cleaning and care</i>	21
<b>18.</b>	<b>STORAGE OF THE STROLLER</b>	<b>21</b>
18.1.	<i>Preparation of the stroller for long-term storage</i>	22
18.2.	<i>Corrosion protection</i>	22
<b>19.</b>	<b>RE-USE</b>	<b>22</b>
<b>20.</b>	<b>SERVICE AND PERIODIC INSPECTIONS</b>	<b>23</b>
<b>21.</b>	<b>DISPOSAL</b>	<b>24</b>
<b>22.</b>	<b>COMPLIANCE WITH REQUIREMENTS FOR MEDICAL DEVICES</b>	<b>24</b>
<b>23.</b>	<b>NOTES</b>	<b>25</b>
<b>24.</b>	<b>LIST OF AUTHORIZED SERVICE POINTS</b>	<b>26</b>
<b>25.</b>	<b>WARRANTY CARD</b>	<b>27</b>

## Introduction

Thank you for purchasing **Mouse** special stroller for young people. We hope that it will be the cause of many wonderful moments with your child and will be helpful every day. We make every effort to ensure that our products meet your expectations and the highest standards. We encourage you to contact our employees and follow us on our website [www.mywam.eu](http://www.mywam.eu) and on social media. There you can find out many interesting things about our products.

### 1. Purpose and indications for the use of the Mouse stroller

"Mouse" is a multifunctional stroller for disabled children. The stroller is based on a frame made of aluminum and steel elements. It allows for quick removal of wheels, headrest, etc. The product is equipped with a number of accessories improving the functionality and comfort of use, such as: table, safety belts, wedge, hood, mosquito net, sleeping bag, shopping bag, etc. Users cannot drive the stroller alone, but only with the help of a guardian. The stroller can be used indoors and outdoors, on various surfaces (asphalt, concrete, stone, gravel) and in various weather conditions. Any other use of the product than described is forbidden. The stroller must be used together with the instructions for use only by adults caring for the user. Young people and children as well as any other persons without proper training should not use the product.

The "Mouse" stroller is a medical device class I indicated for use by disabled children who, due to their health condition, cannot walk or sit down on their own. The product is intended for the transport and movement of the above persons in a sitting or lying position and for posture correction and this is its intended use.

Clinical indications include cerebral palsy, meningeal hernia, dystrophies, paresis, etc.

The product is intended for users aged 2 to 6 years. The product should be used in accordance with the recommendations of the attending physician or physiotherapist.

#### **Contraindications:**

In case of contraindication to the sitting position and/or passive upright positioning, the adjustment of the stroller should be consulted with a physiotherapist or attending physician.

### 2. Basic information about the Mouse stroller

The Mouse youth stroller has many options, e.g.:

- |   |                          |
|---|--------------------------|
| • backrest angle adjustment                   | see point 10.3           |
| • backrest and seat angle adjustment "Cradle" | see point 10.2           |
| • seat depth adjustment - insert              | see point 12.4           |
| • footrest adjustment                         | see point 9.11           |
| • torso stabilization                         | see point 12.11 and 12.2 |
| • optional equipment                          | see point 12             |

### 3. Quick start

To quickly use the stroller it must be unpacked and prepared for use according to item 15 Assembly and disassembly of the stroller and item 16 Preparation for use. The following is a pictorial diagram, but it will not replace careful reading of the complete manual.





---

## Safety

### 4. Important instructions and safety rules

- It is forbidden to use the stroller in a way for which the stroller is not dedicated.
- The product contains small parts that can be dangerous if swallowed by a child!
- Due to the large adjustment possibilities of the stroller, during the angular adjustment of the seat and angular adjustment of the backrest, it is possible to set the position in which the child's head is below the pelvic line. It is forbidden to use this position as it could lead to the risk of the child choking.
- It is forbidden to put additional weight on the stroller, especially when it comes to placing an additional weight on the back of the backrest, e.g. a backpack or a bag. This may cause the Mouse stroller to tip over and injure the user! It does not apply to original additional accessories dedicated to the Mouse stroller that do not exceed the permissible load. All additional elements attached to the stroller affect its stability.
- ATTENTION! Depending on the adjustments made for the Mouse stroller - it concerns the setting of the stroller's push handle, backrest height adjustment, backrest angular adjustment, seat depth adjustment, seat angle adjustment, footrest angle adjustment and footrest adjustment - this stroller may exceed the recommended dimensions.
- ATTENTION! Each time the stroller is used, the steps listed below must be performed and the stroller must not be used in the event of any irregularities. In such a situation, it is recommended to urgently contact the Seller and/or the Manufacturer for advice.
- ATTENTION! When driving through doors or similar restrictions restricting the direction of travel of the stroller, make sure that there is adequate space on both sides of the stroller to prevent pinching of hands or other body parts and to prevent damage to the stroller.
- ATTENTION! The stroller's temperature may increase when the stroller is exposed to sunlight.

---

#### ATTENTION!



- Never leave your child unattended in a stroller, as there is a risk of injury and there is a risk of accident!
- This product is not suitable for running, skating or roller-skating
- This stroller is not intended as a seat in a motor vehicle
- Beware of gaps and moving parts – risk of jamming
- Do not carry the stroller with the child inside it
- The stroller is packed in a plastic bag. The plastic bag is not a toy and should be kept away from infants and young children as it may cause suffocation when put on a head.

---

ATTENTION! It is forbidden to use the stroller if the technical condition of the stroller or a single element of the stroller is bad or the user has any doubts about it, the stroller is damaged or there is a suspicion that damage has occurred or if the operation of the stroller is not correct. In such a situation, it is recommended to urgently contact the Seller and/or the Manufacturer for advice.

## 5. Nameplate

- The nameplate contains basic information about the product: name, model, serial number (ID no.), maximum load, manufacture date, CE mark and also certain warning in accordance with the Operation and Maintenance Manual, the medical device mark and manufacturer's data (see Fig. 5.1).
- The nameplate is located on the crossbeam at the rear of the stroller frame, on its left side, glued from above
- Plate with UDI code  
Plate with UDI code includes the following information: name, model, UDI code, details of manufacturer (see Fig. 5.2)



Fig. 5.1



Fig. 5.2

## 6. Warning labels

- The warning label contains important warnings for the user (see Fig. 6.1).
- The warning label is located on the rear/side part of the frame on both sides.












Fig. 6.1

## 7. Symbols designation

	Product serial number		The manufacturer's name
	Permitted total user weight		The product is a medical device
	Date of manufacture		DESCRIPTION " <b>WARNING!</b> " and the symbol in this manual are used to enhance the reader's attention to the content marked with this sign. Non-compliance with content marked with this description may endanger the health and safety of the user and accompanying persons, and may also damage the product.
	Conformity assessment mark according to applicable European regulations and directives		
	Warning about proceeding in accordance with the product Operation and Maintenance Manual		



## 8. General construction of the stroller

The specialist stroller Mouse consists of two main parts: frame and therapeutic seat with backrest. Mouse offers a wide range of accessories that help to meet individual needs of the users.

The steal frame (Fig.4.1. ) with aluminium elements includes:

Adjustable, ergonomic handlebar	(Fig.8.1.-nr 2)
Front swivel wheels	(Fig.8.1.-nr12)
Front wheel direction lock	(Fig.8.1.-nr15)
Front wheel fork shock absorption system	(Fig.8.1.-nr14)
Front wheel folding system	(Fig.8.1.-nr13)
Rear wheels	(Fig.8.1.-nr 8)
2-step rear wheel shock absorption system	(Fig.8.1.-nr24)
Rear wheel folding	(Fig.8.1.-nr23)
Central breake	(Fig.8.1.-nr 9)
Shopping basket	(Fig.8.1.-nr11)
hooks fixing seat to frame	(Fig.8.1.-nr 6)
Transport lock	(Fig.8.1.-nr 3)
Frame folding mechanism (upper and lower)	(Fig.4.1.-nr 4)

Reclining seat includes:

The seat is made of synthetic material

Side pad / reducer	(Fig.8.1.-nr18)
Adjustable backrest	(Fig.8.1.-nr5)
Safety barrier fixing	(Fig.8.1.-nr7)
Safety barrier	(Fig.8.1.-nr19)
5-point safty belts	(Fig.8.1.-nr20)
Wedge	(Fig.8.1.-nr17)
Adjustable head pads	(Fig.8.1.-nr21)
Adjustable canopy with window	(Fig.8.1.-nr1)
Adjustable and folding footrest	(Fig.8.1.-nr16)
footrest extension	(Fig.8.1.-nr22)
Optional inclination bar	(Fig.8.1.-nr10)



Fig. 8.1  
Main construction of the stroller

## 9. Configuration and use – stroller frame

### 9.1 Frame folding mechanism

The frame folding mechanism is located inside the frame and on the frame in two places in its right part and in three places in its left part. Describing these places from above, from the side of the handle to push the stroller, these are:

- right and left bend - an element that allows folding the upper frame to reduce its dimensions
- mechanism of folding the frame with a lock

### 9.2 Unfolding the stroller frame. Transport lock.

To unfold (open) the stroller frame: lay the stroller frame horizontally on the ground and pull off the transport lock (Fig. 9.2.1). Then gently pull the backrest to the upright position (Fig. 9.2.2.). To set the stroller, pull the handle carefully by lifting it slightly up while bringing the stroller to the desired position until the frame folding lock engages. The frame of the stroller will not unfold properly until the frame folding mechanism has been properly closed. This is evidenced by the clear "click" sound coming from the well-closed mechanism of folding the frame.




Fig. 9.2.1  
Transport lock



Fig. 9.1.  
Folding points



Fig. 9.2.2  
Unfolding the frame

### 9.3 Folding the frame

- Description of the lower folding mechanism of the frame.

To fold the stroller frame: Lift the lever below the push handle on the right side of the frame and simultaneously pull up both handles on the top of the folding mechanism. Then fold the stroller frame and attach the transport lock. Photo 9.3.1.

- Description of the upper folding mechanism of the frame.

To fold the stroller frame: Lift the lever below the push handle on the right side of the frame and simultaneously pull up both handles on the top of the folding mechanism. Photo 9.3.2. Then fold the stroller frame and attach the transport lock.



Fig. 9.3.1  
Lower mechanism



Fig. 9.3.2  
Upper mechanism



Fig 9.4.1  
Quick release button



Fig 9.4.2  
Direction lock knob

#### 9.4 Swivel front wheels, solid. Locking the direction of travel of the front wheels

The used, durable wheels guarantee the highest quality of stroller driving, driving pleasure and elegant appearance. The tires are tubeless, filled with a special gel that has shock absorbing properties very similar to those of inflatable wheels. The wheels are factory-fitted to the front wheel forks, which are finished with top steering axles from above.

Both front wheels can be set for straight-ahead or unlocked for turning, which makes it much easier to maneuver the stroller. Thanks to the locking of the direction of travel of the front wheels and the possibility of putting it on or taking off, there is clearly greater comfort for both the child and the assistant when driving the stroller.

It is recommended to remove the lock and allow the front wheels to turn while driving the stroller on a completely smooth surface, e.g. in shopping centers, school.

- Assembly. To attach the wheels, insert the wheel axle into the folding system slot until you hear a 'click'. The audible clicking sound is the indication that the assembly has been done properly. For dismantling the front wheels press the button on the front wheel folding system (Fig. 9.4.1 ) and eject the wheel axle from the front wheel folding system slot.

- Direction lock. In the event of uneven surfaces, moving on sloping surfaces, when crossing thresholds or curbs, it is recommended to install the lock and block the turning of the front wheels. To block the front wheels, turn the lock knob left. Turn the lock knob right, to unlock. (Fig.9.4.2)

#### 9.5 Rear wheels, solid. Installation and removal of rear wheels on the axle

The used, durable wheels guarantee the highest quality of stroller driving, driving pleasure and elegant appearance. The tires are tubeless, filled with a special gel that has shock absorbing properties very similar to those of inflatable wheels.

Thanks to the possibility of dismounting both the wheels and their axle, the stroller is smaller after folding, which facilitates its transport and storage.

The wheel assembly mechanism is mounted inside the rear wheels and ensures quick assembly and disassembly of the rear wheels. It is started by pressing the quick assembly clamp (Fig. 9.6.1) towards the center of the wheel, after which the rear wheel can be mounted on the rear wheel axle or the rear wheel can be removed from the rear wheel axle.

- Assembly. When installing the rear wheel on the rear wheel axle, press the clamp in the center of the wheel and put the wheel on the beginning of its axle. Then stop pressing the clamp and continue sliding the wheel onto its axle until it clicks.

The correct wheel placement is demonstrated by the clear "click" sound coming from the quick axle rear axle mounting mechanism.

- Dismantling the wheels from the axle. In order to remove the rear wheel from the rear wheel axle, press the clamp in the middle of the wheel and slide the wheel out of its axle.

– Fig 9.5.2.



Fig 9.5.1  
Quick mount buckle



Fig 9.5.2  
Wheel disassemble

### 9.6 Rear wheel amortization system

The stroller is equipped with spring-loaded cushioning of the frame and rear wheels, located at the ends of the rear legs of the frame, when in contact with the crossbeam. Amortization works automatically. The adjustment is done by activating the up-down lever located from the bottom of the shock absorber. In order to facilitate the switchover of settings by the lever, simultaneous pressure on the shock absorber spring is allowed. When the lever is shifted, the shock absorber spring will be compressed (option for hard ground) or extended (soft ground).




Fig 9.6  
Wheel amortization system

### 9.7 Ergonomic handle for pushing the stroller. Mechanism and angle adjustment of the handle

The handle is used to push the stroller, providing a firm grip in various driving conditions. It has the ability to adjust the height of the handle for the child's assistant who drives the stroller - Fig. 9.7

- Adjustment mechanism. These are two buttons on the outside of the frame, located on the joints of the handle to push the stroller.
- Angle adjustment of the stroller push handle. To change the angle of the stroller push handle, push both buttons of the angle adjustment mechanism of the push handle and change its angle. After setting the handle in the correct position, stop pressing the buttons, which will automatically lock the push handle in the desired position.



Fig 9.7  
Handle angular adjustment

### 9.8 Central brake. Applying and deactivating the brake

The central brake pedal (Fig. 9.8), which when pressed blocks the two rear wheels of the stroller directly, is located in the middle of the crossbeam, at the rear of the stroller frame.

- Central brake on and off

To apply the central brake, press the brake pedal down towards the ground until it clicks. The correct activation of the brake is indicated by a clear "click" sound from the central brake. To disengage the central brake, pull the brake pedal up.



Fig 9.8  
Central brake



Fig 9.9  
Tilting foot



### 9.9 Tilting foot

The tilting foot (Fig. 9.9) is located on the right side of the crossbeam, at the rear of the stroller frame. In order to use it, press the tilting foot and at the same time pull the handle to push the stroller, which in turn will cause the front of the stroller to rise up together with the front wheels. Using the tilting foot makes it easier for the stroller to climb onto small steps or curbs.



Fig 9.10  
Transport handle

### 9.10 The transport handle

The transport handle is located on the shopping basket upper fastening (Fig. 9.10)

### 9.11 Footrest. Adjustment and use.

- W In order to change the inclination of the footrest, press the buttons marked on the picture on both sides of the footrest, then set the appropriate position and release the buttons (Fig. 9.11.1)
- In order to change the height of the footrest, unscrew the bolts fixing the footrest and pull the footrest from the angle controllers. In place of the footrest, 2 tubes should be installed and bolts tightened. The footrest should be inserted into the tube so that an additional 3 heights of the footrest are obtained (Fig. 9.11.2). The footrest is equipped with straps that fasten the child's feet. It is recommended to put them on while using the stroller.



Fig. 9.11.1  
Tilt buttons



Fig.. 9.11.2.  
Screws and extensions

## 10. Assembly and disassembly of the seat to the frame

### 10.1. Assembly and disassembly of the seat to the frame

The two frame hooks that are attached to the seat are very important parts of the stroller. They serve the purpose of attaching the seat to the stroller's frame. The proper use of the hooks is essential for the safety of the child as well as that of other people. Make sure there aren't any cracks and damage in the stroller's frame. Also, check if it works properly. Do not use the stroller in case of any suspicion as to the proper operation of the stroller. Keep the safety of your child in mind at all times.

• **Seat assembly.** The seat can face the front or the back. To install the seat, place it on the stroller's frame hooks. Press until you hear a 'click' twice. Fig 10.1.1.

• **Dismantling** the seat. To dismantle the seat, press two colorful buttons situated on the lower part of the stroller's frame hooks. Next, eject the seat from the frame hooks. Fig 10.1.2.



Fig. 10.1.1  
Hooks for mounting the seat



Fig. 10.1.2  
Seat release button

### 10.2. Adjustable seat, „reclining back“

The stroller has a seat that can be adjusted to four different positions including a bassinet position. To change the position, pull back the steel handle at the back of the seat and slowly change the position of the seat. After securing the desired position, let go of the handle.

Caution! Each change of the seat position alters the general center of gravity, which impacts the stability of the stroller.

**Caution! Never leave your child in the stroller while making adjustments! It can put your child in danger.**



Fig. 10.2  
Adjustable seat

### 10.3. Adjustable backrest

The stroller's backrest can be adjusted to various positions varying from a seating position to a lying position. To change the position, move the backrest towards a seating position first. Then, pull the handle located at the back of the backrest and change the position slowly. After reaching the desired position, pull the handle down.



Fig 10.3  
Reclining back handle

**Caution!** Each change of the backrest position requires an adjustment of the 5-point belts.

**Caution!** Never leave your child in the stroller while making adjustments! It can put your child in danger.

**Caution!** Adjustments must be made with the brake engaged.

## 11. Basic equipment

### 11.1. Five-point belts

The Mouse stroller is standard equipped with professional five-point seat belts

- Five-point belts adjustment. Adjust the five-point belts to your child before attaching them. To adjust the length of the belts, move the adjusters on the braces up or down and move the protectors accordingly. To adjust the height of the braces to be attached to the backrest, move the braces from low to high, or vice versa.

- Use. In order to put on the child five-point belts, put the braces on the right and left of the child and then attach both braces to the belt buckle, which is fixed on the seat. In order to remove five-point belts from a child, press the button on the buckles and remove both braces from the buckle and slowly remove these braces from the right and left side of the child.

Attention! Each change in the angle of the backrest requires a readjustment of the five-point belt length.



Fig 11.1.  
Five-point belts



### 11.2. Safety barrier

The Mouse stroller can be equipped with a safety barrier that prevents the sitting child from falling forward in certain indications..

- The safety barrier is mounted with the same fixings as the therapeutic table, so you cannot use both of these elements at the same time. To fastening the safety barrier to the frame slide it into the attachment tubes and listen for the 'click'. The clearly audible 'click' is an indication that the barrier has been securely attached to the seat. To disassemble the barrier pull both barrier latches on the frame and pull out the barrier at the same time.



Fig 11.2  
Safety barrier

### 11.3. Adjustable canopy with window



Fig 11.3  
Adjustable canopy

The Mouse stroller can be equipped with an adjustable canopy with a window that helps protect the child from the sun or rain. The elements of the canopy are: canopy frame (front and rear arch), canopy opening adjustment mechanism (right and left), two snaps with red buttons fastening the canopy to the frame, canopy upholstery, and a mosquito net with a zipper.

- Assembly. The fastening of the adjustable canopy with the window are two plastic elements located on the outer sides of the upper part of the backrest. These attachments are used to mount the canopy with the window.

To install the adjustable canopy with the window, place the backrest upright, and then slide both canopy latches into the backrest mounts in parallel. Mount the canopy in such a way that the upholstery at the

back/bottom of the canopy overlaps the back of the stroller. In addition, make sure that the canopy adjustment element is secured with upholstery inside. This can be done by stretching the upholstery and placing the loop on the button in the corner of the canopy.

- Disassembly. To disassemble the adjustable canopy with window, press the two red buttons on the roof fixing latches and slide the latches out of the backrest mounts.
- Use. To change the size of the canopy, hold the front arch and the rear arch of the canopy frame and bring the arches closer or further apart to set the desired position. The canopy surface adjustment mechanism latches automatically. A special feature of the canopy is a window with a mosquito net for ventilation. The window opens and closes with a zipper.

### 11.4. Shopping basket

The Mouse stroller can be equipped with a shopping basket which is used to transport small, ordinary items of everyday life, resistant to shocks, impacts, bruises, soaking, flooding and dirt.

It is mounted at the bottom of the stroller by wrapping it on the frame of the stroller and fastening it with a zipper. The maximum load of the shopping basket is 3 kg.



Fig 11.4  
Shopping basket

### 11.5. Abduction and stabilizing belts



Fig 11.5

Abduction and stabilizing belts

Abduction and stabilizing belts are used for additional pelvic stability and prevent the child from slipping out of the stroller. The correct way to use them is to fasten the child's thighs so that his pelvis adheres as close as possible to the backrest. This will ensure a healthier and more comfortable position of the child in the stroller.

- **Adjustment.** To adjust the length of the Abduction and stabilizing belts, you must move the adjusters on the front belt straps up or down.
- **Use.** In order to put the child on the Abduction and stabilizing belts, put the child on the installed belts, then wrap the child's leg and fasten both belt buckles to both front straps. In order to remove the Abduction and stabilizing belts from the child, unfasten both belt buckles from the front straps.

### 11.6. Wedge

The Mouse stroller can be equipped with a wedge which in the case of specific indications provides discharge (inclination to the sides) of the child's lower limbs from his body. The wedge has a two-stage adjustment level. Fig.11.6



Fig 11.6

Wedge

### 11.7. Lining pads



Fig 11.7

Lining pads

The Mouse stroller is equipped with lining pads, securing the velcro on the stroller's backrest, if no additional equipment is installed. Fig 11.7.

### 11.8. Upholstery

The Mouse stroller is equipped with upholstery, which is completely removed for cleaning. To do this, detach all the elements fastening it to the stroller frame. How to clean the upholstery is described in point 17.8. Cleaning and care.

## 12. Accessories

### 12.1. Adjustable head supports - ↔ and ↕

To adjust the head support, draw aside the seat upholstery on the front and back of the seat. Next, take off the head support knobs, which are situated on the back of the seat and adjust the head supports at the front of the seat. That accomplished, screw the knobs back in the appropriate place at the back of the seat. Finally, draw the upholstery at front and back of the seat. Fig 12.1

### 12.2. Adjustable chest and loin pads - ↔ and ↕

To adjust the chest and loin pads, draw aside the seat upholstery on the front (Fig 12.2.1) and back of the seat (Fig 12.2.2). Next, take off the chest



Fig 12.1

Side pads

and loin pads knobs, which are situated on the back of the seat (Fig 12.2.3) and adjust the chest and loin pads at the front of the seat. That accomplished, screw the knobs back in the appropriate place at the back of the seat. Finally, draw the upholstery at front and back of the seat. (Fig 12.2.4)



Fig 12.2.1




Fig 12.2.2



Fig 12.2.3



Fig 12.2.4

*Unfasten the upholstery at both sides / Release the knobs / Fasten the upholstery*



Fig 12.3

**Therapeutic table**

### 12.3. Therapeutic table

The Mouse stroller can be equipped with a non-slip therapeutic table with a maximum load of 3kg. – Fig.12.3.

The therapeutic table is mounted with the same fixings as the safety barrier so you cannot use both of these elements at the same time. Slide the table into the attachment tubes and listen for the 'click'. To assemble the table, slide them onto the frame, latch. The table will not be properly installed until the latch engages the frame and snaps into place. This is evidenced by the clear "click" sound coming from a well-connected therapeutic table mount. To disassemble the table, pull both latches on the frame and pull out the table at the same time.

Caution! The horizontal arrangement of the table may change with the adjustment of the cradle and thus items placed on the table must be secured.

### 12.4. Seat minimizing kit

The Mouse stroller is equipped with an insert (Fig. 12.4). It is used to minimize the inside dimensions of the seat width and depth. The upholstery of the seat minimizing kit is removable thanks to the zipper sewn in, which allows it to be replaced or cleaned.

- Assembly. The insert is installed anywhere in the stroller, at any height adapting to the patient's indications. Installation is done with Velcro on the back side of the insert.
- Disassembly is done by removing the element - pulling the velcro from the stroller upholstery.



Fig 12.4

*Seat minimizing kit*



### 12.5. Travel bag

The large and elegant travel bag is mounted on two handles on the inside of the handle to push the stroller (see Fig. 12.5). The maximum bag load is 2kg.

### 12.6. Sun umbrella

The Mouse stroller can be equipped with a Sun umbrella (see Fig. 12.6) that helps protect the child from the sun.

### 12.7. Winter gloves

Two winter gloves that increase the comfort of a child's assistant during colder weather can be put on the handle to push the stroller (see Fig.12.7). First you need to put on two material clamps and fasten them with velcro on the handle to push the stroller, then gloves are put on these material clamps and fastened with snaps.

### 12.8. Winter sleeping bag

The ewa stroller can be equipped with a Winter sleeping bag, increasing the child's comfort during colder weather (Fig. 12.8). This optional equipment requires installation to other optional equipment - a barrier. To mount the Winter sleeping bag, place it on the footrest from below, and transfer it through the safety barrier from above and fasten it with snaps to the frame of the stroller.

### 12.9. Mosquito net

The Mouse stroller can be equipped with a mosquito net (see Fig. 12.9), which prevents insects from coming into direct contact with the child. This additional accessory requires fitting to other optional equipment. In order to mount the mosquito net, put it from above on the entire adjustable canopy with window, and then cover the mosquito net seat and footrest.

### 12.10. Raincover

The Mouse stroller can be equipped with a rain cover that helps protect the child from moisture. This optional equipment requires mounting to other optional equipment. To mount the rain cover, put it on top of the entire adjustable canopy with window and then cover the seat and footrest with the cover (see Fig. 12.10).



Fig 12.5  
Travel bag



Fig 12.6  
Sun umbrella



Fig 12.7  
Winter gloves



Fig 12.8  
Winter sleeping bag



Fig 12.9  
Mosquito net



Fig 12.10  
Raincover

### 13. Technical specification of the Mouse stroller

length of the stroller	95-124 cm
the width of the stroller	60 cm
the height of the stroller	122 - 134 cm
total weight	17.2 – 23.2kg
effective seat depth	22* / 26 cm
seat width	24* / 32 cm
backrest height	52 cm
backrest width	30 cm
distance between footrest and the seat / number of adjustments	14.5-23 cm / 4 positions
length of the stroller, folded up for transport	frame 84 cm /seat 72 cm
width of the stroller, folded up for transport	frame 58 cm /seat 42 cm
height of the stroller, folded up for transport	frame 25 cm /seat 39 cm
capacity of the shopping basket	5 kg
carrying capacity of a tourist bag	2 kg
max. patient weight	30 kg
patient height	80 - 110 cm
seat plane angle / number of adjustments	8°-23° / 4 positions
backrest angle / number of adjustments	16°-88° / 4 positions + lying position
footrest angle / number of adjustment	ok 84°-173° / 6 positions
effectiveness of parking brake of wheelchair positioned forwards down the slope	20°
effectiveness of parking brake of wheelchair positioned backwards down the slope	13°

#### COMMENTS

\*) values shown with inserts

Measurement value tolerance +/- 20mm, +/- 1.5kg, +/- 1 °

The manufacturer reserves the right to introduce changes to the technical parameters

#### 13.1. Material composition of the stroller

Powder coated steel and aluminium, fabrics 100% polyester and plastic (PP, PCV, ABS, PA, PE, PU).

### 14. Delivery and tools

The Mouse stroller is delivered in a cardboard box and pre-assembled. The original packaging should contain the following items:

- Folded stroller frame with seat, backrest and wheels
- Optional equipment as ordered
- User manual with warranty card

Before use, check and make sure that the carton contains all the elements and none of them has been damaged, e.g. during transport. In the event that it is found that in the carton received is missing or damaged any component, contact the Seller and/or the Manufacturer in order to receive the missing elements or replace the product with a new, fully complete and functional.

The multifunctionality and extensive equipment of the Mouse stroller allow you to tailor the stroller to your personal requirements. To do this, read the entire manual carefully before use.

No tools are required to assemble and disassemble the Mouse stroller.

A Phillips screwdriver is required to adjust the height of the footrest. An Allen key is needed to adjust the backrest height and seat depth.

## 15. Assembly and disassembly of the stroller

- **Stroller assembly.** The product pre-packed in the original carton is already pre-assembled. Complete the installation by completing the following steps:
  - a) mount the rear wheels on the axle
  - b) unfold the stroller frame with a footrest adjusting the parameters:
    - enable/disable lock of the direction of travel of the front wheels
    - adjust the height and angle of the stroller push handle
    - adjust the seat tilt angle
    - set the tilt angle of the backrest
    - adjust the length of the five-point belts
    - adjust the height of the footrest and the tilt angle of the footrest
  - c) install appropriate accessories and adjust parameters, for example:
    - adjust the usable width and depth of the seat (set minimizing kit)
    - install a safety barrier or a therapeutic table
    - adjust the deflection and stabilization belts
    - adjust the canopy with the window
    - install and adjust the wedge
    - put on a travel bag
    - put on a shopping basket
    - put on selected items of equipment dedicated to the seasons: winter sleeping bag, winter gloves, mosquito net, raincover, umbrella
- **Dismantling the stroller.** In order to dismantle the stroller, it is necessary to:
  - a) remove the optional equipment
  - b) unclip the seat from the frame of the stroller
  - c) lock the front wheels, fold the footrest and stroller frame
  - d) unclip the rear and front wheels
  - e) fold in the stroller closing it with the transport lock.

## 16. Preparation for use

Before using the stroller:

- Check the good technical condition of the stroller and all its elements, and if the technical condition of the stroller or its single element is bad or the user has any doubts about it, then the stroller must not be used. Please contact the Seller and/or the Manufacturer.
- Check that the stroller frame has been folded out correctly and that the folding mechanism has been properly closed (locked).
- Check that the seat is correctly attached to the stroller frame
- Check that all four wheels are properly installed and that the tread tires are in good condition on all four wheels.
- Check that the central brake is working properly and apply the brake.
- Correctly make all regulations and adjust the stroller to the child's needs, including setting the appropriate depth of the therapeutic seat and the height of the backrest, appropriate inclination of the therapeutic seat and backrest inclination, set the appropriate inclination of the footrest, adjust the appropriate tension of the belts on the backrest and on the seat, check the length of the five-point belts. Check that all adjustments are properly seated and secured.



- Check if the stroller upholstery is dry.
- Secure the child with all available stabilizing elements, such as: five-point belts, narrowing and shallowing inserts, abduction-stabilizing belts, wedge, footrest, footrest tapes, a safety barrier or a therapeutic table. This is to prevent the child from falling out and sliding out of the stroller and to prevent the child's feet from falling off the footrest. Make sure that all adjustments made earlier are suitable and safe for the child.
- Make sure that the child is not placed in the stroller without clothing that protects the body at the points in contact with the stroller's upholstery. This is to prevent chafing of the child's body and to keep the child's hygiene in mind.
- In the case of children who have problems with incontinence, a special oilcloth should be placed on the stroller seat, which can be purchased in pharmacies or medical stores.
- In the case of using any additional aids, such as vents, drains, bags, catheters, etc., make sure that they are well protected against falling out, pulling out, tightening and not interfering with the use of the stroller.

## **17. Using the stroller**

### **17.1. Carrying the stroller**

The best way to carry a stroller is to use it with its wheels. In cases where this is not possible (e.g. to transport the stroller in a car), for safety reasons, the stroller can be lifted and moved after dismantling the additional equipment and the front and rear wheels, in order to reduce the load. When carrying or lifting the stroller frame with the wheels mounted, fold the stroller frame and then hold the stroller frame only by the frame's fixed elements. For safety reasons, do not lift the entire stroller, i.e. the frame with the seat and backrest, wheels, footrest and mounted accessories, as this may damage your health. Do not lift the stroller by its footrest or safety barrier, therapy table or any other accessories, as this may damage the product and endanger the safety of others.

### **17.2. Getting on and off the stroller**

When getting on and off the stroller, place the stroller as close as possible to a stable and safe place where the user will change. After making sure that the central brake is on and the footrest is raised to avoid standing on it, the user can change.

Do not stand on the footplate, as there is a risk of tipping over the stroller, injuring the user and damaging the product. Make sure that when getting in or out of the stroller, the child and any other person should not stand on the footrest, as this may damage it and injure the user!

### **17.3. Correct position in a stroller**

For safety and health reasons, it is recommended that the user have the pelvis and back as close as possible to the backrest as much as possible to support the spine evenly (optimal support for all parts of the spine so that there is no free space).

The thighs should be supported on the largest possible surface, along their entire length, which is used to adjust the depth of the stroller seat, and the feet should be supported on the footrest plate which prevents slipping, maintains optimal leg muscle tension and carries some of the weight from the upper body.

Due to the large adjustment possibilities of the stroller, during the angular adjustment of the seat and angular adjustment of the backrest, it is possible to set the position in which the child's head is below the pelvic line. It is forbidden to use this position as it could lead to the risk of the child choking.

#### 17.4. Riding in a stroller

When using the stroller, always hold the handle to push the stroller to prevent the wheel from rolling slowly. Release the central brake while holding the push handle of the stroller and push the stroller in a safe direction.

#### 17.5. Moving on slopes

Always position the stroller forwards when climbing or driveways. Limit the speed on sloping surfaces and move as slowly as possible. The stroller assistant (parent, guardian) must have sufficient strength to control the stability of the stroller.

Always remember that the stroller's center of gravity changes when it overcomes unevenness or slopes on the ground. In this case, be extremely careful and secure the stroller. When using ramps, make sure that the ramp is suitable, that the wheels of the stroller are high enough to overcome the ramp, that the frame of the stroller will not touch the ramp, and that a second person should be provided as there is a risk of the stroller tipping over, damaging the product and injuring the user.

#### 17.6. Overcoming thresholds or curbs

When overcoming small obstacles (small thresholds, small curbs, etc.), block the tilting foot with the foot and the handle to push the stroller downwards and towards you to raise the front wheels of the stroller and carefully drive the front wheels over the obstacle. Then slowly pull the push handle upwards away from the ground to raise the rear wheels of the stroller and carefully drive the rear wheels onto the obstacle.

#### 17.7. Transport in a car and other vehicles

To transport the Mouse stroller in the car, remove all accessories from the stroller and store them in a safe place. Then dismantle the stroller (see chapter 15) and place it in the trunk.

If the dismantled stroller is not in the trunk of the car separated from the passenger compartment, then the stroller frame must be securely attached to the vehicle. For this purpose, unused seat belts may be used. In order to transport a stroller in other vehicles, such as trains, trams, buses, planes, ships, etc., please ask the persons responsible for transport in these means of transport about the conditions of safe travel to the destination with a stroller and the appropriate requirements.

In order to transport the stroller on other vehicles, remove all accessories from the stroller and store them in a safe place. Then dismantle the stroller (see chapter 15) and place it in the baggage storage area of the given vehicle. If the dismantled stroller is not in the trunk of the vehicle separated from the passenger compartment, the stroller frame must be securely attached to the vehicle.

**The Mouse stroller should not be used as a seat to transport a child in a vehicle.**

Mouse stroller labeling (see Fig. 17.7. NOT intended for use).

The label is located on the frame near the front wheel.

Not AS/NZS 3696.19 compliant

Fig. 17.7. NOT intended for use

#### 17.8. Cleaning and care

The upholstery should be hand washed only (recommended water temperature – 30 degrees). Washing in a washing machine is not allowed. All stains should be cleaned immediately after.

Attention! Do not clean the upholstery with the use of any chemical compounds composed of solvents or bleaches. Manufacturer does not take any responsibility for damages caused by the use of inappropriate cleaning products.

Do not rub off – this may cause the stain to spread on a larger surface but also it may distort the texture and damage the surface.

Always be sure that the product is dry before using.



Hand wash only



Do not bleach



Do not dry in a  
tumble-drier



Dry in a horizontal  
position



Do not iron

In case of liquid stains like food, drink or mineral plant oil stains:

- dry the whole surface of the stain gently pressing it with a white paper towel, tissue or a clean, absorbent wipe or cloth

- in order not to enlarge the stain, moisten a small cloth with a light detergent solution and warm water (or cold water in case of blood stains) then rinse the stain removing it from outside in. Note that an excess amount of water may create circles on the dry and clean material afterwards.

Moreover it may enlarge the stain or even damage the upholstery and the foam inside of it.

Depending on the type of the stain and its size, the described procedure may have to be repeated.

Continue this procedure till complete stain removal however always remember to remove the excess moisture before another rinse.

Next rinse the material with lukewarm water (without detergent) and a clean sponge or clean, white cloth. Dry out.

## 18. Storage of the stroller

The stroller should be stored in a clean, dry and well ventilated room, preferably at 0-35°C.

### 18.1. Preparation of the stroller for long-term storage

In order to prepare the Mouse stroller for long-term storage, secure the product and perform all service operations in the "Periodic Inspections" chapter as if they had occurred.

The stroller should be stored in a clean, dry and well ventilated room.

### 18.2. Corrosion protection

To avoid corrosive damage to the stroller, it should not be used or left in wet and damp rooms.

When the truck is wet, wipe all wet parts with a soft cloth and dry with a light oil preservative. This does not apply to the upholstery of the stroller, the provision of cleanliness and maintenance of which is described in detail in the Cleaning and care chapter.

As a standard, the stroller is covered with the highest quality coatings. In the event of damage to the varnish, protect the damage by applying a touch-up compound using the varnish proposed by the varnish Manufacturer.

Regular cleaning and lubrication of the stroller mechanisms with light oil guarantees long and trouble-free use of the stroller.

## 19. Re-use

The product is suitable for reuse.

All reusable products are subject to a particular burden.

Therefore, their features and performance may deviate from the new product to such an extent that the safety of patients and third parties during their use is at risk.

The product with proper care and maintenance retains its quality for a period of 3 years and this is the expected period of its use.

For secondary use, please:

- Remove all upholstered parts and replace them with new ones.
- The stroller frame with wheels and steering elements should be thoroughly cleaned and disinfected with commercially available cleaning and disinfecting agents.
- The product must then be checked by an Authorized Service Center (see Chapter 24) or the Manufacturer for wear and damage.
- All worn parts and components that do not match or are unsuitable for the user should be replaced with new ones.
- The product prepared in this way can be put back into use.

## 20. Service and periodic inspections

For safe and trouble-free use of the Mouse stroller, it must be regularly and thoroughly inspected and serviced as indicated in the table below.

In a situation where the stroller user is not able to carry out the activities himself, he should return the stroller to a specialized stroller service point or directly to the Manufacturer, who will properly review and will perform necessary service activities. These activities may be charged extra.

All repairs should be carried out in specialized, authorized Services (see point 24) or directly at the Manufacturer.

Only original spare parts that are supplied or recommended by the Manufacturer may be used in the stroller.

### Service should be performed:

- **Each time before using the stroller**

- Check the efficiency of the central brake - After activating the central brake, both rear wheels of the stroller should be locked and cannot rotate.
- Check the attachment of structural elements - It is necessary to check all mechanisms of the stroller and its movable elements, in particular, pay attention to: the correct mounting of the therapeutic seat with the backrest to the stroller frame, the correct mounting of all four wheels, the correct functioning of the five-point seat belts, properly locked seat angle adjustment mechanism, properly locked backrest angle adjustment mechanism and footrest correct mounting.
- Check the condition of the tires - Check that the tires on all four wheels are not damaged or worn.
- Inspection of solid frame parts, seat with backrest and footrest frame - Visually inspect the condition of the structural profiles. These elements must not be damaged or cracked.

- **Every one month**

- Cleaning the stroller - Clean any moving parts, bearings, adjustment mechanisms with a soft cloth and then lubricate them with a few drops of light oil. Be careful not to soil the upholstery of the stroller or your own clothes.
- Check all bolted connections for looseness - Carefully check that all bolts and knobs in the stroller are properly tightened.

- **Every six months**

- Check the cleanliness of the upholstery and inspect the general condition of the stroller - Do not allow the stroller to be neglected as described in the chapter „Cleaning and care“.
- Checking the general condition of the stroller. The stroller must not be neglected

**21. Disposal**

Disposal and recycling: the used stroller should be returned to a recycling center in accordance with the regulations in force in the given country.

**22. Compliance with requirements for medical devices**

Class I medical device in accordance with the Regulation of the European Parliament and of the Council (EU) 2017/745 of April 5, 2017 on medical devices. The manufacturer issued an appropriate declaration of conformity for the product, which is available at [www.mywam.eu](http://www.mywam.eu)

## 23. Notes



## 24. List of authorized service points

### POLAND

MyWam Kupiec, Bartold, Angres sp. j.  
ul. Szczecińska 10, 41-516 Chorzów, PL  
tel: +48 32 733 11 31  
www.mywam.eu  
mywam@mywam.pl

### AUSTRALIA

Dejay Medical & Scientific Pty Ltd  
Unit 2, 6 Money Close, Rouse Hill NSW 2155  
Phone: (02) 9838 8869  
Website: www.dejay.com.au  
E-mail address: info@dejay.com.au

### AUSTRALIA

Motum Pty Ltd  
Unit 1, 7 King Edward Road, Osborne Park, 6017  
Phone: (08) 9381 2087  
Website: www.motum.com.au  
E-mail address: hello@motum.com.au

**25. Warranty Card**

Series and number: .....

Date of purchase:..... - ..... - .....

Seller stamp and signature

**Warranty conditions:**

- The manufacturer grants the buyer a warranty for the purchased stroller, for a period of 24 months from the date of purchase. Defects revealed during this period will be removed free of charge.
- The Buyer submitting a complaint on the stroller should deliver the stroller to the Seller or to the Manufacturer in the safety packaging, together with the elements of the stroller protected against damage.
- Damage caused by insufficient securing of the stroller shipment and its components is not subject to warranty repairs.
- The complainant may not demand gratification if he improperly marks or incorrectly insures the shipment of the stroller and its components.
- The cost of shipping the complained stroller to the Seller or to the Manufacturer shall be borne by the buyer complaining the stroller.
- The warranty period is extended by the time needed to repair the stroller.
- The activities provided for in the manuals, which the user is obliged to perform on their own and at their own expense, cannot be performed as part of a warranty repair.
- The buyer has the right to replace the complained stroller with a new stroller if the same damages or defects that have already been removed twice and which prevent the stroller from being used are still present.
- The stroller should be used in accordance with these operating manual.
- Unauthorized use of the stroller, inconsistent with these operating manual or inconsistent with the principles of correct operation will void the warranty.
- Making any structural changes to the stroller will void the warranty.
- Damage to the stroller caused by: improper operation, improper maintenance, incorrect storage of the stroller, mechanical damage as well as natural abrasion and wear of parts (e.g. upholstery, wheels) are not covered by the warranty. They are repaired at the user's expense.
- The manufacturer warns against performing warranty repairs at unauthorized service facilities. Finding such repair or unauthorized alteration of the stroller will void the warranty.
- The warranty card is valid provided that it has: the date of sale entered, confirmed by the Seller's company seal, the Seller's signature, type and serial number of the stroller.
- The warranty card is valid as long as it does not contain deletions and corrections.
- The stroller manufacturer runs its paid post-warranty service.
- In matters not covered by this warranty card, the provisions of the Civil Code shall apply.
- All matters related to the warranty and complaint of the stroller should be directed to the Seller or directly to the Manufacturer