



**ALLIED
MEDICAL**
LIMITED


USER MANUAL

AML Gutter Frame Walker

MP61040070



Parts of the product:



1	Lever	7	Steel net
2	Braking cable	8	Rear wheel
3	Plastic tray	9	Screw
4	Steel basket	10	Handle tube
5	Front wheel	11	Arm support
6	Seat	12	Handgrip

Warranty:

1 year for main frame only, not including any consuming parts(wheels, brakes, and so on...)


Intended purposes

Walking Aid: including Walkers and Rollators, an all- or partially-wheeled mobility aid in the form of a waist-high framework with handgrips and four or three legs designed to provide a stable support to a person with a disability, a geriatric or an infirm person when standing or ambulating; it may also be used during rehabilitation. It is constructed of lightweight metal, plastic, and rubber materials in the shape of a mobile support structure which the user holds and pushes in front of them as they walk. It is non-foldable and includes a braking mechanism; it may be width-adjustable and some types may include a seat on which the user can rest.

Intended users

A rollator consists of a frame with two, three or four large wheels, handlebars, and a built-in seat. If your condition requires you to stop and rest often, a rollator could be a better choice for you. Whereas a walker has no seat, a rollator with a seat and a crossbar for back support allow the user to stop and rest when needed. A rollator is excellent for elderly with mild to moderate issues with balance, strength or endurance. Elderly often start using them after they have had a fall but I encourage the elderly to use them before they have a fall as a way to prevent injury. The last thing anyone wants is a hip fracture. Rollators also have accessories that can be purchased, such as baskets to hold a shopping bag.

Labels and symbols on the product:



A	Reference Code
B	Product description
C	Manufacturer's address
D	EC Representative
E	Serial number
F	UDI
G	Date of manufacture
H	Maximum user weight
I	EAN 13

Limitation of liability:

Maxhealth Corporation accepts no liability for damage arising from:

- Non-compliance with the User Manual
- Incorrect use
- Natural wear and tear
- Incorrect assembly or set-up by the purchaser or a third party
- Technical modifications
- The usage of unapproved 3rd party accessories instead of a tested and approved accessories
Unauthorized modifications and/or use of unsuitable spare parts

A written authorization of Maxx Corporation must be obtained before installing additional adaptations on this wheelchair. Otherwise no liability claims can be made.

Standards and Regulations

Quality is fundamental to the company's operation, working within the disciplines of ISO 13485.

Maxhealth Corporation is continuously working towards ensuring that the company's impact on the environment, locally and globally, is reduced to a minimum.

Important note

To ensure your safety and to prolong the lifetime of your product, please take good care of it and have it checked and serviced on a regular basis.

Safety instructions

- Read these instructions for use carefully.
- All components should be checked for damage and a secure fit prior to use.
- Use the rollator only as a walking aid.
- Use the rollator only on an even and solid surface.
- Rollator may not be used as a seat in vehicles, such as a bus or train!
- It is prohibited to stand on the rollator or to transport a person or animal on the rollator.
- Take note of the permitted maximum weight of the user. This information can be found in the technical data in these instructions for use as well as on the identification plate
- Use the tray only for small, lightweight objects which are not sharp or hot. Max. load: 5.0 kg

- Use the basket only to transport objects and do not overload it! Max. load: 5.0 kg
- Use on escalators and the like is expressly prohibited!
- The push handles must always be adjusted parallel to the direction of movement – risk of tipping!
- Do not hang any bags on the handle – risk of tipping!
- On steep surfaces, use the service brake (pull the brake lever towards you).

Unfolding



Place the rollator on its wheels. You can unfold the chassis frame by pulling one of the rear wheel tubes.

Push the cross brace straight down until it easily clicks into place.

Rotate the movable folding guard as far as it will go so that it is vertical to the cross brace.

Assembly of accessories

Hang the basket on the small hooks in the front underneath the seat and place the tray down with the small feet in the small holes in the seat section.



Height adjustment of handles



Insert the tube of the adjustable handgrip into the frame at least to the “MAX” marking.

Adjusting the forearm supports

Loosen the locking lever beneath the armrest by turning it anticlockwise. The brake handle can now be pushed until the ideal setting is achieved. The brake handle must be pushed in at least to the “Max” mark. Then tighten the locking lever again by turning it clockwise.



Angle adjustment of the forearm supports



Loosen the star screw on the side until you can move the forearm supports outwards out of the lock. Now you can rotate the forearm supports in the desired angle and push back into the lock. Tighten the star screw again.

Adjusting the brakes

If braking is poor or too hard, the brake can be adjusted either at the brake housing above or using the adjusting screw on the Bowden cable below) to achieve the desired braking effect.



The brake is correctly adjusted when the brake lever as well as the brake wires is free of play and the wheels can be rotated without grinding noises when the brake is released.

Braking | Parking

To brake when walking with the rollator, pull both brake handles upwards with your hands. Always brake slowly while applying equal pressure with both hands.

To apply the brakes/park the rollator, press both brake handles downwards until they audibly click into place.



To release the brake, pull both brake handles upwards once again until they are released (brake handles are in the middle position).

The braking effect can be impaired by wear and tear on the tires. For this reason, regularly check the braking effect.

Sitting

Place the rollator on an even and solid surface. Then apply the brakes. The rollator must be stable and may not slip. Only sit upright; do not lean backwards
Seated position permitted only as shown and against the direction of movement.

Folding



To fold up the rollator, open the folding guard, pull the seat and rollator upwards!
Remove accessories, if necessary.

Specification Sheet

Product	FOREARM SUPPORT STEEL ROLLATOR
Item No.	HE61040070
Weight capacity	300lbs
Product dimension L×W×H	660*592*(1000-1190)mm
Product folding dimension L×W×H	592*240*(1060-1250)mm
Product N.W.	12.9KG
Handle height adjustment min-max	1000-1190mm
Seat dimension L x W	386*162mm
Seat material	PU
Front & Rear wheel dimension	φ191
Frame material	Q195 Steel
Frame color(powder coating/liquid painting)	Powder coating grey
Foldable yes/no?	yes
Basket yes/no/optional?	yes
Cane holder yes/no/optional?	yes
Package	1PC/CTN
Product G.W.	14.2KG
Carton dimension L x W x H	590*255*715mm
1x40ft HQ loading Q'ty	660PCS
FOB Port	XINGANG TIANJIN,CHINA