ENGLISH

MD Medical Device

INTENDED USE

The device is intended to provide gross immobilization to the cervical spine. The device must be fitted and adjusted by a healthcare professional.

Indications for Use

Situations requiring gross immobilization of the cervical spine, these may include:

- Pre and post c-spine surgery
- C-spine precaution for trauma patients

Contraindications

None known.

Warnings and Precautions

- Do not rely solely on extrication collar for sufficient immobilization.
- Do not use improperly sized collars
- To be used only as directed by local Emergency Medical Service Authorities

GENERAL SAFETY INSTRUCTIONS

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

The device is for single patient – multiple use.

FITTING INSTRUCTIONS

SIZING

Position your lowest finger on the sizing line. Select the size closest to the height of your fingers (**Figure 1**). Position the collar with the selected size (**Figure 2**). Press the lock into the collar until secured. Verify by attempting to readjust (**Figure 3**).

APPLICATION

Supine Position

Maintain neutral alignment. Select appropriate size. Slide collar under patient's neck (**Figure 4**).

Seated Position

Maintain neutral alignment. Select appropriate size. Slide collar behind patient's neck. Position chin support.

Secure Collar

Position chin support. Tighten collar and attach hook and loop strap (**Figure 5**).

USAGE

Cleaning and Care

- Wipe the plastic collar shell clean with mild soap and water.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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©Copyright Össur 2021