



Check heel is underneath and not twisted to one side.



## **INSTRUCTIONS:**

- 1. Gather stocking to the toe with your thumbs inside.
- 2. Pull stocking onto plastic blade.
- 3. The stocking must be pulled on tightly to the toe checking that the heel is underneath and not twisted to one side.
- 4. Sitting securely on your bed or chair, lower the stocking aid to the floor while holding the straps in your hands.
- 5. Slide your foot into the plastic blade and wiggle your foot into the toe of the stocking.
- 6. At this point lift heel off the floor, straighten leg and point toes. Pull straps alternately with a quick wiggling action and stocking will be pulled up the leg until you can reach it.