- Note:** Long hair should be placed outside of the back piece.
- Apply the front piece with the chin secured in the recess and center it to secure neutral alignment. The front piece overlaps the back piece to ensure effective immobilization and comfort (**Fig. 6**). The front arrow should point upward.
- Apply the Hook and Loop strap on both sides. Tighten straps alternately to an equal length on both sides (**Fig. 7**). Excess may be trimmed.
- The Plastazote material can be trimmed or heat moulded as necessary depending on the patient's anatomy. Ensure that all edges are rounded off smoothly.

Note: Device application must be secure to ensure proper fit, maintain alignment, and prevent the chin from slipping inside.

Device removal
Ensure that correct spinal protocol is followed. Undo the straps, remove the front piece and slide the back piece out from behind the head.

USAGE
Cleaning and care

- Wipe with a damp cloth and mild soap.
- Air dry.

Note: Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener. **Note:** Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

DISPOSAL
The device and packaging must be disposed of in accordance with respective local national environmental regulations.

LIABILITY
Ossur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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Patient Type	Neck Height	Circumference	Circumference Adjustment	
Stout	2 (2,25" / 5,7 cm)	Large (16" - 19" / 41 - 48 cm)	No tear off	
Average Male	3 (3,25" / 8,3 cm)	Medium (13" - 16" / 33 - 41 cm)	Tear off 1 section on each side	
Thin Female	4 (4,25" / 10,8 cm)	Small (10" - 13" / 25 - 33 cm)	Tear off 2 sections on each side	

Supine Application
Ensure that correct spinal protocol is followed.

- Slide the back piece behind the patient’s neck and center it, noting arrow on panel pointing up (**Fig. 5**).

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