ENGLISH

SYMBOLS

MD	Medical Device
MR	Magnetic Resonance (MR) safe

DESCRIPTION

These Instructions for Use (IFU) provide additional information that enables the usage of an accessory with a particular medical device. Refer to the applicable IFU for the usage of the medical device.

INTENDED USE

The device is intended to enhance gross immobilization of the cervical spine.
The device must be fitted and adjusted by a healthcare professional.
X-ray and CT lucent.

Indications for use

For conditions requiring gross immobilization of the low cervical and/or high-thoracic spine:

- Low C-Spine and high T-Spine precaution for trauma patients.
 - Immobilization for pre and post low C-spine and high T-spine surgery.
 - Other conditions requiring gross immobilization of the low-cervical and/or high-thoracic spine.
- No known contraindications.

Warnings and Cautions:

Regular cleaning of the device and the skin beneath along with inspection of the skin for any signs of irritation is required to reduce the risk of skin ulceration. Bedridden patients are at an increased risk for skin ulceration.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.




Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities. The patient should immediately contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
 - If any pain, skin irritation, or unusual reaction occurs with the use of the device.
- The device is for single patient – multiple use.

FITTING INSTRUCTIONS

Device Application

Proper application of the device is important for immobilization and patient comfort.

Device Removal
Before removing the Stabilizer and collar, mark the position of the hook-and-loop strap location to ensure reapplication maintains the proper fit. With the Stabilizer and collar secured on the patient, loosen and unbuckle the Stabilizer strap on the right side and then on the left side. To remove the collar, unfasten the hook-and-loop strap on the right side and then the left side. Remove Stabilizer and collar from the patient. Plugs can be removed only when the Stabilizer is no longer necessary and only the collar is required as instructed by a qualified healthcare professional.

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Usage

Cleaning and care

- Wipe with a damp cloth and mild soap.
- Air dry.

Note: Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

Note: Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:
• Device not maintained as instructed by the instructions for use.

- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

The two-piece design consists of a front piece and a back piece which are packaged as a set. Always maintain the patient's head in neutral alignment.

Ensure that correct spinal protocol is followed.

1. Insert the plastic reinforcement of the collar front piece into the top of the Stabilizer front piece (**Fig. 2**). Slide the pieces together until the holes are in alignment.
2. After the Stabilizer and collar holes are in alignment, insert the locking plugs (**Fig. 3**). Detach the nylon fastener after the plugs are securely locked into position and the two pieces cannot be separated.
3. Repeat previous steps for the collar back piece and Stabilizer back piece.

4. After properly applying the collar, insert the hook-and-loop straps located on the Stabilizer back piece through the buckles located on the Stabilizer front piece (**Fig. 4**).

5. Pull both straps simultaneously and with equal pressure toward the Stabilizer back piece. Verify that all components are centered on the patient and that the patient is in neutral alignment. Fasten the strap securely and remove excess material, as necessary (**Fig. 5**).

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