



Instructions for Use

PHILADELPHIA[®]
TRACHEOTOMY COLLAR

		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	7
FR	Notice d'utilisation	10
ES	Instrucciones para el uso	13
IT	Istruzioni per l'uso	16
NO	Bruksanvisning	19
DA	Brugsanvisning	22
SV	Bruksanvisning	25
FI	Käyttöohjeet	28
NL	Gebruiksaanwijzing	31
PT	Instruções de Utilização	34
PL	Instrukcja użytkownika	37
RU	Инструкция по использованию	40
JA	取扱説明書	43
ZH	中文说明书	46
KO	사용 설명서	49

SYMBOLS

Medical Device

Magnetic Resonance (MR) safe

INTENDED USE

The device is intended to provide gross immobilization to the cervical spine

The device must be fitted and adjusted by a healthcare professional. X-ray and CT lucent.

Indications for use

Conditions requiring gross immobilization of the cervical spine. This may include:

- C-Spine precaution for trauma patients
- Immobilization for pre and post c-spine surgery

No known contraindications.

Warnings and Cautions:

These instructions do not supersede the specific orders of your healthcare professional.

The use of a cervical collar may increase intracranial pressure (ICP) through jugular venous compression.

If an unstable fracture is suspected or unknown, with or without a sustained trauma, ensure additional spinal precautions are implemented to immobilize the spine.

Regular cleaning of the device and the skin beneath along with inspection of the skin for any signs of irritation is required to reduce the risk of skin ulceration. Bedridden patients are at an increased risk for skin ulceration.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should immediately contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

FITTING INSTRUCTIONS

Size Selection

Proper sizing is critical to achieve comfort and proper immobilization. Ensure that correct spinal protocol is followed.

Height: Measure the vertical distance from the top of the patient's shoulder to the tip of the chin (**Fig. 1-a**). The measurement corresponds to the height of the device. If the patient's measurement falls between two consecutive sizes, apply the smaller size first.

Circumference: Measure around the patient's neck (**Fig. 2-b**). The measurement corresponds to the size circumference of the device.

Measurement A - Neck Height				
Infant	1.25" (3.2 cm)			
Pediatric	1.75" (4.4 cm)			
Adult	2.25" (5.7 cm)	3.25" (8.3 cm)	4.25" (10.8 cm)	5.25" (13.3 cm)

Measurement B - Neck Circumference	
Size	Circumference
Infant	6" - 8" (15.2 - 20.3 cm)
Pediatric	8" - 11" (20.3 - 27.9 cm)
Adult Small	10" - 13" (25.4 - 33.0 cm)
Adult Medium	13" - 16" (33.0 - 40.6 cm)
Adult Large	16" - 19" (40.6 - 48.3 cm)
Adult Extra Large	19"+ (48.3 cm+)

Supine Application

- Slide the back piece behind the patient's neck and center it, noting arrow on panel pointing up (**Fig. 3**).
Note: Long hair should be placed outside of the back piece.
- Apply the front piece with the chin secured in the recess and center it to secure neutral alignment. The front piece overlaps the back piece to ensure effective immobilization and comfort (**Fig. 4**). The front arrow should point upward.
- Apply the Hook and Loop strap on both sides. Tighten straps alternately to an equal length on both sides (**Fig. 5**). Excess may be trimmed.
- The Plastazote material can be trimmed or heat moulded as necessary depending on the patient's anatomy. Ensure that all edges are rounded off smoothly.

Note: Device application must be secure to ensure proper fit, maintain alignment, and prevent the chin from slipping inside.

Device removal

Ensure that correct spinal protocol is followed.

Undo the straps, remove the front piece and slide the back piece out from behind the head.

Accessories and Replacement Parts

Please refer to the Össur catalog for a list of available replacement parts or accessories.

USAGE

Cleaning and care

- Wipe with a damp cloth and mild soap.
- Air dry.

Note: Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

Note: Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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