



Ride CAM[®] Wedge Kit Instructions for Ride Forward[™] Cushion

Symmetrical Use

Use Ride CAM Wedges **symmetrically** to promote more upright posture and decrease ischial pressure.

Maximum 2 wedges/side
(1 Ride CAM wedge kit/side).



Ride Forward Cushion with one pair of Ride CAM wedges.

Asymmetrical Use

Use Ride CAM Wedges **asymmetrically** to help level a flexible pelvic obliquity.

Maximum 2 wedges/side (1 Ride CAM Wedge kit/side).



Uncorrected pelvic obliquity on Ride Forward Cushion.



Partially corrected pelvic obliquity on Ride Forward cushion with use of 1 Ride CAM Wedge on low side.



Corrected pelvic obliquity on Ride Forward cushion with use of 2 Ride CAM Wedges on low side.
