

Paediatric Standing Solutions

	PRONE STANDERS							SUPINE / PRONE STANDERS					
	TOTSTANDER	FREESTANDER			PRONESTANDER			SQUIGGLES+	MYGO STANDER		HORIZON STANDER		
SIZES	1	1	2	3	1	2	3	One Size	1	2	1	2	3
AGE (APPROX.)	1 to 5	1 to 5	4 to 9.5	9 to 17	1 to 5	4 to 9.5	9 to 17	1 to 5	4 to 10	8 to 14	4 to 10	10 to 18	14 to Adult
MAX USER WEIGHT	20 kgs	60 kgs	80 kgs	90 kgs	60 kgs	80 kgs	90 kgs	22 kgs	50 kgs	60 kgs	50 kg	80 kgs	100 kgs
UPRIGHT	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
PRONE (ANGLE ADJUSTABLE)	No	No	No	No	0-40°	0-40°	0-40°	90-160°	90° to vert	90° to vert	90° to vert	90° to vert	90° to vert
SUPINE	No	No	No	No	No	No	No	Yes	Yes	Yes	Yes	Yes	Yes
WHEELED	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
HEIGHT ADJUSTABLE	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
CHEST PLATE HEIGHT - MIN	470 mm	470 mm	680 mm	975 mm	470 mm	680 mm	975 mm	550 mm	760 mm	980 mm	1000 mm	1400 mm	1550 mm
CHEST PLATE HEIGHT - MAX	640 mm	680 mm	975 mm	1350 mm	680 mm	975 mm	1350 mm	820 mm	1080 mm	1300 mm	1400 mm	1800 mm	1800 mm
PRODUCT DIMENSIONS													
BASE LENGTH	450 mm	650 mm	750 mm	750 mm	820 mm	1070 mm	1270 mm	950 mm	978 mm	1130 mm	1135 mm	1280 mm	1280 mm
BASE WIDTH	450 mm	650 mm	750 mm	750 mm	600 mm	650 mm	700 mm	590 mm	630 mm	650 mm	580 mm	640 mm	740 mm
ACTIVITY TRAY DEPTH	No Tray	400 mm	450 mm	500 mm	400 mm	450 mm	500 mm	460 mm	340 mm	340 mm	420 mm	520 mm	520 mm
ACTIVITY TRAY WIDTH	No Tray	525 mm	575 mm	625 mm	525 mm	575 mm	625 mm	560 mm	580 mm	580 mm	440 mm	590 mm	590 mm
KEY FEATURES	Chest band Multi-adjustable pelvic band, Height & angle adjustable kneeplate, Multi-adjustable sandals. Portable.	Angle-adjustable tray with activity bowl Chest support with laterals Belt, knee cups & sandals in various sizes available. Mobile base with lockable castors			Angle-adjustable tray with activity bowl Chest support with laterals . Belt, knee cups & sandals in various sizes available. Mobile base with lockable castors			Hip abduction of up to 60° in prone, upright and supine standing Fully horizontal loading in supine Height-adjustable, foldable push handle for ergonomic use New wrap around knee supports Individually adjustable footplates	Headrest options for supine set up (Whitmeyer compatible). Padded clear tray for prone or supine Independent angle adjustment in knees up to 25° -10 +10° plantar/dorsiflexion Mobile base with lockable castors Manual angle adjustment Belt, knee cups & sandals in a various sizes available		Headrest options for supine set up (Whitmeyer compatible). Separate supine & prone trays Mobile base with lockable castors Powered angle adjustment Belt, knee cups & sandals in a various sizes available. Extremely robust.		
WIPEABLE PADDING	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
REMOVABLE PADDING	No	No	No	No	No	No	No	Yes	Yes Removable covers	Yes Removable covers	Yes	Yes	Yes
COLOUR OPTIONS	Meadow	Meadow	Mondeo Blue	Mondeo Blue	Meadow	Mondeo Blue	Mondeo Blue	Purple	Blue	Blue	Blue	Blue	Blue/Black

Information correct as of August 2023. From time to time, product specifications are changed by manufacturers without prior notice
Transit tie down is available on all chairs for transport purposes.

Paediatric Walking Solutions

OTTOBOCK NURMI NEO



LECKEY MYWAY



MEYWALK 4



SIZES	1	2	3	1	2	3	1	2	3	4
AGE (APPROX.)	1 to 3	3 to 5	5 to 8	1 to 5	4 to 11	9 to 16	1 to 5	4 to 9.5	9 to Adult	9 to Adult
MAX USER WEIGHT	25 kg	40 kg	55 kg	25 kg	50 kg	80 kg	50 kg	100 kg	130 kg	130 kg
UNIT LENGTH	650mm	760mm	830mm	720mm	850mm	1100mm	800mm	850mm	930mm	1100mm
UNIT WIDTH	630mm	680mm	720mm	570mm	670mm	720mm	560mm	620mm	620mm	680mm
MAX. INSIDE LEG HEIGHT	300mm	450mm	550mm	440mm	610mm	810mm	490mm	640mm	810mm	1010mm
MIN. INSIDE LEG HEIGHT	300mm	450mm	550mm	220mm	400mm	550mm	250mm	270mm	570mm	770mm
MAX. THORACIC SUPPORT HEIGHT	640mm	790mm	840mm	920mm	1190mm	1490mm	740mm	810mm	1300mm	1490mm
MIN. THORACIC SUPPORT HEIGHT	530mm	660mm	710mm	635mm	870mm	870mm	520mm	710mm	980mm	1170mm
CHEST PLATE HEIGHT - MAX	Height adjustable grip bar, solid locking mechanism, lightweight and foldable, reverse roll locks on rear wheels, puncture proof low friction castors, stable base, swivelling front castors			Harness to provide core stability, supports the child in an upright posture enabling eye-level contact with standing peers, upper support unit with adjustable height, foldable frame for easy storage			Body posture of user can be set from vertical to 15° forward tilt, suitable for indoor and outdoor use, suspended spring system			

Information correct as of August 2023. From time to time, product specifications are changed by manufacturers without prior notice
 Transit tie down is available on all chairs for transport purposes.