

Ki Mobility Wheelchair Solutions

| | FIXED FRAME | | | | FOLDING FRAME | TILT-IN-SPACE | | |
|----------------------------------|--|---|--|---|--|---|---|--|
| | ETHOS | ROGUE | ROGUE XP | CLIK | CATALYST 5 | LIBERTY FT | FOCUS CR | FLIP |
| DESCRIPTION | Ethos is a highly adjustable frame, which still maintains rigidity in all the right places. Incorporating Iso Tech Suspension', yet still lighter than many frames that lack suspension. | At 4.81kg the Rogue, along with it's lightweight accessory options, is a wheelchair design that requires minimal effort to set up precisely to maximise the distance achieved with each push. | With all the ease of adjustment of the Rogue, the Rogue XP adds in the benefit of width growth, with only a few grams more in total wheelchair weight. | Easy click adjustment and ability to change the set up for kids requiring more stability to transfer from the chair to the floor independantly, the Clik offers more than most fixed frame chairs can. Brilliant for smaller children requiring exact set ups to propel well. | With close on rigid frame performance and weight, the Catalyst 5 is the chair for those needing a model that is highly configurable but still folds for transport. | True portability, low seat to floor heights, built in tilt in space, adjustable COG and a frame that once unfolded locks rigid, make the Liberty an ideal chair for those foot punting and self-propelling. | 5° of anterior tilt for partial stand transfer, full 50° of posterior tilt. Internally routed cabling requires almost no adjustment, narrow footprint, featherlight tilt mechanism with 8 precision sealed bearings. Different wheel sizes and COG adjustment to get the set up just right. | From its flat folding design easily fitting into car boots, to its use of Taper Lok Technology, the Flip is simply the most advanced paediatric tilt-in-space available today. |
| FRAME STYLE (PAEDIATRIC IN BLUE) | Fixed Frame Self-propelling | Fixed Frame Self-propelling | Fixed Frame Self-propelling | Paediatric Fixed Frame Self-propelling | Folding Frame Self-propelling | Tilt-in-space | Tilt-in-space | Paediatric Tilt-in-space |
| SEAT WIDTHS | 30 - 51cm (12" - 20") | 30 - 51cm (12" - 20") | 25 - 46cm (10" - 18") | 20 - 46cm (8" - 18") | 36 - 56cm (14" - 22") | 36 - 56cm (14" - 22") | 36 - 56cm (14" - 22") | 25 - 46cm (10" - 18") |
| SEAT DEPTHS | 36 - 51cm (14" - 20") | 36 - 51cm (14" - 20") | 25 - 50cm (10" - 20") | 20 - 46cm (8" - 18") | 36 - 56cm (14" - 22") | 36 - 51cm (14" - 20") | 36 - 56cm (14" - 22") | 30 - 51cm (12" - 20") |
| FRONT FRAME BENDS | 70, 75, 80, 85 and 90° | 70, 75, 80, 85 and 90° | 70, 75, 80, 85 and 90° | 80, 85 and 90° | 60, 70, 80 and 90° | 60, 70, 80 and 90° | 60, 70, 80 and 90° | 60, 70, 80 and 90° |
| REAR WHEEL SIZE | 22", 24", 25", 26" | 20", 22", 24", 25", 26" | 18", 20" 22", 24", 25" | 18", 20", 22", 24" | 20", 22", 24", 25", 26" | 12", 16", 20", 22", 24" | 12", 16", 20", 22", 24" | 12", 16", 20", 22", 24" |
| FRONT SEAT HEIGHT | 1" - 2" Taper: 42 - 52cm (16.5" - 20.5") 0" Taper: 34 - 52cm (13.5" - 20.5") | 1" - 2" Taper: 42 - 52cm (16.5" - 20.5") 0" Taper: 34 - 52cm (13.5" - 20.5") | 1" Taper: 42 - 52cm (16.5" - 20.5") 0" Taper: 34 - 52cm (13.5" - 20.5") | 29 - 51cm (11.5" - 20") | 34 - 53cm (13.5" to 21") | 33 - 51cm (13" - 20") | 33 - 52cm (13.5" - 20.5") | 34 - 51cm (13.5" - 20") |
| REAR SEAT HEIGHT | 34 - 52cm (13.5" - 20") | 32 - 51cm (12.5" - 20") | 29 - 50.8cm (11.5" - 20") | 29 - 52cm (11.5" - 20.5") | 29 - 50cm (11.5" - 19.5") | 33 - 51cm (13" - 20") | 33 - 52cm (13" - 20.5") | 34 - 51cm (13.5" - 20") |
| TRANSFER WEIGHT | 5.26 kg | 4.81 kg | 5.26 kg | 5.67 kg | 5.80 kg | 11.7 kg | 14.97 kg | 9.53 kg |
| USER WEIGHT LIMIT | 124.7 kg | 124.7 kg | 90.72 kg | 75 kg | 136 kg HD option 158 kg | 113 kg | 136.08 kg HD option 181.48 kg | 75 kg |
| OPTIONS | / | One Arm Drive | One Arm Drive, Drum Brake | One Arm Drive, Drum Brake | One Arm Drive, Drum Brake | Drum Brake | One Arm Drive, Drum Brake | Drum Brake |

Information correct as of February 2024. From time to time, product specifications are changed by manufacturers without prior notice.

Transit tie down is available on all chairs for transport purposes.

*22" depth is achieved with a 20" depth adjustable seat and is only available with T-Arm or Flip-Up Armrests and is not available with back upholstery.

