

Standard Powerchair Solutions

		MID-WHEEL DRIVE		REAR-WHEEL DRIVE					
		CTM HS-2850	CTM HS-2850 REHAB	CTM HS-1500	KARMA KP25.2	CTM HS-5600	WHILL Model C2	WHILL Model F	CTM HS-2750
DRIVE	Seat Type	Captain Style	Rehab Seat	Captain Style	Slung Upholstry	Captain Style	Captain Style	Captain Style	Rehab Seat
	Ground Clearance	2.4" (60mm)	2.4" (60mm)	1.6" (40 mm)	3.0" (76mm)	3.4" (85mm)	2" (51mm)	1.4" (36mm)	3.1" (80mm)
	Turning Radius	24.4" (620mm)	24.4" (620mm)	27.4" (695 mm)	34.0" (864mm)	31.1" (790mm)	29.9" (759mm)	30.7" (780mm)	37.4" (950mm)
FUNCTIONS	Folding	No	No	No	Yes	No	No	Yes	Yes
	Recline Function	No	No	No	No	Yes (Manual)	No	No	100 degree standard (not adjustable)
	Elevate Function	No	No	No	No	No	No	No	No
ELECTRONICS	Motor Type	420W 4900RPM	420W 4900RPM	150W 2800RPM	2 Pole 3600RPM	420W 4600RPM	150W x 2	150W x 2	250W 3250RPM
	Battery Size	2 x 36Ah	2 x 36Ah	2 x 12Ah	2x20Ah (14"W) 2x40Ah (16"+ W)	2 x 50Ah	10.6Ah	10Ah	30.15Ah
DIMENSIONS & WEIGHTS	Seat Width	19.0" (490mm)	16.5" - 18.5" (420 - 470mm)	18.5" (470 mm)	14" - 20" (355mm - 508mm)	20.0" (500mm)	16" - 20" (406mm - 508mm)	18" (457mm)	18.1" (460mm)
	Seat Depth	19.0" (490mm)	17.7" - 18.7" - 19.6" (450 - 475 - 500mm)	18.5" (470 mm)	16.0", 18.0" (406, 457mm)	20.0" (500mm)	16", 18" (406mm, 457mm)	16" (406mm)	18.5" (470mm)
	Seat to Floor Height	20.4" (540mm)	20.4" (540mm)	20.4" (540mm)	19", 483mm (Optional 21")	20.4" (540mm)	15.9" (404mm)	17.4 - 19.5" (443 - 495mm)	21.6" (550mm)
	Overall Length	39.0" (995mm)	39.0" (995mm)	31.9" (810 mm)	39" - 44" (991mm - 1118mm)	44.7" (1135mm)	38.8" (986mm)	36.8" (935mm)	42" (1055mm)
	Overall Width	24.0" (610mm)	24.0" (610mm)	23.6" (600 mm)	Seat width + 7" (177mm)	25.0" (635mm)	21.8" (554mm)	21.9" (556mm)	25" (635mm)
	User Weight Limit	136.0 kg	136.0 kg	115 kg	135.0 kg	181.0 kg	136.1 kg	113.4 kg	150 kg

Note: Measurements are rounded up to one decimal place.

