

Specialty Controls Script



Customer Details					
Client					
Assessor					
Organisation					
	ACC		MOH		Other

Alternative Drive Controls

For assistance with Specialty Controls, please get in touch with your Technical Product Specialist. If you are unsure of who this is, please call us on 0800 31 61 81, to be directed to the right person.

Stealth i-Drive Head Array with Bluetooth (requires EDU or SCIM)



SPIDHBT300-U
Full Head Array

Optional Accessories:

-  **Specs Switch**
AB100-SPBK
-  **Egg Switch**
 - Black** ASL320B
 - Green** ASL320G*
 - Red** ASL320R
 - Yellow** ASL320Y*
-  **Egg Switch Mount**
SPESM
-  **Egg Switch Extension Cable**
PRIDELE155901
-  **9" Lateral Rod with Ring***
SPSU260
-  **13" Lateral Rod with Ring**
SPSU264
-  **18" Lateral Rod with Ring***
SPSU265
-  **SCIM**
PRIDECTL143988
-  **Figure 8**
SPF8675-RING
-  **Headrest Mount**
17mm Ball
SPUHB4D33200FDM



EDU (includes Gooseneck Mount)
.PKEDUQL3
FRMASMB16062



Figure 8 & 13" Rod
SPF8675-RING
SPSU264
SPESM




Headrest Mount
SPUHB4D33200


Full Head Array Pack

- Includes:**
i-Drive Head Array
EDU
Gooseneck mount
Headrest mount
Armrest clamp
Figure 8 & 13" rod
Egg switch & mount
Extension lead (for egg switch)

Note: On a recline, substitute the headrest mount and armrest clamp to these.



Swing-Away Clamp Assembly
PRIDEMEC1809294






Headrest Mount
SPTWB485-17MM

Stealth i-Drive Tri-Array Head Support with Bluetooth (requires EDU or SCIM)

SPIDHBT200-U
Full Tri-Array



Optional Accessories:

-  **Specs Switch**
AB100-SPBK
-  **Egg Switch**
 - Black ASL320B**
 - Green ASL320G***
 - Red ASL320R**
 - Yellow ASL320Y***
-  **Egg Switch Mount**
SPESM
-  **Egg Switch Extension Cable**
PRIDELE155901
-  **SCIM**
PRIDECTL143988
-  **Headrest Mount**
17mm Ball
SPUHB4D33200FDM



EDU (includes Gooseneck Mount)
PKEDUQL3
FRMASMB16062



Figure 8 & 13" Rod
SPF8675-RING
SPSU264
SPESM




Headrest Mount
SPUHB4D33200


Full Tri-Array Pack

- Includes:**
i-Drive Head Array (SPIDHBT200-U)
EDU
Gooseneck mount
Headrest mount
Armrest clamp
Figure 8 & 13" rod
Egg switch & mount
Extension lead (for egg switch)

Note: On a recline, substitute the headrest mount and armrest clamp to the below.



Swing-Away Clamp Assembly
PRIDEMEC1809294









Headrest Mount
SPTWB485-17MM

Stealth i-Drive Tri-Array Head Support with Bluetooth (requires EDU or SCIM)

SPIDHBT310-1
Full Head Array Sip & Puff



Optional Accessories:

-  **Specs Switch**
AB100-SPBK
-  **Egg Switch**
 - Black** ASL320B
 - Green** ASL320G*
 - Red** ASL320R
 - Yellow** ASL320Y*
-  **Egg Switch Mount**
SPESM
-  **Egg Switch Extension Cable**
PRIDELE155901
-  **SCIM**
PRIDECTL143988
-  **Replacement Straw**
SP625-S
-  **Replacement Filter**
SP625-F
-  **Headrest Mount**
17mm Ball
SPUHB4D33200FDM



EDU (includes Gooseneck Mount)
PKEDUQL3
FRMASMB16062



Figure 8 & 13" Rod
SPF8675-RING
SPSU264
SPESM




Headrest Mount
SPUHB4D33200


Full Head Array Sip and Puff Pack

Includes:
i-Drive Head Array Sip and Puff
EDU
Gooseneck mount
Headrest mount
Armrest clamp
Figure 8 & 13" rod
Egg switch & mount
Extension lead (for egg switch)

Note: On a recline, substitute the headrest mount and armrest clamp to the below.



Swing-Away Clamp Assembly
PRIDEMEC1809294



Headrest Mount
SPTWB485-17MM

Switch Driving

.SCSWITCH
Full Switch Driving



Optional Accessories:

-  **Specs Switch**
AB100-SPBK
-  **Buddy Button Black**
DYNSWTA5700
-  **Buddy Button Red**
DYNSWTA5710
-  **Buddy Button Green***
DYNSWTA5720
-  **Buddy Button Yellow***
DYNSWTA5750
-  **5 Way Splitter***
DYNGSM9976
-  **SCIM**
PRIDECTL143988
-  **Tray**
TRAYPRIDExxD1
-  **Attendant Control**
PRIDECTL143981
PRIDEDWR1010D298
PRIDEFRMASMB10455
PRIDEFRMASMB17115
-  **Attendant Control**
PRIDECTL143981
PRIDEDWR1010D298
PRIDEFRMASMB10455
PRIDEFRMASMB17115



EDU (includes Gooseneck Mount)
.PKEDUQL3
FRMASMB16062



Figure 8 & 13" Rod
SPF8675-RING
SPSU264
SPESM



Full Switch Driving Pack

Includes (as pictured above):
5x Buddy buttons
5 way splitter
Tray
EDU

Headrest Mounted


.SCACCDMPHR
Full Chin Control

Optional Accessories:

-  **IDM-Micro-9
Micro Joystick**
-  **Headrest Mount
17mm Ball
SPUHB4D33200FDM**



**EDU (includes
Gooseneck Mount)**
.PKEDUQL3
FRMASMB16062




Headrest Mount
SPUHB4D33200


Note: On a recline, substitute the headrest mount and armrest clamp to the below.



Metric Headpad
SPACP150M



**Swing-Away
Clamp Assembly**
PRIDEMEC1809294



Headrest Mount
SPTWB485-17MM

Full Chin Control Pack

- Includes:
 TWB480-1 FDM - Flipdown
 Headrest Hardware
 CP250 Comfort Plus Headrest
 SM600 Stealth Mini Proportional
 Joystick Mount
 CPS9-R-HW - Stealth Comfort Plus
 R/H Side Swingaway System for 1"
 Ball

Tray Mounted Control Options

Chair width required _____"

Mount type required Drop In
 Slide On
 Bayonet

Control required _____

Optional Accessories:

-  **EDU (includes
Gooseneck Mount)**
.PKEDUQL3
FRMASMB16062
-  **EDU (includes
Gooseneck Mount)**
.PKEDUQL3
FRMASMB16062

Note: For trial a wooden blank tray will be provided which can be drilled to allow for correct position of controller for trial.

Joystick Controller Options

Stand Alone Joystick Pack



Codes Required:
 PRIDECTL143980
 PRIDEMEC142426*
 .PKEDUQL3
 FRMASMB16062

Mini Proportional Joystick Pack



Codes Required:
 SPIDPMP-9
 .PKEDUQL3
 FRMASMB16062

Stand Alone Joystick Pack



Codes Required:
 PRIDECTL143981
 PRIDEDWR1010D298
 PRIDEFRMASMB10455
 PRIDEFRMASMB17115

Codes Required for Recline:

PRIDECTL143981
 PRIDEDWR1010D298
 PRIDEMEC122827

iAccess Enhanced Controls

Buttons



Codes Required:
 PRIDECTL144009
 PRIDEMEC136981 for TB3 Arms-Mount
 .PKMULTIQL3

Toggle



Codes Required:
 PRIDECTL144010
 PRIDEMEC136981 for TB3 Arms-Mount
 .PKMULTIQL3

Bracket for Round Tube Arms PRIDEMEC1906530*

iAccess Enhanced Controls



The Q-Logic 3 enhanced display unit features a full colour LCD display, an industry standard D-9 connector to suit the integration with any aftermarket joystick or specialty control device, two 3.5 mm jack plugs allowing the fitment of auxiliary switches.

Connectivity is fully catered for with the Q-Logic 3 Enhanced Display Unit.

Bluetooth access to computer, tablet and phone functions including iOS and Android
 Infrared remote control to access devices such as TV, Stereo, Sky TV, home automation

Bib Mounted Chin Control (requires EDU or SCIM)

SPIDCH-M



Stealth iDrive Control Harness Medium
This is a chin control harness and requires a mini proportional joystick.



An orbital compression mount and Stealth's proven highly adjustable, link style hardware, along with the addition of the track system allows for unmatched targeted positioning. This will implement the most optimal driving experience for the user.



With infinite positioning possibilities as well as a variety of switch mounting options that will attach to the collar, the Stealth i Drive Control Harness is the RIGHT choice for any lip, cheek, or chin driver!

IDPMP-9
Stealth Min Proportional Joystick



The quick release, swing away, padded collar allows for the ability to remove the device from the right or the left side. The removable collar design also allows height adjustability on the collar frame independent of the adjustable link hardware.



The double padded collar consists of Closed Cell and SunMate® foams. SunMate foam was originally designed for the NASA space program and provides a gentle touch for extended periods of time.

Stealth iDrive Control Harness Medium
This is a chin control harness and requires a mini proportional joystick.

Alternative Joysticks



Includes QL3 Hand Control w/ Lighter Spring

Includes:
PRIDEKIT1902931*
PRIDEDGN2124307
PRIDECTL143976

Lighter spring 120gm joystick force instead of std 240gm

mo-Vis All Round Lite Joystick



STEALTH MO-VIS ALL ROUND JOYSTICK LITE
SPIDM-ARLITE-9

Lighter spring 120gm joystick force instead of std 240gm



STEALTH O-VIS SATELLITE TWISTER RED

SPIDM-17*
Switches for use with SPIDM-ARLITE-9*



STEALTH MO-VIS SATELLITE TWISTER BLUE

SPIDM-20*
Switches for use with SPIDM-ARLITE-9



mo-Vis Heavy Duty Joystick



STEALTH mo-Vis HEAVY DUTY J/STICK
SPIDM-HD*

Requires an enhanced display or SCIM to function

mo-Vis Heavy Duty Joystick



Stealth Mini Proportional Joystick
IDPMP-9

Requires an enhanced display or SCIM to function




STEALTH IDPMP TOPPER KIT

SPIDPMPJ-TOP
Requires an EDU or SCIM to function

Quantum Q-Logic 3E Hand Control (Does not require EDU or SCIM)

Q-Logic 3e Hand Control



GW to add in STD joystick and outline lighter spring option

Q-Logic 3 Handcontrol e



FEATURES:

- Works with Q-logic 3 system
- Operates up to 5 actuator functions
- Simple LCD display (icon based / Min text)
- Bluetooth and cable programming
- Interactive Assist capability
- Left and Right side charger option
- USB-C charger port
- Light functions operate off soft keys
- Programmable soft keys to offer I-access type functionality.
- Day / Night screen via Page button long press

POWER MODULES:

- 120 Amp Q-logic 3

Stealth Alternate Joysticks – Force Required

GW to add in STD joystick and outline lighter spring option

Part Number	Joystick	Grams of Force	Newtons
IDM-MICRO-9	Micro Joystick	≈8.5gf (0.02lbf)	≈8.5gf (0.08N)
IMPMP-9	PMPJ Mini Proportional Joystick	43gf (0.09lbf)	43gf (0.42N)
IDM-MULTI-9*	Multi Joystick	49.89gf (0.11lbf)	49.89gf (0.49N)
IDM-ARLITE-9	All-Round Joystick Lite	120gf (0.26lbf)	120gf (1.18N)
MJM260 MJM266*	Stealth Mushroom Joystick	227gf (0.50lbf)	227gf (2.23N)
IDM-ALLROUND-9*	All-Round Joystick	250gf (0.55lbf)	250gf (2.45N)
IDM-HD*	All-Round HD Joystick	650gf (1.43lbf)	650gf (6.38N)