

Quantum Edge 3 Powerchair Script



Customer Details

Client			
Assessor			
Organisation			
ACC		MOH	Other
User Weight			

Edge 3 Powerchair Script

For assistance with scripting please contact your **Technical Product Specialist** or call **0800 31 61 81** to be put directly through to them.

Edge 3 Powerchair Script

Expandable QL3 Electronics	Switch fender lights
USB charger (plug into controller)	Smooth Ride Suspension (SRS)
iLevel is standard on a REAC Lift and Tilt system	

Seat Size Width x Depth

Width		Depth	
12	19	12	19
13	20	13	20
14	21	14	21
15	22	15	22
16	23	16	23
17	24	17	24
18		18	

Rigidiser Bar

<input type="checkbox"/> Yes	<input type="checkbox"/> No
------------------------------	-----------------------------

Tyres

<input type="checkbox"/> Pneumatic (Standard)	<input type="checkbox"/> Solid
<input type="checkbox"/> Pneumatic All-Terrain (Knobby)	

Charger

<input type="checkbox"/> Yes	<input type="checkbox"/> No
------------------------------	-----------------------------

Spring Type

TRU-Balance 3	
<input type="checkbox"/> Light <small>(client weight: up to 90kg)</small>	<input type="checkbox"/> Standard <small>(client weight: up to 136kg)</small>
TRU-Balance 4	
<input type="checkbox"/> Standard <small>(client weight: up to 113kg)</small>	<input type="checkbox"/> Lift <small>(client weight: up to 136kg)</small>

Powered Seating

<input type="checkbox"/> Non-Tilt	<input type="checkbox"/> Lift & Tilt
<input type="checkbox"/> Tilt	<input type="checkbox"/> Recline



Anterior Tilt with Articulating Platform* (TB4)

<input type="checkbox"/> 10 deg	
<input type="checkbox"/> 30 deg	
Anterior Tilt required:	<input type="text"/> deg

Colour

<input type="checkbox"/> Back in Black	<input type="checkbox"/> Blue Jean <small>(special order)</small>
--	--



<input type="checkbox"/> Green Machine <small>(special order)</small>
--

<input type="checkbox"/> Evolution Orange <small>(special order)</small>



* If Anterior Tilt with Articulating Platform is selected, the PWC will be fitted with the TRU Balance 4 powered positioning system which includes all Powered Seating Functions. Please note, ADI Backs are best suited for TRU Balance 4 seating.

Seat to Floor Height

<input type="checkbox"/>	17.5"	<input type="checkbox"/>	19.5"	<input type="checkbox"/>	
<input type="checkbox"/>	18.5"	<input type="checkbox"/>	Other	<input type="checkbox"/>	

Controller

Controller Type

<input type="checkbox"/>	Standard
<input type="checkbox"/>	Attendant Only
<input type="checkbox"/>	Other (please specify in notes)

Controller Side

<input type="checkbox"/>	Left hand side
<input type="checkbox"/>	Right hand side
<input type="checkbox"/>	Midline
<input type="checkbox"/>	N/A

Controller Mount

<input type="checkbox"/>	Fixed in-line
<input type="checkbox"/>	Swing away
<input type="checkbox"/>	Retract4
<input type="checkbox"/>	N/A

Legrest Type

<input type="checkbox"/>	Centre Mount
<input type="checkbox"/>	Powered Centre Mount
<input type="checkbox"/>	70 Degree Swing Away
<input type="checkbox"/>	60-90 Degree Swing Away
<input type="checkbox"/>	Manual Elevating Swing Away
<input type="checkbox"/>	Manual Elevating Centre Mount
<input type="checkbox"/>	Powered Swing Away
<input type="checkbox"/>	Paediatric

Stump Support

<input type="checkbox"/>	Left
<input type="checkbox"/>	Right
<input type="checkbox"/>	Both
<input type="checkbox"/>	None

Footplate Type & Size

<input type="checkbox"/>	Standard
--------------------------	----------

Angle Adjustable

<input type="checkbox"/>	6 x 8
<input type="checkbox"/>	8 x 11

Centre Mount

<input type="checkbox"/>	5 x 8
<input type="checkbox"/>	5 x 11

<input type="checkbox"/>	Paediatric
--------------------------	------------

<input type="checkbox"/>	Not required
--------------------------	--------------

Seat Pan to Footplate Length

Left

<input type="checkbox"/>	14
<input type="checkbox"/>	15
<input type="checkbox"/>	16
<input type="checkbox"/>	17
<input type="checkbox"/>	18
<input type="checkbox"/>	19
<input type="checkbox"/>	20
<input type="checkbox"/>	Not required
<input type="checkbox"/>	Other (please specify:)
<input type="checkbox"/>	

Right

<input type="checkbox"/>	14
<input type="checkbox"/>	15
<input type="checkbox"/>	16
<input type="checkbox"/>	17
<input type="checkbox"/>	18
<input type="checkbox"/>	19
<input type="checkbox"/>	20
<input type="checkbox"/>	Not required
<input type="checkbox"/>	Other (please specify:)
<input type="checkbox"/>	

Heel Loops

<input type="checkbox"/>	Yes	<input type="checkbox"/>	No
--------------------------	-----	--------------------------	----

Armrest Type

<input type="checkbox"/>	Two post flip back
<input type="checkbox"/>	Flip-up armrests (not compatible with recline)
<input type="checkbox"/>	Quick height adjust drop-in
<input type="checkbox"/>	Heavy duty
<input type="checkbox"/>	Paediatric
<input type="checkbox"/>	Recline armrests (only compatible with recline)
<input type="checkbox"/>	None

Armrest Length & Type

Length

<input type="checkbox"/>	Desk
<input type="checkbox"/>	Full

Type

<input type="checkbox"/>	Standard
<input type="checkbox"/>	Gel
<input type="checkbox"/>	Channel arms

Hip Belt Type & Size

Size

<input type="checkbox"/>	Small	<input type="checkbox"/>	Medium	<input type="checkbox"/>	Large
--------------------------	-------	--------------------------	--------	--------------------------	-------

Type

<input type="checkbox"/>	Standard
<input type="checkbox"/>	2 Point
<input type="checkbox"/>	4 Point
<input type="checkbox"/>	Other (please specify:)
<input type="checkbox"/>	

Headrest

Stealth

<input type="checkbox"/>	Not required	<input type="checkbox"/>	10" Non-flip down mount
<input type="checkbox"/>		<input type="checkbox"/>	10" Flip down mount
<input type="checkbox"/>		<input type="checkbox"/>	14" Non-flip down mount
<input type="checkbox"/>		<input type="checkbox"/>	14" Flip down mount
<input type="checkbox"/>		<input type="checkbox"/>	Other (please specify:)
<input type="checkbox"/>		<input type="checkbox"/>	

Seating - Cushion

<input type="checkbox"/>	Not Required
<input type="checkbox"/>	Varilite Stratus
<input type="checkbox"/>	Varilite Evolution
<input type="checkbox"/>	Ride Java
<input type="checkbox"/>	Ride Forward
<input type="checkbox"/>	Ki Mobility Axiom P
<input type="checkbox"/>	Ki Mobility Axiom SP Visco
<input type="checkbox"/>	Ki Mobility Axiom SP Fluid
<input type="checkbox"/>	Ottobock Cloud
<input type="checkbox"/>	Other (please specify:)
<input type="checkbox"/>	

Seating - Back Support

<input type="checkbox"/>	Not required
<input type="checkbox"/>	Varilite Icon Mid
<input type="checkbox"/>	Varilite Icon Tall
<input type="checkbox"/>	Varilite Icon Deep
<input type="checkbox"/>	Ride Designs Java Back Regular
<input type="checkbox"/>	Ride Designs Java Back Tall
<input type="checkbox"/>	Ride Designs Java Back Grande
<input type="checkbox"/>	Ride Designs Java Back Grande Deep
<input type="checkbox"/>	Axiom Back
<input type="checkbox"/>	ADI Aluminium Back Active Contour (inc height in notes)
<input type="checkbox"/>	ADI Aluminium Back Deep Contour (inc height in notes)
<input type="checkbox"/>	ADI Q-Back Back Active Contour (inc height in notes)
<input type="checkbox"/>	ADI Q-Back Back Deep Contour (inc height in notes)
<input type="checkbox"/>	ADI Flat Top Back Active Contour (inc height in notes)
<input type="checkbox"/>	ADI Flat Top Back Deep Contour (inc height in notes)
<input type="checkbox"/>	ADI Carbon Fibre Back Active Contour
<input type="checkbox"/>	Other (please specify:)

Seating - Laterals (Please Specify)

<input type="checkbox"/>	Not required	<input type="checkbox"/>	Thoracic	<input type="checkbox"/>	Pelvic
--------------------------	--------------	--------------------------	----------	--------------------------	--------

--

Tray

<input type="checkbox"/>	Swing Away	<input type="checkbox"/>	Custom
<input type="checkbox"/>	Drop In	<input type="checkbox"/>	None

Accessories

--

Set-up Instructions & Notes

--