

# karma FLEXX

Specifically designed for functionality with multiple adjustments to suit individual needs and body conditions



ALLIED  
MEDICAL  
LIMITED



Extremely flexible to suit varying needs

## + FLEXIBLE, ROBUST AND STYLISH

The Karma Flexx is designed with plenty of adjustments to create flexible sitting to fit individual needs. With ellipse-shaped frame design, the Flexx is not only stylish but also robust and durable.

The configurable options of the Flexx include variable seat depths, seat to floor heights, backrest heights and backrest angles, castor fork angles and centre of gravity. It can easily be set up to provide optimal comfort and performance.

The bright colour of key adjustment points are self-explanatory; users can easily operate the following functions independently: swing away footrests, flip-back and adjustable armrests.

## + FEATURES

- Configurable dimensions
- Height adjustable flip-back Ergo armrests
- Swing away footrests with flip-up footplates
- Padded height and angle adjustable backrest
- 25 possible rear wheel positions
- Reinforced double crossbrace
- Quick release pneumatic rear wheels
- Push to Lock Brakes with Brake Extension
- Available as a transit wheelchair



## + ACCESSORIES

- Swing-away height adjustable anti-tippers
- Positioning hip belt
- Lap tray
- Cane holder



## + SPECIFICATIONS



**Seat widths**  
14" (35.5 cm) to  
24" (60.9 cm)




**Seat depth**  
16" (40 cm) to  
20" (50.8 cm)




**Wheels**  
6" - 8" (front)  
14", 22"-24" (rear)



**Seat to floor height**  
15.5" - 21"  
(39 - 53 cm)



**Overall width**  
23" (58 cm) to  
28" (71 cm)



**Folded width**  
13.5" (34 cm)



**Back height**  
15" (38 cm) to  
21" (53 cm)



**Overall length**  
40" (102 cm) to  
42" (107 cm)



**Overall height**  
30" (76 cm) to  
42" (106 cm)



**Total weight**  
13.6 kg (14") -  
14.1 kg (24")



**User weight limit**  
130 kg  
(170kg for HD)



**ALLIED  
MEDICAL  
LIMITED**

**0800 31 61 81**

[www.alliedmedical.co.nz](http://www.alliedmedical.co.nz)