

# Karma Flexx Wheelchair Script



Customer Details							
Client							
Assessor							
Organisation							
User Weight	ACC			MOH		Other	

## Karma Flexx Wheelchair

Please note: Scripts are not required for Karma Flexx wheelchairs with the standard specifications (listed below). This script should only be used for clients with alternative requirements.

### Standard Specifications:

- 7" Castor Wheels
- 19" Seat to Floor Height
- 1" Seat Dump
- Height Adjustable Ergo Flipback Arms
- 80° Swing-away Leg supports
- Flip-up Foot Supports
- Push to Lock Wheel Locks
- Swing Away Anti-Tippers
- Calf Strap
- Standard Hip Belt
- Seat widths 14-18": 130kg max user weight
- Seat widths 20"+: 170kg max user weight

## Wheelchair Type

Self Propelling

- Quick Release 24" Pneumatic Spoked Rear Wheel
- Ergo Handrim

Transit

- Fixed 14" Solid Tyre Plastic Rim Rear Wheels

## Seat Size (Width x Depth)

14 x 16	20 x 18 <sup>^</sup>
16 x 16	20 x 20 <sup>^</sup>
16 x 18	22 x 18 <sup>^</sup>
16 x 20	22 x 20 <sup>^</sup>
18 x 16	24 x 18 <sup>^</sup>
18 x 18	24 x 20 <sup>^</sup>
18 x 20	

<sup>^</sup> Heavy Duty Model

## Front Seat Height

		Seat to Floor Height					
Castor Size	1) 5"		16"	17"	18"		
	6"			17"	18"	19"	
	7"			17.5"	18"	19"	
	8"					19"	20" 21"

<sup>1)</sup> 5" Castor not available on models 20" and wider

## Back Cane Height

		Back Height					
Fixed Height	Std			17"			*
Height Adjustable (Adjustable with tools only)	\$	15"	16"	17"	18"	19"	20"

\* Standard on 16 x 20 and 18 x 20 models

## Rear Seat Height

		Rear Seat to Floor Height					
Rear Wheel Size	14"					18"	
	20"	14"	15"	16"	17"	18"	
	22"		15"	16"	17"	18"	19"
	24"			16"	17"	18"	19" 20"

## Anti Tippers

Swing Away Height Adjustable

Front Seat to Floor Height 18-21".....Std
Front Seat to Floor Height 15-17".....\$

Pneumatic tyres       Solid tyres

## Centre of Gravity

Rear       Centre       Forward

**Seat Upholstery**

Standard (Padded) .....	Std
Standard with velcro* .....	\$
Omit Seat Upholstery .....	N/C
Cushion (please specify below) .....	\$
Spare Cushion Cover .....	\$
Cushion	

\*Available on most sizes. Some sizes may be indent only.

**Transit Option**

Transit Tie Down Option .....	\$
-------------------------------	----

\* Not available on HD

**Leg Support**

	Left	Right	
80° Swing-away .....			Std
90° Swing-away .....			\$
Elevating .....			\$
Stump Support .....			\$

**Foot Support**

	Left	Right	
Standard Flip-up .....			Std
Flip-up Angle Adjustable* .....			\$
Flip-up Angle & Depth Adj. Small..			\$
Large (16" or wider chairs only)....			\$

\* Not available on HD

**Calf Strap/Heel Loops**

Calf Strap .....	Std
Heel Loops* .....	\$

\*Standard on HD with Flip-up Foot Supports only, otherwise an upcharge.

**Other Parts Required**

**Back Support**

Standard (Padded) .....	Std
Tension Adjustable (not available on height adj back canes) .....	\$
Omit Back Upholstery .....	N/C
Alternative Back Support (please specify below) .....	\$
Alternative Back Support	

**Arm Support**

Flip-back Ergo (desk length) Height Adjustable .....	Std
Flip-back Boat (full length) Height Adjustable .....	\$

**Wheel Locks**

Push to Lock with Brake extension .....	Std
Pull to Lock* .....	\$

\* Available only on 20, 22 and 24" Rear wheels

**Accessories**

Clear Spoke Guards .....	\$
Oxygen Tank Holder .....	\$
Push Handle Rigidiser Bar* .....	\$
IV Pole .....	\$

\* Std on height adj back canes 18" and wider and HD chairs

**Head Support**

Foldable Medium (16-18" widths) .....	\$
Foldable Large (20" widths) .....	\$
Other .....	\$

**Pelvic Positioning Belt**

Lap Belt .....	Std
2 Point .....	\$
4 Point .....	\$

**Tech Notes**