

# AXIOM

for Kids





# Inspired by Kids

Axiom's pediatric seating line is designed for the unique needs of our smallest clients. Growing bodies have different positioning needs, come in various shapes, and simply weigh less. This is precisely what inspired the **Axiom for Kids** line - to provide kids with better seating products built specifically for them.



## Kids come in all shapes and sizes


Axiom for Kids cushions deliver a better fit for optimal comfort and outcomes by using varying foam softness and shapes by cushion size.

Available in two models: **PXP** ( Positioning ) and **PXSP** ( Skin Protection & Positioning ).



## Accidents happen

Made from high-quality materials, our covers are durable, wipeable, breathable, and easy to clean. Just wipe and continue with your day! Extra cover included with each cushion order.



## Covers in color

Choose from 3 color options to complement the mobility base and the personal style and preferences of the smallest clients. Available in Blue, Purple and Black.

# Let's grow together

Designed with growing kids in mind, the Axiom for Kids back offering provides depth and width adjustment through the back shell, laterals, and mounting hardware. With contour shapes that change with size, our backs provide the confidence of a perfect fit as your child grows. The vapor permeable cover promotes air flow to ensure added comfort.

Photo features the PL Back



AL Back with contoured laterals

## Get the best fit

Setup and growth is easy with intuitive depth and angle adjustments, along with the ability to fit chair widths +/- 1". With three models and three lateral options - Axiom for Kids is the most versatile range of pediatric back solutions on the market.



DEX for Kids

## Hardware that multitasks

DEX for Kids hardware works with all tube diameters and utilizes the DEX 1" locking clamp, providing reliable clamp force, and minimizing back cane interference so you can easily mount the back exactly where you need it.

DEX for Kids comes with a standard lockout feature so curious hands, classmates or backpacks don't accidentally disengage the release.

## Cushion Specs

Description  
 HCPCS Code ( US Only )  
 Size Offering  
 Extra Cover Included  
 Wipeable Cover  
 Color Availability  
 Max User Weight

<b>PXP</b> Positioning		<b>PXSP</b> Skin Protection & Positioning	
Single density contoured foam		Dual density contoured foam	
E2605		E2607	
Width: 8"-16" ( 20 - 40.5cm ) Depth: 8" - 18" ( 20 - 45.5cm )			
Yes			
Yes			
Black, Purple, Blue			
165 lbs. ( 75kg )			

## Back Specs

Description  
 HCPCS Code ( US Only )  
 Size Offering  
 Back Shell Color Availability  
 Attaching Back Hardware  
 Back Adjustability  
 Fits Chair Width  
 Max User Weight  
 Lateral Adjustability  
 Configured Weight  
 Transit  
 Tested Tough

<b>PXPL</b> Posterior Lateral		<b>PXAL</b> Adjustable Lateral		<b>PXALCL</b> Adjustable Contoured Laterals	
Formed back with lateral and lumbar supports.		Planar back with option for fixed or swing-away lateral supports.		Contour back with adjustable contoured lateral supports.	
E2615		E2291		E2293	
Width: 10"-16" ( 23 - 40.5cm )   Depth: 10" - 18" ( 23 - 45.5cm )					
20+ Colors ( see order form for complete offering )					
DEX for Kids ( includes safety lockout ), fits 3/4" to 1-1/8" tube diameter ( no adaptors required )					
Angle: +/- 20 degrees   Depth: 2" standard / 4" with optional extension					
Chair Width +/- 1					
165 lbs. ( 75kg )					
No		+2" / -1" Lateral Adjustment			
12" wide LW Clik wheelchair with Axiom for Kids PL Back - 15.8 lbs. ( 7.2kg )					
Tested according to the dynamic test requirements of ANSI/RESNA WC-4, Section 20					
Tested to ISO-16840 for strength, durability, and reliability					

