

Ki Mobility Rogue XP Express Script



Customer Details							
Client							
Assessor							
Organisation							
	ACC			MOH			Other
User Weight							

Ki Mobility Rogue XP Express Script

If you only require specifications that fall under the Standard Express Package, you do not need to fill in the reverse page.

		Seat Width		
		12"	14"	16"
Seat Depth	12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	13"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	14"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	15"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	16"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	17"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	18"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	19"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Express Codes

- 12W** KIROGUEXP1212STD
20in Maxx Performance Spoke Wheels, 4x1.5 Soft Roll Castors, Height Adjustable T-Arms (7.5-11.5 inch)
- 14W** KIROGUEXP1414STD
22in Maxx Performance Spoke Wheels, 4x1.5 Soft Roll Castors, Height Adjustable T-Arms (8.5-13.5 inch)
- 16W** KIROGUEXP1616STD
22in Maxx Performance Spoke Wheels, 5x1.5 Soft Roll Castors, Height Adjustable T-Arms (8.5-13.5 inch)

Add ROGUEEXPPACKAGE to all orders

Express Package: (Select optional specifications and alternatives on reverse)

- Seat Sling
- Transit Tie Downs
- Depth Adjustment of 3in
- 85 degree front frame bend
- Sling position 1"
- No Taper
- Height adjustable Medium straight back posts with no push handle
- Wheel camber of 2 degrees
- Height adjustable rigidising bar
- Padded Velcro adjust. back upholstery
- Standard fork
- Angle Adjustable Flip Under Footplates
- Street Pneumatic Tyres Black
- Aluminium handrims
- Quick Release axles
- Push to lock wheel locks
- T-Arms with Standard Desk ampads
- Rear Anti Tipppers
- Iridescent Navy Frame
- Double Stroller Handle
- Spoke Guards

Please note: Express Package includes some upcharge products.

Choose Your Specifications (No additional codes are required on the below parts)

Backrest Height						Rear Seat Height					
<input type="checkbox"/> 13"	<input type="checkbox"/> 14"	<input type="checkbox"/> 15"	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	16W: <input type="checkbox"/> 13"	<input type="checkbox"/> 13.5"	<input type="checkbox"/> 14"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
(For 16"-18", .ULROGXPBUP1*)						14W: <input type="checkbox"/> 14"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/>
Backrest Angle						16W: <input type="checkbox"/> 14"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/>
<input type="checkbox"/> Std	<input type="checkbox"/> +2°	<input type="checkbox"/> +4°	<input type="checkbox"/> +6°	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/> +8°	<input type="checkbox"/> +10°	<input type="checkbox"/> -2°	<input type="checkbox"/> -4°	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/> -6°	<input type="checkbox"/> -8°	<input type="checkbox"/> -10°	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Front Seat Height						Seat to Footrest Length					
12W: <input type="checkbox"/> 13.5"	<input type="checkbox"/> 14"	<input type="checkbox"/> 14.5"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	12W: <input type="checkbox"/> 10.5"	<input type="checkbox"/> 11"	<input type="checkbox"/> 11.5"	<input type="checkbox"/> 12"	<input type="checkbox"/>	<input type="checkbox"/>
14W: <input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	14W: <input type="checkbox"/> 11.5"	<input type="checkbox"/> 12"	<input type="checkbox"/> 12.5"	<input type="checkbox"/> 13"	<input type="checkbox"/>	<input type="checkbox"/>
16W: <input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16W: <input type="checkbox"/> 11"	<input type="checkbox"/> 11.5"	<input type="checkbox"/> 12"	<input type="checkbox"/> 12.5"	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> 13"	<input type="checkbox"/> 13.5"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Optional Specifications

Backrest Upholstery

<input type="checkbox"/>	Omit back upholstery						
<input type="checkbox"/>	Tension Adj. (.ULROGXPBUPTA1*)						
<input type="checkbox"/>	13"	<input type="checkbox"/>	14"	<input type="checkbox"/>	15"	<input type="checkbox"/>	16"
<input type="checkbox"/>	17"	<input type="checkbox"/>	18"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Backrest Options

<input type="checkbox"/>	Std Push Handle (.ULROGXPBCKSTD)
<input type="checkbox"/>	Ergo Strollerhandle (.ULROGXPBCKERG****)

Footrest for shorter leg length

<input type="checkbox"/>	High Mount Alum. Ang. Adj. Footplates						
	6"	7"	8"	9"	10"	11"	
12W	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
14W	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
16W	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

(.ULROGXPFPLHM1*)

Camber

<input type="checkbox"/>	0°	<input type="checkbox"/>	4°
--------------------------	----	--------------------------	----

(.ULROGXPCAM0) (.ULROGXPCAM4)

Centre of Gravity

Amputee Positions

<input type="checkbox"/>	-1.5"	<input type="checkbox"/>	-1.25"	<input type="checkbox"/>	-1"	<input type="checkbox"/>	-0.75"
<input type="checkbox"/>	-0.5"	<input type="checkbox"/>	-0.25"	<input type="checkbox"/>	0"	<input type="checkbox"/>	0.25"

Standard Positions

<input type="checkbox"/>	0.5"	<input type="checkbox"/>	0.75"	<input type="checkbox"/>	1"	<input type="checkbox"/>	1.25"
<input type="checkbox"/>	1.5"	<input type="checkbox"/>	1.75"	<input type="checkbox"/>	2"	<input type="checkbox"/>	2.25"
<input type="checkbox"/>	2.5"	<input type="checkbox"/>	2.75"	<input type="checkbox"/>	3"	<input type="checkbox"/>	3.25"
<input type="checkbox"/>	3.5"	<input type="checkbox"/>	3.75"	<input type="checkbox"/>	4"	<input type="checkbox"/>	4.25"

Side Guards (not available with T-Arms)

<input type="checkbox"/>	Composite Clamp on Paed (.ULROGXPSDGCP)
--------------------------	---

Wheel Locks

<input type="checkbox"/>	Low Profile Scissor Locks (.ULROGXPWKLPS)
--------------------------	---

Accessories

<input type="checkbox"/>	Spoke Guard Clear 20 IN w/Velcro (KI501412)
<input type="checkbox"/>	Spoke Guard Clear 22 IN w/Velcro (KI501413)

<input type="checkbox"/>	Spoke Guard Tuxedo Time 20 IN (KI102703)
<input type="checkbox"/>	Spoke Guard Tuxedo Time 22 IN (KI102704)

<input type="checkbox"/>	Spoke Guard Fairy Garden 20 IN (KI102707)
<input type="checkbox"/>	Spoke Guard Fairy Garden 22 IN (KI102708)

Bodypoint Belt Options:

<input type="checkbox"/>	Bodypoint 2pt belt – Centre Pull Flat Mount (12 and 14: BPHB205-M36-A1) (16: BPHB205-M46-A1) + (BPHW200-90-2)
<input type="checkbox"/>	Bodypoint 4pt belt Centre Pull Flat Mount (12 and 14: BPHB405-M36-A1) (16: BPHB405-M46-A1) + (BPHW200-90-2)
<input type="checkbox"/>	Bodypoint Evoflex Pelvic Stabiliser (12 and 14: BPEB235-S38) (16: BPEB205-M46) + (BPHW200-90-2)
<input type="checkbox"/>	Bodypoint Aeromesh Calf Strap (12 and 14: BPSP102S) (16: BPSP102M)

Notes