



















# LECKEY BeMe

Support as Individual as You Are



Customer Details					
Therapist Name					
Therapist Address					
Client Name					
Client Address					
Size Specific					
Small					
Medium					
Large					
BeMe Seat Core					
Code	Description	Quantity	Upholstery	Available for Trial	
150-1601	BeMe Small Seat Core		N/A	✓	 Seat Core
150-2601	BeMe Medium Seat Core		N/A	✓	
150-3601	BeMe Large Seat Core		N/A	✓	
High-Low Chassis Base (Must select one, match seat core size)					
Code	Description	Quantity	Upholstery	Available for Trial	
150-1800-G	Small Manual (Gas Spring) High-Low Chassis		N/A	✓	 Small Chassis
122-2710	Push Handle for Small and Medium Chassis		N/A	✓	
150-2753-G	Medium Manual (Gas Spring) High-Low Chassis		N/A	✓	 Medium Indoor Chassis
150-2766-G	Medium Pump Action (Hydraulic) High-Low Chassis		N/A	✓	
150-2768-G	Medium Pump Action (Hydraulic) Outdoor High-Low Chassis		N/A	✓	
122-2710	Push Handle for Small and Medium Chassis		N/A	✓	 Medium Outdoor Chassis
150-3770-G	Large Pump Action (Hydraulic) Outdoor High-Low Chassis including Push Handle		N/A		
150-3780-G	Large Electric Outdoor High-Low Chassis including Push Handle		N/A		 Large Chassis
150-818	Interface for Small BeMe Seat Core onto Medium Chassis		N/A		
Seat Hardware					
Code	Description	Quantity	Upholstery	Available for Trial	
150-1637	Small Basic Seat		Required	✓	 Basic Seat
150-1638	Small Contoured Seat		Required	✓	
150-2637	Medium Basic Seat		Required	✓	
150-2639	Medium Split Seat		Required	✓	 Split Seat
150-2638	Medium Contoured Seat		Required	✓	
150-4637	Medium-Long Basic Seat		Required		
150-4639	Medium-Long Split Seat		Required		
150-4638	Medium-Long Contoured Seat		Required		
150-3637	Large Basic Seat		Required	✓	 Contoured Seat
150-3639	Large Split Seat		Required	✓	
150-3638	Large Contoured Seat		Required	✓	
150-5637	Large-Long Basic Seat		Required		
150-5639	Large-Long Split Seat		Required		
150-5638	Large-Long Contoured Seat		Required		
Pelvic Belt					
Code	Description	Quantity	Upholstery	Available for Trial	
150-1844	4-point belt Size 1		N/A	✓	 4 Point
150-2844	4-point belt Size 2		N/A	✓	
150-3844	4-point belt Size 3		N/A	✓	
150-4844	4-point belt Size 4		N/A	✓	
117-476	Pelvic Cradle size 2		N/A	✓	

117-477	Pelvic Cradle size 3		N/A	✓		
117-478	Pelvic Cradle size 4		N/A	✓		
117-479	Pelvic Cradle size 5		N/A	✓		
117-480	Pelvic Cradle size 6		N/A	✓		
<b>Femoral Guide (Select up to 2 on Basic seat or 4 on Split Seat)</b>						
<b>Code</b>	<b>Description</b>	<b>Quantity</b>	<b>Upholstery</b>	<b>Available for Trial</b>		
150-1615	Single Femoral Gable - Small		Optional	✓		
150-2615	Single Femoral Gable - Medium		Optional	✓		
150-3615	Single Femoral Gable - Large		Optional	✓		
150-4615	Single Femoral Gable - Extra Large		Optional	✓		
150-1862	Padded Cover for Single Femoral Gable - Small		N/A	✓		
150-2862	Padded Cover for Single Femoral Gable - Medium		N/A	✓		
150-3862	Padded Cover for Single Femoral Gable - Large		N/A	✓		
150-4862	Padded Cover for Single Femoral Gable - Extra Large		N/A	✓		
<b>Pommel (Match seat core size)</b>						
<b>Code</b>	<b>Description</b>	<b>Quantity</b>	<b>Upholstery</b>	<b>Available for Trial</b>		
150-1624	Pommel Small		Optional	✓		
150-1864	Padded Cover - Pommel Small		N/A	✓		
150-2624	Pommel Medium		Optional	✓		
150-2864	Padded Cover - Pommel Medium		N/A	✓		
150-3624	Pommel Large		Optional	✓		
150-3864	Padded Cover - Pommel Large		N/A	✓		
<b>Backrest (does not need to match seat core size)</b>						
<b>Code</b>	<b>Description</b>	<b>Quantity</b>	<b>Upholstery</b>	<b>Available for Trial</b>		
150-1625	Small Basic Backrest		Required	✓		
150-1630	Small Moderate Backrest		Required	✓		
150-2625	Medium Basic Backrest		Required	✓		
150-2630	Medium Moderate Backrest		Required	✓		
150-2633	Shoulder Wing Upgrade for Medium Moderate (Requires Complex 3D Backrest Upholstery)		Required			
150-2631	Medium Complex 3D Backrest		Required			
150-3625	Large Basic Backrest		Required	✓		
150-3630	Large Moderate Backrest		Required	✓		
150-3633	Shoulder Wing Upgrade for Large Moderate (Requires Complex 3D Backrest Upholstery)		Required			
150-3631	Large Complex 3D Backrest		Required			
150-1792	15 degree angle spacer for Moderate shoulder section		N/A			
150-605	Dynamics Upgrade Kit		N/A			
<b>Hip Guides</b>						
<b>Code</b>	<b>Description</b>	<b>Quantity</b>	<b>Upholstery</b>	<b>Available for Trial</b>		
150-1684	Small Hip Guide (Pair)		Included	✓		
150-2684	Medium Hip Guide (Pair)		Included	✓		
150-3684	Large Hip Guide (Pair)		Included	✓		
<b>Hip Guides - Additional Longer Pads (Must Order Hip Guides)</b>						
<b>Code</b>	<b>Description</b>	<b>Quantity</b>	<b>Upholstery</b>	<b>Available for Trial</b>		
150-1653	Hinged Lateral Pad Small (Pair)		Included			
150-3653	Hinged Lateral Pad Large (Pair)		Included			
<b>Chest Harness</b>						
<b>Code</b>	<b>Description</b>	<b>Quantity</b>	<b>Upholstery</b>	<b>Available for Trial</b>		
155-N001-L600	Small Dynamic Butterfly Harness		N/A	✓		
155-N001-L601	Medium Dynamic Butterfly Harness		N/A	✓		
155-N001-L602	Large Dynamic Butterfly Harness		N/A	✓		
155-N001-L603	Extra Large Dynamic Butterfly Harness		N/A	✓		

Chest Lateral Bracket (match backrest size)					
Code	Description	Quantity	Upholstery	Available for Trial	
150-1651	Rigid Lateral Bracket for Size Small/Medium Backrests(Pair)		Included		 Rigid
150-1655	FlipAway Lateral Bracket for Small/Medium Backrests (Pair)		Included	✓	
150-3651	Rigid Lateral Bracket for Large Backrests (Pair)		Included		 Flipaway
150-3655	FlipAway Lateral Bracket for Large Backrests (Pair)		Included	✓	
Chest Lateral Pad					
Code	Description	Quantity	Upholstery	Available for Trial	
150-1675	Size 1 Slimline Lateral Pads (Pair)		Included		 Slimline
150-2675	Size 2 Slimline Lateral Pads (Pair)		Included		
150-3675	Size 3 Slimline Lateral Pads (Pair)		Included		 Wraparound Hugging
150-1260	Small Wraparound Hugging Lateral Pads (Pair)		Included	✓	
150-2260	Medium Wraparound Hugging Lateral Pads (Pair)		Included	✓	 Hinged
150-1653	Hinged Lateral Pad Small (Pair)		Included	✓	
150-3653	Hinged Lateral Pad Large (Pair)		Included	✓	
150-673	Lateral Protraction Pads (Pair)		Included	✓	 Protraction Pads
Armrest (Complete, includes armrest pads)					
Code	Description	Quantity	Upholstery	Available for Trial	
150-1642	Small Armrests (Pair)		N/A	✓	
150-2642	Large Armrests (Pair)		N/A	✓	
Armrest Pads (Extra armrest pads)					
Code	Description	Quantity	Upholstery	Available for Trial	
150-15204	Small Armrest Pad Only (Single)		N/A		
150-25204	Large Armrest Pad Only (Single)		N/A		
Tray (match seat core size)					
Code	Description	Quantity	Upholstery	Available for Trial	
150-1720	Size 1 Tray		N/A	✓	 Tray
150-1888	Size 1 Padded Tray Cover		N/A		
150-1671	Tray Elbow Blockers for Size 1		N/A		 Padded Cover
150-2720	Size 2 Tray		N/A	✓	
150-2888	Size 2 Padded Tray Cover		N/A		
117-641	Tray Elbow Blockers for Size 2		N/A		 Elbow Blockers
150-3720	Size 3 Tray		N/A	✓	
150-3888	Size 3 Padded Tray Cover		N/A		
150-3671	Tray Elbow Blockers for Size 3		N/A		
Leg Support					
Code	Description	Quantity	Upholstery	Available for Trial	
150-1689	Split Footplate Small (Left and Right)		N/A	✓	 Two Piece Footplate
150-2689	Split Footplate Medium (Left and Right)		N/A	✓	
150-3689	Split Footplate Large (Left and Right)		N/A	✓	
117-787	1Inch Sandal Raisers (pair) for split footplate		N/A		 One Piece Footplate
150-1698	One Piece Footplate - Small		Optional	✓	
150-2698	One Piece Footplate - Medium		Optional	✓	
150-3698	One Piece Footplate - Large		Optional	✓	 One Piece Cushion
150-1872	Cushion for Small One-Piece Footplate		N/A	✓	
150-2872	Cushion for Medium One-Piece Footplate		N/A	✓	
150-3872	Cushion for Large One-Piece Footplate		N/A	✓	
Foot Support					
Code	Description	Quantity	Upholstery	Available for Trial	
AH01	Ankle Huggers		N/A		 Ankle Huggers
150-1852	Extra Small Booties		N/A		
152-1600	Sandals Small		N/A	✓	 Extra Small Booties
152-2600	Sandals Medium		N/A	✓	

152-3600	Sandals Large		N/A	✓	
173-2305	Small Ankle Wraps (Only with Sandals)		N/A		
173-3305	Large Ankle Wraps (Only with Sandals)		N/A		
<b>BeMe Headrest</b>					
<b>Code</b>	<b>Description</b>	<b>Quantity</b>	<b>Upholstery</b>	<b>Available for Trial</b>	
150-676	BeMe Headrest		Required	✓	
<b>Universal Headrest Bracket</b>					
<b>Code</b>	<b>Description</b>	<b>Quantity</b>	<b>Upholstery</b>	<b>Available for Trial</b>	
150-781	Universal Headrest Bracket (Whitmyer PRO or Leckey Headrests)		N/A	✓	
<b>Leckey Headrest</b>					
<b>Code</b>	<b>Description</b>	<b>Quantity</b>	<b>Upholstery</b>	<b>Available for Trial</b>	
137-699-05	Leckey Contoured Headrest		Included		
117-699-05	Leckey Contoured Occipital Headrest		Included		
117-851-05	Leckey Flat Headrest		Included	✓	
120-893-05	Leckey Flat Headrest Laterals		Included		
<b>Whitmyer Plush</b>					
<b>Code</b>	<b>Description</b>	<b>Quantity</b>	<b>Upholstery</b>	<b>Available for Trial</b>	
PRO-P1W	Plush Headrest, 10 Inch with 8 Inch PRO Bar		Included	✓	
PRO-P2W	Plush Headrest, 14 Inch with 8 Inch PRO Bar		Included	✓	
PRO-P4W	Plush Headrest, 19 Inch with 8 Inch PRO Bar		Included	✓	
PRO-AP1L	Adjust-A-Plush Narrow Headrest with 8 Inch PRO Bar		Included	✓	
PRO-AP2L	Adjust-A-Plush Standard Headrest with 8 Inch PRO Bar		Included	✓	
<b>Whitmyer Contoured Cradle</b>					
<b>Code</b>	<b>Description</b>	<b>Quantity</b>	<b>Upholstery</b>	<b>Available for Trial</b>	
PRO-CCS2	Small Contoured Cradle Head on 8 Inch PRO Bar		Included	✓	
PRO-CCS3	Medium Contoured Cradle Head on Pro 8 Inch PRO Bar		Included	✓	
PRO-CCS4	Large Contoured Cradle Head on 8 Inch PRO Bar		Included		
<b>Whitmyer Heads Up</b>					
<b>Code</b>	<b>Description</b>	<b>Quantity</b>	<b>Upholstery</b>	<b>Available for Trial</b>	
PRO-HU-ISSCL-O5	Heads Up Short Infant (Size 1 - 8.3 Cm) including chin prompt		Included	✓	
PRO-HU-SSSCL-O6	Heads Up Short Child (Size 2 - 10.1 Cm) including chin prompt		Included	✓	
PRO-HU-MSSCL-O8	Heads Up Short Teenager (Size 3 - 11.6 Cm) including chin prompt		Included	✓	
PRO-HU-LSSCL-O10	Heads Up Short Adult (Size 4 - 14.7 Cm) including chin prompt		Included	✓	
PRO-HU-ITSCL	Heads Up Tall Infant (Size 1 - 8.3 Cm) including chin prompt		Included	✓	
PRO-HU-STSCCL	Heads Up Tall Child (Size 2 - 10.1 Cm) including chin prompt		Included	✓	
PRO-HU-MTSCCL	Heads Up Tall Teenager (Size 3 - 11.6 Cm) including chin prompt		Included	✓	
PRO-HU-LTSCCL	Heads Up Tall Adult (Size 4 - 14.7 Cm) including chin prompt		Included	✓	
<b>Whitmyer Complex</b>					
<b>Code</b>	<b>Description</b>	<b>Quantity</b>	<b>Upholstery</b>	<b>Available for Trial</b>	
PRO-3W	Paediatric S.O.F.T. Three Pad Head Support with 8" PROvertical mounting bar		Included		
PRO-4W	Small Adult S.O.F.T. Three Pad Head Support with 8" PROvertical mounting bar		Included		
PRO-1W	Adult S.O.F.T. Three Pad Head Support with 8" PROvertical mounting bar		Included		
<b>Whitmyer Forehead Support</b>					
<b>Code</b>	<b>Description</b>	<b>Quantity</b>	<b>Upholstery</b>	<b>Available for Trial</b>	
D1010	Whitmyer Bracket for Cradle headrest to attach forehead support		Included		
DFS/3	DYNAMIC FOREHEAD ASSEMBLY - Paediatric		Included		
DFS/1	DYNAMIC FOREHEAD ASSEMBLY - Adult		Included		
D1017	DYNAMIC FOREHEAD STRAP (NEOPRENE AGAINST SKIN)		Included		
D1018	DYNAMIC FOREHEAD STRAP (TOWELLING AGAINST SKIN)		Included		
<b>Whitmyer Bars and Accessories</b>					
<b>Code</b>	<b>Description</b>	<b>Quantity</b>	<b>Upholstery</b>	<b>Available for Trial</b>	
PRC7030	Whit 8" vertical/horizontal PRO mount		Included	✓	
PRC7070	PRO 63mm Offset Bar		Included		
M9450	Dynamic Headrest Kit		Included		
LFC25L	Single Lateral Facial Support for 3 Pad System		Included		
LFC25-CSL	Single Lateral Facial Support for Contoured Cradle		Included		
LFC25-PLL	Single Lateral Facial Support for Plush		Included		
M5200	Swingaway System for headrests		Included		
<b>Miscellaneous</b>					
<b>Code</b>	<b>Description</b>	<b>Quantity</b>	<b>Upholstery</b>	<b>Available for Trial</b>	

150-1645	Sun/Rain Canopy for Small Backrest		N/A	✓	
150-2645	Sun/Rain Canopy for Medium Backrest		N/A	✓	
150-1790	Spacer Pair to decrease minimum seat depth on BeMe Small - 25mm		N/A		
150-1791	Spacer Pair to decrease minimum seat depth on BeMe Small - 50mm		N/A		
150-604	Static Plate to revert dynamic backrest to standard		N/A		
150-699	Van Trolley Set (3 Trolleys + Backpack for Assessments)		N/A		
150-846	Extra Padded Sleeve for Footplate Stem or Lateral Brackets		N/A		
Colours					
Blue - Replace XX with 02					
Purple - Replace XX with 04					
Red - Replace XX with 09					
Grey - Replace XX with 08					
Upholstery					
Code	Description	Quantity	Colour Code	Available for Trial	
150-855-XX	BeMe Headrest Cushion and Cover		02	✓	
150-1840-XX	Small Basic Seat Cushion and Cover		02	✓	
150-1860-XX	Small Contoured Seat Cushion and Cover		02	✓	
150-2840-XX	Medium Basic Seat Cushion and Cover		02	✓	
150-2850-XX	Medium Split Seat Cushion and Cover		02	✓	
150-2860-XX	Medium Contoured Seat Cushion and Cover		02	✓	
150-4840-XX	Medium-Long Basic Seat Cushion and Cover				
150-4850-XX	Medium-Long Split Seat Cushion and Cover				
150-4860-XX	Medium-Long Contoured Seat Cushion and Cover				
150-3840-XX	Large Basic Seat Cushion and Cover		02	✓	
150-3850-XX	Large Split Seat Cushion and Cover		02	✓	
150-3860-XX	Large Contoured Seat Cushion and Cover		02	✓	
150-5840-XX	Large-Long Basic Seat Cushion and Cover				
150-5850-XX	Large-Long Split Seat Cushion and Cover				
150-5860-XX	Large-Long Contoured Seat Cushion and Cover				
150-1870-XX	Small Basic Backrest Cover		02	✓	
150-1800-XX	Small Moderate Backrest Cushions and Covers		02	✓	
150-2870-XX	Medium Basic Backrest Cover		02	✓	
150-2800-XX	Medium Moderate Backrest Cushions and Covers		02	✓	
150-4800-XX	Medium Complex Backrest Cushions and Covers		02	✓	
150-3870-XX	Large Basic Backrest Cover		02	✓	
150-3800-XX	Large Moderate Backrest Cushions and Covers		02	✓	
150-5800-XX	Large Complex Backrest Cushions and Covers		02	✓	

### Technical Measurements

Seat					
	Small	Medium	Large	Medium Long	Large Long
Seat Depth (mm)	170*-325	275-390	350-475	325-440	400-525
Seat Width (mm)	275	325	415	325	415
Max User Weight (kg)**	35	50	70	50	70
Backrest					
	Small	Medium	Large		
Back Height (mm)	360-470	390-520	500-650		
Back Height with Contoured Seat (mm)	370-480	400-530	475-625		
Lateral Width	130-230	200-300	225-350		
Hip Guide Width (mm)	165-230	200-260	250-350		
Backrest Angle	+10 to -25	+10 to -25	+10 to -25		
Max User Weight (kg)**	35	50	70		
Footrest					
	Small	Medium	Large		
Lower Leg length with Basic/Split Seat (mm)	185-285	280-450	290-510		
Lower Leg length with Contoured Seat (mm)	185-285	290-460	300-520		
Maximum Dorsiflexion/Plantarflexion	+/- 10 Degrees	+/- 10 Degrees	+/- 10 Degrees		
Chassis					
	Small	Medium	Large		
Maximum Weight Limit (kg) **	35	50	70		

## Comments

