WE'RE ALWAYS TESTING

Reducing vibration is important for our clients. We tested the Edge 3 and other similar models within our industry using a Vibration Dose Monitor. We used ISO-2631 vibration weighting filters with an ISO-8041 compliance seat sensor.

VIBRATION DOSE VALUE (VDV)

Determined by evaluating the extent of hard impacts.

MAX TRANSIENT VIBE VALUE (MTVV)

Determined by analyzing in-motion vibration and occasional impact.

THE RESULTS ARE IN

31% improvement in VDV* **30**% improvement in MTVV*

* Results are based on the Edge 3 compared with similar models in the industry.



Suspension Advantage SRS

Available Exclusively on the Edge 3, Edge 3 Stretto™ and 4Front®









Providing the smoothest overall ride for driveability and stability in the mid-wheel category.



