

CTM Mobility Scooter Solutions



DESCRIPTION	HS-328	HS-268	HS-520	HS-925	HS-559	HS-745	HS-828	HS-890
DESCRIPTION	Portable Small Range	Automatic Folding	Deluxe Mid Range	3-Wheeled Scooter	Deluxe Mid Range	Mid Range	Upgraded Road Class	Heavy Duty
OVERALL LENGTH	1120 mm	Folded: 775 mm Unfolded: 953 mm	1200 mm	1590 mm	1200 mm	1450 mm	1400 mm	1530 mm
OVERALL Width x Height (mm)	590 x 935 mm	494 x 510mm (folded) 500 x 960mm (un-folded)	640 x 1155 mm	690 x 1200 mm	628 x 1190 mm	685 x 1060 mm	680 x 1325 mm	685 x 1245 mm
FRONT/REAR WHEEL SIZE	9" / 9"	8" / 8"	11" / 11"	15" / 15"	11" / 11"	11" / 13"	13" / 13"	13" / 13"
WEIGHT (incl. batteries)	55 kg	30.6 kg	94 kg	88 kg	101 kg	106 kg	124 kg	114 kg
MAXIMUM SPEED	8 kph	8 kph	12 kph	15 kph	12 kph	12 kph	15 kph	15 kph
USER WEIGHT CAPACITY	136 kg	115 kg	136 kg	136 kg	136 kg	160 kg	160 kg	227 kg
GROUND CLEARANCE	52 mm	43 mm	65 mm	105 mm	80 mm	85 mm	90 mm	85 mm
GRADE CLIMABLE	8 degrees	6 degrees	10 degrees	8 - 10 degree	10 degrees	12 degrees	10 degrees	10 degrees
CURB CLIMBING	50 mm	40 mm	70 mm	100 mm	50 mm	80 mm	100 mm	90 mm
TURNING RADIUS	1450 mm	1120 mm	1400 mm	1800 mm	1500 mm	1460 mm	1660 mm	1730 mm
SUSPENSION	Rear	Rear	Front and Rear	Front and Rear	Front and Rear	Rear	Front and Rear	Front and Rear
SEAT TYPE	Swivel Lightweight Foldable	Padded foldable	Swivel Mid-Back with Headrest and Seat sliding	Scooter Seat with Angle Adjustable Armrests	Premium Upholstery Mid-back with Headrest, Swivel sliding	19" Tall Back Luxury Swivel Seat	Swivel Captain Seat with Sliding Seat and Back Angle Adjustment	Sliding and Swivel Captain Style Seat with Headrest
SEAT WIDTH	460 mm	419 mm	465 mm	456 mm	455 mm	480 mm	515 mm	480 mm
MOTOR SIZE	250W 4000 RPM	250W 4900 RPM	700W 4500 RPM	700W, 4200 RPM	700W 4500 RPM	700W 4100 RPM	850W 5000 RPM	700W 5400 RPM
BATTERY SIZE	2 x 12V 22Ah	1 x 24V 11.6Ah	2 x 12V 50Ah	50Ah battery	2 x 12V 50Ah	2 x 12V 50Ah	2 x 12V 50Ah	2 x 12V 50Ah
BATTERY WEIGHT	15.5 kg	2.8 kg	30.5 kg	30.5 kg	30.5 kg	30.5 kg	30.5 kg	30.5 kg
TRAVEL RANGE*	17.5 km	12 km	39 Km	35 km	42.5 km	36.5 km	27 km	37.5 km
BATTERY CHARGER	3A Off Board	2A Off Board	5A Off Board	5A Off Board	5A Off Board	5A Off Board	5A Off Board	5A Off Board
COLOURS STOCKED	Ice Blue	Sky Blue	Blue, Burgundy	Red, Blue, Silver	Blue, Silver, Champagne, Burgundy	Burgundy	Silver	Burgundy

Information correct as of May 2022. From time to time, product specifications are changed by manufacturers without prior notice
*Refers to transfer weight

CTM Mobility Scooter Solutions



DESCRIPTION	Deluxe Heavy Duty	Hill Climber	Robust	Robust
OVERALL LENGTH	1440 mm	1440 mm	1600 mm	1600 mm
OVERALL Width x Height (mm)	680 x 1360 mm	720 x 1360 mm	720 x 1280 mm	800 x 1620 mm
FRONT/REAR WHEEL SIZE	13" / 13"	13" / 13"	15" / 15"	15" / 15"
WEIGHT (incl. batteries)	117 kg	143.4 kg	178 kg	229 kg
MAXIMUM SPEED	15 kph	10 kph	15 kph	10 kph
USER WEIGHT CAPACITY	181 kg	181 kg	226 kg	200 kg
GROUND CLEARANCE	100 mm	100 mm	110 mm	110 mm
GRADE CLIMABLE	10 degrees	10 degrees	9 degrees	10 degrees
CURB CLIMBING	80 mm	80 mm	90 mm	100 mm
TURNING RADIUS	1600 mm	1600 mm	1860 mm	1860 mm
SUSPENSION	Front and Rear	Front and Rear	Front and Rear	Front and Rear
SEAT TYPE	Sliding and Swivel Captain Style Seat with Reclining Back and Headrest	Sliding and Swivel Captain Style Seat with Reclining Back and Headrest	Sliding and Swivel Captain Style Seat with Reclining Back and Headrest	Swivel Captain Style Seat with Seat Sliding Mechanism
SEAT WIDTH	520 mm	520 mm	559 mm	559 mm
MOTOR SIZE	850W 5000 RPM	850W 3450RPM	800W 4300RPM	1000W 3200RPM
BATTERY SIZE	2 x 12V 50Ah	2 x 12V 75Ah	2 x 12V 100Ah	2 x 12V 100Ah
BATTERY WEIGHT	30.5 kg	55 kg	65.3 kg	65.3 kg
TRAVEL RANGE*	27 km	40 km	50 km	44 km
BATTERY CHARGER	5A Off Board	8A Off Board	12A Off Board	12A Off Board
COLOURS STOCKED	Blue, Burgundy	Blue, Burgundy	Blue, Burgundy	Silver

Information correct as of May 2022. From time to time, product specifications are changed by manufacturers without prior notice
 *Refers to transfer weight

Accessories



Sun Canopy
Scooter Series 5/6/7/8 and 9




Scooter Cover Black
Scooter Series 5/6/7/8



Rear Basket for Seats with Headrest
CTM range of Powerchairs



Seat Bag for Seats with Headrest
Scooter Series 5 and 8 only



Front Basket with Lid
CTM Scooters HS-585/745/890/898



Seat Bag for with no Headrest
Scooter Series 5 and 8 only




CTM Large Rear Bag
Scooter Series 5/6/7/8 and HS-928



Crutch Holder
Scooter Series 5-8 and 9



Cane Holder
Scooter Series 5/6/7 and 8



Oxygen Tank Holder
Scooter Series 5/6/7/8 and 9



Walker Holder
Scooter Series 3/5/6/7 and 8



Safety Flag
Scooter Series 5-8