

Bariatric Manual Wheelchair Solutions



| DESCRIPTION | Low cost self-propelling wheelchair | Low cost transit wheelchair | Cost effective, lightweight and durable alloy wheelchair | Upgraded version of the Flexx which includes heavy duty functions | Robust self-propelling wheelchair specifically designed for larger clients |
|------------------------------|-------------------------------------|-----------------------------------|--|---|--|
| SEAT WIDTH | 50.8 - 60.9 cm (20" - 24") | 55.9 cm (22") | 45.7 - 55.9 cm (18" - 22") | 50.8 - 60.9 cm (20" - 24") | 55.9 - 66.0 cm (22" - 26") |
| SEAT DEPTH | 45.7 cm (18") | 45.7 cm (18") | 40.6 - 50.8 cm (16" - 20") | 40.6 - 50.8 cm (16" - 20") | 45.7 - 50.8 cm (18" - 20") |
| SEATING | Padded PVC upholstery | Slung nylon upholstery | Slung nylon upholstery | Slung nylon upholstery | Slung nylon upholstery |
| TRANSFER WEIGHT | 30 kg | 31.6 kg | 13.5 kg | 16.8 kg | 31.8 kg |
| CLIENT WEIGHT LIMIT | 250 kg | 230 kg | 160 kg | 170 kg | 250 kg |
| FOLDING | Yes | Yes | Yes | Yes | Yes |
| REAR WHEEL SIZE | 24" | 12" | 12", 20", 22" or 24" | 14", 22" or 24" | 20" or 24" |
| CASTOR SIZE | 20.3 cm (8" x 2") | 20.3 cm (8" x 2") | 15.4 or 20.3 cm (6" or 8") | 15.4 or 20.3 cm (6" or 8") | 17.8 cm (7") |
| ARMRESTS | Removable full length | Removable full length | Removable full or desk length, height adjustable | Removable full or desk length, height adjustable | Full length, fixed height, flip back |
| MINIMUM SEAT TO FLOOR HEIGHT | 49.5 cm (19.5") | 49.5 cm (19.5") | 43.2 cm (17") | 45.7 cm (18") | 43.2 cm (17") |
| STANDARD LEGREST STYLE | Swing-away 70°, adjustable length | Swing-away 70°, adjustable length | Swing-away 70°, adjustable length | Swing-in/out 80°, adjustable length | Swing-away 70°, adjustable length |
| MANUFACTURER | AML | AML | Karma | Karma | Karma |

Information correct as of October 2022. From time to time, product specifications are changed by manufacturers without prior notice