

Bariatric Manual Wheelchair Solutions



DESCRIPTION	Low cost self-propelling wheelchair	Low cost transit wheelchair	Cost effective, lightweight and durable alloy wheelchair	Upgraded version of the Flexx which includes heavy duty functions	Robust self-propelling wheelchair specifically designed for larger clients
SEAT WIDTH	50.8 - 60.9 cm (20" - 24")	55.9 cm (22")	45.7 - 55.9 cm (18" - 22")	50.8 - 60.9 cm (20" - 24")	55.9 - 66.0 cm (22" - 26")
SEAT DEPTH	45.7 cm (18")	45.7 cm (18")	40.6 - 50.8 cm (16" - 20")	40.6 - 50.8 cm (16" - 20")	45.7 - 50.8 cm (18" - 20")
SEATING	Padded PVC upholstery	Slung nylon upholstery	Slung nylon upholstery	Slung nylon upholstery	Slung nylon upholstery
TRANSFER WEIGHT	30 kg	31.6 kg	13.5 kg	16.8 kg	31.8 kg
CLIENT WEIGHT LIMIT	250 kg	230 kg	160 kg	170 kg	250 kg
FOLDING	Yes	Yes	Yes	Yes	Yes
REAR WHEEL SIZE	24"	12"	12", 20", 22" or 24"	14", 22" or 24"	20" or 24"
CASTOR SIZE	20.3 cm (8" x 2")	20.3 cm (8" x 2")	15.4 or 20.3 cm (6" or 8")	15.4 or 20.3 cm (6" or 8")	17.8 cm (7")
ARMRESTS	Removable 10" length	Removable 10" length	Removable full or desk length, height adjustable	Removable full or desk length, height adjustable	Full length, fixed height, flip back
MINIMUM SEAT TO FLOOR HEIGHT	49.5 cm (19.5")	49.5 cm (19.5")	43.2 cm (17")	45.7 cm (18")	43.2 cm (17")
STANDARD LEGREST STYLE	Swing-away 70°, adjustable length	Swing-away 70°, adjustable length	Swing-away 70°, adjustable length	Swing-in/out 80°, adjustable length	Swing-away 70°, adjustable length
MANUFACTURER	AML	AML	Karma	Karma	Karma

Information correct as of October 2022. From time to time, product specifications are changed by manufacturers without prior notice