

Paediatric Standing Solutions

	PRONE STANDERS							SUPINE / PRONE STANDERS					
	TOTSTANDER	FREESTANDER			PRONESTANDER			SQUIGGLES+	MYGO STANDER		HORIZON STANDER		
SIZES	1	1	2	3	1	2	3	One Size	1	2	1	2	3
AGE (APPROX.)	1 to 5	1 to 5	4 to 10	9 to 18	1 to 5	4 to 10	9 to 18	1 to 5	4 to 10	8 to 14	4 to 10	10 to 18	14 to Adult
MAX USER WEIGHT	22 kgs	60 kgs	80 kgs	90 kgs	60 kgs	80 kgs	90 kgs	22 kgs	50 kgs	60 kgs	50 kg	80 kgs	100 kgs
UPRIGHT	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
PRONE (ANGLE ADJUSTABLE)	No	No	No	No	0-40°	0-40°	0-40°	90-160°	0-80°	0-80°	90° to vert	90° to vert	90° to vert
SUPINE	No	No	No	No	No	No	No	Yes	Yes	Yes	Yes	Yes	Yes
WHEELED	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
HEIGHT ADJUSTABLE	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
CHEST PLATE HEIGHT - MIN	470 mm	470 mm	680 mm	975 mm	470 mm	680 mm	975 mm	560 mm	760 mm	980 mm	820 mm	1000 mm	1000 mm
CHEST PLATE HEIGHT - MAX	640 mm	680 mm	975 mm	1350 mm	680 mm	975 mm	1350 mm	940 mm	1080 mm	1300 mm	1150 mm	1450 mm	1450 mm
PRODUCT DIMENSIONS													
BASE LENGTH	450 mm	650 mm	750 mm	750 mm	820 mm	1070 mm	1270 mm	950 mm	978 mm	1130 mm	1135 mm	1280 mm	1280 mm
BASE WIDTH	450 mm	650 mm	750 mm	750 mm	600 mm	650 mm	700 mm	590 mm	630 mm	650 mm	580 mm	640 mm	740 mm
ACTIVITY TRAY DEPTH	No Tray	385 mm	410 mm	460 mm	385 mm	410 mm	460 mm	460 mm	340 mm	340 mm	420 mm	520 mm	520 mm
ACTIVITY TRAY WIDTH	No Tray	540 mm	580 mm	645 mm	540 mm	580 mm	645 mm	560 mm	580 mm	580 mm	440 mm	590 mm	590 mm
KEY FEATURES	Chest band Multi-adjustable pelvic band, Height & angle adjustable kneeplate, Multi-adjustable sandals. Portable.	Angle-adjustable tray with activity bowl Chest support with laterals Belt, knee cups & sandals in various sizes available. Mobile base with lockable castors			Angle-adjustable tray with activity bowl Chest support with laterals . Belt, knee cups & sandals in various sizes available. Mobile base with lockable castors			Hip abduction of up to 60° in prone, upright and supine standing Fully horizontal loading in supine Height-adjustable, foldable push handle for ergonomic use New wrap around knee supports Individually adjustable footplates	Headrest options for supine set up (Whitmyer compatible). Padded clear tray for prone or supine Independent angle adjustment in knees up to 25° -10 +10° plantar/dorsiflexion Mobile base with lockable castors Manual angle adjustment Belt, knee cups & sandals in a various sizes available		Headrest options for supine set up (Whitmyer compatible). Separate supine & prone trays Mobile base with lockable castors Powered angle adjustment Belt, knee cups & sandals in a various sizes available. Extremely robust.		
WIPEABLE PADDING	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
REMOVABLE PADDING	No	No	No	No	No	No	No	Yes	Yes Removable covers	Yes Removable covers	Yes	Yes	Yes
COLOUR OPTIONS	Meadow	Meadow	Mondeo Blue	Mondeo Blue	Meadow	Mondeo Blue	Mondeo Blue	Purple	Blue	Blue	Blue	Blue	Blue/Black

Information correct as of June 2024. From time to time, product specifications are changed by manufacturers without prior notice

Paediatric Walking Solutions

OTTOBOCK NURMI NEO



LECKEY MYWAY+



MEYWALK 4



SIZES	1	2	3	1	2	3	1	2	3	4
AGE (APPROX.)	1 to 3	3 to 5	5 to 8	1 to 5	4 to 11	9 to 16	1 to 5	4 to 9.5	9 to Adult	9 to Adult
MAX USER WEIGHT	25 kg	40 kg	55 kg	25 kg	50 kg	80 kg	60 kg	100 kg	130 kg	130 kg
UNIT LENGTH	650 mm	760 mm	830 mm	840 mm	930 mm	1180 mm	800 mm	850 mm	930 mm	1100 mm
UNIT WIDTH	630 mm	680 mm	720 mm	570 mm	670 mm	720 mm	560 mm	620 mm	620 mm	680 mm
MAX. INSIDE LEG HEIGHT	300 mm	450 mm	550 mm	410 mm	610 mm	810 mm	490 mm	640 mm	810 mm	1010 mm
MIN. INSIDE LEG HEIGHT	300 mm	450 mm	550 mm	240 mm	400 mm	550 mm	250 mm	270 mm	570 mm	770 mm
MAX. THORACIC SUPPORT HEIGHT	640 mm	790 mm	840 mm	920 mm	1190 mm	1490 mm	650 mm	999 mm	1330 mm	1520 mm
MIN. THORACIC SUPPORT HEIGHT	530 mm	660 mm	710 mm	635 mm	870 mm	870 mm	520 mm	740 mm	1010 mm	1200 mm
KEY FEATURES	Height adjustable grip bar, solid locking mechanism, lightweight and foldable, reverse roll locks on rear wheels, puncture proof low friction castors, stable base, swivelling front castors			Harness to provide core stability, supports the child in an upright posture enabling eye-level contact with standing peers, upper support unit with adjustable height, foldable frame for easy storage			Body posture of user can be set from vertical to 15 ° forward tilt, suitable for indoor and outdoor use, suspended spring system			

Information correct as of June 2024. From time to time, product specifications are changed by manufacturers without prior notice