

Paediatric Buggy Solutions Chart

	FIREFLY WEGO BUGGY	OTTOBOCK ECO BUGGY	SPECIAL TOMATO EIO PUSH CHAIR	SPECIAL TOMATO JOGGER STROLLER	PATRON PIPER BUGGY	PATRON CORZINO XCOUNTRY BUGGY	MYWAM MOUSE	MYWAM MEWA
DESCRIPTION	Designed to be combined with the Firefly GoTo Seat. Light weight and Portable, It works with both GoTo sizes and allows the child to have the same great head and lateral support while out and about.	Simple Buggy for Children with Mild to Moderate Needs. Thanks to its clever design and lightweight build, the Eco Buggy is a breeze to transport as it folds very easily.	A mild to moderate positioning seat matched with a lightweight stroller designed for children from early intervention to pre-adolescent ages.	Designed for parents and families who want to take long strolls through the city or through rough trails.	The Piper Comfort buggy is used for short term transportation and activities for your child.	The Corzino is a durable, semi-adjustable comfort rehab. buggy, which is suitable for even larger kids up to 17" in seat width.	Rear or forward facing, with tilt and recline as standard, the lightweight Mouse has suspension on all 4 wheels. Making it the perfect solution for smaller kids with supportive needs.	Offering both tilt and recline, as well as adjustable torso and head pads, the MEWA is an all option inclusive solution that doesn't look clinical. Suspension on all wheels and a rear or front facing design.
AGE RANGE (APPROX)	1-6 years	3-8 years	1-10 years	1-10 years	1-12 years	1-15 years	1-6 years	1-10 years
SEAT WIDTH RANGE	35 cm	35 cm	35 cm	35 cm	35 cm	30 / 34 / 38 / 42cm	24 (with insert) / 32 cm	25 (with insert) / 32 cm
SEAT DEPTH RANGE	26 cm	30 cm	30 cm	30 cm	29 cm	28 / 32 / 36 / 38cm	22 (with insert) / 26 cm	24 (with insert) / 28 cm
SEAT BACK HEIGHT	50 cm	55 cm	61 cm	63 cm	68 cm	54 / 67 / 85 / 88cm	52 cm	62 cm
FOLDING	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
RECLINE RANGE	140°	90° - 110°	115° - 165°	100° - 140°	90° - 135°	89 - 98° (depending on size)	16° - 88°	8° - 86°
ADJUSTABLE HANDLEBAR	Yes	No	Yes	Yes	Yes	Yes	Yes	Yes
BRAKE	Foot operated rear brake bar	Foot brake on rear wheels	Rear brake with suspension	Rear parking brake	Foot operated brakes	Foot operated brakes	Central foot operated brake	Central foot operated brake
TOTAL WEIGHT LIMIT	30kg	50kg	41kg	50kg	50kg	45 / 55 / 75 / 75kg	30kg	40kg
STANDARD FEATURES	Easily foldable, wheels suitable for rough terrain, comes with a full accessory pack	Aircraft grade aluminium frame, back angle adjustable, Swivel front castors, Spring mounted front and rear wheels	Lightweight aluminium frame, calf panel and seat back to lay your child flat, padded seat and back laterals, padded head pillow	Lightweight aluminium frame, folding latch for quick and easy one step folding, washable upholstery that is reversible, 5 point harness	Adjustable recline backrest, height adjustable footrest that can be detached, 5 point immobilisation harness with padding	Padded ergonomic seat cover, comfort side protection, adjustable backrest and footrest, 5 point immobilisation harness with padding	Adjustable Footrest, Torso / Head Pads, 5-Point Harness, Seat Minimising Kit, Abduction Belts, Canopy, Adjustable Wedge, Therapy Tray, Shopping Basket, Rain Cover	Adjustable Footrest, Torso / Head Pads, 5-Point Harness, Seat Minimising Kit, Abduction Belts, Canopy, Adjustable Wedge, Therapy Tray, Shopping Basket, Rain Cover
ACCESSORIES AVAILABLE	Standard with Canopy, Rain cover, Foot muff, Mosquito net, Shopping basket, Bottle holder, Wrist strap and compatible with GoTo Seat	Eco Buggy Canopy, 5 Point Harness	Rain cover, Soft-touch liners, Sitters	Rain cover, Head support, Medical necessities bag, Soft-Touch liners, Sitters	Flexipads set, Adjustable headrest, Abduction block, Feet fixation, Classic canopy, Footmuff - Summer, Footmuff - Winter, Raincover, Seat inlay	Basket, Sun hood, 5 point harness system, Tray, Rain cover, Trunk support, Support wing	Wide range of accessories included with buggy package. For the full list, please see our website	Wide range of accessories included with buggy package. For the full list, please see our website
BRAND	Firefly by Leckey	Ottobock	Special Tomato	Special Tomato	Patron	Patron	MyWam	MyWam

Information correct as of August 2024. From time to time, product specifications are changed by manufacturers without prior notice.

