

Tilt-in-Space Wheelchair Solutions

	TILT-IN-SPACE				RECLINING
	VIP2	VIP 515	FOCUS CR	LIBERTY FT	KM5000-18
CONSTRUCTION MATERIAL	Aircraft Grade Aluminium	Aircraft Grade Aluminium	7000 Grade Aluminium	7000 Grade Aluminium	Aircraft Grade Aluminium
SEAT WIDTH	40.6/45.7 cm (16"/18")	40.6/45.7 cm (16"/18")	36 - 56 cm (14" - 22")	36 - 56 cm (14" - 22")	41-46, 51 cm (16", 18", 20")
SEAT DEPTH	40.6/45.7 cm (16"/18")	40.6 cm (17")	36 - 56 cm (14" - 22")	36 - 51 cm (14" - 20")	40/45 cm (16"/18")
FRONT FRAME BENDS	80°	80°	60, 70, 80, and 90°	60, 70, 80, and 90°	70°
FRONT SEAT HEIGHT	46.5/49 cm (18"/19")	51 cm (20")	33 - 51 cm (13.5" to 20.5")	33 - 51 cm (13" - 20")	51 cm (20")
REAR SEAT HEIGHT	46.5/49 cm (18"/19")	51 cm (20")	33 - 51 cm (13.5" to 20.5")	33 - 51 cm (13" - 20")	51 cm (20")
LEGREST OPTIONS	90° Hangers and Elevating Legrest	70° and Elevating Legrest	60, 70, 80, 90° and Elevating Legrest	60, 70, 80, 90° and Elevating Legrest	70° and Elevating Legrest
RESIDUAL LIMB SUPPORT	No	No	Yes	Yes	Yes
TRANSFER WEIGHT	23.4 kg	18.9 kg	14.97 kg	11.7 kg	20 kg
CLIENT WEIGHT LIMIT	115 kg	115 kg	136.08 kg HD option 181.48 kg	113 kg	16" -18" frames - 115 kg 20" frames - 140 kg
TILT RANGE	0 - 35°	0 - 35°	-5° - 50°	0° - 20°	/
RECLINE RANGE	0 - 30°	/	/	/	0° - 180°
REAR WHEEL SIZE	20"	14"	12", 16", 20", 22", 24"	12" - 24"	24"
CASTOR SIZE	17.8 cm (7")	17.8 cm (7")	2.5 - 5.1 cm (5" - 8")	10.2 - 20.3cm (4" - 8")	15.2 - 20.3 cm (6", 8")
ARMREST HEIGHT	22.9 - 27.9 cm (9" - 11")	21 - 29 cm (8" - 11")	19.1 - 36.8 cm (7.5" - 14.5")	21.6 - 34.3cm (8.5" - 13.5")	22.9 - 27.9 cm (9" - 11")
MIN. SEAT TO FLOOR HEIGHT	46.5 or 49 cm (18.3 or 19.3")	50.8 cm (20")	33.0 cm (13")	33.0 cm (13")	50.8 cm (20")
MANUFACTURER	Karma	Karma	Ki Mobility	Ki Mobility	Karma

Information correct as of August 2024. From time to time, product specifications are changed by manufacturers without prior notice
*Refers to transfer weight

