

# Tilt-in-Space Wheelchair Solutions

|                               | ADULT                                                                             |                                                                                   |                                                                                     | PAEDIATRIC                                                                          |                                                                                     |                                                                                     |                                                                                     |
|-------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
|                               | VIP2                                                                              | FOCUS CR                                                                          | LIBERTY FT                                                                          | LITTLE WAVE ARC                                                                     | FLEXX ADAPT                                                                         | KIDEVO PRIME.T                                                                      | KIDEVO MINI.T                                                                       |
|                               |  |  |  |  |  |  |  |
| CONSTRUCTION MATERIAL         | Aircraft Grade Aluminium                                                          | 7000 Grade Aluminium                                                              | 7000 Grade Aluminium                                                                | 7000 Grade Aluminium                                                                | Aircraft Grade Aluminium                                                            | Aluminium                                                                           | Aluminium                                                                           |
| SEAT WIDTH                    | 40.6/45.7 cm<br>(16"/18")                                                         | 35.6 - 55.9 cm<br>(14" - 22")                                                     | 35.6 - 55.9 cm<br>(14" - 22")                                                       | 25.4 - 45.7 cm<br>(10" - 18")                                                       | 25.4 - 35.6 cm<br>(10" - 14")                                                       | 26 - 50 cm<br>(10" - 20")                                                           | 16 - 32 cm<br>(6" - 13")                                                            |
| SEAT DEPTH                    | 40.6/45.7 cm<br>(16"/18")                                                         | 35.6 - 55.9 cm<br>(14" - 22")                                                     | 35.6 - 50.8 cm<br>(14" - 20")                                                       | 30.5 - 50.8 cm<br>(12" - 20")                                                       | 22.9 - 38 cm<br>(9" - 15")                                                          | 20 - 50 cm<br>(8" - 20")                                                            | 16 - 35 cm<br>(6" - 14")                                                            |
| FRAME TYPE                    | Folding                                                                           | Fixed                                                                             | Folding                                                                             | Folding                                                                             | Folding                                                                             | Fixed                                                                               | Fixed                                                                               |
| FRONT SEAT HEIGHT             | 46.5/49 cm<br>(18"/19")                                                           | 33 - 51 cm<br>(13.5" to 20.5")                                                    | 33 - 51 cm<br>(13" - 20")                                                           | 37 - 48 cm<br>(14.5" - 19")                                                         | 43 cm<br>(17")                                                                      | 47 - 57 cm<br>(18.5" - 22.4")                                                       | 36 - 45 cm<br>(14.2" - 17.7")                                                       |
| REAR SEAT HEIGHT              | 46.5/49 cm<br>(18"/19")                                                           | 33 - 51 cm<br>(13.5" to 20.5")                                                    | 33 - 51 cm<br>(13" - 20")                                                           | 37 - 48 cm<br>(14.5" - 19")                                                         | 43cm<br>(17")                                                                       | 47 - 57 cm<br>(18.5" - 22.4")                                                       | 36 - 45 cm<br>(14.2" - 17.7")                                                       |
| GROWABLE SEAT PAN             | No                                                                                | Yes                                                                               | Seat Upholstery (Seat Pan Optional)                                                 | Yes                                                                                 | Yes                                                                                 | Yes                                                                                 | Yes                                                                                 |
| LEGREST OPTIONS               | 90° Hangers and Elevating Legrest                                                 | 60, 70, 80, 90° and Elevating Legrest                                             | 60, 70, 80, 90° and Elevating Legrest                                               | 60, 70, 80, 90° and Elevating Legrest                                               | 80, 90° and Elevating Legrest                                                       | Single Panel, Elevating Leg Rests, Spring Loaded Footplate                          | Single Panel                                                                        |
| RESIDUAL LIMB SUPPORT         | No                                                                                | Yes                                                                               | Yes                                                                                 | Yes                                                                                 | No                                                                                  | /                                                                                   | /                                                                                   |
| TRANSFER WEIGHT               | 23.4 kg                                                                           | 14.97 kg                                                                          | 11.7 kg                                                                             | 9.3 kg                                                                              | 18.6 kg                                                                             | 19.4 kg (Total Weight)                                                              | 9 kg (Total Weight)                                                                 |
| CLIENT WEIGHT LIMIT           | 115 kg                                                                            | 136.08 kg<br>HD Option 181.48 kg                                                  | 113 kg                                                                              | 75 kg                                                                               | 75 kg                                                                               | 75 kg                                                                               | 50 kg                                                                               |
| TILT RANGE                    | 0° - 35°                                                                          | -5° - 50°                                                                         | 0° - 20°                                                                            | 0° - 45°                                                                            | 0° - 45°                                                                            | -10° to 50°                                                                         | -5° to 40°                                                                          |
| ANTERIOR TILT OPTION          | No                                                                                | Up to -5°                                                                         | No                                                                                  | No                                                                                  | No                                                                                  | Up to -10°                                                                          | Up to -5°                                                                           |
| FOOT OPERATED TILT            | No                                                                                | Multiple Mounting Locations                                                       | Tilt Mounted on Push Handle                                                         | Multiple Mounting Locations                                                         | Tilt Mounted on Push Handle                                                         | Yes                                                                                 | Yes                                                                                 |
| RECLINE RANGE                 | 0 - 30°                                                                           | /                                                                                 | /                                                                                   | Optional                                                                            | 0° - 40°                                                                            | 30 °                                                                                | /                                                                                   |
| REAR WHEEL SIZE               | 20"                                                                               | 12", 16", 20", 22", 24"                                                           | 12" - 24"                                                                           | 12" - 16"                                                                           | 16"                                                                                 | 22" - 26"                                                                           | 20" - 24"                                                                           |
| CASTOR SIZE                   | 17.8 cm<br>(7")                                                                   | 12.7 - 20.3 cm<br>(5" - 8")                                                       | 10.2 - 20.3 cm<br>(4" - 8")                                                         | 12.7 - 20.3 cm<br>(5" - 8")                                                         | 15 - 20.3 cm<br>(6" - 8")                                                           | 12.7 cm<br>(5")                                                                     | 10.2 cm, 12.7 cm<br>(4", 5")                                                        |
| ANTITIPPERS                   | Optional                                                                          | Optional                                                                          | Optional                                                                            | Optional                                                                            | Yes                                                                                 | Optional                                                                            | Optional                                                                            |
| ARMREST TYPE                  | Flip Back, Quick Height Adjustable                                                | Multiple Options Including Skirt Guards, Height Adjustable or Flip Back Armrests  | Multiple Options Including Skirt Guards, Height Adjustable or Flip Back Armrests    | Multiple Options Including Skirt Guards, Height Adjustable or Flip Back Armrests    | Height Adjustable                                                                   | Standard, Extended Length, Height Adjustable, with or without Tray Holder           | /                                                                                   |
| ARMREST HEIGHT                | 22.9 - 27.9 cm<br>(9" - 11")                                                      | 19.1 - 36.8 cm<br>(7.5" - 14.5")                                                  | 21.6 - 34.3 cm<br>(8.5" - 13.5")                                                    | 16.5 - 38 cm<br>(6.5" - 15")                                                        | 7.5 - 29.5 cm<br>(3" - 12")                                                         | 18 - 28 cm<br>(7" - 11")                                                            | /                                                                                   |
| MIN. SEAT TO FLOOR HEIGHT     | 46.5 or 49 cm<br>(18.3 or 19.3")                                                  | 33.0 cm<br>(13")                                                                  | 33.0 cm<br>(13")                                                                    | 36.8 cm<br>(14.5")                                                                  | 43 cm<br>(17")                                                                      | 47 cm with a 22" Drive Wheel                                                        | 35.5 cm with a 20" Drive Wheel                                                      |
| HEIGHT ADJUSTABLE PUSH HANDLE | J-Hook, Fixed Height Adjustable                                                   | Multiple Options Including Height Adjustable and Stroller Styles                  | Multiple Options Including Height Adjustable and Stroller Styles                    | Yes                                                                                 | Yes                                                                                 | Yes                                                                                 | Yes                                                                                 |
| FOOT OPERATED BRAKE           | Brakes are Operated from the Push Handle                                          | Yes - Multiple Locations for Mounting                                             | Mounted on Push Handle                                                              | Yes - Multiple Locations for Mounting                                               | Mounted on Push Handle                                                              | Yes - 4 Brake Options                                                               | Yes - 4 Brake Options                                                               |
| TIE DOWNS                     | Yes                                                                               | Yes, Optional                                                                     | Yes, Optional                                                                       | Yes, Optional                                                                       | Yes                                                                                 | Yes                                                                                 | Yes                                                                                 |
| MANUFACTURER                  | Karma                                                                             | Ki Mobility                                                                       | Ki Mobility                                                                         | Ki Mobility                                                                         | Karma                                                                               | Ottobock                                                                            | Ottobock                                                                            |

Information correct as of April 2025. From time to time, product specifications are changed by manufacturers without prior notice