

Tilt-in-Space Wheelchair Solutions

	ADULT			PAEDIATRIC			
	VIP2	FOCUS CR	LIBERTY FT	LITTLE WAVE ARC	FLEXX ADAPT	KIDEVO PRIME.T	KIDEVO MINI.T
							
CONSTRUCTION MATERIAL	Aircraft Grade Aluminium	7000 Grade Aluminium	7000 Grade Aluminium	7000 Grade Aluminium	Aircraft Grade Aluminium	Aluminium	Aluminium
SEAT WIDTH	40.6/45.7 cm (16"/18")	35.6 - 55.9 cm (14" - 22")	35.6 - 55.9 cm (14" - 22")	25.4 - 45.7 cm (10" - 18")	25.4 - 35.6 cm (10" - 14")	26 - 50 cm (10" - 20")	16 - 32 cm (6" - 13")
SEAT DEPTH	40.6/45.7 cm (16"/18")	35.6 - 55.9 cm (14" - 22")	35.6 - 50.8 cm (14" - 20")	30.5 - 50.8 cm (12" - 20")	22.9 - 38 cm (9" - 15")	20 - 50 cm (8" - 20")	16 - 35 cm (6" - 14")
FRAME TYPE	Folding	Fixed	Folding	Folding	Folding	Fixed	Fixed
FRONT SEAT HEIGHT	46.5/49 cm (18"/19")	33 - 51 cm (13.5" to 20.5")	33 - 51 cm (13" - 20")	37 - 48 cm (14.5" - 19")	43 cm (17")	47 - 57 cm (18.5" - 22.4")	36 - 45 cm (14.2" - 17.7")
REAR SEAT HEIGHT	46.5/49 cm (18"/19")	33 - 51 cm (13.5" to 20.5")	33 - 51 cm (13" - 20")	37 - 48 cm (14.5" - 19")	43cm (17")	47 - 57 cm (18.5" - 22.4")	36 - 45 cm (14.2" - 17.7")
GROWABLE SEAT PAN	No	Yes	Seat Upholstery (Seat Pan Optional)	Yes	Yes	Yes	Yes
LEGREST OPTIONS	90° Hangers and Elevating Legrest	60, 70, 80, 90° and Elevating Legrest	60, 70, 80, 90° and Elevating Legrest	60, 70, 80, 90° and Elevating Legrest	80, 90° and Elevating Legrest	Single Panel, Elevating Leg Rests, Spring Loaded Footplate	Single Panel
RESIDUAL LIMB SUPPORT	No	Yes	Yes	Yes	No	/	/
TRANSFER WEIGHT	23.4 kg	14.97 kg	11.7 kg	9.3 kg	18.6 kg	19.4 kg (Total Weight)	9 kg (Total Weight)
CLIENT WEIGHT LIMIT	115 kg	136.08 kg HD Option 181.48 kg	113 kg	75 kg	75 kg	75 kg	50 kg
TILT RANGE	0° - 35°	-5° - 50°	0° - 20°	0° - 45°	0° - 45°	-10° to 50°	-5° to 40°
ANTERIOR TILT OPTION	No	Up to -5°	No	No	No	Up to -10°	Up to -5°
FOOT OPERATED TILT	No	Multiple Mounting Locations	Tilt Mounted on Push Handle	Multiple Mounting Locations	Tilt Mounted on Push Handle	Yes	Yes
RECLINE RANGE	0 - 30°	/	/	Optional	0° - 40°	30 °	/
REAR WHEEL SIZE	20"	12", 16", 20", 22", 24"	12" - 24"	12" - 16"	16"	22" - 26"	20" - 24"
CASTOR SIZE	17.8 cm (7")	12.7 - 20.3 cm (5" - 8")	10.2 - 20.3 cm (4" - 8")	12.7 - 20.3 cm (5" - 8")	15 - 20.3 cm (6" - 8")	12.7 cm (5")	10.2 cm, 12.7 cm (4", 5")
ANTITIPPERS	Optional	Optional	Optional	Optional	Yes	Optional	Optional
ARMREST TYPE	Flip Back, Quick Height Adjustable	Multiple Options Including Skirt Guards, Height Adjustable or Flip Back Armrests	Multiple Options Including Skirt Guards, Height Adjustable or Flip Back Armrests	Multiple Options Including Skirt Guards, Height Adjustable or Flip Back Armrests	Height Adjustable	Standard, Extended Length, Height Adjustable, with or without Tray Holder	/
ARMREST HEIGHT	22.9 - 27.9 cm (9" - 11")	19.1 - 36.8 cm (7.5" - 14.5")	21.6 - 34.3 cm (8.5" - 13.5")	16.5 - 38 cm (6.5" - 15")	7.5 - 29.5 cm (3" - 12")	18 - 28 cm (7" - 11")	/
MIN. SEAT TO FLOOR HEIGHT	46.5 or 49 cm (18.3 or 19.3")	33.0 cm (13")	33.0 cm (13")	36.8 cm (14.5")	43 cm (17")	47 cm with a 22" Drive Wheel	35.5 cm with a 20" Drive Wheel
HEIGHT ADJUSTABLE PUSH HANDLE	J-Hook, Fixed Height Adjustable	Multiple Options Including Height Adjustable and Stroller Styles	Multiple Options Including Height Adjustable and Stroller Styles	Yes	Yes	Yes	Yes
FOOT OPERATED BRAKE	Brakes are Operated from the Push Handle	Yes - Multiple Locations for Mounting	Mounted on Push Handle	Yes - Multiple Locations for Mounting	Mounted on Push Handle	Yes - 4 Brake Options	Yes - 4 Brake Options
TIE DOWNS	Yes	Yes, Optional	Yes, Optional	Yes, Optional	Yes	Yes	Yes
ONE ARM DRIVE	No	Yes	No	No	No	No	No
FOOT PROPELLING	No	Yes	Yes	No	No	No	No
MANUFACTURER	Karma	Ki Mobility	Ki Mobility	Ki Mobility	Karma	Ottobock	Ottobock

Information correct as of May 2025. From time to time, product specifications are changed by manufacturers without prior notice

