

Paediatric Buggy Solutions Chart

	ADAPTIVE STROLLERS					COMPLEX REHABILITATION STROLLERS	
	FIREFLY WEGO BUGGY	SPECIAL TOMATO EIO PUSH CHAIR	SPECIAL TOMATO JOGGER STROLLER	PATRON PIPER BUGGY	PATRON CORZINO XCOUNTRY BUGGY	MYWAM MOUSE	MYWAM MEWA
DESCRIPTION	Designed to be combined with the Firefly GoTo Seat. Light weight and Portable, It works with both GoTo sizes and allows the child to have the same great head and lateral support while out and about.	A mild to moderate positioning seat matched with a lightweight stroller designed for children from early intervention to pre-adolescent ages.	Designed for parents and families who want to take long strolls through the city or through rough trails.	For short term transportation and activities. Included canopy to assist with sensory regulation.	The Corzino is a durable, semi-adjustable comfort rehab. buggy, which is suitable for even larger kids up to 17" in seat width.	Rear or forward facing, with tilt and recline as standard, the lightweight Mouse has suspension on all 4 wheels. Making it the perfect solution for smaller kids with supportive needs.	With all the features of the Mouse but additional seat depth, recline and back height range.
AGE RANGE (APPROX)	1-6 years	1-10 years	1-10years	1-10 years	1-15 years	1-6 years	1-10 years
SEAT WIDTH RANGE	35 cm	35 cm	35 cm	35 cm	30 / 34 / 38 / 42cm	24 (with insert) / 32 cm	25 (with insert) / 32 cm
SEAT DEPTH RANGE	26 cm	30 cm	30 cm	29 cm	28 / 32 / 36 / 38cm	22 (with insert) / 26 cm	24 (with insert) / 28 cm
SEAT BACK HEIGHT	50 cm	61 cm	63 cm	68 cm	54 / 67 / 85 / 88cm	52 cm	62 cm
FOLDING	Yes	Yes	Yes	Yes	Yes	Yes	Yes
RECLINE RANGE	100° - 140° (40° of range)	115° - 165° (50° of range)	100° - 140° (40° of range)	90° - 135° (45° of range)	Backrest angle adjustable (2 positions)	106° - 178° (72° of range)	98° - 176° (78° of range)
TILT 'IN' SPACE	N/A	N/A	N/A	N/A	N/A	8° - 23°	5.5° - 21.5°
ADJUSTABLE HANDLEBAR	Yes	Yes	Yes	Yes	Yes	Yes	Yes
BRAKE	Rear foot brake	Rear foot brake	Rear foot brake	Rear foot brake	Rear mounted push to lock	Rear foot brake	Rear foot brake
USER WEIGHT LIMIT	30kg	41kg	50kg	50kg	45 / 55 / 75 / 75kg	30kg	40kg
PRODUCT WEIGHT*	11 kg	10 kg	12 kg	14.5 kg	11.8 - 14 kg	17.2 kg	20 kg
STANDARD FEATURES	Easily foldable, wheels suitable for rough terrain, comes with a full accessory pack	Easily foldable, built in suspension with Recline. Swivel front castors with lock out feature	Lightweight aluminium frame, folding latch for quick and easy one step folding, washable upholstery that is reversible	Adjustable recline backrest, height adjustable footrest that can be detached	Padded ergonomic seat cover, comfort side protection, adjustable backrest and footrest	Easily foldable, built in suspension with Tilt and Recline. Swivel front castors with lock out feature, reversible seat with Adjustable Footrest	Easily foldable, built in suspension with Tilt and Recline. Swivel front castors with lock out feature, reversible seat
INCLUDED ACCESSORIES	Canopy, Rain cover, Foot muff, Mosquito net, Shopping basket, Bottle holder, Wrist strap	Canopy, 5-Point harness, sunshade, head and thoracic positioning pads	Canopy with window and 5-point harness	5 point immobilisation harness with padding, Canopy	5-Point padded harness	Torso / Head Pads, 5-Point Harness, Seat Minimising Kit, Abduction Belts, Canopy, Pommel, Tray, Grab Bar, Shopping Basket, Rain Cover	Torso / Head Pads, 5-Point Harness, Seat Minimising Kit, Abduction Belts, Canopy, Pommel, Tray, Grab Bar, Shopping Basket, Rain Cover
OPTIONAL ACCESSORIES	GoTo Seat (additional 2kg weight)	Rain Cover, Soft-touch liners, Sitters	Rain cover, necessities bag, Soft-Touch liners, Sitters	Seat inlay reduction kit, Rain cover, Flexipads set, Adjustable headrest and Abduction block	Basket, Canopy, Rain cover, Foot muff, Adjustable headrest with side support	Indoor height adjustable base	Indoor height adjustable base
BRAND	Firefly by Leckey	Special Tomato	Special Tomato	Patron	Patron	MyWam	MyWam
WARRANTY	1 Year warranty	2 Year warranty	2 Year warranty	2 Year warranty	2 Year warranty	2 Year warranty	2 Year warranty

Information correct as of Feb 2026. From time to time, product specifications are changed by manufacturers without prior notice.

*Total weight dependant on accessories being used

