



Selection Back Wide

AL354020015 - XX-Large | AL354020016 - XXX-Large

Designed for people with a more excessive front section that need extra lift and support. The closure on the front is divided into two sections to make it easier to put on and to adjust the height on the front section. The reinforcement in the back divides pressure over a larger surface, yet leaves the spine free from pressure.

It has extra stabilizing straps that "cross-over" in the back. Measurement at Hip: 1250-1400 mm for XX-Large 1400-1600 mm for XXX-Large



Bariatric Perching Stool with Arms & Backrest

PHPAT0911156678

Designed for comfort and stability, this Bariatric Perching Stool features a sturdy frame with an adjustable height to suit users of various needs. The arms and backrest provide added support, making it ideal for those who require extra assistance while sitting or perching. The non-slip feet ensure stability on different floor surfaces, enhancing safety during use. This stool is perfect for use in the kitchen, bathroom, or anywhere a supportive seat is needed.

Suitable for users up to 250 kg



Ossur Wraparound Hinged Knee Support

OS212009 - XX-Large | OS212010 - XXX-Large

A low profile brace with Dual Action Hinges moulded into the brace and adjustable strapping to give extra support. It is made of neoprene lined with Lycra to offer a breathable contact area for comfort and extensive use. Wraparound design makes it very easy to put on and off and also accommodates legs of all shapes and sizes

Colour: Black

Measurement at Mid-Patella:

17" - 18" for XX-Large; 18" - 23" for XXX-Large



Belgrave Bariatric Commode / Shower Chair

DM190

An extremely robust commode with a non-peel frame and a cellular moulded backrest and over seat that can be removed when needed. Comes with a bucket which holds 14l and a lid.

It is a height adjustable commode so the user can achieve the most comfortable position and the arms can also be removed for easy side lateral transfer

Suitable for users up to 280 kg



Under Arm Bariatric Crutches - Pair

SMBE61314B

These heavy-duty steel crutches provide durability and support for users weight up to 249 kg. They are durable crutches for heavier, long term users. The push-button adjusted foot pieces with fitting scale allow for accurate and fast fitting. 1" adjustment increments on both the handle and the height, provide a precise fit which reduces muscle fatigue.

Completely latex free Height from handle for floor: 1575 - 1778 mm



Pride Power Bariatric Lift Chair

PRIDELC-358XLB - Burgundy Vinyl PRIDELC-358XLL - Latte Fabric

An easy-to-operate lift chair that allows users to stand or recline with ease. Engineered for quality and durability, it features a Dual Motor Sync system and high-density REFLEX foam for lasting comfort. The chair includes a battery backup system for power failures.

Available in both fabric and vinyl. Seat Width: 66 cm Weight Capacity: 227 kg



Big John Bariatric Toilet Seat

MP10200S

The Big John Toilet Seat provides extra stability and comfort for larger individuals. It features a contoured design and stabilising synthetic rubber bumpers to fit both round and elongated toilets.

Suitable for users up to 540kg Suitable for toilets with seat fixing centres of 130 – 140mm Made from durable, high-impact ABS plastic Easy to install and clean



Rollator - Bariatric 8" Wheel

MP10200S

he rustproof, robust steel frame is extremely durable. Heavyduty castors can be used both indoors and outdoors. It has a soft padded seat and a removable backrest.

Suitable for users up to 226 kg Anatomic hand grips Folds up for ease of transport Some assembly required



AMLWCTRB22

This sturdy and economical transit wheelchair is ideal for larger individuals needing assistance. It folds easily for storage or transport, making it highly convenient. The removable armrests and swing-away leg rests facilitate easy transfers.

Suitable for users up to 230 kg.

Powder-coated steel frame with double crossbar. Detachable full armrest.

Folds for easy transport. Some assembly required.



Karma Flexx HD Self Propelling

KMFLEXXHDSP2018 - 20" x 18" KMFLEXXHDSP2220 - 22" x 20" KMFLEXXHDSP2020 - 20" x 20" KMFLEXXHDSP2418 - 24" x 18" KMFLEXXHDSP2218 - 22" x 18" KMFLEXXHDSP2420 - 24" x 20"

Designed specifically for larger clients, the Karma BT10 is a versatile and durable bariatric wheelchair. It features an adjustable axle position for customised setup and ease of self-propulsion. The lightweight aluminium frame is easily folded for storage or transport.

Suitable for users up to 170 kg.

Adjustable backrest angle, 80° swing-away leg rests, and flip-back armrests

Folds for easy transport. Some assembly required.



AML Wheelchair Self-Propelling Bariatric

AMLWCSPB20 - 20" x 18" AMLWCSPB22 - 22" x 18" AMLWCSPB24 - 24" x 18"

This robust yet economical self-propelling wheelchair is the ideal choice for larger individuals. It folds easily for storage or transportation, and the removable armrests and swing-away leg rests assist with getting in and out of the chair, providing great benefits for the user and caregivers.

Suitable for users up to 250 kg.

Chrome polish finish steel frame. Detachable full armrest. Folds for easy transport. Some assembly required.



Karma Flexx HD Transit Wheelchair

KMFLEXXHDTR2018 - 20" x 18" KMFLEXXHDTR2220 - 22" x 20" KMFLEXXHDTR2020 - 20" x 20" KMFLEXXHDTR2418 - 24" x 18" KMFLEXXHDTR2420 - 24" x 20"

The Flexx HD is an upgraded version of the standard Flexx, featuring additional heavy-duty functions to assist larger individuals. Its configurable options, including seat depths, backrest angles, and leg rest styles, ensure optimal comfort. The adjustable components are clearly marked in orange for ease of use by both users and carers.

Suitable for users up to 170 kg Ellipse-shaped frame design for robustness. Folds for easy transport. Some assembly required.



Karma BT10 Bariatric Wheelchair

KMBT102220 - 22" x 20" KMBT102420 - 24" x 20" KMBT102620 - 26" x 20"

Designed specifically for larger clients, the Karma BT10 is a versatile and durable bariatric wheelchair. It features an adjustable axle position for customised setup and ease of self-propulsion. The lightweight aluminium frame is easily folded for storage or transport.

Suitable for users up to 250 kg.

Aircraft-grade aluminium frame. Fully detachable flip-back armrests. Swing-away anti-tippers.

Folds for easy transport. Some assembly required.



AML Bariatric Backrest


AMLBB2216 - 22" x 16" AMLBB2816 - 28" x 16" AMLBB2416 - 24" x 16" AMLBB3016 - 30" x 16" AMLBB2616 - 26" x 16"

This heavy-duty back support offers midline support and positioning for bariatric users. The high-strength hardback can support individuals weighing up to 220 kg. The molded foam provides maximum comfort and protection. The outer cover is moisture-resistant and breathable, ensuring durability and ease of maintenance.

Suitable for users up to 220 kg.

Quick-release for easy portability. Adjustable depth, angle, and rotation.

Some assembly required.



AML Heavy Duty Wheelchair Cushion

AMLCHD2220 - 22" x 20" AMLCHD2622 - 26" x 22" AMLCHD2620 - 26" x 20" AMLCHD2820 - 28" x 20"

This cushion is designed for superior positioning, stability, and comfort for larger users. It features an anatomically contoured shape with a waffled ischial relief recess and a reversible outer cover. The durable foam prevents bottoming out, and the moisture-resistant, breathable cover is easy to clean. Helps reduce pressure sores.

Suitable for users up to 266 kg.

Moisture-resistant inner cover. Reversible dartex outer cover



Ki Mobility Axiom G Cushion

KIXG1816 - 18" x 16" KIXG1818 - 18" x 18" KIXG2018 - 20" x 18" KIXG2020 - 20" x 20"

The Axiom G Cushion offers superior comfort and skin protection with its contoured, high-density foam base. It features a copper-infused, anti-microbial cover to reduce odor and improve hygiene. The cushion is designed for durability and effective pressure distribution, making it ideal for users who need long-term support.

*Axiom P Cushion another option for providing positioning and stability

Suitable for users up to 158 kg. Sizes: 14" x 14" to 20" x 20"

Composition: High-density foam Cover: Copper-infused, anti-microbial



Q6 Edge® HD

PRIDEQ6EDGEHDPB

The Q6 Edge HD is a compact mid-wheel heavy duty drive chair. With features like standard 4-pole motors, Mid-Wheel 6 Drive Design and ATX Suspension, the Q6 Edge® HD is engineered to meet the performance needs of the most active user. This chair accepts a complete range of seating and electronics options to deliver exceptional rehab capability at an excellent value.

Weight capacity of 204 kg, with or without power tilt. Heavy-duty seating with seat width from 20" - 32". Standard seating with seat width from 12" - 24" with easy adjustment.



HS-928 Mobility Scooter

HS928BU1

The CTM HS-928 is a robust four wheel scooter for those who want to do the distance! It is a great scooter for larger users or for those wanting to go longer distances, as this heavy duty scooter will deliver. Its sporty design and numerous safety features makes it very comfortable to drive and a really eye catcher.

Colours: red or blue Maximum user weight: 226 kg Maximum speed: 15 kg



Quantum® 1450 Powerchair

PRIDEQ1450PBQ

Designed from the ground up as a bariatric powered mobility solution, the Quantum® 1450 packs the necessary punch and range of appropriate options and sizing to accommodate your requirements. Available in a range of seating configurations from 20" to 32" wide, and with a weight capacity of up to 272kg, the Quantum® 1450 incorporates the most heavy duty of construction with large batteries, high powered controller, and solid wheels all round.

Weight capacity of 272 kg









Please visit our website **www.alliedmedical.co.nz** for our full range of products





0800 31 61 81 helpis@alliedmedical.co.nz

