

Allied Medical

Paediatric Solutions



**ALLIED
MEDICAL**
LIMITED

Product Range

BUGGIES 03

BUGGY SOLUTIONS CHART 05

WHEELCHAIRS 07

POWERCHAIRS 08

WHEELCHAIR SEATING 11

SEATING SOLUTIONS 13

WALKING SOLUTIONS 16

WALKING SOLUTIONS CHART 18

STANDING SOLUTIONS 19

STANDING SOLUTIONS CHART 21

EARLY INTERVENTION 23

BATHING 24

TOILETING 26

24 HOUR POSTURAL CARE 29

OTHER 30



Firefly WeGo Buggy

A lightweight, portable and snazzy-looking pushchair with a few hidden surprises. The WeGo is made for GoTo Seat users, as it slots into the buggy easily.

- Easily foldable and setup in one simple movement
- Wheels suitable for rough terrain
- WeGo reclines back 140° so your child can lie back for a little snooze or for changing
- Comes standard with a full accessory pack including a sun canopy, rain cover, foot muff, mosquito net, shopping basket underneath seat, bottle holder and wrist strap



Special Tomato EIO Push Chair

A stylish, special needs stroller unmatched in features, quality and price. It offers a mild to moderate positioning seat matched with a lightweight stroller designed for children from early intervention to pre-adolescent ages.

- Easy to fold and transport
- Lightweight aluminium frame
- Calf panel and seat back to lay your child flat
- Shopping basket and sun shade
- Significant growth included



Special Tomato Jogger Stroller

This active stroller provides comfort and support to children with special needs while riding on all terrain. It is an easy-to-fold stroller which allows for easy storage wherever you take it.

- Easy to fold and transport
- Lightweight aluminium frame
- Functional rear suspension
- Padded seat can be reclined from 10 degrees to 50 degrees



OttoBock Eco Buggy

Simple Buggy for Children with Mild to Moderate Needs. Thanks to its clever design and lightweight build, the Eco Buggy is a breeze to transport as it folds very easily.

- Small lightweight buggy
- Back angle adjustment 20 degrees (90 - 110 degrees)
- Folding frame, allows it to fold to a compact size to easily fit into most vehicles
- Footrest is height adjustable from 200 - 340 mm (distance from seat edge to footrest)
- Swivel front castors
- Spring mounted front and rear wheels
- Foot brake on back wheels



Patron Corzino XCountry

A semi-adjustable rehabilitation buggy available in four sizes with a wide variety of options available and a high user weight limit of up to 75kg. This buggy is equipped with an optimised buddy-folding system, a new padded ergonomic seat cover and a light frame with 5-point fixation ensuring safety for your child and easy handling.




- Well padded ergonomic seat cover
- Comfortable side protection
- Footrest height adjustable and foldable (angle adj.)
- Backrest angle adjustable (2 positions)
- Backrest height extension (on size 38 and 42 only)



Patron Piper Comfort

This buggy is ideal for short term transportation and activities for your child. It has a fully foldable undercarriage made from lightweight aluminium profiles.

- Extended reclining backrest
- Height adjustable footrest that can be detached
- 5-point immobilisation harness with crotch belt and shoulder padding
- Detachable wheels and forks
- One-piece foldable handlebar

	FIREFLY WEGO BUGGY 	OTTOBOCK ECO BUGGY 	SPECIAL TOMATO EIO PUSH CHAIR 
BRIEF DESCRIPTION	Designed to be combined with the Firefly GoTo Seat. It works with both GoTo sizes and allows the child to have the same great head and lateral support while out and about.	Simple Buggy for Children with Mild to Moderate Needs. Thanks to its clever design and lightweight build, the Eco Buggy is a breeze to transport as it folds very easily.	A mild to moderate positioning seat matched with a lightweight stroller designed for children from early intervention to pre-adolescent ages.
AGE RANGE (APPROX)	1-6 years	3-8 years	1-10 years
SEAT WIDTH RANGE	35 cm	35 cm	35 cm
SEAT DEPTH RANGE	26 cm	30 cm	30 cm
SEAT BACK HEIGHT	50 cm	55 cm	61 cm
FOLDING	Yes	Yes	Yes
RECLINE RANGE	140°	90° - 110°	115° - 165°
ADJUSTABLE HANDLEBAR	Yes	No	Yes
BRAKE	Foot operated rear brake bar	Foot brake on rear wheels	Rear brake with suspension
TOTAL WEIGHT LIMIT	30kg	50kg	41kg
STANDARD FEATURES	Easily foldable, wheels suitable for rough terrain, comes with a full accessory pack	Aircraft grade aluminium frame, back angle adjustable, Swivel front castors Spring mounted front and rear wheels	Lightweight aluminium frame, calf panel and seat back to lay your child flat, padded seat and back laterals, padded head pillow
ACCESSORIES AVAILABLE	Standard with Canopy, Rain cover, Foot muff, Mosquito net, Shopping basket, Bottle holder, Wrist strap and compatible with GoTo Seat	Eco Buggy Canopy, 5 Point Harness	Rain cover, Soft-touch liners, Sitters
BRAND	Firefly by Leckey	OttoBock	Special Tomato

Information correct as of May 2022. From time to time, product specifications are changed by manufacturers without prior notice.

**SPECIAL TOMATO
JOGGER STROLLER**



Designed for parents and families who want to take long strolls through the city or through rough trails.

**PATRON
PIPER BUGGY**



The Piper Comfort buggy is used for short term transportation and activities for your child.

**PATRON CORZINO
XCOUNTRY BUGGY**



The Corzino is a semi-adjustable comfort rehab buggy, which comes in 4 sizes.

1-10 years	1-12 years	1-15 years
35 cm	35 cm	30 - 42 cm
30 cm	29 cm	28 - 39 cm
63 cm	68 cm	54 - 88 cm
Yes	Yes	Yes
100° - 140°	90° - 135°	89 - 98° (depending on size)
Yes	Yes	Yes
Rear parking brake	Foot operated brakes	Foot operated brakes
50kg	50kg	45 - 75kg
Lightweight aluminium frame, folding latch for quick and easy one step folding, washable upholstery that is reversible, 5 point harness	Adjustable recline backrest, height adjustable footrest that can be detached, 5 point immobilisation harness with padding	Padded ergonomic seat cover, comfort side protection, adjustable backrest and footrest, 5 point immobilisation harness with padding
Rain cover, Head support, Medical necessities bag, Soft-Touch liners, Sitters	Flexipads set, Adjustable headrest, Abduction block, Feet fixation, Classic canopy, Footmuff - Summer, Footmuff - Winter, Raincover, Seat inlay	Basket, Sun hood, 5 point harness system, Tray, Rain cover, Trunk support, Support wing
Special Tomato	Patron	Patron



Karma 7520

The Karma lightweight paediatric wheelchair is durable, functional and designed to meet a range of New Zealand kid's needs.

- Ideal as a backup wheelchair for children who use a powerchair, or for those who only need it for shorter periods of time
- Easily adjustable push handle to accommodate carers of all heights
- Quick release, removable rear wheels
- Height adjustable armrests that can be removed for transfers
- Available in self-propelling or transit configuration

Maximum user weight: 100kg | Min seat to floor height: 18.5"
Seat Width & Depth: 11.5" - 13.5"




Karma Flexx Junior

The Karma Flexx Junior wheelchair is an extremely flexible and popular wheelchair solution to suit a range of customer needs. Everything that is orange is adjustable which makes for very easy use.

- Configurable options include seat depth, seat to floor heights, backrest heights, castor fork angles and centre of gravity
- Swing away footrests with flip-up footplates
- Can be set up to provide optimal comfort and performance
- Adjustable backrest angle

Maximum user weight: 80kg | Min seat to floor height: 15.5"
Seat Width & Depth: 10" - 15"



Ki Mobility Rogue XP

Rogue XP is the ultimate performance chair; encouraging you to embrace your personal style at the age when you care the most. It is a teenage chair that gives you the necessary built-in growth; no need to compromise fit, function, ride or style.

- Adjustable fixed frame for growth - up to 3"
- Creative freedom - ability to customise your chair
- Maximum adjustment and growth ability
- Integrated swing-away armrests
- Combination of rigidity and reduced vibration dampening for an ultra-responsive ride
- Adjustable dynamic 5th wheel

Maximum user weight: 90kg | Min seat to floor height: 13.5"
Seat Width & Depth: 10" - 18" or 10" - 20"



Karma 7520 Recliner

This reclining children's wheelchair has been designed especially for the New Zealand market - designed to fit with the needs of children who might be progressing from a standard buggy-type stroller, or have the need to rest in their wheelchair from time to time. It is also ideally suited for use after surgery or during recovery from an accident.

- Easily adjustable push handle - to accommodate carers of all heights
- Quick release, removable rear wheels for ease of transport in vehicle
- Fully removable, height adjustable armrests, that allow for easier transfers
- Reinforced cross bracing - strong but easily folded
- Multiple-position rear wheel and castor mounts to optimise centre of gravity and seat-to-floor height adjustments

Maximum user weight: 100kg | Min seat to floor height: 18.5" | Seat Width & Depth: 11.5" - 13.5"



Ottobock Skippi

Designed specifically for children, this fun rear wheel drive powerchair is perfect for active kids. It is compact and colourful, but still retains large size drive wheels and castors so kids can explore their environment with the push of a joystick.

- Robust paediatric powerchair that will perform well both indoors and outdoors with more active users (has relatively high ground clearance)
- Seat frame can accommodate the width and depth growth

Maximum user weight: 50kg | Min seat to floor height: 17"
Seat Width & Depth: 10" - 15" & 11" - 15"



Jazzy 623 2.0

With highly reliable, efficient motors that deliver a higher range and improved responsiveness than the previous model, the Jazzy® 623 2.0 is engineered to meet the performance needs of the most active user.

- Compact, mid-wheel drive powerchair
- TRU-Balance 3 adjustable seating
- VR2 controller to deliver smooth, nimble manoeuvrability with good indoor and outdoor performance and a tight turning radius

Maximum user weight: 136kg | Min seat to floor height: 17.5"
with tilt | Seat Width & Depth: 12" - 22" or 12" - 22"



Quantum Edge 3 Stretto

The Quantum Edge 3 Stretto is the most capable compact power base in New Zealand with an extremely narrow footprint for outstanding manoeuvrability, making it ideal for children, teenagers and small adults.

- Independent SRS (Smooth Ride Suspension) making for an incredibly easy ride
- The Stretto is available with 12.5" or 14" drive wheels and optional 12" of iLevel power adjustable seat height

Maximum user weight: 136kg | Seat Width & Depth: 12" - 22"
or 12" - 22"



Ottobock Skippi Powerchair

Varilite Junior Back Supports



Lightweight and easy to use, the Varilite Junior Back Support is perfect for children and families on the go. The compact, minimal design directs focus to the child and not the equipment. Available in two styles - Mid and Deep (which offers lateral support).

Specification	Mid	Deep
Back height	9" (23.5cm)	9" (23.5cm)
Outside width of shell	8.5" (21.8cm)	10" (25.1cm)
Inside width of support surface	7" (19cm)	9" (22.6cm)
Fits wheelchair sizes	9" - 11" (23 - 28cm)	9" - 11" (23 - 28cm)
Weight	1.1kg	1.3kg

Ride Designs Java Decaf Back



The Flip XP is designed for the smallest of riders with some of the biggest advancements in the industry. From its flat folding design to its use of Taper Lok Technology, the Flip is simply the most advanced paediatric tilt-in-space available today.

Specification	Java Decaf
Back height	10"- 15" (25-38cm)
Fits wheelchair sizes	10"-15" (25-36cm)
Weight	Up to 1.6kg



Paediatric Curve Cushion

A paediatric cushion that utilises moulded, high-density seat design to properly position the legs, hips, and spine to provide optimal support in a wheelchair. We recommend purchasing a spare cover with all cushion purchases.

Maximum user weight: 50kg | Min seat to floor height: 17"
Seat Width & Depth: 10" - 15" & 11" - 15"



Ki Mobility Axiom P Cushion

A Positioning cushion offering the Axiom Difference design elements with an emphasis on medial and lateral stability. With superior positioning features and surface tension reduction in the well, Axiom P proves positioning cushions can still be comfortable while providing positioning and stability.

Sizing:
Widths: 14" - 22" | Depths: 14" - 22"



Varilite Junior Cushion

The Quantum Edge 3 Stretto is the most capable compact power base in New Zealand with an extremely narrow footprint for outstanding manoeuvrability, making it ideal for children, teenagers and small adults.

Sizes:
Widths: 10" - 12" | Depths: 10" - 12"
Adult sizes available - Reflex cushion



Ride Designs Forward Cushion

The Forward cushion is all about contour and helps prevent sliding while providing superior support and protection. Ride's approach redistributes the forces of pressure and shear from high risk anatomic areas while selectively applying the forces to areas more tolerant. The contours are specifically optimised for a paediatric pelvis.

Sizing:
Widths: 10" - 12" | Depths: 10" - 16"



Special Tomato Soft-Touch Liner Seat and Back

The Soft-Touch Liner offers an affordable alternative to other seating systems. It eliminates the need to transport heavy seating systems. The contoured seat and back cushions promote good posture while sitting in a standard seating system for children with mild to moderate positioning needs.

- Made of a soft-to-touch, easy clean material that is latex-free and anti-microbial
- Can be attached to most seats including but not limited to the following: strollers, highchairs, dining room chairs, restaurant chairs and classroom chairs
- Easy storage as the cushions can be folded together



Special Tomato Soft-Touch Sitter

The Special Tomato Soft-Touch Sitter provides a comfortable and functional seating solution. It has a strong exterior that will stand up, whilst ensuring a comfortable seat for your child.

- Impermeable to fluids
- Peel and tear resistant
- Latex-free material with built in anti-microbial protection
- Soft contoured seat with anti-thrust seating surface
- Built-in abductor
- Contoured head and lateral supports
- Attachment straps for securing to either the Floor Sitter or the Mobile Base



Special Tomato CatchUp Chair

The CatchUp Chair is a stable, highly adjustable, mobile seat. It is easily configured without tools to provide the well fitting seat at the perfect height for a variety of tabletop surfaces. It is appropriate for individuals who require mild to moderate postural support. The CatchUp Chair fits children and adults weighing up to 113kgs.

- Compact design makes it an excellent choice for schools, therapy clinics, and hospitals as well as homes
- Simple, tool-free adjustments
- Pelvic Belt and Chest Harness stabilise a well-aligned posture



Firefly GoTo Seat

The GoTo Seat is a portable support seat for children with extra postural needs. It can be taken anywhere, which makes getting out and about much easier. The GoTo Seat fits easily into the WeGo Buggy, making it that much easier to go places.

- Two sizes available
- Suitable for attachment to most standard chairs in the home, restaurant or classroom
- Adjustable headrest and laterals that adjust for growth
- Two headrest choices: Standard or Advanced
- Optional floor sitter accessory has 26 recline positions, from upright to right back



Leckey PAL Seat

The PAL Seat is designed to give improved stability, which helps reduce fatigue and allows the child longer periods of concentration and fine motor activity. The simplicity of the PAL Seat enables the therapist to assess the child's needs instantly. It is a really simple, fast way to help kids join in with learning.

- Two point lap belt and side pads for lateral stability
- Height adjustable armrests, legs and footplate
- Two sizes accommodating kids from 1 - 12 years



Leckey BeMe

The infinitely configurable BeMe from Leckey, is a revolutionary seating system that puts a child's development and independence at the heart of its design. BeMe is a highly personalised seating system available in three sizes and spec'd to perfection.

- Growable contoured cushion with high-stretch foam and fabric, to maximise pelvic stability
- Interchangeable parts (including seats, back supports, medial, lateral, proximal and distal supports), click and switch with total ease
- Dynamic backrest is air pressure controlled which regulates resistance as the child extends and is speed adjustable




Buckle Guard

A simple solution for the little (and not so little) Houdini's who are able to escape from car safety harnesses and may endanger themselves or others in a vehicle. This strong and difficult to break cover slides over the seatbelt buckle to prevent the wearer from releasing the belt themselves. A caregiver can then safely and easily release the belt by inserting a key or coin into a slot when the destination is reached.

- Designed to go over most car seat buckles - either side release or top release.
- Dimensions: 64 mm W x 64 mm H x 42 mm D
- Comes complete with a seatbelt cutter and an instruction sheet.
- Manufactured from clear and durable ABS plastic.
- Includes a wedge spacer to accommodate for many smaller sized buckles.
- Webbing cutter is provided in case of an emergency.
- Australian design





Firefly Upsee

The Firefly Upsee is an upright mobility device designed to let children stand, walk - and dance - with the help of an adult.

- Sizes available: Small, Medium or Large
- Age: 2-4 / 4-6 / 6-8 years
- Maximum user weight: 15kg / 20kg / 25kg



Ottobock Nurmi Neo Posterior Walker

The Ottobock Nurmi Neo is a highly manoeuvrable kids posterior walker designed to promote an upright posture and natural walking pattern.

- Height-adjustable grip bar
- Solid locking mechanism to stop slipping
- Lightweight and foldable
- Reverse roll locks on the rear wheels
- Puncture proof low friction castors
- Stable base
- Swivelling front castors



Leckey MyWay

MyWay is an upright mobility device that is all about providing children with maximum opportunities for interaction, exploration and participation in a supported, upright position.

- Revolutionary patented harness support
- Adult Connection belt offers an adult to attach to the device and control movements or complete therapeutic exercises
- Upper support unit with adjustable gross height, fine height, prone angle and shoulder height
- Open front design for interaction & exploration and open back design to allow the therapist to work with the user
- MyWay Pedal attachment available to purchase separately for use with MyWay size 2 or 3 frame only

Sizes available: 1-3 | Age: 1-16 years

Maximum user weight: 25-80kg



Meywalk 4

The Meywalk 4 is a high quality, Danish designed walker that offers a suspended spring system that stimulates the muscles and movements to encourage a more natural gait pattern.

- The body posture of the user can be set in any position from vertical to 15° forward tilt
- Suitable option for both indoor and outdoor use
- Suitable for both adolescents and adults, with two larger size options available

Sizes available: Mini - Large | Maximum user weight: 80-130kg



OTTOBOCK NURMI NEO



LECKEY MYWAY



MEYWALK 4

SIZES	1	2	3	1	2	3	1	2	3	4
AGE (APPROX.)	1 to 3	3 to 5	5 to 8	1 to 5	4 to 11	9 to 16	1 to 5	4 to 9.5	9 to adult	9 to adult
MAX USER WEIGHT	25kg	40kg	55kg	25kg	50kg	80kg	50kg	100kg	130kg	130kg
UNIT LENGTH	650mm	760mm	830mm	720mm	850mm	1100mm	800mm	850mm	930mm	1100mm
UNIT WIDTH	630mm	680mm	720mm	570mm	670mm	720mm	560mm	620mm	620mm	680mm
MAX. INSIDE LEG HEIGHT	300mm	450mm	550mm	440mm	610mm	810mm	490mm	640mm	810mm	1010mm
MIN. INSIDE LEG HEIGHT	300mm	450mm	550mm	220mm	400mm	550mm	250mm	270mm	570mm	770mm
MAX. THORACIC SUPPORT HEIGHT	640mm	790mm	840mm	920mm	1190mm	1490mm	740mm	810mm	1300mm	1490mm
MIN. THORACIC SUPPORT HEIGHT	530mm	660mm	710mm	635mm	870mm	1070mm	520mm	710mm	980mm	1170mm
FEATURES	Height adjustable grip bar, solid locking mechanism, lightweight and foldable, reverse roll locks on rear wheels, puncture proof low friction castors, stable base, swivelling front castors			Harness to provide core stability, supports the child in an upright posture enabling eye-level contact with standing peers, upper support unit with adjustable height, foldable frame for easy storage			Body posture of user can be set from vertical to 15° forward tilt, suitable for indoor and outdoor use, suspended spring system			

Information correct as of September 2022. From time to time, product specifications are changed by manufacturers without prior notice. Transit tie down is available on all chairs for transport purposes.



Leckey Totstander

A pre-school standing frame that is adjustable and extremely portable. It is a junior version of the Freestander and is designed to support children in the early stages of standing.

- Simple construction, allowing adjustment and portability
- Allows the child to be positioned close to work tables or benches, whilst keeping them at a similar height to their peers
- Therapists can use the Totstander as an assessment tool to determine a child's ability to stand and bear weight
- Side tubes accommodate growth
- Adjustable pelvic band, knee plate, padded knee cups and sandals that velcro in place

Age: 1-5 years | Maximum user weight: 20kg



Leckey Freestander

An upright stander for children and teens, created to support children with mild to moderate levels of physical involvement. The Freestander is the next size up from the Totstander.

- A lower baseboard allowing for quick and easy transfers from chair to stander
- Pelvic band and optional padded pelvic plate are designed to help promote mid-line symmetry
- Adjustable flexible laterals, padded chest plate, knee plate and padded knee cups, sandals, pelvic band and activity tray

Sizes available: 1-3 | Age: 1-5, 4-9, 9-17 years
Maximum User Weight: 60, 80 or 90kg



Leckey Pronestander

The Leckey Prostander assists in the promotion of weight bearing and mid-line symmetry. The user can be positioned in a fully upright position or at any one of four prone angles.

- Four prone angles are designed to allow weight bearing to be gradually increased as the child's tolerance increases
- Hinged footplate permits wheelchair access and assisted transfer

Sizes available: 1-3 | Age: 1-5, 4-10, 9-18 years
Maximum user weight: 60, 80 or 90kg



Leckey Squiggles+ Standing Frame

The Squiggles+ Stander is an extremely versatile 3-in-1 stander, offering prone, upright and supine standing in one product.

- 3 in 1 design offering prone, upright and supine standing
- Up to 60 degrees of bilateral hip abduction
- Advanced knee support design
- Independent footplates to accommodate leg length discrepancy
- Fully horizontal loading in supine
- Activity tray with height adjustability
- Height-adjustable push handle

Age: 1-5 years | Maximum user weight: 22kg



Leckey Mygo Stander

As the child grows out of the Squiggles Stander, they can move into the Mygo Stander to continue their standing therapy. It encourages natural standing posture and has a range of functions for the growing child.

- This 3-in-1 stander offers prone, upright and supine standing in one product
- It can accommodate up to 25 degree fixed hip and knee contractures
- Supine and prone angles can be changed from 0 to 80 degrees
- Tool free adjustments

Sizes available: 1-2 | Age: 4-10 or 8-14 years
Maximum user weight: 50 or 80kg



Leckey Horizon Stander

This powered stander provides gradual, appropriate weight bearing support in either prone, supine or upright positions. Robust and durable for users aged 4 to adult.

- Included conversion kit offering 90 degrees of adjustment from horizontal to freestanding position in both prone and supine positions
- Hand held, easy to use controls for the powered tilt mechanism

Sizes available: 1-3 | Age: 4-10, 8-14 or 14+ years
Maximum user weight: 50, 80 or 100kg

PRONE STANDERS

TOTSTANDER



FREESTANDER



PRONESTANDER



	TOTSTANDER			FREESTANDER			PRONESTANDER			
SIZES	1	1	2	3	1	2	3	1	2	3
AGE (APPROX.)	1 to 5	1 to 5	4 to 9.5	9 to 17	1 to 5	4 to 9.5	9 to 17	1 to 5	4 to 9.5	9 to 17
MAX USER WEIGHT	20 kgs	60 kgs	80 kgs	90 kgs	60 kgs	80 kgs	90 kgs	60 kgs	80 kgs	90 kgs
UPRIGHT	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
PRONE (ANGLE ADJUSTABLE)	No	No	No	No	0-40°	0-40°	0-40°	0-40°	0-40°	0-40°
SUPINE	No	No	No	No	No	No	No	No	No	No
WHEELED	No	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
HEIGHT ADJUSTABLE	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
CHEST PLATE HEIGHT - MIN	470 mm	470 mm	680 mm	975 mm	470 mm	680 mm	975 mm	470 mm	680 mm	975 mm
CHEST PLATE HEIGHT - MAX	640 mm	680 mm	975 mm	1350 mm	680 mm	975 mm	1350 mm	680 mm	975 mm	1350 mm
PRODUCT DIMENSIONS										
BASE LENGTH	450 mm	650 mm	750 mm	750 mm	820 mm	1070 mm	1270 mm	820 mm	1070 mm	1270 mm
BASE WIDTH	450 mm	650 mm	750 mm	750 mm	600 mm	650 mm	700 mm	600 mm	650 mm	700 mm
ACTIVITY TRAY DEPTH	No Tray	400 mm	450 mm	500 mm	400 mm	450 mm	500 mm	400 mm	450 mm	500 mm
ACTIVITY TRAY WIDTH	No Tray	525 mm	575 mm	625 mm	525 mm	575 mm	625 mm	525 mm	575 mm	625 mm
KEY FEATURES	Chest band Multi-adjustable pelvic band, Height & angle adjustable kneeplate, Multi-adjustable sandals Portable	Angle-adjustable tray with activity bowl Chest support with laterals Belt, knee cups & sandals in various sizes available Mobile base with lockable castors			Angle-adjustable tray with activity bowl Chest support with laterals Belt, knee cups & sandals in various sizes available Mobile base with lockable castors					
WIPEABLE PADDING	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
REMOVABLE PADDING	No	No	No	No	No	No	No	No	No	No
COLOUR OPTIONS	Meadow	Meadow	Mondeo Blue	Mondeo Blue	Meadow	Mondeo Blue	Mondeo Blue	Meadow	Mondeo Blue	Mondeo Blue

Information correct as of September 2022. From time to time, product specifications are changed by manufacturers without prior notice
Transit tie down is available on all chairs for transport purposes.

SUPINE/PRONE STANDERS

SQUIGGLES+



MYGO STANDER



HORIZON STANDER



	SQUIGGLES+	MYGO STANDER	HORIZON STANDER
SIZES	One size	1	2
AGE (APPROX.)	1 to 5	4 to 10	8 to 14
MAX USER WEIGHT	22 kgs	50 kgs	60 kgs
UPRIGHT	Yes	Yes	Yes
PRONE (ANGLE ADJUSTABLE)	90-160°	90° to vert	90° to vert
SUPINE	Yes	Yes	Yes
WHEELED	Yes	Yes	Yes
HEIGHT ADJUSTABLE	Yes	Yes	Yes
CHEST PLATE HEIGHT - MIN	550 mm	760 mm	980 mm
CHEST PLATE HEIGHT - MAX	820 mm	1080 mm	1300 mm
PRODUCT DIMENSIONS			
BASE LENGTH	950 mm	978 mm	1130 mm
BASE WIDTH	590 mm	630 mm	650 mm
ACTIVITY TRAY DEPTH	460 mm	340 mm	340 mm
ACTIVITY TRAY WIDTH	560 mm	580 mm	580 mm
KEY FEATURES	<p>Hip abduction of up to 60° in prone, upright and supine standing</p> <p>Fully horizontal loading in supine</p> <p>Height-adjustable, foldable push handle for ergonomic use</p> <p>New wrap around knee supports</p> <p>Individually adjustable footplates</p>	<p>Headrest options for supine set up (Whitmeyer compatible)</p> <p>Padded clear tray for prone or supine</p> <p>Independent angle adjustment in knees up to 25°</p> <p>-10 +10° plantar/dorsiflexion</p> <p>Mobile base with lockable castors</p> <p>Manual angle adjustment</p> <p>Belt, knee cups & sandals in a various sizes available</p>	<p>Headrest options for supine set up (Whitmeyer compatible)</p> <p>Separate supine & prone trays</p> <p>Mobile base with lockable castors</p> <p>Powered angle adjustment</p> <p>Belt, knee cups & sandals in a various sizes available</p> <p>Extremely robust</p>
WIPEABLE PADDING	Yes	Yes	Yes
REMOVABLE PADDING	Yes	Yes Removable covers	Yes Removable covers
COLOUR OPTIONS	Purple	Blue	Blue



Firefly PlayPak

Playpak is a portable activity kit packed with activities for fun therapy at home or away. Once unzipped, the lightweight bag opens to reveal a colourful floor mat along with rolls, supports and wedges that can be put together in multiple ways to help children with various developmental postures.



Leckey Early Activity System

A colourful, comfortable, supportive activity system for children up to the age of 4 years. The Early Activity System, has been developed to encourage a range of activities based around five key positions: back lying (supine), tummy lying (prone), side lying, floor sitting (long sitting), and hands and knees (4 point kneeling).




Firefly Scoot 4-in-1

Scoot is designed to help children with special needs play, explore and participate in a way that best suits their physical and cognitive ability. Scoot has four configurations: Scoot, Ride, Pooosh and Crawl. Starting with Crawl (allowing kids to explore on their tummy), moving through Scoot (lets kids sit up and use their feet), Ride (lets kids propel themselves) and finally onto Pooosh (lets somebody else do the hard work).




Leckey Therapy Bench with Activity Cradle

Developed to support children in a neutral posture, this therapy bench with innovative activity cradles, facilitates child active therapy by freeing up the hands. Available in one size suitable for children aged 1-14 the Leckey Therapy Bench facilitates anterior/posterior tilt to promote activity and acts as an aid to therapists in the Segmental Assessment of Trunk control.



Daytot Joey

The Daytot Joey is an early intervention crawler designed to encourage the development of motor skills for babies at risk of delay. For use during a baby's first year at home, Daytot Joey helps enhance early motor skills, by encouraging independent movement on their tummy. Daytot Joey supports the baby, allowing them to move freely, while building the all-important foundations for their future physical and cognitive development.



Firefly Splashy

The Splashy is ideal for bathing or getting involved with messy games - as you can take it outside. Little bath buddies get more space to play together as it is really compact and sits nice and low in the bath, so it uses less space and less water.

- 5-point configurable harness to keep kids secure, that can change to 3-point
- The base unit lets you switch easily between 26 recline positions, from upright (106°) to full recline (140°)
- The foam bumpers can be used as a pommel or a support for the head, sides, or bottom as an anti-slide shelf
- Lightweight, detachable parts make it easy to clean
- The Splashy comes with a floorsitter as standard



Firefly SplashyBIG

With all of the features and benefits of Splashy, SplashyBIG opens up many more opportunities for 7-14 year olds.

- 5-point configurable harness to keep kids secure, that can change to 3-point
- SplashyBIG can help promote and encourage independence at bathtime
- The Splashy comes with a floorsitter as standard
- Suitable for 7-14 year old children with a maximum user weight of 55kg



Ottobock Robby Bath Chair

The Robby sits completely flat in the bath, meaning less water. The head, back, seat and lower leg support are separately adjustable - allowing for optimum positioning for children with more complex needs.

- Adjusts to multiple sitting and lying positions
- Extremely lightweight
- Contoured back is designed for comfort and allows for washing of hair
- Trunk belt and lower leg positioning straps are standard



Seahorse Plus Toileting & Shower Chair

This adjustable and supportive toileting & shower chair can be used over a toilet or with the removable potty. Easy tilt operation facilitates showering and bathing. Comes with height adjustable flip-up footrests with foot-straps. The small version is suitable for children aged approx. 5-10 years up to 64kg and the medium version is suitable for children aged approx. 10-15 years up to 102kg (age range is guidance only).

- Gas Assisted Tilt in space
- Up to 40 degree tilt
- 100mm Height Adjustment
- Height and depth adjustable headrest
- Lap Belt included
- Will grow with the child
- Removable Armrests
- Flip-up height adjustable footrests
- Removable Pommel
- Removable Potty



Firefly GottaGo Portable Toilet Seat

The Firefly GottaGo is a lightweight, sturdy, portable toilet seat that makes toileting easier at home and away. Super easy to set up and put away - in just 3 simple steps. The design facilitates a supported squat posture, making toileting easier.

- Suitable for children aged 2-9 years with a maximum user weight of 30kg (66lbs).
- GottaGo frame has 3 height settings: potty height, standard toilet height and raised toilet height.
- Warm-touch, wipe clean materials guarantee a comfortable and snug fit.
- Easy to set up and to put away in just 3 simple steps.
- The entire product weighs less than 10kg



GottaGo - Extra Toilet Seat Only

Seat only for the Firefly GottaGo Portable Toilet Seat. Please note this is only the seat insert shown in the image.

- Seat only for the Firefly GottaGo Portable Toilet. Both sizes of seats fit on the one sized frame from the Firefly GottaGo Toilet Seat Set.



Nuvo Petite Children's Toilet Platform

This platform will enable you to sit comfortably on a range of different toilet heights, with your feet supported. The platform has a smooth and level step up with a slip resistant surface. Smooth handles and non-slip feet offer maximum security.

- Height adjustable platform with slip resistant mat
- Lightweight, powder coated steel frame makes it sturdy yet easy to move
- Non-slip rubber ferrules on the feet ensure additional safety
- Easy assembly without tools



Columbia Toilet Support - Lo Back

Wrap around support for independent toileting. This supportive range of backs allows children, especially those with limited trunk stability, to be safely positioned for any length of time while toileting.

- Easy to attach - just mount the base to the toilet using the same bolts that hold the toilet seat to the bowl
- Removes easily - simply by lifting the back from the base
- Universal base frame fits onto most toilets and fits all sized backs



Columbia Toilet Support Base

Mounting base for the wrap around toilet support for children, ensuring stability and reducing the risk of falling off and to be safely positioned for any length of time while toileting.

- Easy to attach - just mount the base to the toilet using the same bolts that hold the toilet seat to the bowl
- Universal base frame fits onto most toilets and fits all sized backs



Special Tomato Portable Potty Seat - Round

Flexible, durable and soft elongated portable potty seat

- Made from Soft-Touch material that is impermeable to fluids and is peel and tear resistant
- Conforms to the body and doesn't leave a "red mark" after extended periods of sitting
- Comes with its own plastic lined carry bag



Special Tomato Portable Potty Seat - Round

Flexible, durable and soft elongated portable potty seat

- Made from Soft-Touch material that is impermeable to fluids and is peel and tear resistant
- Conforms to the body and doesn't leave a "red mark" after extended periods of sitting
- Comes with its own plastic lined carry bag



Columbia Splash Guard

Toilet Seat Deflector. This soft and flexible product can be installed quickly on almost any toilet seat - all you need to do is attach it to the inside of the toilet seat and do up the two straps.

- Easy to clean
- Fits both open and closed front toilet seats
- Soft padding increases comfort
- Improves cleanliness while toileting



Columbia Padded Reducer Ring

The Columbia Padded Reducer Ring is ideal for children who may require extended toileting time. Available in two sizes, it comes complete with a padded, removable abductor / splash guard.

- Soft yet durable and offers comfort and non-slip security. The soft material conforms to body and sensitive bony protuberances. The abductor/splash guard is removable.
- Sanitary, seamless and urine-proof; the reducer ring can be used anywhere.


Leckey Sleepform

Sleepform is a comfortable, easy to use night-time postural management system. Portable and easily adjusted, each Sleepform system is made of several modular parts making it truly customisable without the requirement of cutting or complex shaping.

Use as much or as little of the system as you like to create an outcome that is just right. The Sleepform starts with a mouldable positioning cushion topped with a soft mattress to ensure airflow is achieved, thereby reducing hot spots. Cushioned chest, hip, neck and knee guides are then added to the mix. Finish with knee rolls and airflow pillows where required. This is a true recipe for success.



- Accommodates supine, prone and side sleeping positions
- Sizes available for kids of all ages. Sizes 1-2 are suitable for a cot, sizes 3-4 are suitable for a single bed
- Sleeping positions can be changed at any given time
- Sleepform mattress moulds to the client, offering a great sleep and excellent support
- Temperature Control Sheet ensures an even temperature is maintained eliminating hotspots for a better sleep (optional accessory)
- Grid sheet allows the product to be precisely repositioned time and time again
- Lightweight and easily transported
- Hand washable and easily cleaned
- Fully supported by instructional DVD, selection guides and sleep diary



PETRA RaceRunner Bike

The RaceRunner bike is perfect for recreational leisure as well as being a fun physical challenge. The bike is fun to take outdoors for activities such as basic mobility, family rides, sports or just for playing! The RaceRunner provides a challenge for the user that will stimulate their senses, and encourage them to develop their navigation and coordination skills.

- Physical activity stimulates the senses and provides endorphins
- Encourages stepping and weight bearing for those who may spend a large amount of time seated
- Encourages better posture due to bone density and muscle bulk increases
- Helps to develop general endurance abilities as well as control of trunk and limbs
- Ideal for stable surfaces in the outdoors such as streets, paths, or a track field
- Fun but challenging exercises such as the RaceRunner can have positive mental impacts on the user.
- Even wheelchair users and severely physically restricted users can enjoy cardiovascular training. As wheelchair users become runners – they literally conquer land transforming their perspective of life
- Suitable for users aged three and above





Please visit our website
www.alliedmedical.co.nz
for our full range of products



0800 31 61 81
helpis@alliedmedical.co.nz

June 2023. This information is correct at the time of printing.
From time to time, product specifications are changed by manufacturers without prior notice.



**ALLIED
MEDICAL**
LIMITED