

Allied Medical

Paediatric Solutions



**ALLIED
MEDICAL**
LIMITED

Product Range

BUGGIES 03

WHEELCHAIRS 06

POWERCHAIRS 9

WHEELCHAIR SEATING 11

SEATING SOLUTIONS 14

WALKING SOLUTIONS 17

STANDING SOLUTIONS 19

EARLY INTERVENTION 21

BATHING 22

TOILETING 24

OTHER 27



Firefly WeGo Buggy

A lightweight, portable and snazzy-looking pushchair with a few hidden surprises. The WeGo is made for GoTo Seat users, as it slots into the buggy easily.

- Easily foldable and setup in one simple movement
- Wheels suitable for rough terrain
- WeGo reclines back 140° so your child can lie back for a little snooze, for changing or feeding
- Comes standard with a full accessory pack including a sun canopy, rain cover, foot muff, mosquito net, shopping basket underneath seat, bottle holder and wrist strap
- Compatible with the Firefly GoTo Seat



Special Tomato EIO Push Chair

A stylish, special needs stroller unmatched in features, quality and price. It offers a mild to moderate positioning seat matched with a lightweight stroller designed for children from early intervention to pre-adolescent ages.

- Easy to fold and transport
- Lightweight aluminium frame
- Calf panel and seat back to lay your child flat
- Shopping basket and sun shade
- Significant growth included
- Compatible with Special Tomato Soft-Touch Liners, size 1 and 2



Special Tomato Jogger Stroller

This active stroller provides comfort and support to children with special needs while riding on all terrain. It is an easy-to-fold stroller which allows for easy storage wherever you take it.

- Easy to fold and transport
- Lightweight aluminium frame
- Functional rear suspension
- Padded seat can be reclined from 10 degrees to 50 degrees
- Compatible with Special Tomato Soft-Touch Liners, size 1 and 2



Ottobock Eco Buggy

Simple buggy for children with mild to moderate needs. Thanks to its clever design and lightweight build, the Eco Buggy is a breeze to transport as it folds very easily.

- Small lightweight buggy
- Back angle adjustment 20 degrees (90 - 110 degrees)
- Folding frame, allows it to fold to a compact size to easily fit into most vehicles
- Footrest is height adjustable from 200 - 340 mm (distance from seat edge to footrest)
- Swivel front castors
- Spring mounted front and rear wheels
- Foot brake on back wheels
- 14" x 12" Seat Size



Patron Corzino XCountry

A semi-adjustable rehabilitation buggy available in four sizes, with a wide variety of options, and a high user weight limit of up to 75kg. This buggy is equipped with an optimised folding system, a padded ergonomic seat cover and a light frame with 5-point fixation ensuring safety for your child and easy handling.

- Well padded ergonomic seat cover
- Comfortable side protection
- Footrest height adjustable and foldable (angle adj.)
- Backrest angle adjustable (2 positions)
- Backrest height extension (on size 38 and 42 only)
- Available in 4 sizes to accommodate children aged 1-15 years



Patron Piper Comfort

This buggy is ideal for short term transportation and activities for your child. It has a fully foldable undercarriage made from lightweight aluminium profiles.

- Extended reclining backrest
- Height adjustable footrest that can be detached
- 5-point immobilisation harness with crotch belt and shoulder padding
- Detachable wheels and forks
- One-piece foldable handlebar



MyWam Mouse

The **MOUSE** stroller is designed for the smallest children who are unable to walk or sit in a correct position as a result of injuries, diseases or motor dysfunctions. The shock-absorbing suspension on all four, large, gel-filled wheels together with the tilt-foot facilitate travel on uneven terrain.

- Backrest tilt and recline feature as well as adjustable torso and head pads so you may easily adjust MOUSE to the required position
- Suitable for children at a height of 80-110 cm and a weight of 30 kg
- The Reducer allows you to quickly and easily extend the width and depth of the seat, allowing the Mouse to grow with the child



MyWam Mewa

The **MEWA** stroller is designed for children who are unable to walk or sit in a correct position as a result of injuries, diseases or motor dysfunctions. The shock-absorbing suspension on all four, large, gel-filled wheels, together with the tilt-foot, facilitate travel on uneven terrain.

- Modern and stylish design
- Backrest tilt and recline feature as well as adjustable torso and head pads so you can easily adjust MEWA to the required position
- Suitable for children at a height of 90-125 cm and a weight of up to 40 kg
- Reversible seat and an array of optional accessories



Karma Flexx Junior

The Karma Flexx Junior wheelchair is an extremely flexible and popular wheelchair solution to suit a range of customer needs. Cost effective and feature packed, with everything that is orange being adjustable, making it easy to use.

- Configurable options include seat depth, seat to floor heights, backrest heights, castor fork angles and centre of gravity. Available in Transit or Self-Propelling setup
- Swing away footrests with flip-up footplates
- Can be set up to provide optimal comfort and performance
- Adjustable backrest angle

Maximum user weight: 80kg | Min seat to floor height: 15.5"
Seat Width & Depth: 10" - 14" | 11" - 15"



Ki Mobility Little Wave Klik XP

A wheelchair designed for the next generation. A system so easy to adjust it promotes proper positioning and wheel access so your child gains the active independence and mobility they want. The Little Wave Klik XP has an adjustment system designed to provide accurate, easy growth options.

- An intuitive 'click based' adjustment system designed to provide accurate, easy growth options
- Adjustable dynamic 5th wheel or anti-tippers
- Reversible front castor arms for added stability
- Seat belt clamp for secure positioning

Maximum user weight: 75kg | Min seat to floor height: 11.5"
Seat Width & Depth: 8" - 18"





Ki Mobility Rogue XP

Rogue XP is the ultimate performance chair; encouraging kids to embrace their personal style at the age when they care the most. This chair has the necessary built-in growth, with no need to compromise on fit, function, ride or style.

- Adjustable fixed frame for growth - up to 3"
- Creative freedom - ability to customise your chair
- Maximum adjustment and growth ability
- Integrated swing-away armrests
- Combination of rigidity and reduced vibration dampening for an ultra-responsive ride
- Adjustable dynamic 5th wheel or anti-tippers
- Suitable for little ones and teenagers alike

Maximum user weight: 90kg | Min seat to floor height: 13.5"
Seat Width & Depth: 10" - 18" | 10" - 20"



Ki Mobility Little Wave Flip XP

The Little Wave Flip XP is a wheelchair designed to be as unique as your child. It is designed for the smallest of riders with some of the biggest advancements in the industry. From its flat folding design to its use of Taper Lok Technology, the Flip is simply the most advanced paediatric tilt-in-space available today.

- Integrates design elements from strollers into the durability of a wheelchair
- Seat height, seat depth, and seat width have never been so easy to change
- Hands-free attendant wheel lock
- Multi-customisable options - refer to script for full details

Maximum user weight: 75kg | Min seat to floor height: 13.5"
Seat Width & Depth: 10" - 18" | 12" - 20"



Ki Mobility Little Wave Arc

The Little Wave Arc combines functionality, comfort, and style for paediatric users. Its ultra-lightweight frame ensures easy mobility and adjustability for growth, without sacrificing durability. With extensive customisation options, it offers independence and comfort for young users.

- Combines the lightweight ease of a stroller with the robustness of a wheelchair
- Easily adjustable seat height, depth, and width for growing children
- Convenient safety features including hands-free attendant wheel lock and flip-up anti-tippers
- New weight-shift tilt mechanism for added stability

Maximum user weight: 75kg | Min seat to floor height: 14.5"
Seat Width & Depth: 10" - 18" | 12" - 20"



Ottobock Kidevo Mini & Mini.t

The Kidevo Mini and Mini.t wheelchairs are designed for children's active and independent mobility, offering customisable seating and adjustable features to grow with them. These wheelchairs provide optimal support and comfort with options for robust foot plates, child-friendly brakes, and functional seat tilts.

- Encourages independent movement and optimal driving behavior from an early age
- Fully adjustable in three dimensions to accommodate growth and changing needs
- Robust Foot Plate: Ensures optimal foot positioning with a sturdy, flat one-piece foot plate and additional positioning option

Maximum user weight: 50kg | Tilt Range (Mini.t Only) -5 to +40 degrees | Seat Width & Depth: 6" - 13" | 6" - 14"



Ottobock Kidevo Prime & Prime.t

The Kidevo Prime and Prime.t wheelchairs are tailored for older children's active and independent mobility, offering customisable seating and adjustable features to grow with them. These wheelchairs provide optimal support and comfort with options for robust foot plates, child-friendly brakes, and functional seat tilts.

- Designed to encourage independent movement and optimal driving behavior
- Fully adjustable in three dimensions to accommodate growth and changing needs
- Robust Foot Plate: Ensures optimal foot positioning with a sturdy, flat one-piece foot plate and additional positioning options

Maximum user weight: 75kg | Tilt Range (Prime.t Only) -10° to +50° | Seat Width & Depth: 10" - 20" | 8" - 20"



Ottobock Kidevo Adapt

The Kidevo Adapt is a versatile mobility base designed for children, offering customisable features to support active and independent movement. With adjustable seat tilt, reliable anti-tippers, and a practical luggage bin, it ensures comfort and safety for everyday use.

- Designed to enhance the child's active and independent movement
- Reliable Anti-Tippers: Ensures safe driving and maneuvering
- Options for drum brake or combi-drum brake systems
- Angle-Adjustable Push Bar: Provides comfortable handling for caregivers.

Maximum user weight: 75kg | Tilt Range -10 to +50 degrees | Seat Width Range: 12" - 24"



Ottobock Skippi

Designed specifically for children, this fun rear wheel drive powerchair is perfect for active kids. It is compact and colourful, but still retains large size drive wheels and castors so kids can explore their environment with the push of a joystick.

- Robust paediatric powerchair that will perform well both indoors and outdoors with more active users (has relatively high ground clearance)
- Seat frame can accommodate the width and depth growth
- Can be disassembled in just a few steps to fit easily into the boot of a small car for transport

Maximum user weight: 50kg | Min seat to floor height: 17"
Seat Width & Depth: 10" - 15" & 11" - 15"



Jazzy 623 2.0

With highly reliable, efficient motors that deliver a higher range and improved responsiveness than the previous model, the Jazzy 623 2.0 is engineered to meet the performance needs of active users.

- Compact, mid-wheel drive powerchair
- TRU-Balance 3 power positioning system
- VR2 controller to deliver smooth, nimble manoeuvrability with good indoor and outdoor performance and a tight turning radius
- Simple to operate hand control

Maximum user weight: 136kg | Min seat to floor height: 17.5"
with tilt | Seat Width & Depth: 12" - 22" or 12" - 22"



Quantum Edge 3 Stretto

The Edge 3 Stretto, with Q-Logic 3 expandable electronics, is the most capable compact power base in New Zealand. Its extremely narrow footprint makes for outstanding manoeuvrability, which is ideal for children, teenagers and small adults.

- Independent SRS (Smooth Ride Suspension) making for an incredibly easy ride
- The Stretto is available with 12.5" or 14" drive wheels, and optional 12" of iLevel power adjustable seat height
- Available with either the TRU-Balance 3 or 4 power positioning system

Maximum user weight: 136kg | Seat Width & Depth: 12" - 22"
or 12" - 22"



Ottobock Skippi Powerchair

Varilite Junior Back Supports



Mid



Deep

Lightweight and easy to use, the Varilite Junior Back Support is perfect for children and families on the go. The compact, minimal design directs focus to the child and not the equipment. Available in two styles - Mid and Deep (which offers lateral support).

Specification	Mid	Deep
Back Height	9" (23.5cm)	9" (23.5cm)
Outside Width of Shell	8.5" (21.8cm)	10" (25.1cm)
Inside Width of Support Surface	7" (19cm)	9" (22.6cm)
Fits Wheelchair Sizes	9" - 11" (23 - 28cm)	9" - 11" (23 - 28cm)
Weight	1.1kg	1.3kg

Ride Designs Java Decaf Back



The Java Back's patented biomechanical design provides the flexibility to easily and independently optimise pelvic support and balance the trunk, supporting the normal rhythm and movement of the body.

Specification	Java Decaf
Back Height	10"-15" (25 -38cm)
Fits Wheelchair Sizes	10" - 15" (25-36cm)
Weight	Up to 1.6kg

Ki Mobility Axiom for Kids Posterior Lateral Back



Unparalleled support and comfort, tailored for the unique needs of young wheelchair users. Its adjustable design and superior materials offer optimal posture support, enhancing mobility and daily independence.

Specification	Posterior Lateral Back
Total Product Width	10"-16" (23 - 40.5cm)
Total Product Depth	10" - 18" (23 - 45.5cm)
User Weight Limit	75 kg

Ki Mobility Axiom for Kids Adjustable Lateral Back



This innovative backrest offers a custom fit that grows with your child, ensuring their wheelchair evolves with their needs. It has a planar back with option for fixed or swing-away lateral supports.

Specification	Adjustable Lateral Back
Total Product Width	10"-16" (23 - 40.5cm)
Total Product Depth	10" - 18" (23 - 45.5cm)
User Weight Limit	75 kg

Ki Mobility Axiom for Kids Adjustable Contoured Lateral Back



This backrest provides customised contoured support, ensuring optimal comfort and posture alignment as your child grows. It contours the back with adjustable contoured lateral supports.

Specification	Contoured Lateral Back
Total Product Width	10"-16" (23 - 40.5cm)
Total Product Depth	10" - 18" (23 - 45.5cm)
User Weight Limit	75 kg

Paediatric Curve Cushion



A paediatric cushion that utilises moulded, high-density seat design to properly position the legs, hips, and spine to provide optimal support in a wheelchair. We recommend purchasing a spare cover with all cushion purchases.

Maximum user weight: 50kg | Min seat to floor height: 17"
Seat Width & Depth: 10" - 15" & 11" - 15"

Ki Mobility Axiom for Kids P Cushion



This Positioning cushion is specifically designed for children's unique positioning needs, adapting to their varying shapes and sizes. This paediatric cushion promotes better clinical outcomes through its tailored design, featuring different foam softness and cushion shapes based on size. It aims to enhance comfort, skin integrity, and stability.

Widths: 8" - 16" | Depths: 8" - 18"

Ki Mobility Axiom for Kids SP Cushion



This Skin Protection and Positioning cushion is specifically designed for children's unique positioning needs, adapting to their varying shapes and sizes. This paediatric cushion promotes better clinical outcomes through its tailored design, featuring different foam softness and cushion shapes based on size. It aims to enhance comfort, skin integrity, and stability.

Widths: 8" - 16" | Depths: 8" - 18"

Varilite Junior Cushion



Designed to provide a comfortable foundation for exploration and adventure. Lightweight multi-stiffness foam promotes immersion and positioning while the washable, water-resistant cover lets kids be kids.

Widths: 10" - 12" | Depths: 10" - 12"
Adult sizes available - Reflex cushion

Ride Designs Forward Cushion



The Forward cushion is all about contour and helps prevent sliding while providing superior support and protection. Ride's approach redistributes the forces of pressure and shear from high risk anatomic areas while selectively applying the forces to areas more tolerant. The contours are specifically optimised for a paediatric pelvis.


Widths: 10" - 12" | Depths: 10" - 16"



Special Tomato Soft-Touch Liner Seat and Back

The Soft-Touch Liner offers an affordable alternative to other seating systems. It eliminates the need to transport heavy seating systems. The contoured seat and back cushions promote good posture while sitting in a standard seating system for children with mild to moderate positioning needs.

- Made of a soft-to-touch, easy clean material that is latex-free and anti-microbial
- Can be attached to most seats including but not limited to the following: strollers, highchairs, dining room chairs, restaurant chairs and classroom chairs
- Easy storage as the cushions can be folded together



Special Tomato Soft-Touch Sitter

The Special Tomato Soft-Touch Sitter provides a comfortable and functional seating solution. It has a strong exterior that will stand up, whilst ensuring a comfortable seat for your child.

- Impermeable to fluids
- Peel and tear resistant
- Latex-free material with built in anti-microbial protection
- Soft contoured seat with anti-thrust seating surface
- Built-in abductor
- Contoured head and lateral supports
- Attachment straps for securing to either the Floor Sitter or the Mobile Base



Special Tomato CatchUp Chair

The CatchUp Chair is a stable, highly adjustable, mobile seat. It is easily configured without tools to provide the well fitting seat at the perfect height for a variety of tabletop surfaces. It is appropriate for individuals who require mild to moderate postural support. The CatchUp Chair fits children and adults weighing up to 113kgs.

- Compact design makes it an excellent choice for schools, therapy clinics, and hospitals as well as homes
- Simple, tool-free adjustments
- Pelvic Belt and Chest Harness stabilise a well-aligned posture



Firefly GoTo Seat

The GoTo Seat is a portable support seat for children with extra postural needs. It can be taken anywhere, which makes getting out and about much easier. The GoTo Seat fits easily into the WeGo Buggy, making it that much easier to go places.

- Two sizes available
- Suitable for attachment to most standard chairs in the home, restaurant or classroom
- Adjustable headrest and laterals that adjust for growth
- Two headrest choices: Standard or Advanced
- Optional floor sitter accessory has 26 recline positions, from upright to right back



Leckey PAL Seat

The PAL Seat is designed to give improved stability, which helps reduce fatigue and allows the child longer periods of concentration and fine motor activity. The simplicity of the PAL Seat enables the therapist to assess the child's needs instantly. It is a really simple, fast way to help kids join in with learning.

- Accessories include two point lap belt, side pads for lateral stability, tray, footplate and lockable ski castors
- Height adjustable armrests, legs and footplate
- Four sizes accommodating kids from 1 - 12 years
- Angle adjustable backrest



Leckey BeMe

The infinitely configurable BeMe from Leckey, is a revolutionary seating system that puts a child's development and independence at the heart of its design. BeMe is a highly personalised seating system available in three sizes and spec'd to perfection.

- Growable contoured cushion with high-stretch foam and fabric, to maximise pelvic stability
- Interchangeable parts (including seats, back supports, medial, lateral, proximal and distal supports), click and switch with total ease
- Dynamic backrest is air pressure controlled which regulates resistance as the child extends and is speed adjustable




Buckle Guard

A simple solution for the little (and not so little) Houdini's who are able to escape from car safety harnesses and may endanger themselves or others in a vehicle. This strong and difficult to break cover slides over the seatbelt buckle to prevent the wearer from releasing the belt themselves. A caregiver can then safely and easily release the belt by inserting a key or coin into a slot when the destination is reached.

- Designed to go over most car seat buckles - either side release or top release
- Dimensions: 64 mm W x 64 mm H x 42 mm D
- Comes complete with a seatbelt cutter and an instruction sheet.
- Manufactured from clear and durable ABS plastic.
- Includes a wedge spacer to accommodate for many smaller sized buckles
- Webbing cutter is provided in case of an emergency.
- Australian design





Firefly Upsee

The Firefly Upsee is an upright mobility device designed to let children stand, walk - and dance - with the help of an adult.

- Sizes available: Small, Medium or Large
- Age: 2-4 / 4-6 / 6-8 years
- Maximum user weight: 15kg / 20kg / 25kg



Ottobock Nurmi Neo Posterior Walker

The Ottobock Nurmi Neo is a highly manoeuvrable kids posterior walker designed to promote an upright posture and natural walking pattern.

- Height-adjustable grip bar
- Solid locking mechanism to stop slipping
- Lightweight and foldable
- Reverse roll locks on the rear wheels
- Puncture proof low friction castors
- Stable base
- Swivelling front castors
- Flip up seat, hip laterals, anti-tippers and gutter arms available as optional accessories



Leckey MyWay+

MyWay+ is an open frame walker which cleverly supports children in an upright position, improving stepping and maximising opportunities for exploration, interaction and participation. A bespoke frame combined with a wraparound hoistable harness, and a wide range of accessories, ensures the MyWay+ can be tailored to meet every need and ability.

- MyWay+'s improved larger castors indoor/outdoor enable children to be active in lots of different environments
- The revolutionary MyWay+ harness has been designed to give 360° wraparound support to promote optimal stepping pattern
- Hoistable For Easier Transitions
- Setup for children with simple or complex needs
- Pedal attachment available to purchase separately

Sizes available: 1-3 | Age: 1-16 years

Maximum user weight: 25-80kg



Meywalk 4

The Meywalk 4 is a high quality, Danish designed walker that offers a suspended spring system that stimulates the muscles and movements to encourage a more natural gait pattern. Plenty of different adjustment options available.

- The body posture of the user can be set in any position from vertical to 15° forward tilt
- Suitable option for both indoor and outdoor use
- Suitable for both adolescents and adults, with two larger size options available
- Differing spring resistance options
- Chest block available in multiple size options

Sizes available: Mini - Large | Maximum user weight: 80-130kg



Leckey Totstander

A pre-school standing frame that is adjustable and extremely portable. It is a junior version of the Freestander and is designed to support children in the early stages of standing.

- Simple construction, allowing adjustment and portability
- Allows the child to be positioned close to work tables or benches, whilst keeping them at a similar height to their peers
- Therapists can use the Totstander as an assessment tool to determine a child's ability to stand and bear weight
- Side tubes accommodate growth
- Adjustable pelvic band, knee plate, padded knee cups and sandals that velcro in place

Age: 1-5 years | Maximum user weight: 20kg



Leckey Freestander

An upright stander for children and teens, created to support children with mild to moderate levels of physical involvement. The Freestander is the next size up from the Totstander.

- A lower baseboard allowing for quick and easy transfers from chair to stander
- Pelvic band and optional padded pelvic plate are designed to help promote mid-line symmetry
- Adjustable flexible laterals, padded chest plate, knee plate and padded knee cups, sandals, pelvic band and activity tray

Sizes available: 1-3 | Age: 1-5, 4-9, 9-17 years

Maximum User Weight: 60, 80 or 90kg



Leckey Pronestander

The Leckey Prostander assists in the promotion of weight bearing and mid-line symmetry. The user can be positioned in a fully upright position or at any one of four prone angles.

- Four prone angles are designed to allow weight bearing to be gradually increased as the child's tolerance increases
- Hinged footplate permits wheelchair access and assisted transfer

Sizes available: 1-3 | Age: 1-5, 4-10, 9-18 years

Maximum user weight: 60, 80 or 90kg



Leckey Squiggles+ Standing Frame

The Squiggles+ Stander is an extremely versatile 3-in-1 stander, offering prone, upright and supine standing in one product.

- 3 in 1 design offering prone, upright and supine standing
- Up to 60 degrees of bilateral hip abduction
- Advanced knee support design
- Independent footplates to accommodate leg length discrepancy
- Fully horizontal loading in supine
- Activity tray with height adjustability
- Height-adjustable push handle

Age: 1-5 years | Maximum user weight: 22kg



Leckey Mygo Stander

As the child grows out of the Squiggles Stander, they can move into the Mygo Stander to continue their standing therapy. It encourages natural standing posture and has a range of functions for the growing child.

- This 3-in-1 stander offers prone, upright and supine standing in one product
- It can accommodate up to 25 degree fixed hip and knee contractures
- Supine and prone angles can be changed from 0 to 80 degrees
- Easily adjustable

Sizes available: 1-2 | Age: 4-10 or 8-14 years

Maximum user weight: 50 or 80kg



Leckey Horizon Stander

This powered stander provides gradual, appropriate weight bearing support in either prone, supine or upright positions. Robust and durable for users aged 4 to adult.

- Offers 90 degrees of adjustment from horizontal to freestanding position in both prone and supine positions dependent on setup
- Hand held, easy to use controls for the powered tilt mechanism

Sizes available: 1-3 | Age: 4-10, 8-14 or 14+ years

Maximum user weight: 50, 80 or 100kg



Firefly PlayPak

Playpak is a portable activity kit packed with activities for fun therapy at home or away. Once unzipped, the lightweight bag opens to reveal a colourful floor mat along with rolls, supports and wedges that can be put together in multiple ways to help children with various developmental postures.



Leckey Early Activity System

A colourful, comfortable, supportive activity system for children up to the age of 4 years. The Early Activity System, has been developed to encourage a range of activities based around five key positions: back lying (supine), tummy lying (prone), side lying, floor sitting (long sitting), and hands and knees (4 point kneeling).



Firefly Scoot 4-in-1

Scoot is designed to help children with special needs play, explore and participate in a way that best suits their physical and cognitive ability. Scoot has four configurations: Scoot, Ride, Pooosh and Crawl. Starting with Crawl (allowing kids to explore on their tummy), moving through Scoot (lets kids sit up and use their feet), Ride (lets kids propel themselves) and finally onto Pooosh (lets somebody else do the hard work).



Leckey Therapy Bench with Activity Cradle

Developed to support children in a neutral posture, this therapy bench with innovative activity cradles, facilitates child active therapy by freeing up the hands. Available in one size suitable for children aged 1-14 the Leckey Therapy Bench facilitates anterior/posterior tilt to promote activity and acts as an aid to therapists in the Segmental Assessment of Trunk control.



Daytot Joey

The Daytot Joey is an early intervention crawler designed to encourage the development of motor skills for babies at risk of delay. For use during a baby's first year at home, Daytot Joey helps enhance early motor skills, by encouraging independent movement on their tummy. Daytot Joey supports the baby, allowing them to move freely, while building the all-important foundations for their future physical and cognitive development.



Firefly Splashy

The Splashy is ideal for bathing or getting involved with messy games - as you can take it outside. Little bath buddies get more space to play together as it is really compact and sits nice and low in the bath, so it uses less space and less water.

- 5-point configurable harness to keep kids secure, that can change to 3-point
- The base unit lets you switch easily between 26 recline positions, from upright (106°) to full recline (140°)
- The foam bumpers can be used as a pommel or a support for the head, sides, or bottom as an anti-slide shelf
- Lightweight, detachable parts make it easy to clean
- The Splashy comes with a floorsitter as standard



Firefly SplashyBIG

With all of the features and benefits of Splashy, SplashyBIG opens up many more opportunities for 7-14 year olds.

- 5-point configurable harness to keep kids secure, that can change to 3-point
- SplashyBIG can help promote and encourage independence at bathtime
- The Splashy comes with a floorsitter as standard
- Suitable for 7-14 year old children with a maximum user weight of 55kg



Ottobock Robby Bath Chair

The Robby sits completely flat in the bath, meaning less water is required. The head, back, seat and lower leg support are separately adjustable - allowing for optimum positioning for children with more complex needs.

- Adjusts to multiple sitting and lying positions
- Extremely lightweight
- Contoured back is designed for comfort and allows for washing of hair
- Trunk belt and lower leg positioning straps are standard



Seahorse Plus Toileting & Shower Chair

This adjustable and supportive toileting & shower chair can be used over a toilet or with the removable potty. Easy tilt operation facilitates showering and bathing. Comes with height adjustable flip-up footrests with foot-straps. The small version is suitable for children aged approx. 5-10 years up to 64kg and the medium version is suitable for children aged approx. 10-15 years up to 102kg (age range is guidance only).

- Gas Assisted Tilt in space
- Up to 40 degree tilt
- 100mm Height Adjustment
- Height and depth adjustable headrest
- Lap Belt included
- Will grow with the child
- Removable Armrests
- Flip-up height adjustable footrests
- Removable Pommel
- Removable Potty



Firefly GottaGo Portable Toilet Seat

The Firefly GottaGo is a lightweight, sturdy, portable toilet seat that makes toileting easier at home and away. Super easy to set up and put away - in just 3 simple steps. The design facilitates a supported squat posture, making toileting easier.

- Available in two sizes, suitable for children aged 2 - 5 and 4 - 9 years, with a maximum user weight of 30kg (66lbs)
- GottaGo frame has 3 height settings: potty height, standard toilet height and raised toilet height
- Warm-touch, wipe clean materials guarantee a comfortable and snug fit
- Easy to set up and to put away in just 3 simple steps
- The entire product weighs less than 10kg

GottaGo - Extra Toilet Seat Only

Seat only for the Firefly GottaGo Portable Toilet Seat. Please note this is only the seat insert shown in the image.

- Seat only for the Firefly GottaGo Portable Toilet. Both sizes of seats fit on the one sized frame from the Firefly GottaGo Toilet Seat Set

Nuvo Petite Children's Commode

Ideal for a child who may have difficulty accessing or using a standard toilet.

- Height and width adjustable
- Easy to clean seat and backrest
- Removable seat for thorough clean
- Discreet removable potty





Nuvo Petite Children's Toilet Platform

This platform will enable you to sit comfortably on a range of different toilet heights, with your feet supported. The platform has a smooth and level step up with a slip resistant surface. Smooth handles and non-slip feet offer maximum security.

- Height adjustable platform with slip resistant mat
- Lightweight, powder coated steel frame makes it sturdy yet easy to move
- Non-slip rubber ferrules on the feet ensure additional safety
- Easy assembly without tools



Columbia Toilet Support - Lo Back

Wrap around support for independent toileting. This supportive range of backs allows children, especially those with limited trunk stability, to be safely positioned for any length of time while toileting. Available in small, medium and large sizes

- Easy to attach - just mount the base to the toilet using the same bolts that hold the toilet seat to the bowl
- Removes easily - simply by lifting the back from the base
- Universal base frame fits onto most toilets and fits all sized backs
- Commonly used with optional padded insert



Columbia Toilet Support Base

Mounting base for the wrap around toilet support for children, ensuring stability and reducing the risk of falling off and to be safely positioned for any length of time while toileting.

- Easy to attach - just mount the base to the toilet using the same bolts that hold the toilet seat to the bowl
- Universal base frame fits onto most toilets and fits all sized backs



Special Tomato Portable Potty Seat - Round

Flexible, durable and soft elongated portable potty seat

- Made from Soft-Touch material that is impermeable to fluids and is peel and tear resistant
- Conforms to the body and doesn't leave a "red mark" after extended periods of sitting
- Comes with its own plastic lined carry bag



Special Tomato Portable Potty Seat - Round

Flexible, durable and soft elongated portable potty seat

- Made from Soft-Touch material that is impermeable to fluids and is peel and tear resistant
- Conforms to the body and doesn't leave a "red mark" after extended periods of sitting
- Comes with its own plastic lined carry bag



Columbia Splash Guard

Toilet Seat Deflector. This soft and flexible product can be installed quickly on almost any toilet seat - all you need to do is attach it to the inside of the toilet seat and do up the two straps.

- Easy to clean
- Fits both open and closed front toilet seats
- Soft padding increases comfort
- Improves cleanliness while toileting



Columbia Padded Reducer Ring

The Columbia Padded Reducer Ring is ideal for children who may require extended toileting time. Available in two sizes, it comes complete with a padded, removable abductor / splash guard.

- Soft yet durable and offers comfort and non-slip security. The soft material conforms to body and sensitive bony protrusions. The abductor/splash guard is removable.
- Sanitary, seamless and urine-proof; the reducer ring can be used anywhere

PETRA RaceRunner Bike

The RaceRunner bike is perfect for recreational leisure as well as being a fun physical challenge. The bike is fun to take outdoors for activities such as basic mobility, family rides, sports or just for playing! The RaceRunner provides a challenge for the user that will stimulate their senses, and encourage them to develop their navigation and coordination skills.

- Physical activity stimulates the senses and provides endorphins
- Encourages stepping and weight bearing for those who may spend a large amount of time seated
- Encourages better posture due to bone density and muscle bulk increases
- Helps to develop general endurance abilities as well as control of trunk and limbs
- Ideal for stable surfaces in the outdoors such as streets, paths, or a track field
- Fun but challenging exercises such as the RaceRunner can have positive mental impacts on the user.
- Even wheelchair users and severely physically restricted users can enjoy cardiovascular training. As wheelchair users become runners – they literally conquer land transforming their perspective of life
- Suitable for users aged three and above





Please visit our website
www.alliedmedical.co.nz
for our full range of products



0800 31 61 81
helpis@alliedmedical.co.nz

August 2024. This information is correct at the time of printing.
From time to time, product specifications are changed by manufacturers without prior notice.

