



Allied Medical

Short Term Loan and

Rehabilitation Solutions



**ALLIED
MEDICAL**
LIMITED

Product Range

MANUAL WHEELCHAIRS 03

ADULT WHEELCHAIRS 03

PAEDIATRIC WHEELCHAIRS 06

BARIATRIC WHEELCHAIRS 07

WHEELCHAIR CUSHIONS 08

WHEELCHAIR ACCESSORIES 10

POWERCHAIRS 13

HOME HEALTHCARE 14

BATHING 14

TOILET AIDS 16

LIVING AIDS 18

WALKING AIDS 25

ROLLATORS 25

CANES 27

CRUTCHES 30

BRACING & SUPPORTS 32

EXERCISE PRODUCTS 51

PAIN MANAGEMENT 52



Karma Agile Self Propelling Wheelchair

KMAG1616P (16" x 16")
KMAG1816P (18" x 16")

KMAG1818P (18" x 18")
KMAG2018P (20" x 18")

- Standard flip-back armrests
- 80° removable footrests to ensure safe transfer in and out of the wheelchair
- Compact and lightweight for easy transport
- Quick-release rear wheels for lighter carrying weight and smaller storage space
- Built with aircraft-grade aluminium
- Elevating legrests and anti-tippers available on request
- Transit kit available on request
- Maximum unit weight 8.8 kg
- Maximum user weight 100 kg



Karma S-Ergo Self Propelling

KM125Q2416 (16" x 17")

KM125Q2418 (18" x 17")

- Flip-back, ergonomically shaped armrests designed to offer comfort and support
- Extended armrests allow for wider seat area and minimise the overall width
- Ergo push rim that aids wheelchair propulsion
- Swing-away footrests swing both inwards & outwards
- Better side support of legs from "tube-in-centre" footplate design
- Aegis upholstery anti-bacterial shielding
- Optional anti-tippers
- Quick release 24" rear wheels
- Maximum unit weight 13.2 kg
- Maximum user weight 115 kg



Karma Lightweight Wheelchair Self Propelling or Transit

KM8020QDXX (Self Propelling)

KM8520FXX (Transit)

- Versatile and easily adjustable - available in 7 sizes
- Lightweight aircraft grade aluminium frame
- Folding frame with reinforced cross bracing
- Removable rear wheels on self propelling model
- Flip back, height adjustable armrests: desk or full length arm pad
- Height adjustable backrest
- Push to lock brakes
- Pelvic positioning belt
- 70-degree legrests with calf strap
- Multiple-position rear wheel mount
- Unit weight 15.5 kg
- Maximum user weight 130 kg (14" - 18" width frames), 160 kg (20" - 22" width frames)





Karma Flexx Wheelchair

KMFLEXXSPXX (Self Propelling)

KMFLEXXSTRXX (Transit)

Seat Width from 14" to 18" and Depth 16" to 18"

- Fixed angle adjustment in back canes -4 to +16° in set increments
- 15.5" minimum seat to floor heights (with 22" wheels and 5" castors)
- Rise and fall in the front castors and 25 rear axle positions
- 80° swing in / swing out legrest hangers and flip back height adjustable armrests
- Options including elevating legrests, 90° hangers, stump supports and attendant brakes available
- Maximum user weight 130kg



Karma Flexx HD Wheelchair

KMFLEXXHDSPXX (Self Propelling)

KMFLEXXHDTRXX (Transit)

Seat Width from 20" to 24" and Depth 18" to 20"

- Fixed angle adjustment in back canes -4 to +16° in set increments
- Rise and fall in the front castors and 25 rear axle positions
- 80° swing in / swing out legrest hangers with aluminium footplates
- Flip back height adjustable armrests
- Options including elevating legrests, 90° hangers, stump supports and attendant brakes available
- Maximum user weight 170 kg
- Tie downs are not available



Karma Standard Lightweight Self-Propelling Reclining Wheelchair

KM5000-16 (16"x 16")

KM5000-20 (20"x 18")

KM5000-18 (18"x 18")

- Reclining backrest
- Headrest with movable cushion to support the head
- Quick release, removable rear wheels for transport
- Fully removable full length armrests
- High mounted, push-to-lock brakes
- Elevating swing-away leg rests, with calf pad
- Reinforced cross bracing – strong but easily folded
- Multiple position rear wheel and castor mounts to optimise centre of gravity and seat to floor height adjustments
- Maximum user weight 100 kg



Karma Ergo Lite Transit Ultra-Light Wheelchair

KM2501-1616 (16" x 16")

KM2501-1816 (18" x 16")

- Folding backrest for ease of storing in vehicles
- Tall backrest for added support and ease for caregiver to push wheelchair
- Extremely lightweight alloy frame, maximum unit weight 8.8 kg
- High mounted, push to lock brakes
- 70 degree leg rests, with calf strap
- Adjustable length leg rests
- Low seat to floor height for ease of transfer
- Modern, high quality anti-bacterial upholstery
- Fixed armrests and leg rests – not able to be misplaced
- Desk armrests allowing wheelchair to easily fit under tables and desks
- Double cross bracing – strong but easily folded
- Maximum user weight 100 kg



AML Self Propelling Wheelchair

AMLWCSP16 (16" x 16")

AMLWCSP18 (18" x 16")

- Powder coated steel frame
- Detachable full armrest
- Swing-away footrest
- Steel chrome handle rim
- Plastic black footplate
- Plastic front fork
- Folds for easy storage and transport
- PU armpads for extra comfort
- Nylon upholstery
- 24" real spoke wheels with PU tires
- 8" x 1.25" front castors with PU tires
- Unit weight 19.5 kg
- Maximum user weight 135 kg



AML Transit Wheelchair

AMLWCTR18 (18" x 16")

- Powder coated steel frame
- Detachable full armrest
- Swing-away footrest
- Steel chrome handle rim
- Plastic black footplate
- Plastic front fork
- Folds for easy storage and transport
- PU armpads for extra comfort
- Nylon upholstery
- 12.5" PU mag wheels
- 8" x 1.25" front castors with PU tires
- Unit weight 19.5 kg
- Maximum user weight 150 kg



Karma Flexx Junior Self Propelling Wheelchair

KMFLEXXSPJR1011 (10" x 11")

KMFLEXXSPJR1212 (12" x 12")

KMFLEXXSPJR1414 (14" x 14")

- Configurable dimensions
- Height adjustable flip-back armrests
- Swing in and out 90° legrest, adjustable length
- Padded height and angle adjustable backrest
- Reinforced double crossbrace
- Quick release pneumatic rear wheels
- Push to lock brakes with brake extension
- 2 point lap belt
- Heel loops
- Maximum user weight limit of up to 75 kg



Karma Flexx Junior Transit Wheelchair

KMFLEXXTRJR1011 (10" x 11")

KMFLEXXTRJR1212 (12" x 12")

KMFLEXXTRJR1414 (14" x 14")

- Configurable dimensions
- Height adjustable flip-back armrests
- Swing in and out 90° legrest, adjustable length
- Padded height and angle adjustable backrest
- Reinforced double crossbrace
- Quick release pneumatic rear wheels
- Push to lock brakes with brake extension
- 2 point lap belt
- Heel loops
- Maximum user weight limit of up to 75 kg



Karma Lightweight Paediatric Self Propelling Wheelchair

KM752011 (11.5" x 11.5")

KM752013 (13.5" x 13.5")

- Easily adjustable push handle
- Quick release, removable rear wheels for ease of storing in vehicle
- Fully removable, height adjustable armrests, that allow for easier transfers
- High mounted, push to lock brakes, that can be set to numerous positions
- Swingaway 90° legrests with calf strap
- Reinforced cross bracing - strong but easily folded
- Durable upholstery with hidden screws
- Multiple-position rear wheel and castor mounts Pump and tool kit for ongoing standard maintenance



Karma BT10 Bariatric Wheelchair

KMBT102220 (22" x 20")

KMBT102420 (24" x 20")

KMBT102620 (26" x 20")

- Angle adjustable backrest (from -5 degrees to +15 degrees)
- Fully detachable flip-back armrests
- Press down, swing-away detachable leg rests with angle adjustable footplates
- Rigidising bar to add strength to frame
- Swing-away anti-tippers
- Hip belt
- Calf strap
- Dual brake system – attendant brake as well as manual brakes
- Pump and tool kit for ongoing standard maintenance
- Maximum user weight 250 kg



AML Self Propelling Bariatric Wheelchair

AMLWCSPB20 (20" x 18")

AMLWCSPB22 (22" x 18")

AMLWCSPB24 (24" x 18")

- Chrome polish finish steel frame
- Steel black side-guards
- Detachable full armrest
- Padded armpads for extra comfort
- Aluminium swing-away footrest
- Folds for easy storage
- Padded PVC upholstery
- 24" mag wheels with PU tires
- 8" x 2" front castors with PU tires
- Maximum user weight 250 kg



AML Transit Bariatric Wheelchair

AMLWCTRB22 (22" x 18")

- Powder-coated steel frame with double crossbar
- Detachable full armrest with PU armpads for extra comfort
- Swing-away footrest
- Plastic black footplate
- Aluminium front fork
- Folds for easy storage and transport
- Padded nylon upholstery
- 12" rear wheels with PU tyres
- 8" x 2" front castors with PU tyres
- Maximum user weight 230 kg



Varilite Reflex Wheelchair Cushion

VL73500114162 14"x 16"

VL73500116162 16"x 16"

VL73500116182 16"x 18"

VL73500118162 18"x 16"

VL73500118182 18"x 18"

VL73500118202 18"x 20"

VL73500120182 20"x 18"

VL73500120202 20"x 20"

- Non-adjustable skin protection cushion featuring VARILITE's unique Air-Release Device
- Achieves adjustment through a self-regulating air chamber within the cushion that adjusts the contour of the cushion to the individual sitting on it
- Ultra-lightweight construction
- Removable incontinent cover that can be machine-washed when the need arises
- Maximum user weight 295 kg



Varilite Stratus Wheelchair Cushion

VL56640 (16X16 INCONTINENT)

VL56840 (16X18 INCONTINENT)

VL58640 (18X16 INCONTINENT)

VL58840 (18X18 INCONTINENT)

- For moderate to low risk of tissue breakdown and moderate positioning needs
- Unique self-inflating valve, which allows air in and out of the cushion so it can contour to the user's shape
- Two types of foam create the functional areas of the cushion: medium foam for the thigh troughs and pelvic bucketing; and firm foam for the perimeter and the medial thigh separator
- User's weight is distributed evenly over the surface of the cushion
- Maximum user weight 320 kg



Varilite Evolution Wheelchair Cushion

VL74620 (14" x 16")

VL74820 (14" x 18")

VL76620 (16" x 16")

VL76820 (16" x 18")

VL76020 (16" x 20")

VL78620 (18" x 16")

VL78820 (18" x 18")

VL78020 (18" x 20")

- Constructed for users at high risk of tissue breakdown and who have moderate positioning needs
- For motor and neurological dysfunction due to spinal cord injury, stroke, multiple sclerosis, cerebral palsy, and traumatic brain injury are the most frequent, though not the only, users of the Evolution cushion
- We recommend purchasing a spare cover with all cushion purchases



Ki Mobility Axiom G Cushion

KIXG1414 (14" x 14")

KIXG1616 (16" x 16")

KIXG1618 (16" x 18")

KIXG1816 (18" x 16")

KIXG1818 (18" x 18")

KIXG2018 (20" x 18")

- Mild contoured well supports the trochanters, transitioning the femurs and enveloping the ischial tuberosities and coccyx for support and stability
- Pre-Contoured Shape and Enhanced Positioning
- Surface Tension Reduction
- High-Quality Material Feature Rich Cover



AML Wheelchair Seat Cushion

AMLWSC1616 (16" x 16")

AMLWSC1816 (18" x 16")

AMLWSC1818 (18" x 18")

AMLWSC2018 (20" x 18")

- Ergonomically designed PU foam contoured cushion for wheelchair seats



AML Wheelchair Bariatric Anti-tippers

AMLWCATB


- Replacement pair of Anti-Tippers for the AML Bariatric wheelchairs
- Sold as a pair



Karma 8000 Series Anti-tippers

KMAT

- Replacement pair of Anti-Tippers for the Karma 8000 series wheelchairs
- Sold as a pair



AML Wheelchair Elevating Legrest

AMLWCELRL (LHS)

AMLWCELRR (RHS)

- Padded Elevating Legrest for the AML wheelchairs
- Ideal for short term use while leg is in plaster or needs elevation
- Sold individually (left or right)



Karma Elevating Legrest

8000 Series

KMELL (LHS)

KMELR (RHS)

7520 Paediatric series

KMEL7520L (LHS)

KMEL7520R (RHS)

- Ideal for short term use while leg is in plaster or needs elevation
- The elevating legrest range is from 60° to 0° (8000 Series)
- The elevating legrest range is from 80° to 0° (7520 Series)



Karma Residual Limb Support

KMSSL (LHS)

KMSSR (RHS)

- Residual limb support for Karma wheelchairs



Karma I.V. Pole Holder

KMIVPOLE

- Ideal for users who frequently visit the hospital for injection treatments
- Step-less length adjustments
- Note: the wheelchair must be equipped with anti-tippers when using this holder for the user's safety



Ottobock Modular Hand Pad - Horn Type

OB436P3=R2-7 (Medium RHS)

OB436P3=R1-7 (Large RHS)

OB436P3=L2-7 (Medium LHS)

OB436P3=L1-7 (Large LHS)

- Can be mounted in connection with the Ottobock modular Channel Forearm Pads
- For people with hand contractures as it can assist in preventing the hand from completely clenching



Ottobock Modular Hand Pad - Cone Type

OB436P4=R2-7 (Medium RHS)

OB436P4=R1-7 (Large RHS)

OB436P4=L2-7 (Medium LHS)

OB436P4=L1-7 (Large LHS)

- Can be mounted in connection with the Ottobock modular Channel Forearm Pads
- For people with hand contractures as it can assist in preventing the hand from completely clenching
- Constructed of polyurethane foam



Ottobock Modular Hand Pad - Extensor Concave

OB436P2=R2-7 (Medium RHS)

OB436P2=R1-7 (Large RHS)

OB436P2=L2-7 (Medium LHS)

OB436P2=L1-7 (Large LHS)

- Ideal for keeping the hand in a relatively neutral position
- Support surface has a built-up palmar arch for mild to moderate positioning needs
- Constructed of polyurethane foam



Ottobock Modular Hand Pad - Flat

OB436P5=2-7 (Medium)

OB436P5=1-7 (Large)

- If no special type of hand pad is required, the flat pad may be a possible solution
- These can also be used for the basis of a custom fabrication
- Constructed of polyurethane foam



Gel Oventions Arm Pad

GO310G (Size 3.5"x10")

GO316G (Size 3.5"x16")

GO312G (Size 3.5"x12")

GO212G (Size 2"x12")

GO314G (Size 3.5"x14")

GO214G (Size 2"x14")

- Easily bolts onto existing armrest T sections
- They consist of a solid wood base gel pad with a neoprene nylon cover and plastic base guard
- Synergel Gel does not absorb liquids and remains consistently soft regardless of ambient temperature.



Ottobock Channel Forearm Pads

OB436A2=3-7 (Small)

OB436A2=2-7 (Medium)

OB436A2=1-7 (Large)

- Arm gutter to help position arm on wheelchair armrest
- Provides functional support for upper extremity weakness or paralysis
- Mounted directly onto the armrest hardware
- Fits either left or right side



Ottobock Channel Forearm Pad

OB436A3=2-7

- One piece channel forearm pad
- Provide functional support if you have upper extremity weakness or paralysis
- It can be mounted directly onto the armrest hardware and can also be used with a Swivel Unit or Offset Plate to get the correct position

POWERCHAIRS



Karma KP25.2 Powerchair

KP2521616L (16" x 16")

KP2521818L (18" x 18")

KP2522018L (20" x 18")

- The new Karma KP25.2 is the updated version of the very durable KP25
- It is a cost effective rear-wheel drive powerchair that provides the user with independence and mobility
- The frame is lightweight, compact and easy to assemble. It can also be folded down to fit into a car boot
- The KP25.2 features a VR2 controller



Jazzy 623 2.0

PRIDEJ6232PB

- Compact full-size mid-wheel drive powerchair with six wheels on the ground for maximum stability
- VR2 controller allowing up to two powered functions from Power tilt, Power elevate and Power elevating leg rests, swing-away or centre mounted
- High efficiency and mileage with Linux 2 pole motor
- Compatible with TRU-Balance Power Positioning
- Quick Height Adjustable (QHA) armrests
- Lower seat to floor heights (down to 17.5" with tilt)
- Maximum user weight 136 kg
- Attendant control



AML Corner Shower Stool

AMLCSS

- Anti-slip suction cups secure the stool to the floor
- Three height adjustable legs
- Maximum user weight 136 kg



Height Adjustable Shower Stool with Arms

AML20107


- Height adjustable shower stool with injection moulded plastic seat
- Powder coated frame
- Armrests assist safe transfer



Height Adjustable Shower Chair With Arms

AML20108

- Height adjustable shower chair with injection moulded plastic seat
- Features a curved plastic seat for comfort and drainage holes for excess water
- Armrests assist safe transfer



Bath Board

MP20301 (Medium)

AML20303 (Large)

- Designed to provide the user with a comfortable and sturdy seat for independent bathing
- One piece injection moulded bath board
- Has four adjustment knobs which enable the bath board to be fitted into the curvature of the bath
- Maximum user weight 250 kg



Perching Stool with Arms and Backrest

AML30305

- Comfortable padded seat
- Curved and removable armrests which help in assisting with safe transfers
- Rubber ferrules provide a sturdy grip on the floor
- Height adjustable legs
- Maximum user weight 113 kg





4" Raised Toilet Seat With Armrests

MP20402

- A moulded plastic 4" Raised Toilet Seat with padded arms that allow for security and independence when getting on and off the toilet
- Ideal for use after hip or knee replacement surgery
- It easily fits to the toilet bowl via adjustable brackets
- Maximum user weight 190 kg



4" Raised Toilet Seat Without Armrests

MP20404

- A one-piece injection moulded 4" high raised toilet seat that fits easily to the toilet bowl with three adjustable brackets
- Maximum user weight 190 kg



6" Raised Toilet Seat

DM12066

- Designed to be used on an existing toilet, it easily attaches to the toilet bowl and locks into place with two rear locks
- There are two hygiene cut-outs for ease of use and with no lid attached to the seat, it is very easy to clean
- No tools are required for installation
- Maximum user weight 200 kg



Bedside Commode

MP20502

- Padded seat & back, moulded armrests
- Powder coated blue steel frame
- Height adjustable legs
- Maximum user weight 125 kg



3 in 1 Commode

AML20505

- Multifunctional - can be used as bedside commode, toilet safety frame and raised toilet seat.
- Removable backrest
- Comfortable, stable and lightweight
- Powder coated steel frame
- Comes complete with bucket, lid and splash guard



Belgrave Bariatric Commode

DM190

- A robust height adjustable commode with a non-peel frame and a cellular moulded backrest and over seat that can be removed when needed
- Bucket holds 14 litres and comes with a lid
- Arms can also be removed for easy side lateral transfer
- Maximum user weight 280 kg



Toilet Safety Frame

MP20703

- Powder Coated Steel frame with plastic arm rests that fits around all three sides of a toilet
- Ideal for use during recovery from hip or knee surgery to ensure extra stability is provided
- The height adjustable legs have rubber tips that ensure it does not slide when being used
- Maximum user weight 130 kg



Good Grips Twisting Jar Opener

GG11332200

- Opens jars with minimal effort and reduces hand strain
- Large handle is easy to grip; perfect for arthritic hands
- Base pad creates a non-slip surface for a better grip
- Sharp teeth grip jars for easy opening
- Twisting design adjusts to all-sized lids



AML Kettle Pourer

SLHA4270

- The frame helps prevent spillage
- Made of wire frame with a cradle
- Holds most kettles and pivots on a raised wire frame to allow gentle pressure to be applied to the handle to pour into the cup
- The kettle sits in the cradle and is secured by a strap



AML Bed & Chair Raiser

SLHA4507

- Moulded plastic cone raisers to make it easier for users to sit/lie down and get up
- Can raise a standard bed or chair up to 150mm
- Steady when in use
- Set of four
- Black only
- Socket is sloped



Over Bed/Chair Table

MP30203

- Adjustable tilting table to be used by someone in bed or sitting in a chair
- Walnut grained melamine laminate for easy care
- Has a solid rim to prevent things from sliding off
- Fitted with four 50 mm smooth glide castors



Bariatric Perching Stool with Arms and Backrest

PHPAT0911156678

- Height-adjustable legs with pin clips
- Non-slip rubber feet for safety
- Depth-adjustable backrest for cushioned support
- Tubular armrests support upright posture
- Maximum user weight 255 kg



Perching Stool

MP30304

- Designed to allow the user to rest and then stand again with ease, which works well for those who tire easily or recovering from surgery or an injury
- Rigid steel, powder coated frame that is able to be used in the kitchen or bathroom
- Maximum user weight 125 kg



AML Bed Cradle

AMLBCHA

- White powder coated steel frame bed cradle helping to eliminate the weight of blankets on feet and ankles
- 5 height adjustments
- Lower frame slides underneath
- Light weight



KCare Patient Assist Rail

KC9077

- Rail can be used with any sized bed, simply slide it under the mattress for easy application
- 45 mm gap between the shield and the bar, allowing for a safe and comfortable grip
- Enclosed return to minimise injury risk
- Inner-shield preventing entrapment risk
- Powder-coated tubular steel frame material



Parsons Comfort Grip Cutlery Individual Cutlery

PA16T001RD (Knife)

PA16T002RD (Fork)

PA16T003RD (Teaspoon)

PA16T004RD (Soup spoon)

- Get maximum user control with minimum effort as the soft vinyl coating has finger indentations for a secure grip
- Also available as full set



Parsons Comfort Grip Weighted Individual Cutlery

PA16T400 (Knife)

PA16T401 (Fork)

PA16T403 (Soup spoon)

- Weighted Cutlery for persons with Parkinson's or other impairments resulting in hand tremors
- Latex-free
- Home dishwasher safe
- Straight handled universal option only
- Cutlery each sold individually
- Right and left-handed versions are available for purchase by indent order only



AML Utensil Strap

SLHA4196

- Ideal for people with a limited grip
- Simply attach the strap to the utensil and slip a hand through the loop
- It fits almost any hand. This item does not include the spoon
- This utensil strap is best suited for the AML utensils with large grips and looped ends



AML Large Scoop Bowl

SLHA4248

- The large mould side to this scoop bowl accommodates individuals who have difficulty when eating particularly with one hand
- It has a non-slip suction cup to keep the bowl in place
- Microwave and household dishwasher safe up to 70° C except for the suction base
- White only



AML Scoop Dish

SLHA4249

- The high back on the scoop dish accommodates individuals who have difficulty when eating particularly with one hand
- It has a non-slip suction cup to keep the dish in place
- Microwave and household dishwasher safe up to 70° C except for the suction base
- White only



AML Plastic Plate Guard

SLHA4250

- An easy to fit and remove plastic plate guard that fits tightly around the rim of any circular plate between 200-300 mm in diameter
- The edge is curved to prevent food from falling off the plate, particularly if eating with one hand
- Microwave and household dishwasher safe up to 70° C
- White only



Parsons Inner Lip Plate

PA16T118 (Sandstone)

PA16T118BL (Blue)

- Holds food on the plate while the user brings the fork or spoon to the edge of the plate making it ideally suited for people with use of only one hand.
- Break proof
- Microwavable
- Available in blue or sandstone



Parsons Double Handled Mug With Lid

PA16T123 (White)

PA16T123RD (Red)


- The handles will accommodate adult fingers
- The transparent design instantly shows the level of contents
- The wide base prevents tipping
- Dishwasher safe up to 90° Celsius
- Microwaveable
- Available with white or red lid



AML Transparent Drinking Cup

SLHA4262

- Transparent cup with a drinking lid to help prevent spillage
- Two lids of varying sizes (large hole and small hole), allowing the user to choose the lid that makes drinking easier
- Microwave and household dishwasher safe up to 70°C



AML Dressing Stick

SLDA5201

- A lightweight device for positioning and locating clothing featuring a smooth wooden shaft hook with a large plastic coated "C" hook at one end and a specially shaped "push/pull" hook at the other
- Length is 685 mm



Sock and Stocking Aid

MP30805

- Lightweight, terry covered plastic cover sock aid makes it easy for socks to be put on, even if you can bend down
- The terry on the outside helps hold the sock in place and the slits in the plastic allow for it to bend around the heel
- The pulling strap makes it possible to use with just one hand
- Fully washable



Parsons Standard Stocking Hose Aid

PA16D005

- A flexible sock aid with nylon webbing straps
- The curved shape holds the stocking on the aid while pulling it up your leg
- Easy to use



AML Pantyhose Aid

SLDA5305

- Lightweight and easy to put on
- Flexible plastic trough can be bent as necessary to fit any size stockings
- Slide feet through and pull the two straps up until the pantyhose can be pulled up within reach




Aluminium Reacher

MP30801 (660 mm)

MP30802 (815 mm)

- A great gadget to help pick up objects from the floor or out of reach
- Lightweight anodized aluminium body
- Plastic handle
- Magnet in jaw for picking up light metal objects



Nipper Reacher

MP30803 (600 mm)

MP30804 (750 mm)

- Multi-Use Reacher with a lightweight anodized aluminium frame and plastic handle
- Magnet in jaw for picking up light metal objects and hook on upper jaw retrieves articles of clothing
- Slip resistant surface on jaws for secure hold
- Jaws rotate to avoid excessive wrist movement



Adjustable Back Rest

MP30101

- A supportive backrest designed to make sitting in bed easier
- Made of a white powder coated frame with a breathable nylon mesh fabric
- Comes complete with a padded, removable head rest
- Has five adjustments of angle which allow it to go from a low lying position to an upright sitting position
- For storage, it folds flat



AML Coccyx Orthopaedic Cushion

AMLCOC


- Coccyx orthopaedic cushion relieves tailbone pressure and encourages a correct sitting position
- Ergonomically designed from PU foam
- Unique cut section for coccyx
- Removable 100% polyester cover
- Easy maintenance



AML Tablet Crusher

SLHA4132

- A grinding mill for medicines.
- Turn the screw top to crush the tablets into powder
- Spare tablets can be stored in a compartment in the lid
- Suitable for carrying in a pocket and easy to clean
- Dimensions are 55 x 44 mm



Long Plastic Shoe Horn

PA16D019

- Strong flexible and latex free, this plastic shoe horn comes with a handy hang up hole
- Length is 18"



Rollator - Petite 6" Wheel

MP10103G

- A smaller rollator suitable for people 5'2" or shorter
- The cable loop brakes are easy to use
- They also have a locking feature which allows the brakes to be "locked on", necessary when using the wheeled walker's built-in seat
- A small bag under the seat allows for personal items to be safely and discreetly carried
- Maximum User Weight 136 kg



Rollator - Bariatric 8" Wheel

MP10200S

- A higher weight capacity than a regular rollator, it offers stability and mobility for the larger adult
- The rustproof, robust aluminium frame is extremely durable
- The heavy duty castors can be used both indoors and outdoors
- It has a soft padded seat and a removable backrest
- Folds up for ease of transport
- Some assembly required
- Maximum user weight 226 kg



AML Gutter Frame Rollator

MP61040070

- Forearm Rollator, which is great for people with weaker arms and shoulders
- Black and grey steel frame
- Folding walker designed to support the user in an upright position with padded armrests for comfort.
- Features a narrow seat and cable brakes on the rear wheels
- Gutter-style arm supports are height adjustable.
- Tray that can be slotted onto the seat for carrying small items
- Removable basket



AML Four Wheeled Rollator - 6\"/>

AMLMHROLBLU6

- The smaller wheels make it more ideal for indoor use, though it can be used outdoors over relatively smooth surfaces such as footpaths
- Rollators' folding mechanism allows for easy storage and transportation
- Comes in blue
- High strength aluminium tubing
- Ergonomic Hand Grips
- Padded seat
- Under-seat storage
- Locking looped hand brakes
- Folds up for ease of transport
- Detachable backrest
- Some assembly required
- Recommended maximum user weight 136 kg



AML Four Wheeled Rollator - 7.3\"/>

AMLMHROLBLU7

- Larger sized castors allow the rollator to be used both indoors and outdoors
- Rollators' folding mechanism allows for easy storage and transportation
- Comes in blue
- High strength aluminium tubing
- Ergonomic Hand Grips
- Padded seat
- Under-seat storage
- Locking looped hand brakes
- Folds up for ease of transport
- Detachable backrest
- Some assembly required
- Recommended maximum user weight 136 kg



Freestyle 3 Wheel Rollator

NRSP64077


- Stylish, contemporary design and purple colour
- Has a clever self-locking mechanism that makes it simple to store and transport the rollator when not in use
- Easy clip on/clip off bag with discreet pockets and zip fastener
- Tried and tested by Occupational Therapists
- Sturdy 190 mm (7.5") wheels enable easy maneuverability
- Maximum User Weight 115 kg



Crook Handle Cane

MP10601K

- Made of aluminium, this cane is lightweight and easy to use
- With a padded foam grip, it is comfortable to hold and its design makes it easy to hook onto a table or chair when not in use
- Black only
- Maximum user weight 125 kg



Freestyle Grab and Go Walking Stick

NRSS08847 (32")

NRSS08859 (34")

NRSS08860 (36")

- Two-in-one walking stick and grabber
- Comfortable, ergonomic rubber grip handle
- Makes reaching small items easier
- Durable walking aid
- Attractive 2-tone colour design
- Maximum user weight 135 kg



Fischer Handle Cane

MP10605K (Black LHS)

MP10606B (Bronze RHS)

MP10606K (Black RHS)

- Lightweight cane designed with arthritic grip
- The shaft is made of lightweight aluminium and it is easily height adjustable



T-Handled Cane

MP10603K (Black)

- Lightweight aluminium, height adjustable straight cane which features an ergonomic T shaped plastic handle
- Maximum user weight 125 kg



Floral Foldable Cane

MP10604V (Violet)

MP10604P (Peony)

MP10604L (Lilac)

- Handy foldable cane ready for when you need it
- The aluminium shaft allows this cane to be lightweight, durable and adjustable
- Comes complete with its own pouch for easy storage



Foldable Cane

MP10604K (Black)

- Handy foldable T-handle cane ready for when you need it
- The aluminium shaft allows this cane to be lightweight, durable and adjustable
- Maximum user weight 125 kg




Quad Cane

MP10501K (Small Base)

MP10502K (Large Base)

- This cane has been adapted for increased stability; through providing a rectangular base for support
- Made of lightweight aluminium, this cane has an offset foam coated handle
- Easily height adjustable



Swan Handled Cane

MP10602B

- Made of lightweight aluminium materials
- Easily height adjustable
- Cane has a soft foam handle, as well as a convenient carrying strap
- Maximum user weight 125 kg



Walking Stick/Cane Wrist Strap

PA14A022

- Attaches to your walking stick or cane with velcro
- Reduces the risk of dropping your cane
- Can be hung or latched around a hook when not in use



Cane Tip - Black

Quad Cane

MP60912 (Small)

MP60913 (Large)

MP60492

- Small tip is 13 mm diameter for use with MP10501K quad cane
- Large tip is 16 mm diameter for use with MP10502K quad cane
- Rubber cane tip
- 19 mm diameter
- Ideal for use with MP10601, MP10602, MP10603 & MP10604 canes



AML Deluxe Rubber Cane Tip - Black

AMLDRCTBLK

- Can tip for walking stick to give added stability and allows the cane to stand alone when needed
- Rubber cane tip
- Offers the user more stability than a standard single cane tip, especially in uneven terrain
- Approx 18mm diameter



Cane Holders

MP10614

- This clip enables cane to suspend from table edge or push clip up cane shaft until it meets underside of table
- Makes it easier than having to reach down and pick the walking stick up from the floor
- Black only



Children's Elbow Crutches

SMBE3640 (4 - 7 Years)

SMBE1070 (6 - 10 Years)

- Strong but lightweight, these aluminium crutches have double height adjustment and are finished in a cheerful buttercup yellow
- They are completely latex free - including ferrule
- Sold as a pair
- Maximum user weight 100 kg



Forearm Crutches

MP10401 (Short/Medium)

MP10402 (Medium/Large)

- Made of lightweight aluminium with a moulded ergonomically contoured arm cuff for comfort
- Double height adjustment with an extra-large rubber tip that provides added stability
- Sold as a pair
- Maximum user weight 115 kg



Forearm Crutches – Comfy Handle

MP10403

- Designed for long term users
- Double height adjustment crutch made from anodised aluminium
- Features a deluxe contoured "Ortho-Grip" Comfy handle that ensures the weight is spread through the palm
- Sold as a pair
- Maximum user weight 115 kg



Forearm Crutches – Heavy Duty

MP10404

- A sturdy elbow crutch for a heavier adult
- Made of lightweight aluminium with a moulded ergonomically contoured arm cuff for comfort
- Double height adjustment with an extra-large rubber tip that provides added stability
- Sold as a pair
- Maximum user weight 160 kg



Elbow Crutches - Adult Tall

SMBE3099

- Sturdy crutches for taller people
- Made of lightweight aluminium with a vinyl coated ergonomically contoured arm cuff for comfort
- Double height adjustment with an extra-large rubber tip that provides added stability
- Sold as a pair
- Maximum user weight 180 kg



Crutch Tip - Round

MP60490

- Grey round shaped rubber crutch tip
- 22 mm diameter
- For use with MP10401 and MP10402





CTS Wrist Brace

AML87162 (LHS XSmall)	AML87152 (RHS XSmall)
AML87163 (LHS Small)	AML87153 (RHS Small)
AML87165 (LHS Medium)	AML87155 (RHS Medium)
AML87167 (LHS Large)	AML87157 (RHS Large)
AML87168 (LHS XLarge)	AML87158 (RHS XLarge)

- Our top selling wrist brace - it's available from X-Small to X-Large in either Left or Right side
- Cotton elastic with circumferential strap and Velcro closure
- Preformed aluminium stay
- An ideal support for Carpal Tunnel Syndrome, as well as post fracture and wrist sprains and strains



Allard S.O.T. Resting Splint

AL287101008 (LHS XXXSmall)	AL287102008 (RHS XXXSmall)
AL287101009 (LHS XXSmall)	AL287102009 (RHS XXSmall)
AL287101010 (LHS XSmall)	AL287102010 (RHS XSmall)
AL287101011 (LHS Small)	AL287102011 (RHS Small)
AL287101012 (LHS Medium)	AL287102012 (RHS Medium)
AL287101013 (LHS Large)	AL287102013 (RHS Large)

- Smooth and lightweight
- Allows adjustment to the desired position
- Aluminium core is embedded into polyethylene foam, covered with fabric
- Offers a resting position which can deliver pain relief, relaxation to the hand and provides a good biomechanical position that may reduce the risk of flexor shortening of the wrist and fingers



Allard Selection Thumb Rigid Support

AL352021011 (LHS Small)	AL352022011 (RHS Small)
AL352021012 (LHS Medium)	AL352022012 (RHS Medium)
AL352021013 (LHS Large)	AL352022013 (RHS Large)

- Keeps thumb in a relaxed and pain free position
- Mild to moderate arthritis in the first CMC-joint and protects the CMC-joint and the first MCP-joint
- Incorporates two elastic materials laminated together with a core of a more rigid material to offer the required stability



Allard Selection Wrist Support For Children

AL353051010 (LHS XSmall)	AL353052010 (RHS XSmall)
AL353051011 (LHS Small)	AL353052011 (RHS Small)
AL353051012 (LHS Medium)	AL353052012 (RHS Medium)
AL353051013 (LHS Large)	AL353052013 (RHS Large)

- Low cut shape in the palm area to allow for the bending zone
- Well cut around the thumb muscles to allow full movement of the thumb
- It comes with a simple aluminium stay
- Helps guide the hand and wrist into a functional position
- Designed to reduce spasticity in hand and wrist in neurological conditions



Allard Selection Thumb Support For Children

AL353001010 (LHS XSmall)	AL353002010 (RHS XSmall)
AL353001011 (LHS Small)	AL353002011 (RHS Small)
AL353001012 (LHS Medium)	AL353002012 (RHS Medium)
AL353001013 (LHS Large)	AL353002013 (RHS Large)

- Elastic all the way around
- Designed to keep thumb in a relaxed and pain free position
- Designed to guide the thumb and hand into a more functional position to facilitate the gripping ability



Ossur Wrist Brace Exoform Carpal Tunnel

OS517083 (LHS Small)	OS517073 (RHS Small)
OS517085 (LHS Medium)	OS517075 (RHS Medium)
OS517087- (LHS Large)	OS517077 (RHS Large)

- Anatomically correct fit
- Pre-shaped into a neutral position for comfortable limited motion/flexibility
- Soft touch, contoured interior ensures a comfortable fit
- Excellent for more active users - especially those involved in more extreme sports
- Ideal for carpal tunnel syndrome, tendonitis and soft tissue injuries



Functional/Night Resting Splint

AML35360 (LHS Small)

AML35361 (LHS Medium)

AML35362 (LHS Large)

AML35365 (RHS Small)

AML35366 (RHS Medium)

AML35367 (RHS Large)

- Forearm brace, ideal for use post surgical, CVA, TBI and carpal tunnel
- Aluminium stays throughout the foam, the thumb position adjustable
- Brace can be bent as required
- Thick padding over wrist ensures comfort



Allard Lacer Wrist Brace

AL514271011 (LHS Small)

AL514271012 (LHS Medium)

AL514271013 (LHS Large)

AL514271014 (LHS X-Large)

AL514272011 (RHS Small)

AL514272012 (RHS Medium)

AL514272013 (RHS Large)

AL514272014 (RHS X-Large)

- Easy-to-adjust lacing system for a secure, customised fit
- Breathable Material: High-quality fabric for comfort and durability, keeping the skin dry
- Aluminum Stay Support: Sturdy support promoting proper wrist alignment and motion control
- Ideal for recovery, pain management, and support during various activities
- Simple application and removal process for convenience
- Long design offers added stability to the wrist and lower forearm.



Scott Universal Wrist Brace

SC1385

- Simple design, with cotton sleeve that allows additional comfort and easy one handed application
- Designed to fit left or right wrist
- Cotton-elastic construction with heavy-duty vinyl strip protecting wrist
- 50 mm wide elastic wrist strap for greater wrist stability, along with a removable aluminium stay
- One size fits most



Scott Unifoam Wrist Support

SC4005LR (LHS Regular)
SC4005LXL (LHS X-Large)

SC4005RR (RHS Regular)
SC4005RXL (RHS X-Large)

- Extremely comfortable, soft foam wrist brace
- Ideal for night resting and for carpal tunnel syndrome
- Removable aluminium stay and extra padding around key joints



Scott Elastic Wrist Brace With Palm Stay

SC4039LS (LHS Sml)
SC4039LM (LHS Med)
SC4039LL (LHS Lge)
SC4039LXL (LHS X-Lge)

SC4039RS (RHS Sml)
SC4039RM (RHS Med)
SC4039RL (RHS Lge)
SC4039RXL (RHS X-Lge)

- Immobilises the wrist and alleviates pain
- Ideal for carpal tunnel syndrome
- Constructed of 175 mm beige duo-tension elastic for varied tension around the palm and wrist area



Foam Limb Restraint

PR90110

- Contact closure with double "D" ring cinch lock
- Foam Limb Restraint, Adult, Universal Size



Allard Selection Ulnar Deviation Splint

AL352181011 (LHS Sml)

AL352181012 (LHS Med)

AL352181013 (LHS Lge)

AL352182011 (RHS Sml)

AL352182012 (RHS Med)

AL352182013 (RHS Lge)

- Splint to improve finger alignment and enhance hand function, and functional grip
- Designed to meet the contours of the hand, the splint supports the arches of the hand whilst providing re-alignment and stabilisation
- Helps to manage pain in the hand



LMB Finger Splints

LMB501AA (X-Small Size 1)

LMB501A (Small Size 2)

LMB501B (Medium Size 3)

LMB501C (Large Size 4)

- Bendable wire in closed-cell polyethylene foam
- Wire allows for easy adjustment to meet individual splinting needs and to accommodate edema
- Large curved dorsal pad helps distribute pressure
- Measure finger from palmar crease to distal finger crease



Finger Cot - 12 Pack

PR71903(Small)

PR71905 (Medium)

PR71907 (Large)

- Lightweight
- Padded aluminium
- Easy application and removal
- Sold in 12 pack



Staxx Finger Splint

PR72241 (Size 1)

PR72242 (Size 2)

PR72243 (Size 3)

PR72244 (Size 4)

PR72245 (Size 5)

PR72246 (Size 5.5)

PR72247 (Size 6)

PR72248 (Size 7)

- Protect and support the distal joint of the finger in extension
- Permitting unrestricted movement of the proximal interphalangeal joint
- Sold individually and comes in skin tone



Procare Contoured Wrist Support

PR87007

- Perforated suede construction with soft flannel lining and cotton stockinette for patient comfort
- D-rings with loop/lock closure provide a secure fit
- Aluminium palmer stay can be contoured to fit each patient
- Designed to prevent brace from extending beyond palm crease
- Ideal for Carpal Tunnel Syndrome and immobilisation of the wrist



Procare Abducted Thumb Spica

PR87100

- Volara foam laminate construction with contact closure
- Provides support to the thumb and wrist
- Can be used for soft tissue injuries to the thumb, tendonitis, arthritis or deQuervain's Syndrome
- One size fits most



Scott Wrist Wrap

SC1307

- Made of breathable elastic with Velcro closures
- Wrap around strap has double locking system for desired tension and comfort
- 3" wide, beige in colour, one size fits most



Scott Universal Tennis Elbow Support

SC1969

- All round tennis and golfers elbow support
- 50 mm wide Royal blue plush fabric on both sides, accepts Velcro closure anywhere
- Easy to adjust with loop lock ring



AML Elbow Ranger

AML94210 (Medium)

AML94220 (Large)

- For post-operative immobilisation and protected range of motion
- Elbow Ranger is ideal when a neutral position is desired
- For non-operative application, collateral ligament strains and ruptures and ulnar nerve transposition
- May be used following cast removal to stabilise fractures and for post-surgical immobilisation
- Removable handgrip assists in the prevention of pronation/supination
- Aluminium R.O.M. bilateral hinges allow for adjustments in 30 degree increments with support by the formed Kydex shells



AML Shoulder Abduction Sling

AML31201 (Small)

AML31202 (Medium)

AML31203 (Large)

AML31204 (X-Large)

- Fits left or right
- Allows shoulder and arm to remain in a neutral position
- Pillow allows auxiliary air exchange to reduce the risk of secondary infections
- Breathable Fabric and Padding Sling for comfort



AML Clinic Shoulder Immobiliser

AMLCSIXSML (X-Small)

AMLCSISML (Small)

AMLCSIMED (Medium)

AMLCSILGE (Large)

AMLCSIXLG (X-Large)

- Simple to use, multipurpose shoulder immobiliser
- Sized canvas sling and soft foam waist strap immobilise the shoulder while supporting the elbow, wrist and hand



Ottobock Omo Neurexa Plus

OB5065N=R-XS-7 (RHS X-Small) OB5065N=L-XS-7 (LHS X-Small)
 OB5065N=R-S-7 (RHS Small) OB5065N=L-S-7 (LHS Small)
 OB5065N=R-M-7 (RHS Medium) OB5065N=L-M-7 (LHS Medium)
 OB5065N=R-L-7 (RHS Large) OB5065N=L-L-7 (LHS Large)
 OB5065N=R-XL-7 (RHS X-Large) OB5065N=L-XL-7 (LHS X-Large)

- Will correct the alignment of the humerus head within the joint socket to stabilise and position the shoulder correctly
- The design of the Omo Neurexa Plus makes it easy to put on with one hand
- Made of comfortable, long lasting and machine washable materials
- The material is soft and will not rub against the skin
- Conditions this brace may be used for include: dislocation, stroke, hemiplegie and paralysis in the shoulder



Scott Four Way Clavicle Strap

SC1594M (Medium)

SC1594L (Large)

SC1594XL (X-Large)

- Support for recovery from a collarbone fracture
- Features a strap of orthopaedic polyester felt, covered with stockinette
- The closure in the back uses a D-ring and pressure sensitive velcro over felt pad



Procure Clavicle Splint

PR85003 (Small)

PR85007 (Large)

PR85005 (Medium)

PR85008 (X-Large)

- Adjustable prong buckles ensure correct fit and prevent slippage
- Stockinette-covered felt shoulder straps
- Ideal for clavicular fractures and postural problems
- Washable



SAE - Spinal Anterior Extension

AMLSAE42215

- Thoraco Lumbar Sacral Orthoses (TLSO)
- Comfortable, lightweight, universal sized
- Anatomically shaped hyper-extension orthosis
- Low profile and very easy to fit
- Constructed of anodized aluminium. Kydex plastic pads lined with soft foam for easy modification with a heat gun
- Hinged sternum and pubic assembly
- Ideal for compression fractures, kyphosis, osteoporosis and other conditions that require hyper-extension of the spine



AML Multifit Back Support

AMLMULTIFITB (Black)

- Basic back support with double elastic pulls on each side
- Ideal for short-term use providing a very economic solution
- Fits waist 800 - 1125 mm



Scott Universal 6" Sacral Support

SC903

- Support of the lower back and spine
- 6" wide elastic and Durafoam, which combine to provide support to the sacral area
- Adds pressure to the lower spine to help stabilise the spine and relieve back pain
- 3" wide side tension strap for added compression
- Universal size - fits hip circumference of 700 - 1200 mm



Scott Elastic Duo-Adjustable Back Support

SC3245S (Small)

SC3245L (Large)

SC3245M (Medium)

SC3245XL (X-Large)

- Elastic back support provides excellent abdominal and lumbosacral support
- 225 mm wide elastic back panel with 75 mm wide elastic duo-tension straps
- Ideal for women and long term use



Scott Scrotal Support

SCSA1503S (Small)

SCSA1503M (Medium)

SCSA1503L (Large)

- Scrotal athletic support
- The fabric is comfortable and breathable
- Measure around waist for size
- Hand wash with a mild detergent in cold water



Procure Pavlik Harness

PR90512 (XSmall)

PR90515 (Medium)

PR90513 (Small)

PR90517 (Large)

- For treatment of hip dysplasia and dislocation
- Chest band and shoulder straps made with heavy webbing with Kodel Polyester pile backing for durability and comfort
- Follow a qualified healthcare provider's instructions regarding usage and infection control




Headmaster Collar

SDHMC03-G (Small)

SDHMC02-G (Medium)

SDHMC01-G (Large)

- Improves breathing and swallowing by opening the trachea and esophagus
- Cool and comfortable brace, with low heat retention
- Provides firm support for the head
- Quick release velcro makes it easy to put on and take off



Ossur Miami J Collar Set

OSMJ-300 (Short)

OSMJR-300 (SML W/REPL PAD)

OSMJR-400 (MED W/REPL PAD)


- Offers superior immobilisation that simplifies size selection by having 6 unique front and back components
- Engineered to work optimally together for various precautions



Miami J Collar Thoracic Extension

OSMTR900 (W/REPL PADS)

- Fits all of the Miami J collars and is ideally used whenever greater stability of the cervical junctions is required and the control of the collar is needed to be extended
- Fitting is tool free and it is metal free to allow for ease of imaging
- Please note that the collar needs to be ordered separately



S. O. M. I. Brace

SOMI

- Sternal Occipital Mandibular Immobilizer, designed to provide flexion, neutral or hypertension positioning and immobilization to the neck and spine
- Very adjustable and has a removable chin plate which helps the patient with their hygiene and also eating
- It is made of vinyl-coated aluminium so it is lightweight and can also be contoured to the user



Philadelphia Adjustable Collar

OSPHP-AA

- Immobilisation and comfort
- Anterior height adjustability accommodates neck height from 2¼ to 4¼ inches
- Tear-off tabs accommodate neck circumference from 10 to 19 inches
- Restricts cervical spine flexion, extension and rotation
- Total cervical arch support maintains neutral alignment of the c-spine
- Water-resistant for bathing and aqua-therapy



Philadelphia CoolMax Neck Liners

OSPHPLCS (Small)

OSPHPLCL (Large)

OSPHPLCM (Medium)

- Sold separately



Philadelphia Tracheotomy Rigid Collar

OSPHP-TXX

(available a a broad range of sizes from Infant to Adult)

- Provides strong, open, anterior support
- Constructed of laminated foam plastic
- Recommended for single patient use



Patriot Extrication Collar

OSPHP-PP

- Either a supine or seated position
- Restricts cervical spine flexion, extension and rotation
- Adjustable bilateral Velcro® hook and loop fasteners
- Polyethylene Front and Back pieces for improved rigidity
- Large trachea opening for quick access to carotid pulse monitoring and emergency tracheotomies



AML R.O.M. Walker

AMLROMWS (Small)

AMLROMWL (Large)

AMLROMWM (Medium)

- Low profile brace with Dual Action Hinges
- Adjustable strapping to give extra support.
- Accordion moulding in the popliteal area eliminates bunching and increases comfort
- Neoprene lined with Lycra to offer a breathable contact area for comfort and extensive use



AML Lower Leg Walker

AMLLLWS (Small)

AMLLLWM (Medium)

AMLLLWL (Large)

- Low profile walker that provides stability for lower leg injuries
- Simple and cost-effective solution that provides the support you need to rehabilitate your injury
- It is made from a breathable fibre liner and a non-skid rocker sole
- Comes in three sizes and is universal to fit your left or right foot



Ossur Formfit Night Splint With Tread

OS50023 (Small)

OS50025 (Medium)

OS50027 (Large)

- 3-dimensional padding throughout the shell
- Adjustable flexion straps that control the amount of stretch provided for effective treatment and provide a gentle stretching of the plantar fascia and Achilles tendon
- Is ideal for treatment of Plantar Fasciitis, Achilles Tendonitis and other lower extremity overuse injuries



Ossur Walker Equaliser

OSW0600BLK

- Unique shock absorbing sole minimises the impact on the heel while walking, preventing further irritation of injuries
- Rocker bottom walking base on this moonboot is specially designed to be low profile and wide to promote a natural and stable gait
- It is ideal for soft tissue injuries, grade 2 and 3 sprains, or stable fractures as well as post-operative, trauma, or rehab use



Ossur Airform Night Splint

OS50045 (Medium)

OS50047 (Large)

- Designed to maintain the ankle in a neutral position (right angle, as if standing)
- Provides a gentle stretch on the Plantar Fascia and Achilles Tendon
- Lightweight, low profile slim design makes it comfortable to wear at night which means better compliance
- Effective for treating plantar fasciitis, achilles tendonitis, achilles strains, achilles bursitis, drop foot and post-static pain



Ossur Formfit Ankle Brace (Speed Brace)

OSB212000001 (X-Small)

OSB212000004 (Large)

OSB212000002 (Small)

OSB212000005 (X-Large)

OSB212000003 (Medium)

- Provides the compression and comfort of a soft ankle support with 20 - 35% more protection than traditional ankle stirrups
- Exoform technology forms a plastic exoskeleton directly onto soft breathable fabric
- Speed laces allows the user to lace the brace up in seconds, saving time and encouraging the brace's use



Allard ToeOFF 2.0

AL289221011 (LHS Small)

AL289222011 (RHS Small)

AL289221012 (LHS Medium)

AL289222012 (RHS Medium)

AL289221013 (LHS Large)

AL289222013 (RHS Large)

- Carbon composite dynamic Ankle Foot Orthosis solution designed for users with mild to severe foot drop
- Organic shape, with short side extensions and more rounded shape
- Carbon and Kevlar design absorbs energy at the heel strike and returns it at toe-off, restoring a more normal and dynamic gait



Allard Blue Rocker 2.0

AL289401012 (LHS Medium) AL289402011 (RHS Small)
 AL289401013 (LHS Large) AL289402012 (RHS Medium)
 AL289402013 (RHS Large)
 AL289402014 (RHS X-Large)

- Carbon composite dynamic Ankle Foot Orthosis solution designed for users with mild to severe foot drop
- Organic shape, with short side extensions and more rounded shape
- Carbon and Kevlar design absorbs energy at the heel strike and returns it at toe-off, restoring a more normal and dynamic gait



Allard SoftKIT 2.0

AL287500011 (Small)
 AL287500012 (Medium)
 AL287500013 (Large)

- For ToeOFF 2.0 and Blue Rocker 2.0 AFO
- Soft interface between the orthosis and the tibial crest
- Padding system that includes two narrow neoprene pads to create a pressure-free channel
- Helps relieve any potential discomfort from bony prominences



Knee Immobiliser

PR80023 (Small)
 PR80025 (Medium)
 PR80027 (Large)
 PR80028 (X-Large)

- Comfortable, sturdy, foam/mesh construction with alternating web straps and loop/lock contact closure
- All sizes are 20" / 50.8 cm in length
- Comfortable and sturdy
- Loop/lock contact closure for easy fitting
- Combination of stays for proper immobilisation
- HI-beam stays
- Foam/mesh construction with alternating web straps



Knee Ranger

PR94240

- Ideal for providing protected range of motion for postoperative immobilization
- The dual-axis, polycentric hinge allows for convenient dial-in adjustment with flexion and extension stops at 0°, 15°, 30°, 60° and 90°.
- Trimmable foam leg wraps for a customised fit
- Medial/lateral plastic shells riveted to hinge bars provide additional support following post-operative procedures
- Accommodates Ankle/Shoe Insert for added rotational control



Ottobock Genu Neurexa Knee Brace

OB8165=S-7 (Small)

OB8165=L-7 (Large)

OB8165=M-7 (Medium)

OB8165=XL-7 (X-Large)

- Comfortable to wear and should be worn directly on the skin to take full advantage of the material properties
- Prevents overextension of your knee through the secure cross-over straps at the back of the brace
- The range of motion can be individually adjusted according to the needs of the user
- Temperature-regulating material, so perspiration that may usually occur when wearing a brace is reduced, and the brace can be worn for prolonged periods of time



Ossur Wraparound Hinged Knee Support

OS212003 (Small)

OS212008 (X-Large)

OS212005 (Medium)

OS212009 (XX-Large)

OS212007 (Large)

OS212010 (XXX-Large)

- Low profile brace with Dual Action Hinges
- Adjustable strapping to give extra support
- Accordion moulding in the popliteal area eliminates bunching and increases comfort
- Neoprene lined with Lycra to offer a breathable contact area for comfort and extensive use



AML Supralite Thermoplastic Ankle Foot Orthosis (A.F.O.)

AMLSUPAFOFL (Female LHS) AMLSUPAFOML (Male LHS)
AMLSUPAFOFR (Female RHS) AMLSUPAFOMR (Male RHS)

- Made of durable, orthopaedic grade thermoplastic
- Can be custom trimmed and reshaped using cast scissors and a heat gun
- This product is recommended for single patient use



Lace Up Ankle Brace

PR81313 (Small) PR81317 (Large)
PR81315 (Medium) PR81318 (X-Large)

- Fits easily into any athletic shoe
- Sturdy canvas stirrup brace
- Designed with removable plastic medial / lateral stays
- Provide support for acute ankle ligament sprains and chronic instabilities



Scott Double Strap Ankle Support

SC325S (Small) SC325L (Large)
SC325M (Medium) SC325XL (X-Large)

- Figure of 8 Ankle Brace
- Slip-on ankle support with open heel, made of beige four-way stretch spandex
- Double strap style offers compression to both medial and lateral sides of the ankle
- Provides mild ankle protection to help prevent injury or re-injury for chronic ankle instabilities or strains
- May help control over rotation of the ankle.



Black Dictus Band Shoeless Cuff

ERDS03L (LHS)

ERDS03R (RHS)

- Simple AFO for Drop Foot that can be used without shoes in conjunction with the ankle strap from the Outdoor Dictus Band or Dictus Flex
- Please note: Only lower foot cuff included, Black Dictus Band (sold separately)



Black Dictus Band - Outdoor

ERDS02J (Junior)

ERDS02L (Large)

ERDS02M (Medium)

ERDS02XL (X-Large)

- An orthosis to help people who suffer from Drop Foot
- Drop Foot (or Foot Drop) is a neuromuscular (nerve and muscle) disorder that affects the patient's ability to raise their foot at the ankle
- Included parts: 1 x Ankle band, 2 x Latex Rubber Bands,
- 1 x Monohook, 2 x Steel hooks



Replacement Mono Hook for Dictus Band

ERDS0102

- 'Hook' system that can be used with the outdoor Dictus Band
- It slides into the tongue area of the shoe underneath the laces, with the hook coming through the laces
- Latex Rubber Band then hooks into the hook of the Mono Hook and also the hook on the Dictus Band to assist with drop foot
- Sold individually



Replacement Spring Steel Hooks for Dictus Band - Pair

ERDS0103

- For use with an outdoor Dictus Band to hook through shoe laces and attach Latex Rubber Bands to
- Sold as a pair



Replacement Latex Rubber Band

ERDS01041

- For use with Dictus Bands - either indoor or outdoor
- Fantastic to have on hand for when your current bands wear out over time
- Sold individually



Dictus Hip-Foot

ERD530500-311 (Medium)

ERD530500-321 (Large)

- Designed to assist individuals with weak hip flexors, such as those with Multiple Sclerosis (MS)
- It provides support and enhances mobility, adapting quickly and easily to the user's needs
- Elastic bands assist with hip flexion
- Adjustable level and pull direction for tailored support
- Can correct foot inversion
- Suitable to wear over or under clothes
- Can be used on one or both legs



Dictus Hip-Knee

ERD530500-311 (Medium)

ERD530500-321 (Large)

- Supports individuals with leg weakness, specifically targeting the hip and knee areas
- Offers a comfortable and easy-to-use solution to enhance leg movement
- Provides assistance for hip flexion
- Easily adjustable for customised support
- Simple to change flexion, adduction, and inward rotation pull direction
- Comfortable for wear over or under garments
- Suitable for unilateral or bilateral use



Tulis Classic Gel Heel Cups

TULIS10215 (Regular)

TULIS10216 (Large)

- Tulis Gel is a specially developed polymer that is lighter, softer and more resilient than the material in similar heel cups
- The patented waffle design provide the ultimate shock absorption for use in shoes of all types
- Sold as a pair



Tulis Heavy-Duty Gel Heel Cups

TULIS10219 (Small)

TULIS10217 (Regular)

TULIS10218 (Large)

- Tulis Heavy-Duty Gel heel cups provide the ultimate in shock absorption for high-impact activities
- Sold as a pair



Pedal Exerciser

MP30401

- Created especially for people with limited mobility, who cannot use conventional exercise equipment
- The pedal exerciser helps to strengthen your legs and/or arms and hands while increasing your blood circulation
- Chrome plated with full tension range



AML Folding Pedal Exerciser with Pedometer

AMLFPE

- Full tension range
- Can also exercise upper body if required
- Wrap around foot straps
- Folds down making it easier to slide out of the way when not in use, no tools required
- Recommended maximum user weight 100 kg

Gemstim TENs Machine

GEMSTIM

- A cost effective TENs unit ideal for single patient use
- The unit is a dual channel electric stimulator for active treatment applications
- Comes complete with: Four pre-wired reusable self adhesive electrodes; 9-volt battery; Leads (two channel); Instruction manual; compact carry case



Self Adhesive Electrodes TENS machine


ELECTSAD

- Sold separately

Obstetric TENS Machine

GEMSTIM

- A cost effective TENs unit ideal for single patient use
- The unit is a dual channel electric stimulator for active treatment applications
- Comes complete with: Four pre-wired reusable self adhesive electrodes; 9-volt battery; Leads (two channel); Instruction manual; compact carry case



Self Adhesive Electrodes Obstetric TENS machine

OBSELECT

- Sold separately

Leads For Gemstim TENs Machine

LEADSGS

- Replacement electrode leads to suit GEMSTIM TENS machines
- Sold as a pair





Gel for TENs Machines

GEL-60ML

- Gel for use with TENS Machines that use silicone electrodes
- Can also be used with products such as foetal dopplers
- 60 ml tube





Please visit our website
www.alliedmedical.co.nz
for our full range of products



0800 31 61 81
helpis@alliedmedical.co.nz

April 2024. This information is correct at the time of printing.
From time to time, product specifications are changed by manufacturers without prior notice.



**ALLIED
MEDICAL**
LIMITED