

ACTIVE X™

- LOW PROFILE AND ULTRA LIGHT VENTILATED MATERIAL
- ANATOMICALLY SHAPED PATELLA-OPENING
- STRAPS EASY TO RE-POSITION
- HINGE WITH ADJUSTABLE FLEXION - EXTENSION
- AVAILABLE AS TUBE OR WRAP
- "STAY-UP" SILICONE ON INSIDE TO PREVENT MIGRATION

Another New
Product!

ActiveX Patella



A COMPLETELY NEW LINE OF FUNCTIONAL KNEE BRACES!

TUBE SHORT • WRAP SHORT • TUBE LONG • WRAP LONG • PATELLA STABILIZATION

ACTIVE X™

LONG TUBE & WRAP

This longer version (40 cm) is intended for different kinds of knee pain, such as mild to moderate sideways instability. The orthosis has aluminum uprights with adjustable flexion and extension stops to allow for different range of motion adjustments.

It is available in two versions, a pull-on tube and an open wrap-around model. Both models have circular straps with d-ring. Two on the thigh section and two on the calf section. These can easily be adjusted and re-positioned. The braces are delivered with the joint set on 10° of flexion. A full range of stops are included in the package.

Intended use: Mild to moderate MCL/LCL instability. Minor medial or lateral meniscus injuries. Mild osteoarthritis to the knee.



ACTIVE X™

SHORT TUBE & WRAP

This shorter version (32 cm) is intended for different kinds of knee pain, such as mild to moderate sideways instability. The orthosis has aluminum uprights with adjustable flexion and extension stops to allow for different range of motion adjustments.

It is available in two versions, a pull-on tube and an open wrap-around model. Both models have circular straps with d-ring. One on the thigh section and one on the calf section. These can easily be adjusted and re-positioned. The braces are delivered with the joint set in 10° of flexion. A full range of stops are included in the package.

Intended use: Mild to moderate MCL/LCL instability. Minor medial or lateral meniscus injuries.



Item No.	Description	Sizes	Left/Right	Total brace length
28050	ActiveX Short Tube	Small - XXX-Large	Bilateral	32 cm
28051	ActiveX Short Wrap	Small - XXX-Large	Bilateral	32 cm
28052	ActiveX Long Tube	Small - XXX-Large	Bilateral	40 cm
28053	ActiveX Long Wrap	Small - XXX-Large	Bilateral	40 cm

Sizes	Circumference 15 cm above mid-patella
Small	35 - 43 cm
Medium	43 - 50 cm
Large	50 - 57 cm
X-Large	57 - 64 cm
XX-Large	64 - 71 cm
XXX-Large	71 - 79 cm

PATELLA

This open model is designed to stabilize the patella and prevent luxation and sub-luxation. The built-in pad fits into the groove on the side of the patella and grabs hold to keep it in place.

As an extra support it is combined with a plastic half moon shaped splint that is removable.

ActiveX Patella is made from the ventilated ActiveX-material to give comfort and ventilation.

Intended use: PFPS (Patellar Femoral Pain Syndrom), Chondromalacia patella, luxated or subluxated patella.

Another New Product!



Item No.	Description	Sizes	Left/Right	Knee circumference mid-patella
28055	ActiveX Patella	XX-Small	Bilateral	20 - 25 cm
28055	ActiveX Patella	X-Small	Bilateral	25 - 30 cm
28055	ActiveX Patella	Small	Bilateral	30 - 36 cm
28055	ActiveX Patella	Medium	Bilateral	36 - 42 cm
28055	ActiveX Patella	Large	Bilateral	42 - 48 cm

ACTIVE X™

LATERAL PATELLA

A smooth and nice conforming knee brace with a pad on the lateral side for stabilizing the patella. For extra stabilization two straps cross above and below patella and secures with Velcro on the opposite side.

This brace is made from Latex free Neoprene. On the inside silicone "Stay-Up" material prevents the brace from migrating down during exercise.

Intended use: Non-specific patella pain, subluxation and luxation of patella.



Item No.	Description	Sizes	Left/Right	Total brace length
28054	ActiveX Lateral Patella	X-Small - XX-Large	Left/Right	33 cm

Sizes	Circumference 15 cm above mid-patella
X-Small	28 - 35 cm
Small	35 - 43 cm
Medium	43 - 50 cm
Large	50 - 57 cm
X-Large	57 - 64 cm
XX-Large	64 - 71 cm