

## Lacer Wrist Short

A premium orthosis for rigid support and wrist stabilization, with a conformable stay and wrist strap to aid stabilization. Available in both a long (art.no 51427) and short version, with various sizes to suit all conditions. The lacing closure allows an easy one-handed application, while the palmar pad provides additional comfort and support. Comes with both volar and dorsal stays. Provides rigid wrist support for various injuries and conditions, pre-and post-operations.

### Indications for use:

Arthritis, sprains, strains & instability. Soft tissue injuries and conditions. Post fracture management.



| Art.nr.   | Model             | Size | Left/<br>right | Length | Circumference |
|-----------|-------------------|------|----------------|--------|---------------|
| 514261009 | Lacer Wrist Short | XXS  | Left           | 18 cm  | 14-17 cm      |
| 514261010 | Lacer Wrist Short | XS   | Left           | 20 cm  | 16-19 cm      |
| 514261011 | Lacer Wrist Short | S    | Left           | 21 cm  | 17-20 cm      |
| 514261012 | Lacer Wrist Short | M    | Left           | 21 cm  | 20-23 cm      |
| 514261013 | Lacer Wrist Short | L    | Left           | 21 cm  | 22-27 cm      |
| 514261014 | Lacer Wrist Short | XL   | Left           | 21 cm  | 25-31 cm      |
| 514262009 | Lacer Wrist Short | XXS  | Right          | 18 cm  | 14-17 cm      |
| 514262010 | Lacer Wrist Short | XS   | Right          | 20 cm  | 16-19 cm      |
| 514262011 | Lacer Wrist Short | S    | Right          | 21 cm  | 17-20 cm      |
| 514262012 | Lacer Wrist Short | M    | Right          | 21 cm  | 20-23 cm      |
| 514262013 | Lacer Wrist Short | L    | Right          | 21 cm  | 22-27 cm      |
| 514262014 | Lacer Wrist Short | XL   | Right          | 21 cm  | 25-31 cm      |

