

ToeUP™ - textile foot drop orthosis

A new foot drop orthosis from Allard!

ToeUP is a textile ankle foot orthosis designed for individuals with a mild form of foot drop. The textile surface makes the orthosis comfortable and soft against the ankle while it, together with the elastic cords, provides stability and support for a lift of the forefoot during walking.

- Provides a lift of the forefoot during walking
- The support is regulated by the placement of hooks and elastic cords
- Soft and flexible, for outdoor and indoor use
- The D-ring can hold the elastic bands together for a more discreet silhouette



Support and comfort for mild foot drop

ToeUP™ features a low-profile textile ankle section that fastens above the malleoli. The ankle section is connected using an elastic cord with hooks, which are secured beneath the lacing or within the lacing holes on the shoe. This mechanism provides elevation for the forefoot during ambulation. The level of elevation can be adjusted by altering the placement of the hooks or the strap on the ankle section.

The supplementary D-ring can be used to hold the elastic bands together and allows for individualized positioning.

ToeUP offers an optional foot part component for indoor use, or with unlaced shoes such as sandals.



The D-ring fastens the straps together.



Foot part for indoor use or with unlaced shoes.



Art.no	Model	Size	L/R
28500	ToeUP - package incl. 28502-28506	XS-L	Bilateral
28501	Foot part for indoor use	S-L	L and R
28502	Extra Elastic Cord	XS-S, M-L	-
28503	Extra Hooks	One size, pair	-
28504	Extra D-Ring	One size	-
28505	Extra Under Lace Hook	One size	-
28506	Extra Ankle Section	XS-L	Bilateral



Scan or click for more information!



**ALLIED
MEDICAL**