

Bariatric Wheelchairs



	AML BARIATRIC SELF-PROPELLING	AML BARIATRIC TRANSIT	KARMA 8020	KARMA FLEXX HD	KARMA BT10
DESCRIPTION	Low cost self-propelling wheelchair	Low cost transit wheelchair	Cost effective, lightweight and durable alloy wheelchair	Upgraded version of the Flexx which includes heavy duty functions	Robust self-propelling wheelchair specifically designed for larger clients
SEAT WIDTH	50.8 - 60.9 cm (20" - 24")	55.9 cm (22")	45.7 - 55.9 cm (18" - 22")	50.8 - 60.9 cm (20" - 24")	55.9 - 66.0 cm (22" - 26")
SEAT DEPTH	45.7 cm (18")	45.7 cm (18")	40.6 - 50.8 cm (16" - 20")	45.7 - 50.8 cm (18" - 20")	45.7 - 50.8 cm (18" - 20")
SEATING	Padded PVC upholstery	Slung nylon upholstery	Slung nylon upholstery	Slung nylon upholstery	Slung nylon upholstery
TRANSFER WEIGHT	30 kg	31.6 kg	13.5 kg	16.8 kg	31.8 kg
CLIENT WEIGHT LIMIT	250 kg	230 kg	160 kg	170 kg	250 kg
FOLDING	Yes	Yes	Yes	Yes	Yes
REAR WHEEL SIZE	24"	12"	12", 20", 22" or 24"	14", 20, 22" or 24"	20" or 24"
CASTOR SIZE	20.3 cm (8" x 2")	20.3 cm (8" x 2")	15.4 or 20.3 cm (6" or 8")	15.4 or 20.3 cm (6" or 8")	17.8 cm (7")
ARMRESTS	Removable full length	Removable full length	Removable full or desk length, height adjustable	Removable full or desk length, height adjustable	Full length, fixed height, flip back
MINIMUM SEAT TO FLOOR HEIGHT	49.5 cm (19.5")	49.5 cm (19.5")	43.2 cm (17")	45.7 cm (18")	43.2 cm (17")
STANDARD LEGREST STYLE	Swing-away 70°, adjustable length	Swing-away 70°, adjustable length	Swing-away 70°, adjustable length	Swing-in/out 80°, adjustable length	Swing-away 70°, adjustable length
MANUFACTURER	AML	AML	Karma	Karma	Karma

Information correct as of September 2019. From time to time, product specifications are changed by manufacturers without prior notice



	UNDER ARM BARIATRIC CRUTCHES	BARIATRIC ROLLATOR	BARIATRIC COMMODE	HS-928 MOBILITY SCOOTER
DETAILS	Maximum user weight: 249kg Height from handle to floor: 157-177cm	Castors: 7.5" Maximum user weight: 226kg Seat size: 550 x 380mm	Maximum user weight: 125kg Adjustable height: 400 - 610mm Width between arms: 460mm	Maximum user weight: 226kg Maximum speed: 15kg



	SELECTION BACK WIDE	ABDOMINAL BINDER	OSSUR WRAPAROUND HINGED KNEE SUPPORT
DETAILS	Size XXL: 125 - 140cm Size XXXL: 140 - 160cm (Pelvic circumference)	Fits hip circumference: 750 - 1125mm	Size XL: 16 - 17" Size XXL: 17 - 18" Size XXXL: 18 - 23" (Circumference of knee at mid-patella)