

BROADBAND

STURDY, LIGHTWEIGHT RAMPS



FEATURES AT A GLANCE

- Our new lightweight broad-track series is based on a strong honeycomb design.
- The ramps are made of a newly developed, lightweight yet high-strength material.
- This manufacturing process makes our lightweight construction ramps up to 20% lighter than comparable competing models, without compromising the extremely high load capacity (up to 300 kg).
- The wide and non-slip surface ensures safe handling.
- The sturdy and exceptionally strong design makes the BROADBAND Ramps suitable for most wheelchairs, electric wheelchairs and scooters as well as pushchairs and walkers.
- The ramps are available in different models and lengths to meet different needs.

SPECTRUM

1 pcs. foldable
in length



DEFRAG

Set (2 pcs.) foldable
in length



SPECTRUM

LENGTH	ITEM NO. Set	WIDTH	WEIGHT Set	LOAD CAPACITY max	STEP HEIGHT max
70cm	BS1P700	71 cm	6,3 kg	300 kg	16 cm
90cm	BS1P900	71 cm	7,6 kg	300 kg	19 cm
120 cm	BS1P1200	71 cm	9,5 kg	300 kg	25 cm
150 cm	BS1P1500	71 cm	11,5 kg	300 kg	31 cm

DEFRAG

LENGTH	ITEM NO. Set	WIDTH	WEIGHT Set	LOAD CAPACITY max	STEP HEIGHT max	LENGTH FOLDED
180 cm	BD2P1800	73 cm	14,2 kg	300 kg	37 cm	94 cm
240 cm	BD2P2400	73 cm	21,2 kg	300 kg	50 cm	125 cm
300 cm	BD2P3000	73 cm	24,8 kg	300 kg	60 cm	155 cm