



PORTABLE MINI FOUR-WHEEL MOBILITY SCOOTER

The HS-226 offers a stylish design and user-friendly features for a comfortable driving experience and excellent performance. It features a quick-removable battery pack for easy storage and car boot transportation. With its compact size, advanced features, and extended travel range, it's your ideal companion for deeper adventures.



HS-226

Features



Easy-To-Use Control Panel



Multi-Angle Adjustable Steering



Standard Front Basket



Power Off Safety Mechanism



Tool-Free Quick Release System



LED Headlight



ND Gear Shifting Device



Lithium Battery Information

Weight W/ Battery	39.7 kgs
Weight of Battery (Pack)	4 kgs
Battery	30.15AH
Battery Charger	4A Off Board
Travel Range	40km /24.85 mi

Overall Length	1000 mm / 39.37"	
Overall Width	610 mm / 24" (between Armrests)	
Overall Height	870 mm / 34.2"	
Front Wheels	200 mm / 8"	
Rear Wheels	200 mm / 8"	
Weight W/ Batteries	39.7 kgs / 87.5 lbs	
Weight of Heaviest Piece	15 kgs / 4.97 lbs	
Max. Speed	8 kmph / 5 mph	
Weight Capacity	136 kgs / 300 lbs	
Ground Clearance	35 mm / 1.4"	
Grade Climable	6 degree	
Curb Climable	50 mm / 1.97"	
Turning Radius	1130 mm / 44.49"	
Suspension	Front & Rear	
Brake	Electro-mechanical	
Seat Type	Swivel Padded Foldable	
Seat Width	465 mm / 18.3"	
Motor Size	250W, 4400 r.p.m	
Battery Size	Li 30.15AH	
Weight of Battery Pack	4 kgs	
Travel Range	40 km / 24.85 mi	
Battery Charger	4A Off Board	
Electronics	On / Off Key Switch, Battery Level Indicator, Speed Control Knob	

※ Subject to change without notice.

Colors

Special request only.
MOQ's may apply.

Tuscan Yellow

Pistachio Green

Fjord Blue

Vapour Grey

Special request only.
MOQ's may apply.

Wisteria Pink

