



HALF COVER MOPED STYLE MOBILITY SCOOTER

Taking inspiration from the motorcycle industry, the HS-926 Half Cover Moped Style Mobility Scooter offers a remarkable riding experience with its expansive seating and ergonomic aesthetics. Its sleek and sporty design exudes a stylish charm from every perspective, making heads turn wherever your journey takes you!



HS-926

Features



LED Display



Storage Compartment and Hook




Adj. Tiller



LED Lighting Package & Daytime Running Light



15" Large Wheels



Front and Rear Suspension



Back Sliding Seat


Options



Sun Canopy



Cane and Crutch Holder



Walker Holder



USB Charger



Rear Box

Colors



Titanium Silver

Rose Red

Pearl Blue

~~Fine white~~

Roland Purple

Overall Length	1580 mm / 62.2"
Overall Width	690 mm / 26.8"
Overall Height	1380 mm / 54.3"
Front Wheels	370 mm / 15"
Rear Wheels	370 mm / 15"
Weight W/ Batteries	143 kgs / 315.3 lbs
Max. Speed	12.8 kmph / 8 mph
Weight Capacity	180 kgs / 396.8 lbs
Ground Clearance	80 mm / 3.1"
Grade Climbable	10 degree
Curb Climbable	100 mm / 3.9"
Turning Radius	2200 mm / 86.8"
Suspension	Front & Rear
Brake	Electro-mechanical & Hand Brake
Seat Type	Scooter Seat w/ Angle Adj. Armrests
Seat Width	525 mm / 20.7"
Motor Size	1000W, 5000 r.p.m.
Battery Size	(2) 12V. 75Ah (50Ah Option)
Weight of Battery	51.5 kgs / 113.5 lbs
Travel Range	60 km / 37.3 mi
Battery Charger	8A Off Board
Electronics	On / Off Key Switch, Battery Level Indicator, Speed Control Button

*Subject to change without notice.