

DOORSTEP

RAMPS MADE OF HIGH-STRENGTH ALUMINIUM



ADJUSTABLE THRESHOLD RAMPS

There are 3 different threshold ramps that are practical and very easy to position. You can overcome various thresholds, small individual steps and other obstacles. The corrugated and anodized surface ensures a reliable hold even in wet conditions and is highly resistance to scratches.

The ramps are easy and flexible to use.

FEATURES AT A GLANCE

- Innovative design for easy access and for overcoming small obstacles (indoors and outdoors)
- Ribbed and anodized surface provides additional grip, even in wet conditions
- Made of high-strength, scratch-resistant aluminium
- The mini-threshold ramps have practical carrying handles and can be folded
- Easy to handle


MICRO



MINI



ADJUSTABLE



TECHNICAL DETAILS

DOORSTEP	WIDTH (mm)	MAX HEIGHT FOR STEP (mm)	MAX LOAD (kg)	NET WEIGHT (kg)	ARTICLE NUMBER
Mini	760	30-70	350	3,3	PFD04076
Micro	2 x 250	30-70	350	2 x 1,1	PFD04225
Adjustable	760	30-70	350	1,6	PAW04276