

Quantum Edge 3 Stretto Powerchair Script



Customer Details

Client						
Assessor						
Organisation						
	ACC		DSS		Other	
User Weight						

Quantum Edge 3 Stretto Powerchair Script

For assistance with scripting please contact your **Technical Product Specialist** or call **0800 31 61 81** to be put directly through to them.

Edge 3 Stretto Powerchair Script

- Expandable QL3 Electronics
- Switch fender lights
- USB charger (plug into controller)
- Smooth Ride Suspension (SRS)
- iLevel is standard on a REAC Lift and Tilt system
- Overall width 14" drive wheels: 21.75" (55.2cm)
- Back in Black shrouds

Seat Size Width x Depth

Width		Depth	
12	17	12	19
13	18	13	20
14	19	14	21
15	20	15	22
16		16	23
		17	24
		18	



Wheel Size (Pneumatic Tyre)

3 x 14" (standard)	2 x 12.5"
--------------------	-----------

Charger

<input type="checkbox"/> Yes	<input type="checkbox"/> No
------------------------------	-----------------------------

Powered Seating

<input type="checkbox"/> Non-Tilt	<input type="checkbox"/> Lift & Tilt (iLevel is standard)
<input type="checkbox"/> Tilt	<input type="checkbox"/> Recline

Anterior Tilt with Articulating Platform* (TB4)

<input type="checkbox"/> 10 deg	
<input type="checkbox"/> 30 deg	
Anterior Tilt required:	deg

Seat to Floor Height

<input type="checkbox"/> 17.5"	<input type="checkbox"/> 19.5"
<input type="checkbox"/> 18.5"	<input type="checkbox"/> Other

Controller

Controller Type

<input type="checkbox"/> Standard
<input type="checkbox"/> Attendant Only
<input type="checkbox"/> Q-Logic 3e

Other (please specify in notes)

Controller Side

<input type="checkbox"/> Left hand side
<input type="checkbox"/> Right hand side
<input type="checkbox"/> Midline
<input type="checkbox"/> N/A

Controller Mount

<input type="checkbox"/> Fixed in-line
<input type="checkbox"/> Swing away
<input type="checkbox"/> Retract4
<input type="checkbox"/> N/A

* If Anterior Tilt with Articulating Platform is selected, the PWC will be fitted with the TRU Balance 4 powered positioning system which includes all Powered Seating Functions. Please note, ADI Backs are best suited for TRU Balance 4 seating.

* Special order item. Not available for trial.

Legrest Type	
<input type="checkbox"/>	Centre Mount
<input type="checkbox"/>	Powered Centre Mount
<input type="checkbox"/>	70 Degree Swing Away
<input type="checkbox"/>	60-90 Degree Swing Away
<input type="checkbox"/>	Manual Elevating Swing Away
<input type="checkbox"/>	Manual Elevating Centre Mount
<input type="checkbox"/>	Powered Swing Away
<input type="checkbox"/>	Paediatric

Stump Support	
<input type="checkbox"/>	Left
<input type="checkbox"/>	Right
<input type="checkbox"/>	Both
<input type="checkbox"/>	None

Footplate Type & Size	
<input type="checkbox"/>	Standard

Angle Adjustable		Centre Mount	
<input type="checkbox"/>	6 x 8	<input type="checkbox"/>	5 x 8
<input type="checkbox"/>	8 x 11	<input type="checkbox"/>	5 x 11
<input type="checkbox"/>	Paediatric	<input type="checkbox"/>	Not required

Seat Pan to Footplate Length			
Left		Right	
<input type="checkbox"/>	14	<input type="checkbox"/>	14
<input type="checkbox"/>	15	<input type="checkbox"/>	15
<input type="checkbox"/>	16	<input type="checkbox"/>	16
<input type="checkbox"/>	17	<input type="checkbox"/>	17
<input type="checkbox"/>	18	<input type="checkbox"/>	18
<input type="checkbox"/>	19	<input type="checkbox"/>	19
<input type="checkbox"/>	20	<input type="checkbox"/>	20
<input type="checkbox"/>	Not required	<input type="checkbox"/>	Not required
<input type="checkbox"/>	Other (please specify:)	<input type="checkbox"/>	Other (please specify:)
<input type="checkbox"/>		<input type="checkbox"/>	

Heel Loops	
<input type="checkbox"/>	Yes
<input type="checkbox"/>	No

Armrest Type	
<input type="checkbox"/>	Two post flip back
<input type="checkbox"/>	Flip-up armrests (not compatible with recline)
<input type="checkbox"/>	Quick height adjust drop-in
<input type="checkbox"/>	Heavy duty
<input type="checkbox"/>	Paediatric
<input type="checkbox"/>	Recline armrests (only compatible with recline)
<input type="checkbox"/>	None

Armrest Length & Type			
Length		Type	
<input type="checkbox"/>	Desk	<input type="checkbox"/>	Standard
<input type="checkbox"/>	Full	<input type="checkbox"/>	Gel
		<input type="checkbox"/>	Channel arms

Hip Belt Type & Size			
Size			
<input type="checkbox"/>	Small	<input type="checkbox"/>	Medium
<input type="checkbox"/>		<input type="checkbox"/>	Large
Type			
<input type="checkbox"/>	Standard		
<input type="checkbox"/>	2 Point		
<input type="checkbox"/>	4 Point		
<input type="checkbox"/>	Other (please specify:)		
<input type="checkbox"/>			

Headrest	
<input type="checkbox"/>	Not required
Stealth	
<input type="checkbox"/>	10" Non-flip down mount
<input type="checkbox"/>	10" Flip down mount
<input type="checkbox"/>	14" Non-flip down mount
<input type="checkbox"/>	14" Flip down mount
<input type="checkbox"/>	Other (please specify:)
<input type="checkbox"/>	

Seating - Cushion	
<input type="checkbox"/>	Not Required
<input type="checkbox"/>	Varilite Junior
<input type="checkbox"/>	Varilite Stratus
<input type="checkbox"/>	Varilite Evolution
<input type="checkbox"/>	Ride Java
<input type="checkbox"/>	Ride Forward
<input type="checkbox"/>	Ki Mobility Axiom P
<input type="checkbox"/>	Ki Mobility Axiom SP Visco
<input type="checkbox"/>	Ki Mobility Axiom SP Fluid
<input type="checkbox"/>	Ottobock Cloud
<input type="checkbox"/>	Other (please specify:)
<input type="checkbox"/>	

Seating - Back Support

<input type="checkbox"/>	Not required
<input type="checkbox"/>	Varilite Junior Mid
<input type="checkbox"/>	Varilite Junior Deep
<input type="checkbox"/>	Varilite Icon Mid
<input type="checkbox"/>	Varilite Icon Tall
<input type="checkbox"/>	Varilite Icon Deep
<input type="checkbox"/>	Ride Designs Paediatric Java Decaf Back
<input type="checkbox"/>	Ride Designs Java Back Regular
<input type="checkbox"/>	Ride Designs Java Back Tall
<input type="checkbox"/>	Ride Designs Java Back Tall Deep
<input type="checkbox"/>	Ride Designs Java Back Grande
<input type="checkbox"/>	Ride Designs Java Back Grande Deep
<input type="checkbox"/>	Axiom Back
<input type="checkbox"/>	ADI Aluminium Back Active Contour (inc height in notes)
<input type="checkbox"/>	ADI Aluminium Back Deep Contour (inc height in notes)
<input type="checkbox"/>	ADI Q-Back Back Active Contour (inc height in notes)
<input type="checkbox"/>	ADI Q-Back Back Deep Contour (inc height in notes)
<input type="checkbox"/>	ADI Flat Top Back Active Contour (inc height in notes)
<input type="checkbox"/>	ADI Flat Top Back Deep Contour (inc height in notes)
<input type="checkbox"/>	ADI Carbon Fibre Back Active Contour
<input type="checkbox"/>	Other (please specify:)

Seating - Laterals (Please Specify)

<input type="checkbox"/>	Not required	<input type="checkbox"/>	Thoracic	<input type="checkbox"/>	Pelvic
--------------------------	--------------	--------------------------	----------	--------------------------	--------

--

Tray

<input type="checkbox"/>	Swing Away	<input type="checkbox"/>	Custom
<input type="checkbox"/>	Drop In	<input type="checkbox"/>	None

Accessories

--

Set-up Instructions & Notes

--

*Green, blue and raspberry shrouds available for special orders