

Sled Impact Test

PM 2517  
Pride Mobility

**Frontal Impact of a Edge 4 (TR# 306-25)  
Secured by a Surrogate Six-Point, Strap-Type Tiedown  
and Loaded with a Hybrid III Midsize Male ATD  
Restrained by a Surrogate Vehicle-Anchored Three-Point Belt**

Tested in accordance with Annex A of  
ISO 7176-19 (2022): *Wheeled Mobility Devices for Use in Motor Vehicles*


Test Date: June 24, 2025


Submitted to:  
Pride Mobility  
401 York Avenue  
Duryea, Pennsylvania 18642

The University of Michigan  
Transportation Research Institute  
2901 Baxter Road  
Ann Arbor, Michigan 48109  
Accredited to ISO/IEC 17025:2017



Authorized Signatures

  
Nichole Orton  
Senior  
Research Engineer

  
Miriam Manary  
Lead  
Research Engineer