

Compact in design with a folding frame for easy storage and transportation



The K Care Electric Bath Lift is designed to assist users to safely transfer into and out of a bath.

## Electric Bath Lift

KA610

Using a bath lift can enhance safety and independence in the bathroom, allowing individuals to enjoy a bath whilst minimising assistance from others.

### Specifications

#### KA610

Overall Height	1071mm
Seat Height Range	60-480mm
Overall Width	360mm
Seat Width	350mm
Seat Depth	430mm
Overall Length - Folded	728mm
Overall Length - Reclined	950mm
Backrest Angle	140°
Weight Without Accessories	11kg
Max User Weight	140kg

**Important:** If you or someone you know is considering a bath lift, it's advisable to consult with a healthcare professional or occupational therapist to determine the most suitable option based on individual requirements.



Using a bath lift can enhance safety and independence in the bathroom, allowing individuals to enjoy a bath whilst minimising assistance from others.

## Features

- Easy to assemble and use.
- Compact in design with a folding frame for easy storage and transportation.
- High-quality durable seating pad and backrest.
- Easily positioned in standard-sized bathtubs (non-standard or specially customised bathtub sizes may not be suitable).
- Suction cups for secure positioning.
- The bath lift can be raised and lowered during bathing.
- Fold out side flaps assist with transfer to and from the bathlift.
- The backrest can be locked to prevent accidental tilting.
- The bath lift is waterproof (except for the charger).
- The bath lift is controlled with a hand set – using a rechargeable battery that is integrated into the handset. The handset includes a battery indicator, up/down controls and an emergency stop button.
- Australian Register of Therapeutic Goods – ARTG no 390530.
- Independently tested to International standards for Product Safety, Electrical Safety, and EMC.
- Power adaptor has Australian Electrical Compliance.

