

Electric Toilet Lift Seat

Maintain your dignity and independence with this electric toilet lift seat. A safe, sleek solution to help you sit and stand with ease.



+ A LITTLE LIFT FOR GREATER INDEPENDENCE

Designed to make sitting down and standing up from the toilet safer and easier, the Electric Toilet Lift Seat provides smooth powered assistance at the press of a button.

Ideal for use at home or in supported living environments, it combines a compact bathroom-friendly design with secure support, simple controls and up to 14" (35.56 cm) of powered lift.

+ SIT. LIFT. STAND WITH CONFIDENCE.

Designed for comfort & hygiene

The contoured seat features an integrated chute to guide waste directly into the toilet bowl, while the sleek design and minimal seams make cleaning easier.

Its lower seated position can also support a more natural toileting posture.

Power that lasts

A single full charge provides more than 120 lift cycles and can last up to 30 days, depending on frequency of use.

For maximum battery life, simply switch the wireless remote off when it's not in use.

Simple to install

No major bathroom modifications required. Simply remove the existing toilet seat and secure the lift in place using the included hardware.

Supporting independence at home

Whether operated independently or with assistance from a caregiver, the Electric Toilet Lift Seat can help make an everyday task safer, easier and more dignified.

Recommended weight capacity: 136 kg.

The unit can support up to 150 kg; however, users over 136 kg should have their seated body width checked against the distance between the armrests to ensure a comfortable fit.

+ SPECIFICATIONS

- 14" powered lift to assist with sitting and standing
- Easy-to-use controls on both armrests
- Wireless remote allows operation by the user or a caregiver
- Swing-up armrests for easier access and transfers
- Stable four-foot base for added support
- Fits most toilets with bowl heights from 14" to 18.5"
- Under 25" wide to suit standard NZ bathrooms
- 136 kg recommended user weight
- Rechargeable battery with more than 120 lift cycles per full charge
- Plug-in operation also available



Scan to view product details

Speak with our team about whether this solution is suitable for your home or client.

0800 31 61 81 | helpis@alliedmedical.co.nz | alliedmedical.co.nz

