

Ki Mobility Focus Express Script



Customer Details

Client				
Assessor				
Organisation				
User Weight	ACC		DSS	Other

Ki Mobility Focus Express Script

If you only require specifications that fall under the Standard Express Package, you do not need to fill in the reverse page.

Seat Depth	Seat Width			
	14"	16"	18"	20"
	16"	16"	16"	16"
	18"	18"	18"	18"
	20"	20"	20"	20"

Code for all Express Packages:

KIFOCUSCR14181620STDNAV

OR

KIFOCUSCR16201822STDNAV

Please add Seat Pan/Footplate code found in table to the left, based on your size selection, as well as the Standard Express Package code, above.

Express Package: (Select optional specifications and alternatives on reverse)

- Growing adjustable Aluminium seat pan
- Transit tie downs
- Height adjustable straight back cane with adjustable push handle
- Fold down back (angle adjustable)
- 80° legrest hanger
- Aluminium angle adjustable footplates
- Long extension tube hangers
- 4-Way in line hanger release
- 6" castors - 2" Poly
- 16" Pneumatic with Airless Insert
- Combination Attendant Foot Lock with Push to Lock with Extension Handles
- Quick release axles
- Height adjustable T-arm armrests with full arm pads
- Rear anti-tips
- Iridescent Navy Frame

Please note: Express Package includes some upcharge products.

Choose Your Specifications (No additional codes are required on the below parts)

Backrest Height			
20"	21"	22"	23"
24"	25"	26"	27"
Front Seat Height			
19"	20"	21"	

Seat to Footrest Length			
14"	14.5"	15"	15.5"
16"	16.5"	17"	17.5"
18"			

Other specifications are available on reverse as a non-standard express package

Optional Specifications

Leg Hanger / Font Frame Bend 1

70° (.ULEXPLEG70NAVY)	Omit Hanger
-----------------------	-------------

Medium (please specify):

11.5"	12"	12.5"	13"
-------	-----	-------	-----

13.5"

70° (.ULEXPSTF70M) | 80° (.ULEXPSTF80M)

Armpad 2

Desk arm pad (.ULEXPAPRD)

Adjustable Heel Loops 3

14-15W (.ULLIBHL1415)
16-22W (.ULLIBHL1622)

Tall Armrests 4

Short side shield and standard pad Desk, 14-17 deep (.ULEXP TARSD)
Short side shield and standard pad Full, 14-17 deep (.ULEXP TARSF)
Long side shield and standard pad Desk, 18-20 deep (.ULEXP TARLD)
Long side shield and standard pad Full, 18-20 deep (.ULEXP TARLF)

Rear Wheel Kit - Self propel 5

24" Maxx Spoke Street pneumatic tyres with Aluminium anodised Hand rims
Plastic coated push rim (.ULEXP HRMPC)

*Selecting this option will change Wheel Locks to Push to Lock

Natural Fit:
Thumb grip regular (.ULEXP HRMTGREG)
Ergonomic Fit:
Ergo Lite Regular (.ULEXP HRMELREG)

PRO Elevated Legrest - Angle Adj. FP Long 6

16/17 LHS (.ULEXP ELRLHS16)	16/17 RHS (.ULEXP ELRRHS16)
18-25 LHS (.ULEXP ELRLHS18)	18-25 RHS (.ULEXP ELRRHS18)

Residual Limb Support 7

5" x 8" with smooth cover	7" x 10" with smooth cover
Left (.ULEXP RLS5X8L)	Left (.ULEXP RLS7X10L)
Right (.ULEXP RLS5X8R)	Right (.ULEXP RLS7X10R)

Accessories 8

Black neoprene impact guards (.ULEXP IMGSA)
Clear spoke guards - 24" wheel only (.ULEXP SPGCL24)
Tray

Notes