



Gel Ovations

Solutions

Gel Ovations™



Armrest Pad

GO310G (3.5" x 10")

GO312G (3.5" x 12")

GO314G (3.5" x 14")

GO316G (3.5" x 16")

GO210G (2" x 10")

GO212G (2" x 12")

GO214G (2" x 14")

- These are available in five different sizes to suit all different needs and will easily bolt onto existing armrest T sections
- They consist of a solid wood base gel pad with a neoprene nylon cover and plastic base guard
- Designed for the ultimate in comfort and protection and can aid in healing existing sores and/or abrasions
- Sold individually



Ottobock Trough Liner

GOOTLM (Medium)

GOOTLL (Large)

- These are pads designed to fit inside Ottobock Channel Forearm Pads
- They simply Velcro in place and are long enough to pad up behind the elbow, making it considerably softer for the user
- Medium fits OB436A2=2-7 (13")
- Large fits OB436A2=1-7 (19")



Gel Ovations Waterfall Armpad

GO414GWT (4" x 14")

- GEL waterfall style pads with Nylon / Lycra cover
- Designed to be easily fitted to both manual and power wheelchairs
- Sold as a pair



Sta-Rite Forehead Stabiliser

GOSRAFS

- Headrest and Forehead stabiliser are padded with dimensional gel
- Fully adjustable stabiliser with quick release swing away function
- Mounts to standard ball joint mount



Split Calf Support Panel With Strap

GOGHS2

- Split calf panel allows use of swingaway leg without removal of pad
- Fits 16" to 20" wide seats
- Positioning strap included
- Fits Right or Left leg
- Sold Individually



Leg Protector Wraps

GOLPW-LOGO

- With a velcro closure these are designed to fit around a standard 70° legrest hanger
- Sold as a pair



Footrest Pads

GOGFPC

- 6" x 8" Bolt on gel footplate
- Designed to easily fit both manual and powered wheelchairs
- Wood base gel pad with neoprene cover
- Sold individually



Gel & Foam Footplate Protector

GOGFPCA1

- 6" x 8" Lightweight gel and foam footplate covers
- Designed to stick on with an adhesive back onto standard adult sized footplates
- Sold Individually



Footrest Pads - 1-Piece Footplate

GOGFPCS

- 11" x 9" gel footboard
- Designed to bolt onto centre-mounted footboard, common on powerchairs including the Quantum range - regular size
- It can also be used on some scooters
- Wood base gel pad with neoprene cover
- Sold individually




Gel Knee Protector Buttons

GOKBM (Regular)

GOKBMT (Thin)

- Designed to fit swingaway legrests, these pads come
- complete with mounting hardware to interface with
- popular brands
- Regular is 1 1/4" thick, 4" diameter
- Thin is 7/8" thick, 4" diameter



Sta-Rite Fix Bracket Adductor

GOSRABC (7/8")

GOSRABD (1")

- Quadra-pivot positioning allows pad position adjustment in 3 planes and 6 directions
- Comes standard with 3" x 4" GOAEP pad
- Serrated leg clamps and set screws at every pivot means the pad will stay right where you want it
- Size taken from 7/8" or 1" round tube leg rests
- Sold individually



Sta-Rite Swing-Away Bracket Adductor

GOSRASC (7/8")

GOSRASD (1")

GOSRASF (Flat)

- Quadra-pivot positioning allows pad position adjustment in 3 planes and 6 directions
- Comes standard with 3" x 4" GOAEP pad
- Serrated leg clamps and set screws at every pivot means the pad will stay right where you want it
- 7/8" or 1" fits round tube leg rests
- FM fits flat to square tube seating systems or surfaces
- Sold as a pair



Sta-Rite 6-Way Lateral Thigh Support

GOSRMHG10 (3.5" x 10")

GOSRMHG12 (3.5" x 12")

GOSRMHG14 (3.5" x 14")

- Multi-axis Hip Guide that is adjustable in 6 directions
- They mount easily to the side of virtually any seating system and can be easily adjusted up and down and can be locked in the position of choice
- Made of stainless steel on the primary support part with 2 interlocking crown pivots to allow adjustment in and out. The track on the pad allows continuous forward and back adjustment
- Sold individually



Oval Stick On Pad

GO34SP (3" x 4")

GO48SP (4" x 8")

GO412SP (4.5" x 12")

- Flat stick on gel pad with neoprene nylon cover
- Each is 0.375" in thickness
- Ideal to place on trays, calf pads and any surface that the user's body may come into contact with
- Sold individually



Egg Shaped Pad

GOESPLOGO (2.625" x 4.5")

- Domed stick on gel pad with neoprene nylon cover
- Each is 0.625" in thickness
- Sold individually



Lateral Gel Pad

GO46DGLP1

- 4" x 6" Wheelchair Lateral Pad made with silicone Dimensional gel
- Removable washable Neoprene / Nylon cover
- Standard 1" T-nut spacing fits most lateral support hardware
- Sold individually



Adductor / Elbow Gel Pad

GOAEP

- Flat 3" x 4" gel pad comes with a solid wooden bottom
- Bottom has a plastic cover and 2 T-nut connections with 1" spacing
- Sold individually



DeZyne Gel Sheet Pad

GODG14145 (14" x 14")

GODG16165 (16" x 16")

GODG16185 (16" x 18")

GODG18165 (18" x 16")

GODG18185 (18" x 18")

GODG18205 (18" x 20")

GODG20185 (20" x 18")

GODG20205 (20" x 20")

GODG24205 (24" x 20")

- The unique lightweight cellular design can redistribute pressure, reduce abrasion, absorb light shocks and dampen vibration while minimising heat and moisture build up
- Made from water insoluble proprietary silicone gel so that it is ideal for use in places where the pad might get wet, so it's easy to clean
- Adapts to all surfaces and is long lasting as it will not flow or leak if cut or punctured, and is also odourless, hypoallergenic and latex free
- Each is 0.5" in thickness



DeZyne Gel Pad

GODG65955 (6.5" x 9.5")

GODG895 (8" x 9")

- Pad sized version of the DeZyne Gel Sheet
- The unique lightweight cellular design can redistribute pressure, reduce abrasion, absorb light shocks and dampen vibration while minimising heat and moisture build up
- Made from water insoluble proprietary silicone gel so that it is ideal for use in places where the pad might get wet, so it's easy to clean
- Long lasting, will not flow or leak if cut or punctured, it is also odourless, hypoallergenic and latex free
- Adapts to all surfaces
- Each is 0.5" in thickness

ACCESSORIES



Screw Sets for Gel Ovation Armpads

GORTFSL

- Adapter set used to convert 3/4" round tube to flat surface using 2 bolts
- Most commonly used with flat bottom arm pads
- The 2 adapter pieces accommodate various bolt configurations
- Order 2 sets to adapt a pair of flat bottom arm rests

Sta-Rite™

DeZyne™ GEL



Please visit our website
www.alliedmedical.co.nz
for our full range of products



0800 31 61 81
helpis@alliedmedical.co.nz

July 2024. This information is correct at the time of printing.
From time to time, product specifications are changed by manufacturers without prior notice.



**ALLIED
MEDICAL**
LIMITED