

# **Seating System Armrests IPB, International**

This Illustrated Parts Breakdown (IPB) Catalog has been compiled from the latest product information available at the time of publication. Quantum Rehab® / Pride Mobility Products® Corporation reserves the right to make changes and updates to this manual and to the product(s) it represents. Any changes, updates, or improvements to the product(s) may result in slight differences between this manual and the product that has been purchased.

An IPB is a graphically portrayed parts listing of product components which have been deemed to have serviceable, maintenance significance. When applicable, IPBs utilize “exploded” illustrations that detail serviceable components and a corresponding itemized parts list of reference numbers, part numbers, descriptions, quantities, and any applicable pricing. Most product components are incorporated into assemblies and may not be available for replacement.

Serviceable components will be noted by a callout in the illustration that will correspond to a alphanumeric reference number in the parts list. Items may also be listed in the parts list by parent assemblies. In most cases when a parent assembly is listed, there will be child parts or sub-assemblies listed under that parent.

For example: The parent assembly reference number is A1 — A5. That part contains all of the items listed as and contained within A1, A2, A3, A4, & A5. As long as these subcomponents have a callout associated with them, they can be selected for replacement or quote.

Updates to the contents of an IPB will be made on a regular basis based upon applicable product changes and enhancements. This PDF catalog has been export from the Provider based service site. The most up to date parts look up tool and ordering application is located at: [prideprovider.com](http://prideprovider.com). Please utilize this system to enhance your user experience. It is built with highly interactive materials that are set with relational properties that will facilitate easier part look-up, more accurate parts selection, notification of onhand part availability, and quick links to related documents and products.

We welcome all questions and comments regarding the Quantum Rehab® / Pride Mobility Products® Corporation IPB systems.

Any pricing that is included in Quantum Rehab® /Pride Mobility Products® Corporation IPBs is subject to change without prior notice.

Thank you for choosing products from Quantum Rehab® and Pride Mobility Products® Corporation.




6.3. 10"-11" Height w/ Std. Armpad

10"-11"H, Full Length, Right, Std. - Recline Armrest Assy	1
10"-11"H, Full Length, Left, Std. - Recline Armrest Assy	3
10"-11"H, Desk Length, Right, Std. - Recline Armrest Assy	5
10"-11"H, Desk Length, Left, Std. - Recline Armrest Assy	7
Appendix A:Pages Index	9
Appendix B:Parts Index	10

6.3. 10"-11" Height w/ Std. Armpad


\*\*\* Individual components are not available for replacement. The entire assembly must be purchased. \*\*\*



Item	Part No.	Name	Qty
A1	SETLRAM1116	ARMREST,PAD,LONG,QUANTUM,V2	1
B1	FRMASMB7134	ARMREST,ASSY,RECLINE,10"H-11"H,FULL LENGTH,RIGHT	1
C1	DWR1367C043	END CAP,RECLINE ARMREST,TRU-BALANCE POWER RECLINE,C-1367-043	1
D1	SCRSHFH1158	SCREW, SHFH,BLACK ZINC,10-32 X 0.875",GR 8,W/PATCH	4
E1	ARMASMB7110006	ARMREST,ASSY,RECLINE,FULL LENGTH,RIGHT,W/RECEIVER BLOCK	1
F1	FRMASMB7066	ARMREST,ASSY,RECEIVER BLOCK,TB2 RECLINE,C-1302-027	1
G1	DWR1367C042	CAP,PUSH BUTTON,RECLINE ARMREST QUICK RELEASE	1
H1	DWR1302B003	PUSH PIN,3/8-24 X 1.88"L,QUICK RELEASE,TB2 RECLINE ARMREST	1
I1	NUTNJAM1034	NUT,NYLOCK JAM,ZINC,3/8-24	1
J1	FRMASMB7142	ARMREST,ASSY,PIVOT,10"H-11"H,TRU-BALANCE POWER RECLINE,W/MOUNTING HARDWARE,RIGHT	1

6.3. 10"-11" Height w/ Std. Armpad



\*\*\* Individual components are not available for replacement. The entire assembly must be purchased. \*\*\*



Item	Part No.	Name	Qty
A1	SETLRAM1116	ARMREST,PAD,LONG,QUANTUM,V2	1
B1	FRMASMB7133	ARMREST,ASSY,RECLINE,10"H-11"H,FULL LENGTH,LEFT	1
C1	DWR1367C043	END CAP,RECLINE ARMREST,TRU-BALANCE POWER RECLINE,C-1367-043	1
D1	SCRSHFH1158	SCREW, SHFH,BLACK ZINC,10-32 X 0.875",GR 8,W/PATCH	4
E1	ARMASMB7110005	ARMREST,ASSY,RECLINE,FULL LENGTH,LEFT,W/RECEIVER BLOCK	1
F1	FRMASMB7066	ARMREST,ASSY,RECEIVER BLOCK,TB2 RECLINE,C-1302-027	1
G1	DWR1367C042	CAP,PUSH BUTTON,RECLINE ARMREST QUICK RELEASE	1
H1	DWR1302B003	PUSH PIN,3/8-24 X 1.88"L,QUICK RELEASE,TB2 RECLINE ARMREST	1
I1	NUTNJAM1034	NUT,NYLOCK JAM,ZINC,3/8-24	1
J1	FRMASMB7141	ARMREST,ASSY,PIVOT,10"H-11"H,TRU-BALANCE POWER RECLINE,W/MOUNTING HARDWARE,LEFT	1

6.3. 10"-11" Height w/ Std. Armpad

\*\*\* Individual components are not available for replacement. The entire assembly must be purchased. \*\*\*






Item	Part No.	Name	Qty
A1	SETLRAM1115	ARMREST,PAD,SHORT,QUANTUM,V2	1
B1	FRMASMB7136	ARMREST,ASSY,RECLINE,10"H-11"H,DESK LENGTH,RIGHT	1
C1	DWR1367C043	END CAP,RECLINE ARMREST,TRU-BALANCE POWER RECLINE,C-1367-043	1
D1	SCRSHFH1158	SCREW, SHFH,BLACK ZINC,10-32 X 0.875",GR 8,W/PATCH	4
E1	ARMASMB7110004	ARMREST,ASSY,RECLINE,DESK LENGTH,RIGHT,W/RECEIVER BLOCK	1
F1	FRMASMB7066	ARMREST,ASSY,RECEIVER BLOCK,TB2 RECLINE,C-1302-027	1
G1	DWR1367C042	CAP,PUSH BUTTON,RECLINE ARMREST QUICK RELEASE	1
H1	DWR1302B003	PUSH PIN,3/8-24 X 1.88"L,QUICK RELEASE,TB2 RECLINE ARMREST	1
I1	NUTNJAM1034	NUT,NYLOCK JAM,ZINC,3/8-24	1
J1	FRMASMB7142	ARMREST,ASSY,PIVOT,10"H-11"H,TRU-BALANCE POWER RECLINE,W/MOUNTING HARDWARE,RIGHT	1

6.3. 10"-11" Height w/ Std. Armpad

\*\*\* Individual components are not available for replacement. The entire assembly must be purchased. \*\*\*



Item	Part No.	Name	Qty
A1	SETLRAM1115	ARMREST,PAD,SHORT,QUANTUM,V2	1
B1	FRMASMB7135	ARMREST,ASSY,RECLINE,10"H-11"H,DESK LENGTH,LEFT	1
C1	DWR1367C043	END CAP,RECLINE ARMREST,TRU-BALANCE POWER RECLINE,C-1367-043	1
D1	SCRSHFH1158	SCREW, SHFH,BLACK ZINC,10-32 X 0.875",GR 8,W/PATCH	4
E1	ARMASMB7110003	ARMREST,ASSY,RECLINE,DESK LENGTH,LEFT,W/RECEIVER BLOCK	1
F1	FRMASMB7066	ARMREST,ASSY,RECEIVER BLOCK,TB2 RECLINE,C-1302-027	1
G1	DWR1367C042	CAP,PUSH BUTTON,RECLINE ARMREST QUICK RELEASE	1
H1	DWR1302B003	PUSH PIN,3/8-24 X 1.88"L,QUICK RELEASE,TB2 RECLINE ARMREST	1
I1	NUTNJAM1034	NUT,NYLOCK JAM,ZINC,3/8-24	1
J1	FRMASMB7141	ARMREST,ASSY,PIVOT,10"H-11"H,TRU-BALANCE POWER RECLINE,W/MOUNTING HARDWARE,LEFT	1

Appendix A:Pages Index

10"-11"H, Desk Length, Left, Std. - Recline Armrest Assy . . . . . 7

10"-11"H, Desk Length, Right, Std. - Recline Armrest Assy . . . . . 5

10"-11"H, Full Length, Left, Std. - Recline Armrest Assy . . . . . 3

10"-11"H, Full Length, Right, Std. - Recline Armrest Assy . . . . . 1

## Appendix B:Parts Index

ARMASMB7110003	8
ARMASMB7110004	6
ARMASMB7110005	4
ARMASMB7110006	2
DWR1302B003	2, 4, 6, 8
DWR1367C042	2, 4, 6, 8
DWR1367C043	2, 4, 6, 8
FRMASMB7066	2, 4, 6, 8
FRMASMB7133	4
FRMASMB7134	2
FRMASMB7135	8
FRMASMB7136	6
FRMASMB7141	4, 8
FRMASMB7142	2, 6
NUTNJAM1034	2, 4, 6, 8
SCRSHFH1158	2, 4, 6, 8
SETLRAM1115	6, 8
SETLRAM1116	2, 4

# Illustrated Parts Breakdown

## **USA - Quantum Rehab**

401 York Avenue  
Duryea, Pa 18642  
[www.quantumrehab.com](http://www.quantumrehab.com)

## **USA - Pride Mobility Products Corp.**

182 Susquehanna Avenue  
Exeter, Pa 18643  
[www.pridemobility.com](http://www.pridemobility.com)

## **Canada**

5096 South Service Road  
Beamsville, Ontario L0R 1B3  
[www.pridemobility.ca](http://www.pridemobility.ca)  
[www.quantumrehab.com](http://www.quantumrehab.com)

## **Australia**

20-24 Apollo Drive  
Hallam, Victoria 3803  
[www.pridemobility.com.au](http://www.pridemobility.com.au)  
[www.quantumrehab.com.au](http://www.quantumrehab.com.au)

## **New Zealand**

38 Lansford Crescent  
Avondale, Auckland 0600  
[www.pridemobility.co.nz](http://www.pridemobility.co.nz)

## **UK**

32 Wedgwood Road  
Bicester, Oxfordshire OX26 4UL  
[www.pride-mobility.co.uk](http://www.pride-mobility.co.uk)  
[www.quantumrehab.co.uk](http://www.quantumrehab.co.uk)

## **Netherlands**

### **(Authorised EU Representative)**

De Zwaan 3  
1601 MS Enkhuizen  
[www.pridemobility.eu](http://www.pridemobility.eu)  
[www.quantumrehab.eu](http://www.quantumrehab.eu)

## **Italy**

Via del Progresso, ang. Via del Lavoro  
Loc. Prato della Corte  
00065 Fiano Romano (RM)  
[www.pride-italia.it](http://www.pride-italia.it)  
[www.quantumrehab-italia.it](http://www.quantumrehab-italia.it)

## **France**

26 Rue Monseigneur Ancel  
69800 Saint-Priest  
[www.pridemobility.fr](http://www.pridemobility.fr)  
[www.quantumrehab.fr](http://www.quantumrehab.fr)

## **Spain**

Calle Las Minas Número 67  
Polígono Industrial Urtinsa II, de Alcorcón  
28923 Madrid

## **Germany**

Hövelrieger Str. 28  
33161 Hövelhof  
[www.pridemobility.de](http://www.pridemobility.de)

## **China**

Room 508, Building #4  
TianNa Business Zone  
No. 500 Jianyun Road  
Pudong New Area  
Shanghai 201318

