

Seating System Armrests IPB, International

Illustrated Parts Breakdown

USA - Quantum Rehab

401 York Avenue
Duryea, Pa 18642
www.quantumrehab.com

USA - Pride Mobility Products Corp.

182 Susquehanna Avenue
Exeter, Pa 18643
www.pridemobility.com

Canada

5096 South Service Road
Beamsville, Ontario L0R 1B3
www.pridemobility.ca
www.quantumrehab.com

Australia

20-24 Apollo Drive
Hallam, Victoria 3803
www.pridemobility.com.au
www.quantumrehab.com.au

New Zealand

38 Lansford Crescent
Avondale, Auckland 0600
www.pridemobility.co.nz

UK

32 Wedgwood Road
Bicester, Oxfordshire OX26 4UL
www.pride-mobility.co.uk
www.quantumrehab.co.uk

Netherlands

(Authorised EU Representative)

De Zwaan 3
1601 MS Enkhuizen
www.pridemobility.eu
www.quantumrehab.eu

Italy

Via del Progresso, ang. Via del Lavoro
Loc. Prato della Corte
00065 Fiano Romano (RM)
www.pride-italia.it
www.quantumrehab-italia.it

France

26 Rue Monseigneur Ancel
69800 Saint-Priest
www.pridemobility.fr
www.quantumrehab.fr

Spain

Calle Las Minas Número 67
Polígono Industrial Urtinsa II, de Alcorcón
28923 Madrid

Germany

Hövelrieger Str. 28
33161 Hövelhof

China

Room 508, Building #4
TianNa Business Zone
No. 500 Jianyun Road
Pudong New Area
Shanghai 201318



Standard Armpad Configuration

Desk Lenght, Left, 2 Post Flip Up Armrest 1

Desk Length, Right, 2 Post Flip Up Armrest 3

Full Length, Left, 2 Post Flip Up Armrest 5


Full Length, Right, 2 Post Flip Up Armrest 7

Appendix A:Pages Index 9

Appendix B:Parts Index 10

Standard Armpad Configuration

*** Individual components are not available for replacement. The entire assembly must be purchased. ***






Item	Part No.	Name	Qty
A1	SETLRAM1115	ARMREST,PAD,SHORT,QUANTUM,V2	1
B1	SETLRAM1164	ARMPAD,ASSY,WATERFALL,BLACK,9", DESK LENGTH	1
C1	FRMASMB7060	ARMREST,TWO POST FLIP-UP,DESK LENGTH,LEFT	1
D1	SCRASBH1053	SCREW,BHSCS,BLACK ZINC,10-32 X 1.25"	2
E1	FRMASMB6792	ARMREST,TWO POST FLIP-UP,LOWER,LEFT	1
F1	PROHDLC1062	HAND GRIP,BLACK,FOAM,.75"ID X 10.00"L X .25" THICKNESS	1
G1		END CAP,0.875" DIA,FITS 18 GAUGE TUBE,DOMED	1
H1	FRMASMB6374	ARMREST,LATCH BLOCK,TWO POST FLIP-UP,FRONT,LEFT	1
I1	FRMASMB6967	ARMREST,TWO POST FLIP-UP,DESK LENGTH,UPPER	1

Standard Armpad Configuration

*** Individual components are not available for replacement. The entire assembly must be purchased. ***






Item	Part No.	Name	Qty
A1	SETLRAM1115	ARMREST,PAD,SHORT,QUANTUM,V2	1
B1	SETLRAM1164	ARMPAD,ASSY,WATERFALL,BLACK,9", DESK LENGTH	1
C1	FRMASMB7061	ARMREST,TWO POST FLIP-UP,DESK LENGTH,RIGHT	1
D1	SCRASBH1053	SCREW,BHSCS,BLACK ZINC,10-32 X 1.25"	2
E1	FRMASMB6793	ARMREST,TWO POST FLIP-UP,LOWER,RIGHT	1
F1	PROHDLC1062	HAND GRIP,BLACK,FOAM,.75"ID X 10.00"L X .25" THICKNESS	1
G1		END CAP,0.875" DIA,FITS 18 GAUGE TUBE,DOMED	1
H1	FRMASMB6375	ARMREST,LATCH BLOCK,TWO POST FLIP-UP,FRONT,RIGHT	1
I1	FRMASMB6967	ARMREST,TWO POST FLIP-UP,DESK LENGTH,UPPER	1

Standard Armpad Configuration

*** Individual components are not available for replacement. The entire assembly must be purchased. ***




Full Length, Left, 2 Post Flip Up Armrest



Item	Part No.	Name	Qty
A1	SETLRAM1116	ARMREST,PAD,LONG,QUANTUM,V2	1
B1	SETLRAM1165	ARMPAD,ASSY,WATERFALL,BLACK,14", FULL LENGTH	1
C1	FRMASMB7062	ARMREST,TWO POST FLIP-UP,FULL LENGTH,LEFT	1
D1	SCRASBH1053	SCREW,BHSCS,BLACK ZINC,10-32 X 1.25"	2
E1	FRMASMB6792	ARMREST,TWO POST FLIP-UP,LOWER,LEFT	1
F1	PROHDLC1062	HAND GRIP,BLACK,FOAM,.75"ID X 10.00"L X .25" THICKNESS	1
G1		END CAP,0.875" DIA,FITS 18 GAUGE TUBE,DOMED	1
H1	FRMASMB6374	ARMREST,LATCH BLOCK,TWO POST FLIP-UP,FRONT,LEFT	1
I1	FRMASMB6968	ARMREST,TWO POST FLIP-UP,FULL LENGTH,UPPER	1

Standard Armpad Configuration

*** Individual components are not available for replacement. The entire assembly must be purchased. ***



Full Length, Right, 2 Post Flip Up Armrest



Item	Part No.	Name	Qty
A1	SETLRAM1116	ARMREST,PAD,LONG,QUANTUM,V2	1
B1	SETLRAM1165	ARMPAD,ASSY,WATERFALL,BLACK,14", FULL LENGTH	1
C1	FRMASMB7063	ARMREST,TWO POST FLIP-UP,FULL LENGTH,RIGHT	1
D1	SCRASBH1053	SCREW,BHSCS,BLACK ZINC,10-32 X 1.25"	2
E1	FRMASMB6793	ARMREST,TWO POST FLIP-UP,LOWER,RIGHT	1
F1	PROHDL1062	HAND GRIP,BLACK,FOAM,.75"ID X 10.00"L X .25" THICKNESS	1
G1		END CAP,0.875" DIA,FITS 18 GAUGE TUBE,DOMED	1
H1	FRMASMB6375	ARMREST,LATCH BLOCK,TWO POST FLIP-UP,FRONT,RIGHT	1
I1	FRMASMB6968	ARMREST,TWO POST FLIP-UP,FULL LENGTH,UPPER	1

Appendix A:Pages Index

Desk Lenght, Left, 2 Post Flip Up Armrest 1

Desk Length, Right, 2 Post Flip Up Armrest 3

Full Length, Left, 2 Post Flip Up Armrest 5

Full Length, Right, 2 Post Flip Up Armrest 7

Appendix B:Parts Index

FRMASMB6374	2, 6
FRMASMB6375	4, 8
FRMASMB6792	2, 6
FRMASMB6793	4, 8
FRMASMB6967	2, 4
FRMASMB6968	6, 8
FRMASMB7060	2
FRMASMB7061	4
FRMASMB7062	6
FRMASMB7063	8
PROECAP1054	2, 4, 6, 8
PROHDL1062	2, 4, 6, 8
SCRASBH1053	2, 4, 6, 8
SETLRAM1115	2, 4
SETLRAM1116	6, 8
SETLRAM1164	2, 4
SETLRAM1165	6, 8

Illustrated Parts Breakdown

USA - Quantum Rehab

401 York Avenue
Duryea, Pa 18642
www.quantumrehab.com

USA - Pride Mobility Products Corp.

182 Susquehanna Avenue
Exeter, Pa 18643
www.pridemobility.com

Canada

5096 South Service Road
Beamsville, Ontario L0R 1B3
www.pridemobility.ca
www.quantumrehab.com

Australia

20-24 Apollo Drive
Hallam, Victoria 3803
www.pridemobility.com.au
www.quantumrehab.com.au

New Zealand

38 Lansford Crescent
Avondale, Auckland 0600
www.pridemobility.co.nz

UK

32 Wedgwood Road
Bicester, Oxfordshire OX26 4UL
www.pride-mobility.co.uk
www.quantumrehab.co.uk

Netherlands

(Authorised EU Representative)

De Zwaan 3
1601 MS Enkhuizen
www.pridemobility.eu
www.quantumrehab.eu

Italy

Via del Progresso, ang. Via del Lavoro
Loc. Prato della Corte
00065 Fiano Romano (RM)
www.pride-italia.it
www.quantumrehab-italia.it

France

26 Rue Monseigneur Ancel
69800 Saint-Priest
www.pridemobility.fr
www.quantumrehab.fr

Spain

Calle Las Minas Número 67
Polígono Industrial Urtinsa II, de Alcorcón
28923 Madrid

Germany

Hövelrieger Str. 28
33161 Hövelhof
www.pridemobility.de

China

Room 508, Building #4
TianNa Business Zone
No. 500 Jianyun Road
Pudong New Area
Shanghai 201318

